

June 7, 2023

CS Medical LLC Kendall Ashe Vice President, General Manager 2179 East Lyon Station Rd Creedmoor, North Carolina 27522

Re: K230381

Trade/Device Name: Ethos Automated Ultrasound Probe Cleaner Disinfector with AquaCide

Cleaner/High-Level Disinfectant and QwikCheck Chemical Indicator

Regulation Number: 21 CFR 892.1570

Regulation Name: Diagnostic Ultrasonic Transducer

Regulatory Class: Class II Product Code: PSW

Dated: May 12, 2023 Received: May 12, 2023

#### Dear Kendall Ashe:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <a href="https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm">https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm</a> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

K230381 - Kendall Ashe Page 2

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <a href="https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products">https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products</a>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <a href="https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems">https://www.fda.gov/medical-device-problems</a>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<a href="https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</a>) and CDRH Learn (<a href="https://www.fda.gov/training-and-continuing-education/cdrh-learn">https://www.fda.gov/training-and-continuing-education/cdrh-learn</a>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<a href="https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</a>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

# Bifeng Qian -S

Bifeng Qian, M.D., Ph.D.
Assistant Director
DHT4B: Division of Infection Control
and Plastic Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

# DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

## Indications for Use

Form Approved: OMB No. 0910-0120

Expiration Date: 06/30/2023 See PRA Statement below.

510(k) Number *(if known)* K230381

**Device Name** 

Ethos Automated Ultrasound Probe Cleaner Disinfector with AquaCide Cleaner/High-Level Disinfectant and QwikCheck Chemical Indicator

Indications for Use (Describe)

The Ethos automated ultrasound probe cleaner disinfector is designed to provide cleaning and high-level disinfection of surface and endocavity ultrasound probes. The system uses AquaCide cleaner/disinfectant, which is designed to be used only with the Ethos automated cleaner disinfector. The disinfectant bottles cannot be reused in the system.

AquaCide cleaner/disinfectant is intended for use as a single use cleaner and high-level disinfectant, used exclusively in the Ethos automated cleaner disinfector for cleaning and high-level disinfection of surface and endocavity ultrasound probes.

AquaCide cleaner/disinfectant should be used with the following contact conditions in the Ethos automated cleaner disinfector:

High-level Time Temperature Minimum Recommended Concentration

Disinfectant

AquaCide 3 minutes 47oC 1750 ppm peracetic acid

The QwikCheck Chemical Indicator is for use in Ethos to determine whether the concentration of peracetic acid, the active ingredient in AquaCide, is above or below the Minimum Recommended Concentration (MRC) of 1750 ppm.

The Ethos cleaner disinfector with AquaCide and QwikCheck are intended for use by qualified individuals trained in its use.

Type of Use (Select one or both, as applicable)	
Prescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)

#### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

#### \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

# 510(K) SUMMARY - K230381

510(k) Owner

CS Medical L.L.C. 2179 E. Lyon Station Road Creedmoor, North Carolina 27522 Phone 919-255-9472

Fax 919-528-3400

Contact Name

Kendall Ashe

Vice President and General Manager, CS Medical, L.L.C.

2179 E. Lvon Station Road

Creedmoor, North Carolina 27522

Email: kendallashe@csmedicalllc.com

**Submission Prepared** 

6-June-2023

**CS Medical Trade Name** 

Ethos Automated Ultrasound Probe Cleaner Disinfector with AquaCide Cleaner/High-

Level Disinfectant and

QwikCheck Chemical Indicator

Common Name

High level disinfection reprocessing instrument for ultrasonic transducers.

liquid.

Classification Name

Diagnostic Ultrasound Transducer (21 CFR 892.1570, Product Code

PSW)

## Legally Marketed Predicate Devices

TEEClean Automated TEE Probe Cleaner Disinfector with TEEZyme and TD-5 or TD-8 (K182891) – Primary Predicate Device

#### Reference Device

TD 200 Automated TEE Probe Disinfector with TD-12 and QwikCheck (K192228)

Reference Device

# Description of the CS Medical Ethos Automated Ultrasound Probe Cleaner Disinfector with AquaCide Cleaner/High-Level Disinfectant and QwikCheck **Chemical Indicator**

The Ethos cleaner disinfector provides cleaning and high-level disinfection of surface and endocavity ultrasound probes when used according to the operating instructions and when used with AquaCide cleaner/disinfectant. The Ethos cleaner disinfector is for use only with AquaCide cleaner/disinfectant. The AquaCide cleaner/disinfectant is for use only in the Ethos cleaner disinfector. Thus, the Ethos cleaner disinfector, AquaCide cleaner/disinfectant represent a dedicated system. Each soiled ultrasound probe has the condom/cover removed and is bedside cleaned according to the ultrasound probe manufacturer's instructions before insertion into the Ethos cleaner disinfector. If a

condom/cover is not used the user must manually clean the probe. A fresh, unopened bottle of granular PAA AquaCide cleaner/disinfectant is loaded into the Ethos. The Ethos brings in water, mixes and heats the AquaCide solution to a minimum of 47°. While that is occurring, the Ethos brings in water and the ultrasound probe is pre-rinsed. After temperature is achieved and MRC is confirmed, the Ethos bring the AquaCide solution to the probe and cleans and disinfects the ultrasound probe for at least three minutes. Then the Ethos thoroughly rinses the AquaCide off the ultrasound probe before the cycle is complete. The ultrasound probe is then removed from the Ethos and dried according to the ultrasound probe manufacturer's instructions. The Ethos is ready for a new cycle immediately after the preceding cycle is completed. The Ethos cleaner disinfector incorporates a method for validating the PAA solution through an automatic chemical indicator to ensure each dose of PAA is at or above the MRC. The Ethos cleaner disinfector prints a verification report indicating a successful cleaning and disinfection cycles as well as the time and the average temperature during the cleaning/disinfection. The ultrasound probe is then removed from the Ethos cleaner disinfector and dried according to the ultrasound probe manufacturer's instructions. The Ethos cleaner disinfector is ready for a new cycle immediately after the preceding cycle is completed.

#### **Indications for Use Statement:**

The Ethos automated ultrasound probe cleaner disinfector is designed to provide cleaning and high-level disinfection of surface and endocavity ultrasound probes. The system uses AquaCide cleaner/disinfectant, which is designed to be used only with the Ethos automated cleaner disinfector. The disinfectant bottles cannot be reused in the system.

AquaCide cleaner/disinfectant is intended for use as a single use cleaner and high-level disinfectant, used exclusively in the Ethos automated cleaner disinfector for cleaning and high-level disinfection of surface and endocavity ultrasound probes.

AquaCide cleaner/disinfectant should be used with the following contact conditions in the Ethos automated cleaner disinfector:

High-level	Time	Temperature	Minimum Recommended
Disinfectant			Concentration
AquaCide	3 minutes	47°C	1750 ppm peracetic acid

The QwikCheck Chemical Indicator is for use in Ethos to determine whether the concentration of peracetic acid, the active ingredient in AquaCide, is above or below the Minimum Recommended Concentration (MRC) of 1750 ppm.

The Ethos cleaner disinfector with AquaCide and QwikCheck are intended for use by qualified individuals trained in its use.

Comparison of Proposed Device to Reference Device and Primary Predicate Device

Comp	Proposed Device	Primary Predicate	Reference Device	Comparison
Element	CS Medical Ethos Automated Ultrasound Probe Cleaner Disinfector with AquaCide Cleaner/High-Level Disinfectant and QwikCheck Chemical Indicator (K230381)	TEEClean® Automated TEE Probe Cleaner Disinfector with TEEZyme® Cleaner and TD-5® or TD-8® High-Level Disinfectants (K182891)	TD 200 <sup>®</sup> Transesophageal Probe Disinfector and TD-12 High-level Disinfectant (K192228)	Companson
Classificatio n Name (CFR; Product code)	Diagnostic ultrasonic transducer (21 C.F.R. § 892.1570, Product Code PSW)	Diagnostic ultrasonic transducer (21 C.F.R. § 892.1570, Product Code PSW)	Diagnostic ultrasonic transducer (21 C.F.R. § 892.1570, Product Code PSW)	Proposed and primary predicate same classification
Indications for Use	The Ethos automated ultrasound probe cleaner disinfector is designed to provide cleaning and high-level disinfection of surface and endocavity ultrasound probes. The system uses AquaCide cleaner/disinfectant, which is designed to be used only with the Ethos automated cleaner disinfector. The disinfectant bottles cannot be reused in the system.  AquaCide cleaner/disinfectant is intended for use as a single use cleaner and high-level disinfectant, used exclusively in the Ethos automated cleaner disinfector for cleaning and high-level disinfection of surface and endocavity ultrasound	The TEEClean automated cleaner disinfector is intended to replace manual cleaning of Transesophageal (TEE) ultrasound probes and automate high-level disinfection of TEE probes. The system uses TEEZyme enzymatic cleaner to clean TEE probes as well as TD-5 or TD-8 disinfectant to high level disinfect TEE probes. TEE probes must undergo bedside cleaning prior to insertion into the TEEClean.  The TD-5 or TD-8 disinfectant bottles cannot be reused in the system.  TD-5 disinfectant is intended for use as a single use high-level disinfectant used	The TD 200 Automated TEE Probe Disinfector with TD-12 High Level Disinfectant is designed to provide high-level disinfection of Transesophageal (TEE) ultrasound probes. The system can use TD-12 disinfectant, which is designed to be used only with the TD 200 disinfector. The disinfectant bottles cannot be reused in the system.  TD-12 disinfectant is intended for use as single use high-level disinfectant to be used exclusively in the TD 200 disinfector for high-level disinfection of TEE ultrasound probes.  TD-12 high level disinfectant and TD	Proposed device and primary predicate similar indications for use, reference device provides same disinfectant

probes.

AquaCide cleaner/disinfectant should be used with the following contact conditions in the Ethos automated cleaner disinfector:

High-level Disinfectant AquaCide

Time 3 minutes

Temperature 47°C

Minimum Recommended Concentration

1750 ppm peracetic acid

The QwikCheck
Chemical Indicator is
for use in Ethos to
determine whether the
concentration of
peracetic acid, the
active ingredient in
AquaCide, is above or
below the Minimum
Recommended
Concentration (MRC)
of 1750 ppm.

The Ethos cleaner disinfector with AquaCide and QwikCheck are intended for use by qualified individuals trained in its use.

exclusively in the TEEClean automated cleaner disinfector for high-level disinfection of TEE ultrasound probes. TD-5 disinfectant should be used with the following contact conditions in the TEEClean automated cleaner disinfector: High-level disinfectant TD-5 Time - 5 minutes Temperature - 38° -40°C

40°C
Minimum
Recommended
Concentration
1.7% glutaraldehyde

TD-8 disinfectant is intended for use as a single use high-level disinfectant used exclusively in the TEEClean automated cleaner disinfector for high-level disinfection of TEE ultrasound probes. TD-8 disinfectant should be used with the following contact conditions in TEEClean automated cleaner disinfector: High-level disinfectant TD-8 Time - 5 minutes Temperature - 38° -40°C Minimum Recommended Concentration 0.3% orthophthalaldehyde TEEZyme enzymatic cleaner, TD-5 and TD-8 high level

200 disinfector is intended for use by qualified individuals trained in its use.

TD-12 disinfectant should be used with the following contact conditions in TD 200 disinfector:

High-level Disinfectant TD-12

Time 3 minutes

Temperature 38°C

Minimum Recommended Concentration

1750 ppm peracetic acid

disinfectant, and

		TEEClean automated cleaner disinfector system are intended for use by qualified individuals trained in its use.		
Instrument ation For Automation	Ethos® with AquaCide® disinfector is automated for single use with only AquaCide cleaner disinfectant. The user initiates the automated cleaning and disinfection cycles via touchpad and receives cleaning and disinfection verification ticket and data is stored electronically.	TEEClean® with TEEZyme® enzymatic cleaner and TD-5® or TD-8® disinfectors is automated for single use with only TD-5® or TD-8® disinfectants. The user initiates the automated cleaning and disinfection cycles via touchpad and receives cleaning and disinfection verification ticket and data is stored electronically.	TD 200® with TD-12® disinfector is automated for single use with only TD-12® disinfectant. The user initiates the automated cycle via touchpad and receives disinfection verification ticket.	Proposed device and primary predicate device similar, reference device provides same disinfectant

Comparison of Operational Principles

Ethos (K230381)	TEEClean (K182891)	TD 200 (K192228)
The Ethos cleaner disinfector	The TEEClean cleaner disinfector	The TD 200 disinfector
provides high-level disinfection of	provides high-level disinfection of	provides high-level
surface and endocavity	transesophageal (TEE) ultrasound	disinfection of
ultrasound probes when used	probes when used according to	transesophageal (TEE)
according to the operating	the operating instructions, and	ultrasound probes when
instructions, and when used with	when used with TEEZyme	used according to the
AquaCide cleaner disinfectant.	enzymatic cleaner and either TD-5	operating instructions, and
Each soiled ultrasound probe has	or TD-8 disinfectant. Each soiled	when used with TD-12
the condom/cover removed and	ultrasound probe is pre-cleaned	disinfectant. Each soiled
pre-cleaned before insertion into	manually before insertion into the	ultrasound probe is pre-
the Ethos. A fresh, unopened	TEEClean cleaner disinfector. A	cleaned and manually
bottle of AquaCide cleaner	fresh, unopened bottle of TD-5 or	cleaned before insertion into
disinfectant is loaded into the	TD-8 disinfectant is loaded into	the TD 200 disinfector. A
Ethos. The Ethos heats the	the TEEClean cleaner disinfector.	fresh, unopened bottle of
AquaCide cleaner disinfectant to	The TEEClean cleaner disinfector	TD-12 disinfectant is loaded
the correct temperature and MRC	heats the TEEZyme enzymatic	into the TD 200 disinfector.
while spraying the ultrasound	cleaner to the correct temperature,	The TD 200 disinfector heats
probe with water. After MRC is	soaks the ultrasound probe, and	the TD-12 disinfectant to the
confirmed the AquaCide solution	then thoroughly rinses the	correct
is sprayed in the ultrasound	enzymatic cleaner off the	temperature, soaks the
probe at the correct temperature	ultrasound probe before the cycle	ultrasound probe, and then
and time for disinfection. Then	is complete. The TEEClean	thoroughly rinses the
the Ethos thoroughly rinses the	cleaner disinfector then heats the	disinfectant off the
AquaCide off the	TD5 or TD-8 disinfectant to the	ultrasound probe before the

ultrasound probe before the cycle is complete. The ultrasound probe is then removed from the Ethos and dried according to the ultrasound probe manufacturer's instructions. The Ethos cleaner disinfector is ready for a new cycle immediately after the preceding cycle is completed. A fresh bottle of AquaCide cleaner disinfectant is used with each cycle and mixed inside the Ethos, monitoring of the disinfectant's potency required at the MRC of 1750 ppm and QwikCheck is used. Due to the disinfectant cycling through the entire AquaCide system, the Ethos disinfects itself by the conclusion of the cycle.

correct temperature, soaks the ultrasound probe, and then thoroughly rinses the disinfectant off the ultrasound probe before the cycle is complete. The ultrasound probe is then removed from the TEEClean cleaner disinfector and dried according to the ultrasound probe manufacturer's instructions. The TEEClean cleaner disinfector is ready for a new cycle immediately after the preceding cycle is completed. Because a fresh bottle of TD-5 or TD-8 disinfectant is used with each disinfection cycle, no monitoring of the disinfectant's potency is required, nor is there any requirement for daily testing of the disinfectant solution. Due to the disinfectant cycling through entire TEEClean system, the TEEClean disinfects itself at the conclusion of the disinfection cycle.

cycle is complete. The ultrasound probe is then removed from the TD 200 disinfector and dried according to the TEE probe manufacturer's instructions. The TD 200 disinfector is ready for a new cycle immediately after the preceding cycle is completed. A fresh bottle of TD-12 disinfectant is used with each cycle and mixed inside the TD 200 monitoring of the disinfectant's potency required at the MRC of 1750 ppm and QwikCheck is used. Due to the disinfectant cycling through the entire TD 200 system, the TD 200 disinfects itself by the conclusion of the disinfection cycle.

Comparison of Critical Design Features, Process Monitors, and Process Parameters

Characteristic	Ultrasound Probe Cleaner Disinfector with AquaCide and QwikCheck system	TEE Probe Cleaner	TD 200/TD-12 disinfector/disinfecta nt system	Compare
510(k) number	K230381	K182891	K192228	
	high-level disinfection of	Automated cleaning and high-level disinfection of TEE ultrasound probes	disinfection of TEE	Similar
Cleaner for use with the device	AquaCide (Peracetic Acid) 1750 ppm	TEEZyme < 1% Subtilisins	N/A	Similar
Dedicated Disinfectants for use with the device	Acid)	TD-5 (Glutaraldehyde); TD-8 (Ortho- phthalaldehyde)	TD-12 (Peracetic Acid)	Similar
	AquaCide (Peracetic Acid) 1750 ppm	TD-5 (1.7% Glutaraldehyde); TD-8 (0.3% Ortho- phthalaldehyde)	TD-12 (Peracetic Acid) 1750 ppm	Similar
Disinfectant Buffer System	Sodium Carbonate and Sulfamic Acid	Phosphates	Sodium Carbonate and Sulfamic Acid	Similar

Disinfectant pH	8.5 – 9.0	7.45 – 7.55	8.5 – 9.0	Similar
Operating principles	Peracetic acid	Aldehyde sterilization	Peracetic acid	Similar
	sterilization	,	sterilization	
Process monitors	Digital display screen,	Digital display screen,	Digital display screen,	Same
	printout	printout	printout	
Process parameters	Cleaning: 7 min contact		Cleaning: N/A	Similar
	38 - 54°C	at least 45°C.		
	Disinfection: 3 min	Disinfection: 5 min	Disinfection: 3 min	
	contact at least 47°C.	contact at 38 – 40°C.	contact at least 38°C.	
Software/firmware	Yes	Yes	Yes	Same
control				
Performance claims	Cleaning and High-level disinfection	Cleaning and High-level disinfection	High-level disinfection	Same
Drovidos high lovel	Yes	Yes	Yes	Same
Provides high-level disinfection for heat-	res	res	res	Same
sensitive ultrasound				
probes				
Uses cleaner	Yes	Yes	No	Same
Uses peracetic acid	Yes	No	Yes	Similar
based disinfectant	. 33		. 55	Jiiiiai
	No	Yes	No	Similar
disinfectants				
Single-use disinfectants				Same
only i.e. no open bottle		Yes	Yes	
shelf life of disinfectant				
Single-use disinfectants	18 months	12 months	18 months	Similar
unopened shelf life				
Diagona potibility alsin	Yes	Vaa	Vaa	Como
		Yes	Yes	Same
irritation passing results				
Biocompatibility –	Yes	Yes	Yes	Same
sensitization passing				
results				
Biocompatibility –	Yes	Yes	Yes	Same
cytotoxicity passing				
results				
Toxicology assessment	Yes	Yes	Yes	Same
passing results				
Automated disinfection	Yes	Yes	Yes	Same
cycle	30	1.55	. 30	Janio
Uses 5nm water filter	Yes	Yes	No	Same
for rinse water	. 33	. 33		24110
Disinfection process				Same
user control via	Yes	Yes	Yes	
software functions with				
success/failure print out				
User hazard to				Same
disinfectant	Yes	Yes	Yes	
contact reduced by				

bottle loading system			
User hazard to vapor exposure controlled by vapor management system utilizing air circulation and filtration with no room air	Yes	Yes	Same
circulation required			

Summary of Non-Clinical Studies

Title			Result
Discours of the life of Olding		Acceptance Criteria	Oli sula ti issait a sa t
Biocompatibility Skin	To ensure contact with a	•	Slight irritant
Irritation Test	I •	the study the device extract can be	
		a slight irritant	
	patients		
Biocompatibility	To ensure contact with a	•	Non-sensitizing
Sensitization Test	I •	the study the device extract can be	
		a slight sensitizing	
	patients		
Biocompatibility	To ensure contact with a	ISO 10993-5, under conditions of	Not cytotoxic
Cytotoxicity Test	I -	the study, device extract is not	
		cytotoxic	
	patients		
<b>Toxicology Assessment</b>	To ensure there was no	FDA Submissions for Liquid	Low risk of Toxic
	other toxicological risks	Chemical Sterilants/High Level	
		Disinfectants, evidence/	
		justifications provided device is	
		non-toxic	
Bacillus subtilis Bench	To ensure AquaCide has	(AOAC) Official Method 966.04 for	No growth
Test	sporicidal efficacy for a	5.0 hrs at 45°C have no growth	
	high-level disinfectant	-	
	To ensure AquaCide has	(AOAC) Official Method 966.04 for	No growth
Bench Test	sporicidal efficacy for a	5.0 hrs at 45°C have no growth	
	high-level disinfectant		
Mycobacterium terrae	To ensure AquaCide has	(AOAC) Official Method 965.12 for	No growth
Bench Test	tuberculocidal efficacy for a	3.0 mins at 45°C have at least 6log	
	high-level disinfectant	reduction	
Staphylococcus aureus	To ensure AquaCide has	(AOAC) Official Method 955.15 for	No growth
Bench Test	bactericidal efficacy for a	3.0 mins at 45°C have at least 6log	
	high-level disinfectant	reduction	
Salmonella enterica	To ensure AquaCide has	(AOAC) Official Method 955.14 for	No growth
Bench Test	bactericidal efficacy for a	3.0 mins at 45°C have at least 6log	
	high-level disinfectant	reduction	
Pseudomonas	To ensure AquaCide has	(AOAC) Official Method 964.02 for	No growth
aeruginosa Bench Test	bactericidal efficacy for a	3.0 mins at 45°C have at least 6log	
	high-level disinfectant	reduction	

Trichophyton	To ensure AquaCide has	(AOAC) Official Method 955.17 for	No growth
interdigitale Bench Test	fungicidal efficacy for a	3.0 mins at 45°C have at least 6log	
	high-level disinfectant	reduction	
Herpes Simplex Virus	To ensure AquaCide has	ASTM E1053-20 for 3.0 mins at	Greater than 6log
Type 1 Bench Test	virucidal efficacy for a high-	45°C have greater than 6log	reduction
	level disinfectant	reduction	
Human Influenza Virus	To ensure AquaCide has	ASTM E1053-20 for 3.0 mins at	Greater than 6log
A (H1N1) Bench Test	virucidal efficacy for a high-	45°C have greater than 6log	reduction
	level disinfectant	reduction	
Adenovirus Type 1	To ensure AquaCide has	ASTM E1053-20 for 3.0 mins at	Greater than 6log
Bench Test	virucidal efficacy for a high-	45°C have greater than 6log	reduction
	level disinfectant	reduction	
Mycobacterium Terrae	To ensure when inoculated	Performed in Ethos per FDA	All lots and all
Simulated Use Test	with the most robust	Submissions for Liquid Chemical	probes had greater
	organism, in the most	Sterilants/High Level Disinfectants	than a 6log
	challenging places on	for 3.0 mins at 45°C have greater	reduction
	probes, in worst case	than 6log reduction for all probes	
	conditions, the Ethos and	and all lots	
	AquaCide still perform		
	high-level disinfection		
AquaCide Storage	To ensure that at the end	FDA Submissions for Liquid	18 month shelf life
Stability Test	of the 18 month shelf life	Chemical	
	time the AquaCide can still	Sterilants/High Level Disinfectants	
	meet specifications for	solubilized to meet MRC of 1750	
	high-level disinfection	ppm PAA at 18 months	

#### **Clinical In-Use Testing**

After routine use, soiled surface and endocavity probes were subjected to cleaning and disinfection in the Ethos® cleaner disinfector with AquaCide® and QwikCheck® under standard operating parameters. In all cases there was a measurable complete kill of microorganisms after ultrasound probe processing.

#### Conclusion

Based on the intended use, technological characteristics, non-clinical performance data, and clinical in-use testing, Ethos<sup>®</sup> Automated Ultrasound Probe Cleaner Disinfector with AquaCide<sup>®</sup> Cleaner/High-Level Disinfectant and QwikCheck<sup>®</sup> Chemical Indicator (K230381) is as safe, as effective, and performs as well as or better than the legally marketed predicate device TEEClean<sup>®</sup> Automated TEE Probe Cleaner Disinfector with TEEZyme<sup>®</sup> Cleaner and TD-5<sup>®</sup> or TD-8<sup>®</sup> High-Level Disinfectants (K182891).