

April 13, 2023

Inter Medical Medizintechnik GmbH % Hans-Guenter Osiek Managing Director Daimlerstrasse 34-36 Luebbecke, D-32312 Germany

Re: K230393

Trade/Device Name: UniCam Evo Software Regulation Number: 21 CFR 892.2050

Regulation Name: Medical Image Management And Processing System

Regulatory Class: Class II

Product Code: LLZ Dated: February 14, 2023 Received: February 14, 2023

Dear Hans-Guenter Osiek:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Daniel M. Krainak, Ph.D.

Assistant Director

Magnetic Resonance and Nuclear Medicine Team

Division of Radiological Imaging

and Radiation Therapy Devices

OHT8: Office of Radiological Health

Office of Product Evaluation and Quality

Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

Form Approved: OMB No. 0910-0120 Expiration Date: 06/30/2023

See PRA Statement below.

510(k) Number (if known)
k230393
Device Name
UniCam Evo Software
Indications for Use (Describe)
The UniCam Evo Software package is indicated for nuclear medicine image post-processing software for scintigraphic
and SPECT imaging. UniCam Evo is indicated for processing and quantification of planar and tomographic bone
scintigraphy scans, processing of cerebrovascular scintigraphic scans, processing of planar and tomographic
cardiovascular scintigraphic studies, renal dynamic scans, planar and SPECT lung scintigraphy, planar gastrointestinal
scintigraphic scans, liver scintigraphy, thyroid, parathyroid scintigraphy and displaying multimodal image fusion between
SPECT/CT/MR/PET/ultrasound modalities.
The UniCam Evo Software is intended for nuclear medicine specialists, nuclear medicine radiologists, or trained medicine
technologists. The operator shall have basic computer operation skills.
technologists. The operator shall have basic computer operation skins.
The software loads scintigraphic data sets in Dicom format from the local database or queries the corresponding dicom
data from a PACS archive. The operator may perform analysis on the data, ROI analysis, tomographic reconstruction,
organ dependent procedure steps and the software displays the results in form of images, numerical data or curves. The
operator may print the results, store as an image or dicom secondary capture. Reconstructed transversal, coronal and
sagittal slices can be also stored in Dicom 3.0 format.
Type of Use (Select one or both, as applicable)
Prescription Use (Part 21 CFR 801 Subpart D)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



Inter Medical Medizintechnik GmbH Daimlerstr. 34-36 D-32312 Luebbecke Phone +49 (57 41) 23 59-16 Fax +49 (57 41) 23 59-10 Contact: Mr. H-G. Osiek E Mail Osiek@Intmed.de www.IntMed.de

510(k) Summary

InterMedical Medizintechnik GmbH 510(k) Submitter's Name:

> Daimlerstrasse 34-36 D-32312 Luebbecke

Germany

Contact Person: Hans Guenter Osiek

> Tel: 00495741 2359-0 Fax: 00495741 2359-10 Email: osiek@intmed.de

US Agent: William J. Carroll, Eclipse Systems Inc.

> 422 Briarwood Drive, Guildford, CT 06437

USA

US Agent Contact: William J. Carroll, Managing Director

Eclipse Systems Inc.

14 Commercial st., Unit B

Branford, Connecticut, 06405, USA

Phone: 1- (203) 4830665 Fax: 1-(203) 4830669

Email:wcarroll@eclipsesys.com

Establishment Registration Number:

Date Prepared: 17 July 2022

Device Name:

UniCam Evo Software Trade name:

Evo Software Common Name:

Regulatory Class: Ш

510(k) Summary

Inter Medical Medizintechnik GmbH
Daimlerstr. 34-36
D-32312 Luebbecke
Phone +49 (57 41) 23 59-16
Fax +49 (57 41) 23 59-10
Contact: Mr. H-G. Osiek
E Mail Osiek@Intmed.de

www.IntMed.de

Classification Name: 21CFR 892.2050, Medical image management and

processing system

Panel: Radiology Devices

Product Code: LLZ

Identification of Predicate Device(s):

Manufacturer	Device name	510(k) Number
PHILIPS MEDICAL		
SYSTEMS (CLEVELAND),	ODYSSEY LX, MODEL 211320	K003437
INC. 595 MINER RD.		
CONVERGENT IMAGING	UniSyn Image Fusion	K081987
SOLUTIONS		
49 FIRST AVE. SUITE B		
OTTAWA, ONTARIO, CA		
K1S 2G1		

Device Description:

UniCam Evo is a software package for processing planar, whole body, dynamic planar and SPECT nuclear medicine data sets. The UniCam Evo Software is used for digital images capture, pre-processing, saving, post-processing, multiplanar reconstruction, multimodal image fusion, multimodal image registration, 3D rendering, viewing, printing, archiving and transferring between image viewing workstations. The UniCam Evo Software is intended for nuclear medicine specialists, nuclear medicine radiologists, or trained medicine technologists. The operator shall have basic computer operation skills.

The software loads scintigraphic data sets in Dicom 3.0 format from the local database or queries the corresponding dicom data from a PACS archive. The operator may perform analysis on the data, ROI analysis, tomographic reconstruction, organ dependent procedure steps and the software displays the results in form of images, numerical data or curves. The operator may print the results, store as an image or dicom secondary capture. Reconstructed transversal, coronal and sagittal slices can be also stored in Dicom format.

The data to be processed can be transferred via the DICOM 3.0 Standard to another Nuclear Medicine Workstations of another Manufacturer or from a Gamma Camera / SPECT acquisition workstation.

The software runs under standard Windows operating systems, such as XP, 7, 8, 10, 11. The basic operation principle of the processing program is displayed as flowchart below:



Inter Medical Medizintechnik GmbH
Daimlerstr. 34-36
D-32312 Luebbecke
Phone +49 (57 41) 23 59-16
Fax +49 (57 41) 23 59-10
Contact: Mr. H-G. Osiek
E Mail Osiek@Intmed.de
www.IntMed.de

Intended Use Statement:

The UniCam Evo Software package is indicated for nuclear medicine image post-processing software for scintigraphic and SPECT imaging. UniCam Evo is indicated for processing and quantification of planar and tomographic bone scintigraphy scans, processing of cerebrovascular scintigraphic scans, processing of planar and tomographic cardiovascular scintigraphic studies, renal dynamic scans, planar and SPECT lung scintigraphy, planar gastrointestinal scintigraphic scans, liver scintigraphy, thyroid, parathyroid scintigraphy and displaying multimodal image fusion between SPECT/CT/MR/PET/ultrasound modalities.

The UniCam Evo Software is intended for nuclear medicine specialists, nuclear medicine radiologists, or trained medicine technologists. The operator shall have basic computer operation skills.

The software loads scintigraphic data sets in Dicom format from the local database or queries the corresponding dicom data from a PACS archive. The operator may perform analysis on the data, ROI analysis, tomographic reconstruction, organ dependent procedure steps and the software displays the results in form of images, numerical data or curves. The operator may print the results, store as an image or dicom secondary capture. Reconstructed transversal, coronal and sagittal slices can be also stored in Dicom 3.0 format.

Predicate Device Comparison

The primary legally marketed device is the ODYSSEY LX, MODEL 211320 (K003437) and the second predicate device is the UniSyn Image Fusion (K081987) Software. The intended use of the device and the predicate devices are similar and the device under subject combines features from both predicate devices. The combined features are quantitative analysis of nuclear medicine data processing, reporting, displaying images, storing images in Dicom format, ROI analysis from predicate device 1, and multimodal image fusion, multimodal image registration, displaying, storing images in Dicom format from predicate device 2.

Inter Medical Medizintechnik GmbH
Daimlerstr. 34-36
D-32312 Luebbecke
Phone +49 (57 41) 23 59-16
Fax +49 (57 41) 23 59-10
Contact: Mr. H-G. Osiek
E Mail Osiek@Intmed.de
www.IntMed.de

Device Comparison Chart:

Description	Subject	Primary	Predicate	Significant difference
·	device	Predicate	device 2	, and the second
		device		
Device name and	UniCam Evo Software	ODYSSEY	UniSyn Image Fusion	
510k number		LX, MODEL	(K081987)	
		211320		
		(K003437)		
Intended	The UniCam Evo Software	The Philips	UniSyn is a software	Similar
use/Indications	package is indicated for	Medical	application for image	
for use	nuclear medicine image	Systems	registration and	
	post-processing software for	ODYSSEY	fusion display of	
	scintigraphic and SPECT	LX	scanned image data	
	imaging. UniCam Evo is	computer	from CT, PET, SPECT	
	indicated for processing and	workstatio	and MR scanners. It	
	quantification of planar and	n performs	is to be used by	
	tomographic bone	acquisition	qualified radiology	
	scintigraphy scans,	,	and nuclear medicine	
	processing of	processing,	professionals. UniSyn	
	cerebrovascular	display,	creates multi-planar	
	scintigraphic scans,	archiving,	reformat and	
	processing of planar and	printing	maximum intensity	
	tomographic cardiovascular	and	projection displays of	
	scintigraphic studies, renal	networkin	the data and	
	dynamic scans, planar and	g of	provides	
	SPECT lung scintigraphy,	Nuclear	measurements such	
	planar gastrointestinal	Medicine	as area, volume and	
	scintigraphic scans, liver	data.	Standard Uptake	
	scintigraphy, thyroid,		Values for user	
	parathyroid scintigraphy		defined regions on	
	and displaying multimodal		the image.	
	image fusion between			
	SPECT/CT/MR/PET/ultrasou			
	nd modalities.			
	The UniCam Evo Software is			
	intended for nuclear			
	medicine specialists, nuclear			
	medicine radiologists, or			
	trained medicine			
	technologists. The operator			
	shall have basic computer			
	operation skills.			

Inter Medical Medizintechnik GmbH Daimlerstr. 34-36 D-32312 Luebbecke Phone +49 (57 41) 23 59-16 Fax +49 (57 41) 23 59-10 Contact: Mr. H-G. Osiek E Mail Osiek@Intmed.de www.IntMed.de

				www.IntMed.de
	The software loads			
	scintigraphic data sets in			
	Dicom format from the local			
	database or queries the			
	corresponding dicom data			
	from a PACS archive. The			
	operator may perform			
	1			
	analysis on the data, ROI			
	analysis, tomographic			
	reconstruction, organ			
	dependent procedure steps			
	and the software displays			
	the results in form of			
	images, numerical data or			
	curves. The operator may			
	print the results, store as an			
	image or dicom secondary			
	capture. Reconstructed			
	transversal, coronal and			
	sagittal slices can be also			
	stored in Dicom 3.0 format.			
Device		L		
Where to use	Office settings in clinic or	Office	Office settings in	No significant
	hospital	settings in	clinic or hospital	difference
	nospital	clinic or	ciniic or nospital	direction
		hospital		
Stand-alone	yes	yes	yes	No significant
Software	yes	yes	yes	difference
Web application	no	no	no	No significant
				difference
Mobile medical	no	no	no	No significant
арр				difference
Operating system	Windows 7 64 bit or later	Linux	Microsoft Windows	Odyssey LX and UniSyn
			XP and later	operate on different
				operating system, for
				further details see
				Explanation A) below
User interface	Graphical User Interface	Graphical	Graphical User	No significant
	•	User	Interface	difference
		Interface		
Nuclear Medicine	yes	yes	yes	No significant
data	, 53	,	, 53	difference
quantification				difference
algorithms				
_				
according to the				

Inter Medical Medizintechnik GmbH Daimlerstr. 34-36 D-32312 Luebbecke Phone +49 (57 41) 23 59-16 Fax +49 (57 41) 23 59-10 Contact: Mr. H-G. Osiek
E Mail Osiek@Intmed.de
www.IntMed.de

				www.IntMed.de
SNMMI and				
EANM Procedure				
Standards and				
Guidelines				
ROI drawing,	yes	yes	yes	No significant
reporting, 3D	,	,	,	difference
rendering,				
Multiplanar				
reconstruction,				
Study				
comparison				
	l dicine Scintigraphy Processing	footures		
Thyroid uptake,	1		Not specified	Primary Predicate
1	yes	yes	Not specified	
Parathyroid				device:
suppression,				No significant
Bone two and				difference
three phase,				Predicate device 2 is
Brain				not intended for the
scintigraphy,				full spectrum of
Renal				scintigraphic image
scintigraphy,				processing.
Lung planar,				
Multiple gated				
cardiac blood				
pool imaging,				
First-pass				
radionuclide				
ventriculography,				
HIDA, Gallbladder				
ejection fraction				
calculation,				
Gastric emptying,				
Sentinel node,				
Esophageal				
transit, Salivary				
Gland				
SPECT processing f	eatures	<u> </u>		
Lung, bone, brain	yes	yes	Not specified	Primary Predicate
perfusion,	, 55	,		device:
myocardial				No significant
perfusion				difference
periodicii				uniciciice
				Predicate device 2 is
				not intended for the
	<u> </u>			full spectrum of



Inter Medical Medizintechnik GmbH
Daimlerstr. 34-36
D-32312 Luebbecke
Phone +49 (57 41) 23 59-16
Fax +49 (57 41) 23 59-10
Contact: Mr. H-G. Osiek
E Mail Osiek@Intmed.de

www.IntMed.de

				scintigraphic image
Multimodal image fo	usion			processing.
Image fusion / overlay display, Opacity control, Image registration	yes	Not specified	yes	Primary Predicate device: is not intended for multimodal image fusion. Predicate device 2: No significant difference

Explanation A) The technological characteristics between the UniCam Evo Software and Odyssey LX are different, as the software operates on different operating system. Although Odyssey LX runs under a different operating system with another graphical user interface, both devices are image-processing programs, the different operating systems, do not limit the functionality of the software.

The Unicam Evo Software and the predicate devices are software product that run on PC-based workstations. Image data is input to the devices and used to generate 2D and 3D views, perform image processing, quantitative processing, image fusion and co-registration. Like the two predicate devices, the software has image processing, quantitative analysis, Region Of Interest drawing and image registration abilities, fusion of images from different-modalities, image storage and retrieval, as well as patient information management functions. The UniCam Evo and UniSyn Image Fusion are software products that accept multiple image data types including magnetic resonance, computed tomography, single photon emission computed tomography, or positron emission tomography. Odyssey LX accepts only nuclear medicine, namely, planar and SPECT scintigraphic data sets. The UniCam Evo Software and the predicate devices are capable to display, load, process and store image data in Dicom format.

Nonclinical testing summary

Verification and validation testing confirms that product specifications are met which are equivalent in design and technological characteristics as the predicate devices. The testing results support that the functional testing met for the acceptance of the device. The UniCam Evo Software passed all testing and supports the claims of substantial equivalence to the predicate devices.

The software has been extensively tested using verification and validation test protocols. An emphasis on early testing improves the quality and reliability of software prior to release.



Inter Medical Medizintechnik GmbH
Daimlerstr. 34-36
D-32312 Luebbecke
Phone +49 (57 41) 23 59-16
Fax +49 (57 41) 23 59-10
Contact: Mr. H-G. Osiek
E Mail Osiek@Intmed.de
www.IntMed.de

by comparing - for the selected test studies - the results from the current version with the results from reference values wherever possible the phantom studies are generated and the results from the program are validated with the expected values calculated in an independent way. The results of processing are validated (for a carefully selected series of studies) against the values from the published original reference.

The carried out tests are finally being evaluated and released in an all-embracing risk analysis. The results of the software tests are part of the evaluation of the remaining risks (see Device hazard analysis).

During a review it was checked, if all specifications were implemented in the software as descript in the specification and matches to the requirements.

Each system integration / validation tests contain the following information: the test description, name of the tester, its signature, date of the test, software version number, steps required to perform the test, pass fail criteria for the test step and a checkbox for documenting the passes/failed step.

Summary of the performed tests:

All performed tests passed and have met the acceptable test criteria. There were no test failures.

Summary conclusion:

The device has been designed, verified and validated complying with applicable safety standards for this type of medical equipment. Bench and clinical data demonstrate that processing methods, images and results are equivalent comparing to the predicate devices. No adverse effect has been detected.

Before placing the system on the market and use on human beings, Inter Medical has reviewed all known information and carried out a risk analysis for the software. The comparison table between UniCam Evo Software to predicate devices showed the close similarities and therefore we concluded that it is substantially equivalent to the legally marketed device.