

June 1, 2023

SB-Kawasumi Laboratories, Inc. % Valerie Followell
Official Correspondent for SB-KAWASUMI LABORATORIES, INC.
Regulatory Compliance Associates, Inc. (RCA)
10411 Corporate Drive, Suite 102
Pleasant Prairie, Wisconsin 53158

Re: K230408

Trade/Device Name: K-SHIELD Zen (Model Numbers: PBMH/ PBMA/ PBM)

Regulation Number: 21 CFR 862.1675

Regulation Name: Blood Specimen Collection Device

Regulatory Class: Class II

Product Code: JKA Dated: May 2, 2023 Received: May 3, 2023

Dear Valerie Followell:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/edrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance) and CDRH Learn (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

David Wolloscheck, Ph.D.

Assistant Director

DHT3C: Division of Drug Delivery and

Davil Wallarche 2

General Hospital Devices,

and Human Factors

OHT3: Office of GastroRenal, ObGyn, General Hospital and Urology Devices Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

510(k) Number (if known)

Form Approved: OMB No. 0910-0120 Expiration Date: 06/30/2020

See PRA Statement below.

| K230408 |
|--|
| Device Name K-SHIELD Zen (Model Numbers: PBMH/ PBMA/ PBM) |
| Indications for Use (Describe) K-SHIELD Zen is designed for accessing peripheral veins to collect blood. The Press Button anti-needle stick protector is an integral safety device intended to minimize accidental needlestick injuries. |
| |
| |
| |
| |
| |
| Type of Use (Select one or both, as applicable) |
| Prescription Use (Part 21 CFR 801 Subpart D) Over-The-Counter Use (21 CFR 801 Subpart C) |
| CONTINUE ON A SEPARATE PAGE IF NEEDED. |

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

> Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary -K230408

DATE: 06/01/2023

I. SUBMITTER

Sponsor/Manufacturer

SB-KAWASUMI LABORATORIES, INC.

3-25-4 Tonomachi,

Kawasaki-ku, Kawasaki-shi,

Kanagawa 210-8602, Japan

Tel:+81-44-589-9750

Fax:+81-44-589-9768

US Contact and Official Correspondent for SB-KAWASUMI LABORATORIES, INC.

Valerie Followell

Official Correspondent

Regulatory Compliance Associates®, Inc.

10411 Corporate Drive, Suite 102

Pleasant Prairie, Wisconsin 53158

Phone: (847) 400-6187

Email: v.followell@rcainc.com

II. <u>DEVICE</u>

| Device Trade/Proprietary Name: | K-SHIELD Zen (Model Numbers: PBMH/ PBMA/ PBM) |
|-----------------------------------|---|
| Device Common or Usual Name: | Blood Specimen Collection Device |
| Device Classification Name: | Tubes, Vials, Systems, Serum Separators, Blood Collection |
| Device Regulatory Classification: | Class II |
| Device Classification Regulation: | 21 CFR 862.1675 |
| Product Code: | JKA – Blood Specimen Collection Device |
| | FMI – Needle, Hypodermic, Single Lumen |
| Submission Type: | 510(k) |
| Classification Panel: | General Hospital |

III. PREDICATE DEVICE

Predicate Device:

• K-SHIELD Zen (K220799)

IV. DEVICE DESCRIPTION

The K-SHIELD Zen is a sterile, single use winged blood collection set consisting with a cannula, flexible tube, connection parts and the Press Button anti-needle stick protector designed to minimize accidental needlestick injuries when accessing peripheral veins to collect patient blood. The Press Button anti-needle stick protector stores the needle within the housing when the activation button is pressed after blood collection procedure.

This product has a variety of Needle gauge sizes (21G, 23G, 25G) and a lineup of 3 types of connection parts, Luer connector (PBM), Multiple Sample Luer Adapter (PBMA) and Preattached Holder (PBMH). During blood collection, the Luer connector is connected to a syringe. Multiple Sample Luer Adapter and Pre-attached holder are connected to vacuum blood collection tubes / blood culture bottles to obtain multiple samples. The different configurations (connection options) are chosen by the user based on facility blood collection procedure.

The mechanism of action for the device, all configurations, is blood collection using syringes or vacuum blood collection tubes / blood culture bottles to assist the user with obtaining blood from a patient.

The K-SHIELD Zen is a single use, disposable Ethylene Oxide Sterilized medical device. Sterilization Validation was conducted per ISO 11135: 2014 – Sterilization of Health-Care Products – Ethylene Oxide –Requirements for the Development, Validation, and Routine Control of a Sterilization Process for Medical Devices.

The K-SHIELD Zen is considered blood path, indirect, limited contact (< 24 hours) per *ISO* 10993-1: 2018 - Biological Evaluation of Medical Devices – Part 1: Evaluation and Testing within a Risk Management Process. All device materials comply with ISO 10993-1.

The K-SHIELD Zen is a sterile, single-use device that consists of the following components:

- Winged needle with Press Button
 Anti-needle Stick Protector
- Main Tube
- Clip
- Connector Tube

- Luer Connector
- Locking Cap (optional)
- Multiple Sample Luer Adapter (optional)
- Pre-attached Holder (optional)

The K-SHIELD Zen is intended for use at healthcare facilities or in hospitals.

V. <u>INTENDED USE</u>

K-SHIELD Zen is designed for accessing peripheral veins to collect blood. The Press Button antineedle stick protector is an integral safety device intended to minimize accidental needlestick injuries.

VI. INDICATIONS FOR USE

K-SHIELD Zen is designed for accessing peripheral veins to collect blood. The Press Button antineedle stick protector is an integral safety device intended to minimize accidental needlestick injuries.

VII. COMPARISON OF TECHNOLOGICAL CHARACTERISTICS WITH THE PREDICATE DEVICE

 The comparison of the similarities and differences between the K-SHIELD Zen (subject device) and the K-SHIELD Zen (K220799) (predicate device) is provided in the Substantial Equivalence (SE) Comparison Table below.

| Substantial Equivalence (SE) Comparison Table | | | |
|---|---|---|--------------------------------|
| | K-SHIELD Zen (Subject Device) | K-SHIELD Zen (K220799) (Predicate Device) | Same/Similarities/ Differences |
| Manufacturer | SB-KAWASUMI LABORATORIES, INC | SB-KAWASUMI LABORATORIES, INC | |
| Device Trade or Proprietary Name | K-SHIELD Zen | K-SHIELD Zen | |
| 510(k) Number | K230408 | К220799 | |
| Device Class | Class II | Class II | Same |
| Device Classification Name | Tubes, Vials, Systems, Serum Separators, Blood Collection | Tubes, Vials, Systems, Serum Separators, Blood Collection | Same |
| Device Common Name | Blood specimen collection device | Blood specimen collection device | Same |
| Product Code | JKA and FMI | JKA and FMI | Same |

| Regulation Number | 21 CFR 862.1675 | 21 CFR 862.1675 | Same |
|------------------------------|--|--|---|
| Indications for Use | K-SHIELD Zen is designed for accessing peripheral veins to collect blood. The Press Button anti-needle stick protector is an integral safety device intended to minimize accidental needlestick injuries. | K-SHIELD Zen is designed for accessing peripheral veins to collect blood. The Press Button anti-needle stick protector is an integral safety device intended to minimize accidental needlestick injuries. | Same |
| Intended Use | Same as Indications for Use | Same as Indications for Use | Same |
| Patient Population | Adults and Pediatric | Adults Only | Different (Biocompatibility assessments (DEHP extraction testing and risk assessment supports comparison) |
| Device Description | The K-SHIELD Zen is a sterile, single use winged blood collection set consisting with cannula, flexible tube, connection parts and the Press Button anti-needle stick protector designed to minimize accidental needlestick injuries when accessing peripheral veins to collect patient blood. The Press Button anti-needle stick protector stores the needle within the housing when the activation button is pressed after blood collection procedure. | The K-SHIELD Zen is a sterile, single use winged blood collection set consisting with cannula, flexible tube, connection parts and the Press Button anti-needle stick protector designed to minimize accidental needlestick injuries when accessing peripheral veins to collect patient blood. The Press Button anti-needle stick protector stores the needle within the housing when the activation button is pressed after blood collection procedure. | Same |
| Press Button Feature | Yes | Yes | Same |
| Anti-Needle Stick Prevention | Yes | Yes | Same |
| Color-Coded Wing | Yes | Yes | Same |
| Needle | Stainless Steel | Stainless Steel | Same |

| Noodlo Caugo | 21G, 23G, 25G | 216 226 256 | Same |
|-----------------------------|------------------------------------|------------------------------------|------|
| Needle Gauge Diameter OD | 210, 230, 230 | 21G, 23G, 25G | Same |
| | Lilland Theory Marell | LUkus This MASH | Same |
| Needle | Ultra-Thin Wall | Ultra-Thin Wall | Sume |
| Diameter ID | | | |
| | | | |
| Needle Protector | PP | PP | Same |
| Needle Protector | rr | rr | Same |
| | | | |
| Needle Sheath | Isoprene Rubber | Isoprene Rubber | Same |
| Cover | isopi elle kubbei | isopi elle kubbei | Same |
| | | | |
| Wing | PP | PP | Same |
| Tubing | PVC | PVC | Same |
| Female Luer | ABS | ABS | Same |
| Connector | | | |
| Luer Adaptor | PP | PP | Same |
| Hub | | | |
| Adhesive | Heat Curing Epoxy | Heat Curing Epoxy | Same |
| Packaging | Envelope Pouch | Envelope Pouch | Same |
| Materials Comply | Materials comply with | Materials comply with ISO | Same |
| with | ISO 10993-1 | 10993-1 | |
| ISO 10993-1 | | | |
| Non-Pyrogenic | Non-Pyrogenic | Non-Pyrogenic | Same |
| Non-Toxic | Non-Toxic | Non-Toxic | Same |
| Prescription (Rx) | Prescription Use | Prescription Use | Same |
| or Over- the | | | |
| Counter | | | |
| (OTC) Use | | | |
| Use | Hospital / Clinic | Hospital / Clinic | Same |
| Environment | | | |
| Single-Use | Single-Use Only | Single-Use Only | Same |
| Sterile | Sterile | Sterile | Same |
| Sterility | SAL 10 ⁻⁶ | SAL 10 ⁻⁶ | Same |
| Assurance Level | | | |
| (SAL 10 ⁻ | | | |
| 6) | | | |
| Sterilization Method | Ethylene Oxide Sterilization (EtO) | Ethylene Oxide Sterilization (EtO) | Same |

| Shelf-Life 3 Year Shelf-Life | 3 Year Shelf-Life | Same |
|------------------------------|-------------------|------|
|------------------------------|-------------------|------|

VIII. SUMMARY OF PERFORMANCE DATA AND PERFORMANCE TEST CONCLUSIONS

The following Performance Data were provided in support of the substantial equivalence (SE) determination.

Summary of Non-Clinical Bench Performance Testing

This submission is for the same device but for an expanded patient population and non-clinical performance and biocompatibility data is leveraged/referenced from K220799.

The table below includes the list of the testing results submitted specific to supporting this expanded patient population.

Biocompatibility Testing

Biocompatibility Testing including:

- DEHP Extraction test
- Risk Analysis Report of DEHP Exposure

FDA Recognized Testing Standards:

- ISO 10993-17:2002(First Edition 2002-12-01) Biological Evaluation of Medical Devices – Part 17: Establishment of allowable limits for leachable substances
- ISO 10993-18:2020 (Second edition 2020-01 Amendment 1 2022-05) Biological evaluation of medical devices Part 18: Chemical characterization of
 medical device materials within a risk management process [Including
 Amendment 1 (2022)]

Test Results: PASSED

The results of these Non-Clinical Bench Performance Data were provided in support of the substantial equivalence determination.

Chemical Testing

Chemical Testing

EtO Residuals

FDA Recognized Testing Standards:

ISO 10993-7:2008 (Second Edition 2008-10-15) —
 Biological Evaluation of Medical Devices — Part 7: Ethylene Oxide Sterilization
 Residuals [Including: Technical Corrigendum 1 (2009), AMENDMENT 1: Applicability
 of allowable limits for neonates and infants (2019)]

Test Results: PASSED

The result of the Non-Clinical Bench Performance Data was provided in support of the substantial equivalence determination.

IX. OVERALL CONCLUSIONS

• The K-SHIELD Zen [subject device] is substantially equivalent to the K-SHIELD Zen (K220799) [predicate device] with respect to indications for use, treatment method and technological characteristics.