

April 14, 2023

Abbott Medical Mingzi Deng Associate Director, Regulatory Affairs 4 Robbins Road Westford, Massachusetts 01886

Re: K230411

Trade/Device Name: Dragonfly OpStar<sup>TM</sup> Imaging Catheter

Regulation Number: 21 CFR 870.1200

Regulation Name: Diagnostic intravascular catheter

Regulatory Class: Class II Product Code: DQO Dated: February 14, 2023 Received: February 15, 2023

### Dear Mingzi Deng:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <a href="https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm">https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm</a> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <a href="https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products">https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products</a>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <a href="https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems">https://www.fda.gov/medical-device-problems</a>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<a href="https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</a>) and CDRH Learn (<a href="https://www.fda.gov/training-and-continuing-education/cdrh-learn">https://www.fda.gov/training-and-continuing-education/cdrh-learn</a>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<a href="https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</a>) for more information or contact DICE by email (<a href="DICE@fda.hhs.gov">DICE@fda.hhs.gov</a>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

# Aneesh S. Deoras -S

Aneesh Deoras
Assistant Director
Division of Cardiac Electrophysiology,
Diagnostics and Monitoring Devices
Office of Cardiovascular Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

## DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

### **Indications for Use**

510(k) Number (if known)

Form Approved: OMB No. 0910-0120

Expiration Date: 06/30/2023 See PRA Statement below.

X230411	
Device Name Dragonfly OpStar™ Imaging Catheter	
ndications for Use (Describe) The Dragonfly OpStar <sup>TM</sup> Imaging Catheter with the OCT imaging system is intended for the imaging of coronary arter and is indicated in patients who are candidates for transluminal interventional procedures. The Dragonfly OpStar Imag Catheter is intended for use in vessels 2.0 to 3.5 mm in diameter. The Dragonfly OpStar Imaging Catheter is not intended for use in the left main coronary artery or in a target vessel which has undergone a previous bypass procedure.	ging
Type of Use <i>(Select one or both, as applicable)</i> ☑ Prescription Use (Part 21 CFR 801 Subpart D) ☐ Over-The-Counter Use (21 CFR 801 Subpart C)	

This section applies only to requirements of the Paperwork Reduction Act of 1995.

CONTINUE ON A SEPARATE PAGE IF NEEDED.

#### \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(K) SUMMARY

510(K) SUMMARY						
510(k) Summary Per 21 CFR §807.92						
510(k) Number	K230411					
Date Prepared	February 14, 2023					
	Abbott Medical					
<b>Submitter Name</b>	4 Robbins Road					
& Address	Westford, MA, 01886					
Carrier at Danner	Derek Pike					
Contact Person	978-577-3595					
Alternative	Mingzi Deng					
<b>Contact Person</b>	781-640-4474					
Proprietary /	Dragonfly OpStar <sup>™</sup> Imaging Catheter					
Trade Name						
Common / Usual	Diagnostic Imaging Catheter					
Name						
Product	Product Code: DQO					
Classification						
Product						
Regulation	21 CFR 870.1200					
Number						
Predicate Device	Dragonfly OpStar <sup>TM</sup> Imaging Catheter, AptiVue <sup>TM</sup> Software version E.5.1 (K192019), cleared 11 November 2019					
Device Description	The Dragonfly OpStar Imaging Catheter is a sterile, single-use intravascular catheter consisting of a catheter body external sheath and an internal rotating fiber optic imaging core. The external sheath serves two primary functions: 1) to facilitate placement of the device into the coronary artery and 2) to cover and protect the internal rotating fiber optic imaging core. The inner rotating fiber optic imaging core emits near infrared light to tissues and receives reflected light It is driven by a stainless-steel torque wire visible under fluoroscopy and pulled back through the window tube of the external sheath by the Drive-motor and Optical Controller (DOC). The emitted and returned reflected light are combined and processed by the OPTIS System software to construct an Optical Coherence Tomography (OCT) image. The patient is never exposed to moving parts as the external sheath completely covers the rotating imaging core.					
Indications for Use / Intended Use	The Dragonfly OpStar <sup>™</sup> Imaging Catheter with the OCT imaging system is intended for the imaging of coronary arteries and is indicated in patients who are candidates for transluminal interventional procedures. The Dragonfly OpStar Imaging Catheter is intended for use in vessels 2.0 to 3.5 mm in diameter. The Dragonfly OpStar Imaging Catheter is not intended for use in the left main coronary artery or in a target vessel which has undergone a previous bypass procedure.					

The Dragonfly OpStar Imaging Catheter is substantially equivalent to the predicate Dragonfly OpStar Imaging Catheter (K192019) in terms of intended use, indications for use, operational characteristics, fundamental design, and technological characteristics. There are no technological differences between the predicate device and new device.

	Feature	Predicate Device: Dragonfly OpStar Imaging	Proposed Device: Dragonfly OpStar Imaging		
	reature	Catheter (K192019)	Catheter		
	<b>Intended Use</b>	The Dragonfly OpStar	Same		
		Imaging Catheter with the			
		OCT imaging system is			
		intended for the visualization			
Comparison of		and imaging of coronary			
-		arteries during an			
Subject to		interventional procedure.			
<b>Predicate Device</b>	Indications	Intended for the imaging of	Same		
	for Use	coronary arteries and is			
		indicated in patients who are			
		candidates for transluminal interventional procedures.			
		The Dragonfly OpStar			
		Imaging Catheter is intended			
		for use in vessels 2.0 to 3.5			
		mm in diameter. The			
		Dragonfly OpStar Imaging			
		Catheter is not intended for			
		use in the left main coronary			
		artery or in a target vessel			
		which has undergone a			
		previous bypass procedure.			
	Design verification and validation bench tests were performed on the Dragonfly				
Summary on			ternal design control procedures.		
Non-Clinical	The results demonstrate that the Dragonfly OpStar Imaging Catheter meets the				
Testing	user needs and product specifications and is appropriate for its intended use and				
8	does not raise any new issues of safety and effectiveness.				
Summary of	No clinical testing is provided in this pre-market notification.				
Clinical Testing					
		OpStar Imaging Catheter is subs			
Statement of	predicate Dragonfly OpStar Imaging Catheter (K192019) in terms of intended				
	use, indications for use, operational characteristics, fundamental design, and				
Equivalence	technological characteristics.				