

April 26, 2023

LivsMed Inc. Dong Wook Lee QMR (Quality Management Representative) #304, D-dong, 700, Pangyo-ro, Bundang-gu Seongnam-si, Gyeonggi-do 13516 Korea, South

Re: K230499

Trade/Device Name: ArtiSential Laparoscopic Instruments-Electrodes, Bipolar series (four versions, ABF01series, ABD01 series, ABD02 series and ABD04 series)

Regulation Number: 21 CFR 878.4400 Regulation Name: Electrosurgical cutting and coagulation device and accessories Regulatory Class: Class II Product Code: GEI Dated: February 24, 2023 Received: February 24, 2023

Dear Dong Lee:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Mark Trumbore, Ph.D. Assistant Director DHT4A: Division of General Surgery Devices OHT4: Office of Surgical and Infection Control Devices Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure



ArtiSential Laparoscopic Instruments-Electrodes, Bipolar Series (four versions, ABF01 series, ABD01 series, ABD02 series and ABD04 series) Indications for use Statement

3. Indications for Use Statement

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

510(k) Number *(if known)* K230499

Device Name

ArtiSential Laparoscopic Instruments-Electrodes, Bipolar series (four versions, ABF01 series, ABD01 series, ABD02 series and ABD04 series)

Indications for Use (Describe)

Indications for use include electrosurgical coagulation, dissection, and grasping of tissue during the performance of laparoscopic and general surgical procedures.

Type of Use (Select one or both, as applicable)				
Prescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)			
CONTINUE ON A SEPAR	ATE PAGE IF NEEDED.			
This section applies only to requirements of the Paperwork Reduction Act of 1995. *DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.*				
The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:				
Department of Heal Food and Drug Adm Office of Chief Infor Paperwork Reductio <i>PRAStaff@fda.hhs.</i>	mation Officer on Act (PRA) Staff			

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

1. General Information

Applicant/Submitter:	LivsMed Inc.
Address:	#304, D-dong, 700, Pangyo-ro, Bundang-gu, 13516 Seongnam-si, Gyeonggi-do, Republic of Korea Tel) +82-70-4282-7652 Fax) +82-31-706-3211
Contact Person:	Dong Wook Lee / QMR (Quality Management Representative)
Address:	#304, D-dong, 700, Pangyo-ro, Bundang-gu, 13516 Seongnam-si, Gyeonggi-do, Republic of Korea Tel) +82-70-7709-4993 Fax) +82-31-706-3211 Email) dongwook.livsmed@gmail.com
Preparation Date:	02-24-2023

2. Device Name and Code

Device Trade Name	ArtiSential Laparoscopic Instruments-Electrodes,	
	Bipolar Series (four versions, ABF01 series, ABD01	
	series, ABD02 series and ABD04 series)	
Common Name	Electrosurgical Instruments	
Classification Name	Electrosurgical, cutting & coagulation & accessories	
Product Code	GEI	
Regulation Number	21 CFR 878.4400	
Classification	Class II	
Review Panel	General & Plastic Surgery	

3. Predicate Devices

ArtiSential Laparoscopic Instruments-Electrodes, Bipolar Series(four versions, ABF01 series, ABD01 series, ABD02 series and ABD04 series) are substantially equivalent to the following devices

Table 3.1 Predicate device 1

Applicant	Device Name	510(k) Number
LivsMed Inc.	ArtiSential Bipolar	K200875
	Fenestrated Forceps	

Table 3.2 Predicate device 2

Applicant	Device Name	510(k) Number
LivsMed Inc.	ArtiSential Bipolar	K220384
	Fenestrated Forceps	

4. Device Description

The ArtiSential Laparoscopic Instruments – Electrodes, Bipolar series(four versions, ABF01 series, ABD01 series, ABD02 series and ABD04 series) are sterile, single-use, invasive instruments that used in laparoscopic surgery. There are four versions, ABF01 series, ABD01 series, ABD02 series and ABD04 series. Three versions are same except for jaw shape. This product is a specific component, but not the entire electrosurgical device. The device is not intended to be marketed with multiple components, accessories, and as part of a system.

5. Indications for use

5.1 Indications for use

Indications for use include electrosurgical coagulation, dissection, and grasping of tissue during the performance of laparoscopic and general surgical procedures.

6. Technical Characteristics in Comparison to Predicate Devices

	Proposed device	Predicate Device 1	Predicate Device 2
510(K)	In process	K200875	K220384
Number			
Manufacture	LivsMed, Inc.	LivsMed, Inc.	LivsMed, Inc.
Device Name	ArtiSential Laparoscopic	ArtiSential Laparoscopic	ArtiSential Laparoscopic
	Instruments-Electrodes	Instruments-Electrodes	Instruments-Electrodes
Clearance Date	N/A	05-21-2020	02-24-2022
Classification /	Class 2 / 878.4400	Class 2 / 878.4400	Class 2 / 878.4400
Regulation			
Product Code	GEI	GEI	GEI
Intended for	Prescription Use	Prescription Use	Prescription Use
Indications for	Electrosurgical	Electrosurgical	Electrosurgical
Use	coagulation, dissection,	coagulation, dissection,	coagulation, dissection,
	and grasping of tissue	and grasping of tissue	and grasping of tissue
	during the performance	during the performance	during the performance
	of laparoscopic and	of laparoscopic and	of laparoscopic and
	general surgical	general surgical	general surgical
	procedures.	procedures.	procedures.

Deinsinters			T1 1
Principles of	This product is a single-	This product is a single-	This product is a single-
operation	use instrument used in	use instrument used in	use instrument used in
	electrosurgical units to	electrosurgical units to	electrosurgical units to
	hold soft tissues or	hold soft tissues or	hold soft tissues or
	coagulate and make an	coagulate and make an	coagulate and make an
	incision (tissue	incision (tissue	incision (tissue
	dissection) during	dissection) during	dissection) during
	general laparoscopic	general laparoscopic	general laparoscopic
	surgery, which uses the	surgery, which uses the	surgery, which uses the
	principle of applying	principle of applying	principle of applying
	high-frequency currents	high-frequency currents	high-frequency currents
	from the electrode to the	from the electrode to the	from the electrode to the
	human body to generate	human body to generate	human body to generate
	heat by bioimpedance	heat by bioimpedance	heat by bioimpedance
	when radio frequency	when radio frequency	when radio frequency
	(RF) energy from the	(RF) energy from the	(RF) energy from the
	electrosurgical unit	electrosurgical unit	electrosurgical unit
	applies an electric	applies an electric	applies an electric
	current to the electrode	current to the electrode	current to the electrode
	part, and using the	part, and using the	part, and using the
	generated heat to incise	generated heat to incise	generated heat to incise
	cellular tissues and cause	cellular tissues and cause	cellular tissues and cause
	coagulation.	coagulation.	coagulation.
	It is composed of a jaw,	It is composed of a jaw,	It is composed of a jaw,
	Φ 8 diameter shaft, grip	Φ 8 diameter shaft, grip	Φ 8 diameter shaft, grip
	(including a control	(including a control	(including a control
	ring), and electrosurgical	ring), and electrosurgical	ring), and electrosurgical
	unit connection	unit connection	unit connection
	electrode connector.	electrode connector.	electrode connector.
	During a procedure with	During a procedure with	During a procedure with
	this product, the jaw	this product, the jaw	this product, the jaw
	opens if the control ring	opens if the control ring	opens if the control ring
	opens, and jaw closes if	opens, and jaw closes if	opens, and jaw closes if
	the control ring closes.	the control ring closes.	the control ring closes.
	In addition, the jaw is	In addition, the jaw is	In addition, the jaw is
	also bent up, down, left	also bent up, down, left	also bent up, down, left
	and right within a range	and right within a range	and right within a range
	of $\pm 80^{\circ}$ or more by	of $\pm 80^{\circ}$ or more by	of $\pm 80^{\circ}$ or more by
	moving the grip up,	moving the grip up,	moving the grip up,
	down, left and right, and	down, left and right, and	down, left and right, and
	the jaw can also turn	the jaw can also turn	the jaw can also turn
	360° when rotating the	360° when rotating the	360° when rotating the
	grip.	grip.	grip.
Energy Type	Radiofrequency	Radiofrequency	Radiofrequency
Electrode type	Bipolar	Bipolar	Bipolar
(monopolar or			
bipolar)			
Physical	- Shaft diameter: 8mm	- Shaft diameter: 8mm	- Shaft diameter: 8mm
dimensions	- Shaft Length: 250mm,	- Shaft Length: 250mm,	- Shaft Length: 250mm,
and design	380mm, 450mm	380mm, 450mm	380mm, 450mm
(size, length)			
Rated voltage	200Vp	200Vp	200Vp
Materials	Stainless steel	Stainless steel	Stainless steel
(electrode)			

ArtiSential Laparoscopic Instruments-Electrodes, Bipolar Series (four versions, ABF01 series, ABD01 series, ABD02 series and ABD04 series)

K230499

Materials (insulation)	Polyetherimide	Polyetherimide	Polyetherimide
Materials (Shaft)	Glass fiber	Glass fiber	Glass fiber
Articulating feature	Pitch:±80° or more, Yaw:±80° or more and Open-Close	Pitch:±80° or more, Yaw:±80° or more and Open-Close	Pitch:±80° or more, Yaw:±80° or more and Open-Close
Tip rotation	360°	360°	360°
Sterilization	EO	EO	EO

7. Performance Data

7.1 Biocompatibility

The device has been evaluated for its biological safety according to ISO 10993-1 "Biological Evaluation of Medical Devices – Part 1 : Evaluation and Testing Within a Risk Management Process". Following endpoints have been assessed during the evaluation:

- Cytotoxicity
- Intracutaneous reactivity
- Skin Sensitization
- Acute systemic toxicity
- Pyrogenicity

7.2 Electrical Safety

The ArtiSential Laparoscopic Instruments-Electrodes, Bipolar Series(four versions, ABF01 series, ABD01 series, ABD02 series and ABD04 series) have been tested according to IEC 60601-1, IEC 60601-1-2, IEC 60601-1-6, IEC 60601-2-18 and IEC 60601-2-2. The test setup included:

The device had passed all performed tests.

7.3 Sterilization

ArtiSential Laparoscopic Instruments-Electrodes, Bipolar Series(four versions, ABF01 series, ABD01 series, ABD02 series and ABD04 series) are provided sterile, intended to be single-use. This product is EO-Sterilization in accordance with ISO-11135.

7.4 Shelf life

The proposed expiration date is 3 years from the manufacturing date. The real-time testing will be performed to confirm the shelf-life for 3 years

7.5 Performance test

The device had passed all performed tests.

	m	
Test clause and specification	Test requirement	Results- Remarks
1. Appearance	There should be no defects in	No crack, stain or no
: Implement the visual	the appearance of the product	substances on the surface
inspection for surface of the	and there should be no	of the product
device	problem in use.	
2. Dimension	It shall be within \pm 5% of the	Pass
: Measure by ruler and vernier	indicated value of the	Refer to [Test result] on
calipers	dimensional term.	9-50 page at attachment
		12
3. Operational test	The jaw must be smoothly	The jaw and hub are bent
: Manipulating the grip and	opened and closed and free	up, down, left and right
control ring, and measure the	from jamming, the jaw and	within above 80° and can
angle at bending and rotation	hub can be bent up, down,	rotate 360°.
by goniometer.	left, and right a range of	
	above $\pm 80^{\circ}$ and are capable	
	of 360° rotation.	
4. Tensile strength	The jaw and shaft	No damage to the
: Hold the jaw and shaft	connections shall not be	connection when applying
connections respectively and	damaged from pulling of 20	a force of 20N
apply a force of 20 N using	N.	
Push pull gauge.	14.	
5. Feedthrough test	Electricity should be	The resistance value
: Electrical conduction	transmitted between the	between the electrode tip
		and the connector is less
between the electrode tip and	electrode tip and the	
the connector is tested using a	connector.	than 1Ω
DMM (digital multi meter).		

7.6 Thermal effect

Thermal effects on tissue were also tested. A histological analysis was performed on thermal effect to porcine tissues(liver, kidney and abdominal muscle) through an electrosurgical device.

Based on these performance characteristics, the results demonstrate that the performance requirements were met, the device performs as intended and that the subject device has substantially equivalent performance characteristics to the predicate devices.

8. Substantial Equivalence

ArtiSential Laparoscopic Instruments-Electrodes, Bipolar Series(four versions, ABF01 series, ABD01 series, ABD02 series and ABD04 series) indication for use is same to the predicate device 1 (K200875) and predicate device 2 (K220384). The energy type, electrode type, sterilization as well as physical characteristics are the same. Although there are some minor differences with each product, these differences between the ArtiSential Bipolar Series and the predicate device do not raise new or different questions of safety and efficacy. There is no new technology and no difference that would raise new or different questions of safety or efficacy.

9. Conclusions

In conclusion, the comparison carried out covers all products, models, sizes, and the entire intended purpose of the device under evaluation. The subject device which is the ArtiSential Laparoscopic Instruments-Electrodes, Bipolar Series(four versions, ABF01 series, ABD01 series, ABD02 series and ABD04 series) are same to the predicate device in principles of operation, technological characteristics, as well as performance characteristics. The testing was conducted to evaluate the performance of subject device in comparison to the predicate device. Results of validation and verification activities in design control that include testing/certification to designated standards and performance testing of the devices has demonstrated substantial equivalence of the subject device to the predicate in terms of safety and effectiveness for requested intended use.