

Barco NV % Ms. Julie Vandecandelaere Regulatory Affairs Officer President Kennedypark 35 Kortrijk, W-VL 8500 BELGIUM

March 21, 2023

Re: K230520

Trade/Device Name: Nio Color 2MP (MDNC-2521); Nio Color 3MP (MDNC-3521)

Regulation Number: 21 CFR 892.2050

Regulation Name: Medical image management and processing system

Regulatory Class: Class II

Product Code: PGY

Dated: February 27, 2023 Received: February 27, 2023

Dear Ms. Vandecandelaere:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for

devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/training-and-continuing-education/cdrh-learn) and CDRH Learn (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Jessica Lamb, Ph.D. Assistant Director

Imaging Software Team

DHT8B: Division of Radiological Imaging Devices

and Electronic Products

OHT8: Office of Radiological Health Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

510(k) Number (if known)

Form Approved: OMB No. 0910-0120 Expiration Date: 06/30/2020

See PRA Statement below.

K230520		
Device Name Nio Color 2MP (MDNC-2521) Nio Color 3MP (MDNC-3521)		
Indications for Use (Describe) The display is intended to be used for displaying and viewing digital images (excluding digital mammography) for review and analysis by trained medical practitioners.		
The display is intended to be used for displaying and viewing digital images (excluding digital mammography) for review and analysis by trained medical practitioners. The display may be used in dental applications.		
Type of Use (Select one or both, as applicable)		
Prescription Use (Part 21 CFR 801 Subpart D) Over-The-Counter Use (21 CFR 801 Subpart C)		
CONTINUE ON A SEPARATE PAGE IF NEEDED.		

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary (in accordance with 21 CFR 807.92) K230520		
1. Company	Barco N.V. Healthcare Division 35 President Kennedyr B-8500 Kortrijk BELGIUM	park
2. Contact person	Julie Vandecandelaere Regulatory Affairs Officer Tel: +32 (0)56 26 13 19 julie.vandecandelaere@barco.com	
3. Date of submission	27 February 2023	
4. Device information	Trade name/model: Common name:	Nio Color 2MP (MDNC-2521) Nio Color 3MP (MDNC-3521) MDNC-2521 MDNC-3521
	Classification name: Classification code: Device classification: Regulation number: Regulation name:	System, image processing, Radiological PGY Class 2 892.2050 Medical Image Management and Processing System
5. Predicate device	Common name: Classification name: Classification code: Device classification: Regulation number: Regulation name: MDNC-3521: Nio Color Common name:	MDNC-2221 cleared under 510(K) K133663 MDNC-2221 System, image processing, Radiological LLZ Class 2 892.2050 Medical Image Management and Processing System 3MP (MDNC-3421) cleared under 510(k) K170837 MDNC-3421 System, image processing, Radiological
	Classification name: Classification code:	System, image processing, Radiological PGY



	Device classification: Class 2
	Regulation number: 892.2050
	Regulation name: Medical Image Management and Processing System
6. Device description	Both the Nio Color 2MP (MDNC-2521) and the Nio Color 3MP (MDNC-3521) are medical computer displays designed for general radiology imaging applications. The MDNC-3521 model can also be used in dental applications. The devices can also be used for home reading in radiology.
	The MDNC-2521 is a derivative of the MDNC-2221. The MDNC-3521 is a derivative of the MDNC-3421.
	The modified displays are effectively identical to the respective predicate devices except for the following changes:
	 ✓ Updated LCD panel with same resolution and dimensions compared to the respective predicate devices ✓ New housing, display stand and internal mechanics, with similar functionality and design principle compared to the respective predicate devices ✓ Updated internal electronics boards, with similar functionality and design principle compared to the respective predicate devices ✓ Updated firmware, with similar functionality and design principle compared to the respective predicate devices ✓ New packaging, with similar functionality and design principle compared to the respective predicate devices or any other Barco diagnostic display ✓ Rephrasing of the intended usage environment, because of the evolution towards more home reading in radiology ✓ Small update in the intended user description The modified device has the following similarities compared to the unmodified device: ✓ The same intended use ✓ The same operating principle ✓ The same fundamental technology
	The displays can be used optionally with QAWeb Enterprise software, listed under D332294 as a class 1 device with product code LHO. QAWeb Enterprise is a calibration software that is intended as a quality assurance software for the displays. QAWeb Enterprise software helps to keep the display DICOM compliant.
	The display can be used optionally with Intuitive Workflow Tools, cleared in K191845 as a class 2 device with product code PGY. The Intuitive Workflow Tools are accessories for image enhancement on diagnostic displays:
	 SpotView: The Barco SpotView display feature allows focusing on a region of interest in an image by boosting the display's backlight such that the maximum luminance in provided inside the region of interest. SpotView also enables focused observation during reading by dimming images outside the region of interest and increasing the contrast in this region. Also magnification and inversion of pixels are possible with SpotView.
	 AAM – Application Appearance Manager: This workflow tool allows you to set the luminance as well as the color space for each application that is on the workstation. There are often multiple windows open on a screen, but not all of them need the high brightness of the diagnostic applications.



These applications are addressed by the general term 'Intuitive Workflow Tools'.

The integration of the Intuitive Workflow Tools with the displays have been de-risked, verified and validated to ensure that they do not affect the safety and effectiveness of the displays.

The display can be marketed with or without the Barco MXRT display controller boards. The display controller board is installed in a PACS workstation computer, connected to the display.

7. Intended Use of the Device

MDNC-2521: The display is intended to be used for displaying and viewing digital images (excluding digital mammography) for review and analysis by trained medical practitioners.

MDNC-3521: The display is intended to be used for displaying and viewing digital images (excluding digital mammography) for review and analysis by trained medical practitioners. The display may be used in dental applications.

Note: There are no changes to the indications for use statement from that of the unmodified device.

8. Comparison of technological characteristics

Item	Predicate Device (K133663)	Device for which listing is sought
Device name	Nio Color 2MP (MDNC-2221)	Nio Color 2MP (MDNC-2521)
Screen technology	IPS-Pro	IPS-SFT Color LCD
Active screen size (diagonal)	540 mm (21.3")	541 mm (21.3")
Active screen size (HxV)	432 x 324 mm (17.0 x 12.8")	433 x 325 mm (17.1 x 12.8")
Aspect ratio (H:V)	4:3	4:3
Resolution	2MP (1600 x 1200 pixels)	2MP (1600 x 1200 pixels)
Color imaging	Yes	Yes
Gray imaging	Yes	Yes
Bit depth	30 bit	30 bit
Viewing angle (H, V)	178°	178°
Uniformity Correction	ULT	ULT
SteadyColor Calibration	No	Yes (in MXRT Display Controller), when used as a system with MXRT Display Controller & QAWeb Enterprise
Ambient Light Compensation (ALC)	Yes, reading room selection	Yes, reading room selection
Ambient Light Sensor	Yes	Yes
Front sensor	Yes	Yes
Maximum luminance (panel typical)	800 cd/m ²	1000 cd/m ²



DICOM calibrated luminance	500 cd/m ²	600 Cd/m ²
Contrast ratio (panel typical)	1400:1	2000:1
Response time ((Tr + Tf)/2) (typical)	10 ms	12 ms ^[*] (gray-to-gray average)
Housing color	RAL 9003 / RAL 9004	Black (RAL 9004) / White (RAL 9003)
Video input signals	1x DVI 1x DisplayPort	2x DisplayPort 1.4
USB ports	1x USB 2.0 upstream (endpoint) 2x USB 2.0 downstream	2x USB-B 2.0 upstream (endpoint) 5x USB-A 2.0 downstream (of which 1 charge port)
Power rating	24 VDC, 3.75 A	24 VDC, 4 A
Power consumption	50 W (nominal) < 1 W (hibernate)	37 W (nominal) < 0.35 W (hibernate)
Dimensions with stand (W x H x D)	Portrait: 378 x 525~625 x 235 mm Landscape: 491 x 466~566 x 235 mm	x 225 mm
Dimensions without stand (W x H x D)	Portrait: 378 x 491 x 83 mm Landscape: 491 x 378 x 83 mm	Portrait: 351 x 491 x 64 mm Landscape: 491 x 351 x 64 mm
Dimensions packaged (W x H x D)	655 x 388 x 495 mm	455 x 210 x 770 mm
Net weight with stand	With protective cover: 12.6 kg Without protective cover: 11.3 kg	With protective cover: 8.8 kg Without protective cover: 7.7 kg
Net weight without stand	With protective cover: 7.6 kg Without protective cover: 6.3 kg	With protective cover: 5.8 kg Without protective cover: 4.7 kg
Net weight packaged	With protective cover: 16.8 kg (without optional accessories) Without protective cover: 15.5 kg (without optional accessories)	With protective cover: 12.2 kg (without optional accessories) Without protective cover: 11.2 kg (without optional accessories)
Tilt	-5° to +25°	-10° to +30°
Swivel	-30° to +30°	-30° to +30°
Pivot	90°	90°
Height adjustment range	100 mm	100 mm
Mounting standard	VESA (100 mm)	VESA (100 mm)
Screen protection	Protective, anti-reflective glass cover (optional)	Protective, anti-reflective front glass (optional)
Recommended	All digital images, except	All digital images, except
modalities	digital mammography	digital mammography



Supplied accessories	User guide	User guide
	Documentation disc	Documentation disc
	System sheet	System sheet
	Video cable	Video cables
	Mains cable(s)	USB cables
	USB cable	Mains cables
	External power supply	External power supply
Optional accessories	Graphics board	Display controller
QA software	QAWeb	QAWeb
Warranty	5 years, including 20000 hrs	5 years, including 20000 hrs
	backlight warranty	backlight warranty
Operating temperature	0 °C to 35 °C (15 °C to 30	0 °C to 35 °C (20 °C to 30
	°C within specs)	°C within specs)
Storage temperature	-20 °C to 60 °C	-20 °C to 60 °C
Operating humidity	8 % to 80 % (non-	8% to 80% (non-
	condensing)	condensing)
Storage humidity	5% to 85% (non-	5% to 85% (non-
	condensing)	condensing)
Minimum operating	70 kPa minimum	70 kPa minimum
pressure		
Storage pressure	50 to 106 kPa	50 to 106 kPa

Item	Predicate Device (K170837)	Device for which listing is sought
Device name	Nio Color 3MP (MDNC-3421)	Nio Color 3MP (MDNC-3521)
Screen technology	IPS-TFT Color LCD	IPS-SFT Color LCD
Active screen size (diagonal)	540 mm (21.3")	541 mm (21.3")
Active screen size (HxV)	432 x 324 mm (17.0 x 12.8")	433 x 325 mm (17.1 x 12.8")
Aspect ratio (H:V)	4:3	4:3
Resolution	3MP (2048 x 1536 pixels)	3MP (2048 x 1536 pixels)
Color imaging	Yes	Yes
Gray imaging	Yes	Yes
Bit depth	30 bit	30 bit
Viewing angle (H, V)	178°	178°
Uniformity Correction	ULT	ULT
SteadyColor Calibration	No	Yes (in MXRT Display Controller), when used as a system with MXRT Display Controller & QAWeb Enterprise
Ambient Light Compensation (ALC)	Yes, reading room selection	Yes, reading room selection
Ambient Light Sensor	Yes	Yes
Front sensor	Yes	Yes
Maximum luminance (panel typical)	900 cd/m ²	1050 cd/m ²





11	T===		•
DICOM calibrated luminance	500 cd/m ²	600 Cd/m ²	
Contrast ratio (panel	1400:1	2000:1	
typical)	1100.1	2000.1	
Response time ((Tr +	20 ms	12 ms ^[*] (grey-to-grey	
Tf)/2) (typical)		average)	
Housing color	RAL 9003 / RAL 9004	Black (RAL 9004) / White	
Node a la susta di sucella	DVT D Deval Link Disales Dest	(RAL 9003)	
Video input signals	DVI-D Dual Link DisplayPort	2x DisplayPort 1.4	
USB ports	1x USB 2.0 upstream	2x USB-B 2.0 upstream	
	(endpoint)	(endpoint)	
	3x USB 2.0 downstream	5x USB-A 2.0 downstream	
	24.1/0.0.4.4	(of which 1 charge port)	
Power rating	24 VDC, 4 A	24 VDC, 4 A	
Power consumption	50 W (nominal)	45 W (nominal)	
		< 0.35 W (hibernate)	
Dimensions with stand	Portrait: 378 x 528~628 x	Portrait: 351 x 531~631	
(W x H x D)	235 mm	x 225 mm	
	Landscape: 491 x 472~572 x		
	235 mm	x 225 mm	
Dimensions without	Portrait: 378 x 491 x 84 mm	Portrait: 351 x 491 x 64 mm	
stand (W x H x D)	Landscape: 491 x 378 x 84 mm	Landscape: 491 x 351 x 64 mm	
Dimensions packaged (W x H x D)	500 x 280 x 670 mm	455 x 210 x 770 mm	
Net weight with stand	With protective cover: 11.2	With protective cover: 8.8 kg	
	kg	Without protective cover: 7.7	
	Without protective cover:	kg	
Net weight without	10.7 kg With protective cover: 6.2 kg	With protective cover: 5.8 kg	
stand	Without protective cover: 6.2 kg	With protective cover: 5.8 kg Without protective cover: 4.7	
Staria	kg	kg	
Net weight packaged	With protective cover: 15.7	With protective cover: 12.2	
	kg (without optional	kg (without optional	
	accessories)	accessories)	
	Without protective cover:	Without protective cover:	
	15.2 kg (without optional	11.2 kg (without optional	
	accessories)	accessories)	
Tilt	-10° to +30°	-10° to +30°	
Swivel	-45° to +45°	-30° to +30°	
Pivot	90°	90°	
Height adjustment range	100 mm	100 mm	
Mounting standard	VESA (100 mm)	VESA (100 mm)	
Screen protection	Protective, anti-reflective	Protective, anti-reflective	
	glass cover (optional)	front glass (optional)	
Recommended	All digital images, except	All digital images, except	
modalities	digital mammography	digital mammography	



Supplied accessories	User guide Documentation disc System sheet Video cable Mains cable(s) USB cable External power supply	User guide Documentation disc System sheet Video cables USB cables Mains cables External power supply
Optional accessories	Graphics board	Display controller
QA software	QAWeb	QAWeb
Warranty	5 years, including 20000 hrs backlight warranty	5 years, including 20000 hrs backlight warranty
Operating temperature	0 °C to 40 °C (15 °C to 35 °C within specs)	0 °C to 35 °C (20 °C to 30 °C within specs)
Storage temperature	-20 °C to 60 °C	-20 °C to 60 °C
Operating humidity	8 % to 80 % (non- condensing)	8% to 80% (non- condensing)
Storage humidity	5% to 85% (non- condensing)	5% to 85% (non- condensing)
Minimum operating pressure	70 kPa minimum	70 kPa minimum
Storage pressure	50 to 106 kPa	50 to 106 kPa

[*] The intrinsic response time of the new LCD panel on the modified device is further improved with Barco's RapidFrame technology, a proprietary medical overdrive algorithm which improves the temporal response of the system.

The technological characteristics shown above show that the devices MDNC-2521 and MDNC-3521 are substantially equivalent to their respective predicate devices MDNC-2221 and MDNC-3421 and they do not reveal new issues of safety and performance.

9. Performance testing

The below performance bench tests are performed and corresponding results reported for the modified devices MDNC-2521 and MDNC-3521 in comparison to the already cleared device MDNC-3421 as per the *Physical Laboratory Testing* instructions in "Guidance for Industry and FDA Staff: Display Devices for Diagnostic Radiology", issued in 2022:

- Spatial resolution MTF
- Pixel defects, Artifacts
- Temporal Response
- Maximum and Minimum Luminance
- Luminance response, Conformance to DICOM GSDF
- Angular Dependency of Luminance
- Luminance uniformity
- Reflection coefficient Display Reflectance incl. Specular, Diffuse & Haze coefficients
- Veiling glare or small-spot contrast
- Color tracking

The tests showed that the modified devices MDNC-2521 and MDNC-3521 have similar characteristics compared to already cleared device MDNC-3421 and did not reveal new issues of safety and performance.

Additionally, the modified devices MDNC-2521 and MDNC-3521 are compliant to EMC and Safety standards.

No animal testing or clinical testing has been performed.



10. Conclusion

The Nio Color 2MP (MDNC-2521) and Nio Color 3MP (MDNC-3521) were found to be substantially equivalent to their respective predicate devices MDNC-2221 and MDNC-3421, due to the following reasons:

- a) Device and predicate device have the same intended use respectively
- b) The technological characteristics differences from the predicate device do not affect safety or effectiveness
- c) Bench testing showed that the device has similar characteristics compared to the already cleared device and did not reveal new issues of safety and performance.

