

June 22, 2023

RecensMedical, Inc. Lee Yeonui Regulatory Affairs Manager #507, #908, SK V1 center, 830 Dongtansunhwan-daero Hwaseong-si, Gyeonggi-do 18468 Korea, South

Re: K230599

Trade/Device Name: TargetCoolTM Regulation Number: 21 CFR 878.4350 Regulation Name: Cryosurgical Unit And Accessories Regulatory Class: Class II Product Code: GEH, MLY Dated: May 26, 2023 Received: May 26, 2023

Dear Lee Yeonui:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801; medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,



Mark Trumbore Assistant Director, THT4A1: Robotically-Assisted Surgical Devices Team DHT4A: Division of General Surgery Devices OHT4: Office of Surgical and Infection Control Devices Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number *(if known)* K230599

Device Name TargetCoolTM

Indications for Use (Describe)

TargetCoolTM (Cooling mode / Boosting mode) is indicated for the temporary reduction of pain, swelling, inflammation, and hematoma from minor surgical procedures, minor sprains or other minor sports injuries, and as an adjunct to rehabilitative treatment (e.g., intermittent cold with stretch).

TargetCoolTM (Freezing mode) is indicated for the surgical destruction of target tissue by applying cryogenic gases at extreme low temperatures

- Molluscum Contagiosum
- Skin Tags
- Actinic Keratosis
- Lentigo
- Verruca Plana
- Verruca Vulgaris
- Verruca Lesions
- Genital Lesions
- Seborrheic Keratosis

Type of Use (Select one or both, as applicable)	
Prescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov*

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

1. ADMINISTRATIVE INFORMATION

Manufacturer Name	RecensMedical, Inc. #507, #908, SK V1 center, 830 Dongtansunhwan- daero, Hwaseong-si, Gyeonggi-do, Republic of Korea
Phone:	+82 31 8043 3064
Fax:	+82 31 630 2092
Date of the summary preparation	March 01, 2023
Contact person:	Yeonui Lee / Regulatory Affairs Manager yui.lee@recensmedical.com

2. DEVICE NAME AND CLASSIFICATION

TargetCool TM
Cryosurgical Device
Cryosurgical Unit and Accessories
21 CFR 878.4350 / 21 CFR 878.4810
Class ll
General & Plastic Surgery
GEH, MLY

3. PRIMARY PREDICATE DEVICE

510(k) Number:	K203481
Trade name:	CryoVIVE
Classification name:	Cryosurgical Unit and Accessories
Classification Regulations:	21 CFR 878.4350
Class:	Class ll
Classification Panel:	General & Plastic Surgery
Product code:	GEH, MLY

4. SECOND PRIDICATE DEVICE

510(k) Number: Trade name: Classification name: K221234 TargetCool™ Cryosurgical Unit and Accessories



Classification Regulations:	21 CFR 878.4350
Class:	Class ll
Classification Panel:	General & Plastic Surgery
Product code:	GEH, MLY

5. INDICATIONS FOR USE

TargetCoolTM (Cooling mode, Boosting) is indicated for the temporary reduction of pain, swelling, inflammation, and hematoma from minor surgical procedures, minor sprains or other minor sports injuries, and as an adjunct to rehabilitative treatment (e.g., intermittent cold with stretch).

TargetCoolTM (Freezing mode) is indicated for the surgical destruction of target tissue by applying cryogenic gases at extreme low temperatures

- Molluscum Contagiosum
- Skin Tags
- Actinic Keratosis
- Lentigo
- Verruca Plana
- Verruca Vulgaris
- Verruca Lesions
- Genital Lesions
- Seborrheic Keratosis

6. DEVICE DESCRIPTION

TargetCoolTM is a handheld device that can deliver rapid, precise, and controlled cooling to the skin tissue. The TargetCoolTM device consists of a main device, nozzle, a guard, a filter, and a cartridge.

The main device produces controlled cooling based on thermoelectric cooling, which controls the temperature of the targeted area. TargetCoolTM displays the skin temperature measured in real time, the set cooling temperature and time, and the device status through the LCD display. Also, if the measured temperature is below -1°C and lasts for more than 1 second, the status light blinks in blue with a beep sound.

The principle of pain relief through skin cooling is as follows: As the skin cools by spraying the cryogen onto the skin, the nerve conduction velocity (NCV) of the skin decreases, increasing the pain threshold (PTH) and pain tolerance (PT). Through this mechanism, pain in the skin is relieved. The advances represented by TargetCoolTM lie in the fact that the device executes this process in a precise, controlled, and rapid way.

7. PERFORMANCE DATA

The Company's Performance Data for TargetCool[™] is as follows:

Bench Testing

TargetCoolTM complies with all applicable standards, including ISO 13485:2016, IEC 60601-1 for electrical safety and IEC 60601-1-2 for electromagnetic compatibility.

Biocompatibility (ISO 10993) was also performed to demonstrate conformance with established industry standards. The only patient-contacting material on TargetCoolTM is the guard, which is comprised of polycarbonate. The polycarbonate is the exact same material used in CryoVIVE, which was FDA-cleared under K203481. Therefore, the biocompatibility test is not applicable for the TargetCoolTM.

The device hazard analysis was completed and risk-control implemented to mitigate identified hazards. The testing results support that all the specifications have met the acceptance criteria of each module and interaction of processes.

TargetCoolTM passed all testing and supports the claims of substantial equivalence and safe operation.

Clinical Testing

Clinical publications and clinical data demonstrated the safety and effectiveness of the TargetCoolTM for treatment of pediatric patient.

8. SUBSTANTIAL EQUIVALENCE

The comparison chart below provides evidence to support the equivalence determination between TargetCool[™] and the predicate device (K220674) with respect to intended use, technological characteristics and principles of operation. TargetCool[™] shares the same indications for use, device operation, technical and functional capabilities, and therefore is substantially equivalent to the predicate device.

Product Name	TargetCool [™]	TargetCool™	Comparison
	(Subject Device)	(K221234)	-
Indications for Use / Intended Use	TargetCool [™] is indicated for the	TargetCool [™] is indicated for the	
	temporary reduction of pain, swelling,	temporary reduction of pain, swelling,	
	inflammation, and hematoma from	inflammation, and hematoma from	
	minor surgical procedures, minor	minor surgical procedures, minor	Same
	sprains or other minor sports injuries,	sprains or other minor sports injuries,	Same
	and as an adjunct to rehabilitative	and as an adjunct to rehabilitative	
	treatment (e.g., intermittent cold with	treatment (e.g., intermittent cold with	
	stretch).	stretch).	
	Main system, Control button, LCD,	Main system, Trigger, LCD Display,	
	Cooling nozzle, Freezing nozzle,	Cooling-Nozzle, Guard, Filter, CO ₂	
Component	Cooling guard, filter, CO ₂ cartridge	Cartridge	Difference
	and Boosting accessories (Boosting		
	nozzle, a Boosting guard, and a		

[TargetCoolTM-Cooling mode, Boosting mode]

	Boosting Contain	ner).		
	The unit blows very low-temperature		The unit blows very low-temperature	
Mechanism of	-	er at temperature and	gas at temperature and time settings,	Similar
Action	time settings, on	to the desired	onto the desired treatment area	~~~~~~
	treatment area			
Mode type	Cooling Mod	le, Boosting Mode	Cooling mode	Difference
Cryogen Type	CO ₂ , Cold-wa	ter (Saline solution)	CO_2	Same
Temperature	2-4 °C within 5 sec		2-4 °C within 5 sec	Same
		Continuous (0 sec)		
Treatment	Cooling Mode	PL (2 sec pre-set)	0~60 sec	Similar
Duration		1~60sec	0~00 sec	Similar
	Boosting Mode	3~4 min		
Gas Volume	65g	cartridge	65g cartridge	Same
	Alarm and statu	s light blinking if the	Alarm and status light blinking if the	
Safety feature	temperature of the skin is determined to		temperature of the skin is determined to	Same
	be less than -1 $^{\circ}$ C	C for 1 second.	be less than -1 °C for 1 second.	

The Cooling mode & Boosting mode of TargetCool[™] is equivalent in intended use, principles of operation, and performance temperature to TargetCool[™] (K221234) and raises no new issues of safety or effectiveness.

The differences are the Boosting mode and the Boosting accessories for the Boosting mode.

In Boosting mode, saline cooled by CO_2 gas is sprayed onto the skin. The saline solution is used for the purpose of keeping the skin temperature constant at 2-4°C for 3-4 minutes and for safety purposes to prevent the skin temperature from dropping below freezing.

Cooling mode uses CO_2 gas and Boosting mode uses cold-water, but the temperature of Cooling mode and Boosting mode is the same (2~4°C). Therefore, the difference does not raise new questions of safety and effectiveness

Product Name	TargetCool TM	CryoVIVE	Comparison
FIGURET Name	(Subject Device)	(K203481)	Comparison
	TargetCool [™] are intended for the	The CryoVIVE are intended for the	
	surgical destruction of target tissue by	surgical destruction of target tissue by	
	applying cryogenic gases at extreme	applying cryogenic gases at extreme	
	low temperatures	low temperatures	
	- Molluscum Contagiosum	- Molluscum Contagiosum	
Indications for	- Skin Tags	- Skin Tags	
Use	- Actinic Keratosis	- Actinic Keratosis	Same
/ Intended Use	- Lentigo	- Lentigo	
	- Verruca Plana	- Verruca Plana	
	- Verruca Vulgaris	- Verruca Vulgaris	
	- Verruca Lesions	- Verruca Lesions	
	- Genital Lesions	- Genital Lesions	
	- Seborrheic Keratosis	- Seborrheic Keratosis	
	Main system, Control button, LCD,	Main system, Control button, LCD,	
	Cooling nozzle, Freezing nozzle,	Nozzle (Cooling and Freezing), Guide	
Component	Cooling guard, filter, CO ₂ cartridge	tip, filter, CO ₂ cartridge	Same
	and Boosting accessories (Boosting		
	Nozzle, a Boosting guard, and a		

[TargetCoolTM-Freezing mode]

TargetCool[™]

	Boosting Container).		
Mechanism of Action	Cryogen, CO ₂ is delivered to the treatment site to effect cellular destruction	Cryogen, CO ₂ is delivered to the treatment site to effect cellular destruction	Same
Temperature	CO ₂ (-79°C)	CO ₂ (-79°C)	Same
Gas dispensing rate	0.578 g/sec	0.578 g/sec	Same
Gas Volume	65g cartridge	65g cartridge	Same
Tissue damage	Cell necrosis occurs only inside the ice ball	Cell necrosis occurs only inside the ice ball	Same

Cryoablation is the fundamental technological principle for the freezing mode of the subject device, TargetCoolTM, and the predicate devices (CryoVIVE).

The Freezing mode of TargetCoolTM is equivalent in intended use, principles of operation, and performance temperature to TargetCoolTM (K221234) and raises no new issues of safety or effectiveness. The only difference is that the freezing nozzle's length is shortened and the guard is not inserted during the procedure in the freezing mode. Because of this the device can be placed closer to the lesion, so that it can be applied to smaller lesion sizes. However, these differences do not significantly affect safety and/or effectiveness.

8. CONCLUSION

TargetCoolTM and the legally marketed predicate devices have the same intended use, Indications for Use statement and the technological characteristics. While the technological characteristics differ between the two systems, the differences are minor. Performance testing data established that the TargetCoolTM is safe and effective as the legally marked predicate devices and that the TargetCoolTM does not raise any different questions of safety and effectiveness than the predicate.

On this basis and in accordance with 21 CFR§ 807.100(b), TargetCool[™] is substantially equivalent to the predicate device.