

June 1, 2023

Henan Chaoya Medical Equipment Co., Ltd. % Jarvis Wu Consultant Shanghai Sungo Management Consulting Company Limited 14th Floor, 1500# Central Avenue Shanghai, Shanghai 200122 China

Re: K230617

Trade/Device Name: Surgical Face Mask (L 175*95mm)

Regulation Number: 21 CFR 878.4040 Regulation Name: Surgical Apparel

Regulatory Class: Class II Product Code: FXX Dated: March 6, 2023 Received: March 6, 2023

Dear Jarvis Wu:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

K230617 - Jarvis Wu Page 2

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance) and CDRH Learn (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Brent Showalter -S

Brent Showalter, Ph.D.
Assistant Director
DHT6B: Division of Spinal Devices
OHT6: Office of Orthopedic Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

510(k) Number (if known)

Form Approved: OMB No. 0910-0120
Expiration Date: 06/30/2023

Expiration Date: 06/30/2023 See PRA Statement below.

K230617	
Device Name Surgical Face Mask (L 175*95mm)	
ndications for Use (Describe) The Surgical Face Mask is intended to be worn to protect both the microorganisms, body fluids and particulate material. These face to reduce the potential exposure to blood and body fluids. This is	e masks are intended for use ininfection control practices
Type of Use <i>(Select one or both, as applicable)</i> Prescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary K230617

Document Date Prepared: 2023/5/15

A. Applicant:

Name: Henan Chaoya Medical Equipment Co., Ltd.

Address: No. 19, 13th Street, Economic and Technological Development, Zhengzhou, Henan, China

Contact Person: Ruihua Miao

Tel: +86 18039303229

EMail: 949940319@qq.com

Submission Correspondent:

Primary contact: Mr. Jarvis Wu

Shanghai SUNGO Management Consulting Co., Ltd.

14th floor, 1500# Century Ave., Shanghai 200122, China

Tel: +86-21-58817802

Email: haiyu.wang@sungoglobal.com

Secondary contact: Mr. Raymond Luo

14th floor, 1500# Century Ave., Shanghai 200122, China

Tel: +86-21-68828050

Email: fda.sungo@gmail.com

B. Device:

Proprietary Name: Surgical Face Mask (L 175*95mm)

Common Name: Surgical Face Mask

Model(s): L 175*95mm

Regulatory Information

Classification Name: Surgical Face Mask

Classification: Class II

Product code: FXX

Regulation Number: 878.4040

Review Panel: Surgical Apparel

C. Primary Predicate Device:

510K	Device name	ASTM F2100-19 level	Manufacturer
K221272	Medical Face Mask	Level 2 &3	Shandong Xingyu Gloves Co., Ltd.

(Note: Predicate device has NOT been subject to any Medical Device Recalls, including design-related recall.)

D. Indications use of the device:

The Surgical Face Mask is intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids and particulate material. These face masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single use, disposable device(s), provided non-sterile.

E. Device Description:

The Proposed device(s) are blue color, and flat pleated type mask, utilizing ear loops way for wearing, and they all has Nose clips design for fitting the face mask around the nose.

The proposed device(s) are manufactured with three layers, the inner and outer layers are made of non-woven fabric (polypropylene), and the middle layer is made of melt blown fabric (polypropylene). The model of proposed device, ear loops, is held in place over the users' mouth and nose by two elastic ear loops welded to the face mask. The elastic ear loops are made of Polyester and Spandex. The nose piece contained in the proposed device(s) is in the layers of face mask to allow the user to fit the face mask around their nose, which is made of polypropylene, iron and zinc. The proposed device(s) are sold non-sterile and are intended to be single use, disposable device.

F. Comparison with predicate device

Table 1 General Comparison

Device	Proposed Device	Predicate device	Comparison
	Henan Chaoya Medical Equipment	Shandong Xingyu Gloves Co.,	
Manufacturer	Co., Ltd.	Ltd.	-
510K number	K230617	K221272	-
Classification	Class II Device,	Class II Device,	Same
Classification	FXX (21CFR878.4040)	FXX (21CFR878.4040)	

		The Surgical Face Mask is	The Medical Face Mask is	
		intended to be worn to protect both	intended to be worn to protect	
		the patient and healthcare	both the patient and healthcare	
		personnel from transfer of	personnel from transfer of	
		microorganisms, body fluids and	microorganisms, body fluids and	
Indicatio	ons for use	particulate material. These face	particulate material. The	
inuicatio	ons for use	masks are intended for use in	Disposable Medical Face Mask is	Same
		infection control practices to	intended for use in infection	
		reduce the potential exposure to	control practices to reduce the	
		blood and body fluids. This is a	potential exposure to blood and	
		single use, disposable device(s),	body fluids. This is a single use,	
		provided non-sterile.	disposable device(s), provided	
			non-sterile.	
	Outer	Non-woven fabric(polypropylene)	Non-woven fabric	Same
	layer			
	Middle	Melt blown fabric (polypropylene) Melt-blown Non-woven fabric		Same
Material layer Inner				Same
		Non-woven fabric (polypropylene)	Non-woven fabric	C
layer				Same
Nose		Polypropylene ,iron and zinc	Metal Core Plastic	Different
clip			(Iron wire & Polypropylene)	Different
Ear		Polyester and Spandex	Polyester and Spandex	Same
loops			,	Surre
Color		Blue	Blue, Black	Same
Dimensi	on	175mm ±8mm	175mm+/-5mm	Similar
(length)				
Dimension		95mm±4mm	95mm+/-5mm	Similar
(Width)				
OTC use		Yes	Yes	Same
Sterility		Non-Sterile	Non-Sterile	Same
Use		Single Use, Disposable	Single Use, Disposable	Same
ASTM F2100-19 Lev		Level 2	Level 2 & 3	Same

Henan Chaoya Medical Equipment Co., Ltd. No. 19, 13th Street, Economic and Technological Development, Zhengzhou, Henan, China

level			
	Meet ISO10993 ,proved non-	Meet ISO10993 ,proved non-	Same
Biocompatibility	cytotoxicity, non-irritating and non-	cytotoxicity, non-irritating and non-	
	sensitizing	sensitizing	

From the comparison we found the material of proposed device's nose clip was different from the predicate device. The biocompatibility tests were conducted to both components to ensure their compliance to the ISO10993-5 and ISO10993-10. There is no new risk generated from the difference of the material.

G. Non-Clinical Test Conclusion

Non-clinical tests were conducted to verify that the proposed device met all design specifications as was similar to the predicate device. The test results demonstrated that the proposed device complies with the following standards and the requirements stated in the Guidance for Industry and FDA Staff: Surgical Masks – Premarket Notification [510(k)] Submission issued on March 5, 2004:

- ➤ ISO 10993-5: 2009 Biological Evaluation of Medical Devices -- Part 5: Tests For In Vitro Cytotoxicity
- ➤ ISO 10993-10: 2010 Biological Evaluation of Medical Devices Part 10: Tests For Irritation And Skin Sensitization
- ASTM F2100, Standard Specification for Performance of Materials Used In Medical Face Masks
- ➤ ASTM F1862, Standard Test Method for Resistance of Medical Face Masks To Penetration by Synthetic Blood (Horizontal Projection of Fixed Volume At A Known Velocity);
- ➤ EN 14683, Medical Face Masks—Requirements and Test Methods;
- ➤ ASTM F2101, Standard Test Method for Evaluating the Bacterial Filtration Efficiency (BFE) Of Medical Face Mask Materials, Using A Biological Aerosol of Staphylococcus Aureus;
- ASTM F2299, Stand test method for determining the initial efficiency of materials used in medical face masks to penetration by particulates using latex spheres;
- ➤ 16 CFR 1610, Standard for the Flammability of clothingtextiles;

Table 2 - Performance Testing

Item	Purpose	Proposed device	Acceptance Criteria	Result
Fluid Resistance Performance ASTM F1862	Assess the performance of a mask to resistance to a synthetic blood preparation targeted toward the mask at a set pressure	3 non-consecutive lots tested, using a sample size of 32/lot. 32 out of 32 pass at 120 mmHg	29 out of 32 pass at 120 mmHg for level 2	PASS
Particulate Filtration Efficiency ASTM F2299	Assess the performance of a mask to penetration by sub-micron polystyrene latex particles of 0.1 micron	3 non-consecutive lots tested, using a sample size of 32/lot. Lot1: 98.21% Lot2: 98.21% Lot3: 98.16%	≥ 98%	PASS
Bacterial Filtration Efficiency ASTM F2101	Assess the performance of a mask to penetration by a prepared solution with known concentration of an indicator bacterial organism	3 non-consecutive lots tested, using a sample size of 32/lot. Lot1: 98.21% Lot2: 98.24% Lot3: 98.20%	≥ 98%	PASS
Differential Pressure (Delta P) EN 14683 Annex C	Assess the performance of a mask for resistance to air movement through the materials of the face of the mask	3 non-consecutive lots tested, using a sample size of 32/lot. Lot1:2.52 mmH ₂ O/cm ² Lot2: 2.60 mmH ₂ O/cm ² Lot3: 2.60	< 6.0mmH ₂ O/cm ²	PASS

Henan Chaoya Medical Equipment Co., Ltd. No. 19, 13th Street, Economic and Technological Development, Zhengzhou, Henan, China

		mmH ₂ O/cm ²		
Flammability 16 CFR 1610	Assess the resistance of a mask to ignition	3 non-consecutive lots tested, using a sample size of 32/lot. Class I	Class I	PASS

Table 3 – Biocompatibility Testing

Test Method	Purpose	Acceptance Criteria	Result
	Assess the potential risk		PASS
Cytotoxicity	of Cytotoxicity of mask	Non-Cytotoxic	Under the conditions of the
	material		study, the device is non-
			cytotoxic.
	Assess the potential risk		PASS
Irritation	of Irritation of mask	Non-Irritating	Under the conditions of the
	material		study, the device is non-
			irritating.
	Assess the potential risk		PASS
Sensitization	of Sensitization of mask	Non-Sensitizing	Under the conditions of the
	material		study, the device is non-
			sensitizing

H. Clinical Test Conclusion

No clinical study is included in this submission.

I. Conclusion

The conclusion drawn from the nonclinical tests demonstrates that the subject device in 510(K) submission, the Surgical Face Mask (Model: L 175*95mm) is as safe, as effective, and performs as well as or better than the legally marketed predicate device cleared under K221272.