

April 14, 2023

Radformation, Inc. % Kurt Sysock Co-founder/CEO 335 Madison Avenue, 4th floor NEW YORK NY 10017

Re: K230685

Trade/Device Name: AutoContour Model RADAC V3 Regulation Number: 21 CFR 892.2050 Regulation Name: Medical Image Management And Processing System Regulatory Class: Class II Product Code: QKB Dated: March 9, 2023 Received: March 13, 2023

Dear Kurt Sysock:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part

801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,



Lora D. Weidner, Ph.D. Assistant Director Radiation Therapy Team DHT8C: Division of Radiological Imaging and Radiation Therapy Devices OHT8: Office of Radiological Health Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number *(if known)* K230685

Device Name AutoContour Model RADAC V3

Indications for Use (Describe)

AutoContour is intended to assist radiation treatment planners in contouring and reviewing structures within medical images in preparation for radiation therapy treatment planning.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov*

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

K230685

AutoContour Software

Radformation, Inc.

Special 510(k) Summary

Table of Contents

5.1. Submitter's Information	3
5.2. Device Information	3
5.3. Predicate Device Information	4
5.4. Device Description	4
5.5. Indications for Use	4
5.6. Technological Characteristics	4
5.7. Discussion of differences	10
5.9. Conclusion	25

This 510(k) Summary has been created per the requirements of the Safe Medical Device Act (SMDA) of 1990, and the content is provided in conformance with 21 CFR Part 807.92.

5.1. Submitter's Information

Table 1 : Submitter's Information			
Submitter's Name:	Kurt Sysock		
Company:	Radformation, Inc.		
Address:	335 Madison Avenue, 4th Floor New York, NY 10017		
Contact Person:	Alan Nelson Chief Science Officer, Radformation		
Phone:	518-888-5727		
Fax:			
Email:	anelson@radformation.com		
Date of Summary Preparation	03/09/2023		

5.2. Device Information

Table 2 : Device Information				
Trade Name: AutoContour Model RADAC V3				
Common Name:	AutoContour, AutoContouring, AutoContour Agent, AutoContour Cloud Server			
Classification Name:	Class II			
Classification:	Medical image management and processing system			
Regulation Number:	892.2050			
Product Code:	QKB			
Classification Panel:	Radiology			

5.3. Predicate Device Information

AutoContour Model RADAC V3 (Subject Device) makes use of its prior submissions -AutoContour Model RADAC V2 (K220598) - as the Predicate Device.

5.4. Device Description

As with AutoContour Model RADAC V2, the AutoContour Model RADAC V3 device is software that uses DICOM-compliant image data (CT or MR) as input to: (1) automatically contour various structures of interest for radiation therapy treatment planning using machine learning based contouring. The deep-learning based structure models are trained using imaging datasets consisting of anatomical organs of the head and neck, thorax, abdomen and pelvis for adult male and female patients, (2) allow the user to review and modify the resulting contours, and (3) generate DICOM-compliant structure set data the can be imported into a radiation therapy treatment planning system.

AutoContour Model RADAC V3 consists of 3 main components:

- 1. A .NET client application designed to run on the Windows Operating System allowing the user to load image and structure sets for upload to the cloud-based server for automatic contouring, perform registration with other image sets, as well as review, edit, and export the structure set.
- 2. A local "agent" service designed to run on the Windows Operating System that is configured by the user to monitor a network storage location for new CT and MR datasets that are to be automatically contoured.
- 3. A cloud-based automatic contouring service that produces initial contours based on image sets sent by the user from the .NET client application.

5.5. Indications for Use

AutoContour is intended to assist radiation treatment planners in contouring and reviewing structures within medical images in preparation for radiation therapy treatment planning.

5.6. Technological Characteristics

The Subject Device, AutoContour Model RADAC V3 makes use of AutoContour Model RADAC V2 (K220598) as the Predicate Device for substantial equivalence comparison. The functionality and technical components of this prior submission remain unchanged in AutoContour Model RADAC V3. This submission is intended to build on the technological characteristics of the 510(k) cleared AutoContour Model RADAC V2 pertaining to new structure models for both CT and MRI.

5.6.1. Updates vs. AutoContour (K220598)

The updated submission expands the use of machine-learning based contouring to include additional organs and volumes of Interest found in MR and CT image types.

Table 3: Technological Characteristics AutoContour Model RADAC V3 vs. AutoContour Model RADAC V2 (K220598)				
Characteristic	Subject Device: AutoContour Model RADAC V3Predicate Device: AutoContour M RADAC V2 (K220598)			
Indications for Use	AutoContour is intended to assist radiation treatment planners in contouring and reviewing structures within medical images in preparation for radiation therapy treatment planning.	ng and treatment planners in contouring and reviewing structures within medical images		
Design: Image registration	Manual and Automatic Rigid registration. Automatic Deformable Registration	Manual and Automatic Rigid registration. Automatic Deformable Registration		
Design: Supported modalities	CT or MR input for contouring or registration/fusion. PET/CT input for registration/fusion only. DICOM RTSTRUCT for output	CT or MR input for contouring or registration/fusion. PET/CT input for registration/fusion only. DICOM RTSTRUCT for output		

Regions and	CT or MR input for contouring of	CT or MR input for contouring of		
Volumes of	anatomical regions: Head and Neck,	anatomical regions: Head and Neck,		
interest (ROI)	Thorax, Abdomen and Pelvis.	Thorax, Abdomen and Pelvis.		
	CT Models:	CT Models:		
	A_Aorta	A_Aorta		
	A_Aorta_Asc	A_Aorta_Asc		
	 A_Aorta_Dsc 	 A_Aorta_Dsc 		
	A_LAD	A_LAD		
	A_Pulmonary	Bladder		
	Bladder	Bone_Ilium_L		
	Bladder_F	Bone_Ilium_R		
	Bone_Ilium_L	Bone_Mandible		
	Bone_Ilium_R	Bowel_Bag		
	Bone_Mandible	BrachialPlex_L		
	Bone_Pelvic	BrachialPlex R		
	Bone_Skull	Brain		
	Bone_Sternum	Brainstem		
	Bowel	Breast_L		
	Bowel Bag	Breast R		
	Bowel_Large	Bronchus		
	Bowel_Small	Carina		
	BrachialPlex_L	CaudaEquina		
	BrachialPlex R	Cavity_Oral		
	Brain	Cochlea L		
	Brainstem	Cochlea_R		
	Breast_L			
	Breast R	Ear_Internal_L Ear_Internal_R		
	Bronchus	Esophagus		
	BuccalMucosa	External		
	Carina	• Eye_L		
	CaudaEquina	• Eye_R		
	Cavity_Oral	• Femur_L		
	Cavity_Oral_Ext	 Femul_L Femur_R 		
	Chestwall_L	 Femur_R Femur_RTOG_L 		
	Chestwall_OAR	Femur_RTOG_R		
	Chestwall_R	Glnd_Lacrimal_L		
	Chestwall_RC_L	Gind_Lacrimal_R		
	Chestwall_RC_R	Gind_Lacimia_it Gind_Submand_L		
	Cochlea L	Gind_Submand_R		
	Cochlea R	Glnd_Thyroid		
	Colon_Sigmoid	HDR_Cylinder		
	Cornea_L	Heart		
	Cornea_R	Humerus_L		
	Duodenum	Humerus R		
	Ear_Internal_L	Kidney_L		
	Ear_Internal_R	Kidney_R		
	Esophagus	Kidney_Cuter_L		
	External	Kidney_Outer_R		
	• Eye_L	Larynx		
	• Eye_L • Eye_R	• Lens_L		

Femur_Head_I	_
Femur_Head_I	
Femur_L	LN_Ax_L
Femur_R	LN_Ax_R
Femur_RTOG_	L • LN_IMN_L
Femur_RTOG_	R IN_IMN_R
GallBladder	LN_Neck_IA
 Genitals_F 	 LN_Neck_IB-V_L
 Genitals_M 	 LN_Neck_IB-V_R
 Glnd_Lacrimal_ 	L • LN_Neck_II_L
 Glnd_Lacrimal_ 	R • LN_Neck_II_R
 Glnd_Submand 	_L • LN_Neck_II-IV_L
 Glnd_Submand 	_R
 GInd_Thyroid 	LN_Neck_III_L
HDR_Cylinder	LN_Neck_III_R
Heart	LN_Neck_IV_L
Hippocampus_	L • LN_Neck_IV_R
Hippocampus_	R • LN_Neck_VIA
Humerus_L	LN_Neck_VIIA_L
Humerus_R	 LN_Neck_VIIA_R
 Kidney_L 	 LN_Neck_VIIB_L
 Kidney_R 	 LN_Neck_VIIB_R
 Kidney_Outer_I 	LN_Pelvics
 Kidney_Outer_I 	R • LN_Sclav_L
Larynx	 LN_Sclav_R
Larynx_Glottic	Liver
Larynx_NRG	Lung_L
 Larynx_SG 	 Lung_R
Lens_L	Marrow_Ilium_L
Lens_R	Marrow_Ilium_R
• Lips	Musc_Constrict
Liver	OpticChiasm
 LN_Ax_L 	OpticNrv_L
LN_Ax_L1_L	OpticNrv_R
 LN_Ax_L1_R 	Parotid_L
• LN_Ax_L2_L	Parotid_R
• LN_Ax_L2_L3_	L • PenileBulb
• LN_Ax_L2_L3_	R • Pituitary
• LN_Ax_L2_R	Prostate
• LN_Ax_L3_L	Rectum
• LN_Ax_L3_R	Rib
• LN_Ax_R	SeminalVes
• LN_IMN_L	SpinalCanal
● LN_IMN_R	SpinalCord
• LN_IMN_RC_L	Stomach
• LN_IMN_RC_R	Trachea
● LN_Inguinofen	
LN_Inguinofen	
LN_Neck_IA	MR Models:
LN_Neck_IB-V	L • OpticChiasm
 LN_Neck_IB-V_ 	R • OpticNrv_L

of	25
	of

	LN_Neck_II_L	OpticNrv_R
	 LN_Neck_II_R 	Brainstem
	 LN_Neck_II-IV_L 	Hippocampus_L
	 LN_Neck_II-IV_R 	 Hippocampus_R
	 LN_Neck_II-V_L 	
	 LN_Neck_II-V_R 	
	LN_Neck_III_L	
	LN_Neck_III_R	
	LN_Neck_IV_L	
	LN_Neck_IV_R	
	LN_Neck_V_L	
	 LN_Neck_V_R 	
	LN_Neck_VIA	
	LN_Neck_VIIA_L	
	LN_Neck_VIIA_R	
	LN_Neck_VIIB_L	
	LN_Neck_VIIB_R	
	LN_Paraaortic	
	LN_Pelvics	
	LN_Pelvic_NRG	
	LN_Sclav_L	
	LN_Sclav_R	
	LN_Sclav_RADCOMP_L	
	 LN_Sclav_RADCOMP_R 	
	Lobe_Temporal_L	
	Lobe_Temporal_R	
	Lung_L	
	Lung_R	
	Macula_L	
	Macula_R	
	Marrow_Ilium_L	
	Marrow_Ilium_R	
	Musc_Constrict	
	Nipple_L	
	Nipple_R	
	OpticChiasm	
	OpticNrv_L	
	OpticNrv_R	
	Pancreas	
	Parotid_L	
	Parotid_R	
	PenileBulb	
	Pericardium	
	Pituitary	
	Prostate	
	Rectum	
	Rectum_F	
	Retina_L	
	Retina_R	
	• Rib	
	• Rib_L	
l		

	Rib_R
	SeminalVes
	SpinalCanal
	SpinalCord
	Spleen
	Stomach
	Trachea
	• V_Venacava_I
	V_Venacava_S
	• VB
	• VB_C1
	• VB_C2
	• VB_C3
	• VB_C4
	• VB_C5
	• VB_C6
	• VB_C7
	• VB_L1
	• VB_L2
	• VB_L3
	• VB_L4
	• VB_L5
	• VB_T01
	• VB_T02
	• VB_T03
	• VB_T05
	• VB_T06
	• VB_T07
	• VB_T08
	• VB_T09
	• VB_T10
	• VB_T11
	• VB_T12
Ι.	
	MR Models:
	Brainstem
	Cerebellum
	• Eye_L
	• Eye_R
	GInd_Prostate
	Hippocampus_L
	Hippocampus_R
	● Hypo_True
	Hypothalamus
	OpticChiasm
	OpticNrv_L
	OpticNrv_R
	• OpticTract_L
	• OpticTract_R
	·····

	 Pituitary Prostate SeminalVes 	
Computer platform & Operating System	Windows based .NET front-end application that also serves as agent Uploader supporting Microsoft Windows 10 (64-bit) and Microsoft Windows Server 2016.	Windows based .NET front-end application that also serves as agent Uploader supporting Microsoft Windows 10 (64-bit) and Microsoft Windows Server 2016.
	Cloud-based Server based automatic contouring application compatible with Linux.	Cloud-based Server based automatic contouring application compatible with Linux.
	Windows python-based automatic contouring application supporting Microsoft Windows 10 (64-bit) and Microsoft Windows Server 2016.	Windows python-based automatic contouring application supporting Microsoft Windows 10 (64-bit) and Microsoft Windows Server 2016.

As shown in Table 3, almost all technological characteristics are either substantially equivalent or a subset of the Predicate Device's technological characteristics.

5.7. Discussion of differences

Minor differences

The following minor differences exist, but do not represent any significant additional risks or decreased effectiveness for the device for its intended use:

• New CT Models:

Compared with the Predicate Device, AutoContour Model RADAC V3 supports contouring **90** new models on CT images (the new models are listed below). The addition of these models do not represent a significant deviation from the intended use and operation of AutoContour, nor does it represent a new significant unmitigated risk because:

(a) very similar CNN architecture was used to train these new CT models(b) all new models passed the same DSC test protocol criteria that was applied to the models in the predicate device for similar structure sizes(c) the same risk mitigations that have been applied to the predicate device models have also been applied to all new models

- A_Pulmonary
- Bladder_F
- Bone_Pelvic
- Bone_Skull
- Bone_Sternum
- Bowel
- Bowel_Large
- Bowel_Small

- BuccalMucosa
- Cavity_Oral_Ext
- Chestwall_L
- ChestWall_OAR
- Chestwall_R
- Chestwall_RC_L
- Chestwall_RC_R
- Colon_Sigmoid
- Cornea_L
- Cornea_R
- Duodenum
- Femur_Head_L
- \circ Femur_Head_R
- GallBladder
- Genitals_F
- Genitals_M
- Hippocampus_L
- Hippocampus_R
- Larynx_Glottic
- Larynx_NRG
- Larynx_SG
- LN_Ax_L1_L
- LN_Ax_L1_R
- LN_Ax_L2_L
- LN_Ax_L2_L3_L
- LN_Ax_L2_L3_R
- LN_Ax_L2_R
- LN_Ax_L3_L
- LN_Ax_L3_R
- LN_IMN_RC_L
- LN_IMN_RC_R
- \circ LN_Inguinofem_L
- LN_Inguinofem_R
- $\circ \ LN_Neck_II-V_L$
- $\circ \ LN_Neck_II-V_R$
- \circ LN_Neck_V_L
- \circ LN_Neck_V_R
- LN_Paraaortic
- LN_Pelvics_NRG
- LN_Sclav_RC_L
- LN_Sclav_RC_R
- Lobe_Temporal_L
- Lobe_Temporal_R
- Macula_L

- Macula_R
- Nipple_L
- Nipple_R
- Pancreas
- Pericardium
- Rectum_F
- Retina_L
- Retina_R
- Rib_L
- Rib_R
- \circ Spleen
- UteroCervix
- V_Venacava_I
- VB
- \circ VB_C1
- VB_C2
- VB_C3
- VB_C4
- VB_C5
- VB_C6
- VB_C7
- VB_L1
- VB_L2
- VB_L3VB_L4
- VB_L5
- VB_L3
- VB_T01
- VB_102
- VB_100
- VB T05
- VB_T06
- VB_T07
- VB_T08
- VB_T09
- VB_T10
- VB_T11
- VB_T12
- New MR Models:

Compared with the Predicate Device, AutoContour Model RADAC V3 supports contouring **11** new models on MR images (the new models are listed below). The addition of these models do not represent a significant deviation from the intended use and operation of AutoContour, nor does it represent a new significant unmitigated risk because:

(a) very similar CNN architecture was used to train these new CT models(b) all new models passed the same DSC test protocol criteria that was applied to the models in the predicate device for similar structure sizes(c) the same risk mitigations that have been applied to the predicate device models have also been applied to all new models

- Cerebellum
- Eye_L
- Eye_R
- \circ GInd_Prostate
- Hypo_True
- Hypothalamus
- OpticTract_L
- OpticTract_R
- Pituitary
- Prostate
- \circ SeminalVes

5.8. Performance Data

The following performance data were provided in support of the substantial equivalence determination.

Sterilization & Shelf-life Testing

AutoContour is a pure software device and is not supplied sterile because the device doesn't come in contact with the patient. AutoContour is a pure software device and does not have a Shelf Life.

Biocompatibility

AutoContour is a pure software device and does not come in contact with the patient.

Electrical safety and electromagnetic compatibility (EMC)

AutoContour is a pure software device, hence no Electromagnetic Compatibility and Electrical Safety testing was conducted for the Subject Device.

Software Verification and Validation Testing

<u>Summary</u>

As with the Predicate Device, no clinical trials were performed for AutoContour Model RADAC V3. Non-clinical tests were performed according to Radformation's AutoContour Complete Test Protocol and Report, which demonstrates that AutoContour Model RADAC V3 performs as intended per its indications for use. Further tests were performed on independent datasets from those included in training and validation sets in order to validate the generalizability of the machine learning model.

Description of Changes to Test Protocol

Changes to the testing protocol between AutoContour RADAC V2 and RADAC V3 were made to improve reviewer independence and the validation dataset for the intended population. These changes better demonstrate the ability of the structure model outputs in assisting the user to contour more efficiently as per AutoContour's Indications for Use.

We do not feel that the changes are a significant deviation from the past report as the primary passing criteria is still based on the same minimum mean DSC score and a qualitative review of the structure model output. For RADAC V3 structure models, additional DSC and qualitative review validation was performed on image data that was acquired uniquely from the data used for training. Additionally, independent reviewers (not employed by Radformation) were used to evaluate the clinical appropriateness of structure models as they would be evaluated for the purposes of treatment planning. This external review was performed as a replacement to intraobserver variability testing done with the RADAC V2 structure models as it better quantified the usefulness of the structure model outputs in an unbiased clinical setting.

The RADAC V3 test protocol also adds a section that addresses the validation of any existing structure models that were approved within previous releases of the software (RADAC V2). This regression testing was added as a way to confirm that updates made to the software made for any new releases do not affect the output of any previously approved models. The addition of this test is not significant to the testing of the new structure models as it only confirms that no structure output "drift" has occurred between version releases.

Testing Summary

Mean Dice Similarity Coefficient (DSC) was used to validate the accuracy of structure model outputs when tested on image data sequestered from the original training data population. The test datasets were independent from those used for training and consisted of approximately 10% of the number of training image sets used as input for the model. For CT structure models there were an average of 373 training and 50 testing image sets. Among the patients used for CT training and testing 51.7% were male and 48.3% female. Patient ages range 11-30 : 0.3%, 31-50 : 6.2%, 51-70 : 43.3%, 71-100 : 50.3%. Race 84.0% White, 12.8% Black or African American, 3.2% Other. CT datasets spanned across treatment subgroups most typically found in a radiation therapy treatment clinic with the most common diagnosis being cancers of the Prostate (21%), Breast (21%), Lung (29%), Head and Neck (16%), Other (13%). Images were acquired using a Philips Big Bore CT simulator with the majority of scans having an average slice thickness of 2mm, In-plane resolution between 1-1.2 mm, and acquisition parameters of 120kVp, 674+/-329 average mAs.

Ground truthing of each test data set were generated manually using consensus (NRG/RTOG) guidelines as appropriate by three clinically experienced experts consisting of 2 radiation therapy physicists and 1 radiation dosimetrist.

Structure models were categorized into three size categories as DSC metrics can be sensitive to structure volume. A structure would pass initial validation if the mean DSC exceeded 0.8 for Large volume structures (eg. Bladder, Spleen) 0.65 for Medium volume structures (eg. gallbladder, duodenum) and 0.5 for Small structures (eg Cornea, Retina). For CT Structure models large, medium and small structures resulted in a mean DSC of 0.88+/-0.06, 0.88+/-0.08, and 0.75+/-0.12 respectively. A full summary of the CT structure DSC is available below:

Table 4: CT Training Data Results for AutoContour Model RADAC V3							
CT Structure	Size	Pass Criteria	# of Training Sets	# of Testing Sets	DSC (Avg)	DSC Std Dev	Lower Bound 95% Confidence Interval
A_Pulmonary	Medium	0.65	169	43	0.88	0.03	0.83
Bladder_F	Large	0.8	252	63	0.94	0.03	0.89
Bone_Pelvic	Large	0.8	201	51	0.94	0.01	0.92
Bone_Skull	Large	0.8	80	20	0.92	0.01	0.90
Bone_Sternum	Medium	0.65	80	20	0.9	0.02	0.87
Bowel	Medium	0.65	705	45	0.93	0.08	0.80
Bowel_Large	Medium	0.65	805	52	0.89	0.17	0.61
Bowel_Small	Medium	0.65	705	45	0.93	0.05	0.85
BuccalMucosa	Medium	0.65	392	98	0.7	0.05	0.62
Cavity_Oral_Ext	Medium	0.65	392	98	0.94	0.02	0.91
Chestwall_L	Large	0.8	79	20	0.9	0.03	0.85
Chestwall_R	Large	0.8	79	20	0.9	0.03	0.85
Chestwall_OAR	Large	0.8	118	30	0.9	0.03	0.85
Chestwall_RC_L	Large	0.8	80	20	0.91	0.04	0.84
Chestwall_RC_R	Large	0.8	80	20	0.91	0.04	0.84
Colon_Sigmoid	Medium	0.65	392	98	0.66	0.28	0.20
Cornea_L	Small	0.5	N/A*	N/A*	N/A*	N/A*	N/A
Cornea_R	Small	0.5	N/A*	N/A*	N/A*	N/A*	N/A
Duodenum	Medium	0.65	659	44	0.88	0.16	0.62
Femur_Head_L	Medium	0.65	160	40	0.95	0.04	0.88
Femur_Head_R	Medium	0.65	160	40	0.95	0.04	0.88
Gallbladder	Medium	0.65	512	32	0.96	0.03	0.91
Genitals_F	Large	0.8	233	59	0.92	0.02	0.89
Genitals_M	Large	0.8	173	44	0.93	0.03	0.88
Hippocampus_L	Medium	0.65	226	57	0.67	0.1	0.51
Hippocampus_R	Medium	0.65	226	57	0.67	0.1	0.51

Larynz,NRG Medium 0.65 449 113 0.84 0.04 0.73 Larynz,SG Medium 0.65 449 113 0.84 0.04 0.73 Larynz,SG Medium 0.65 443 110 0.81 0.06 0.71 LN_Ax_L1_R Large 0.8 437 110 0.81 0.06 0.71 LN_Ax_L2_L3_L Medium 0.65 203 51 0.79 0.04 0.72 LN_Ax_L2_R Medium 0.65 203 51 0.79 0.04 0.72 LN_Ax_L2_R Medium 0.65 203 51 0.74 0.07 0.62 LN_Ax_L3_R Medium 0.65 203 51 0.74 0.07 0.62 LN_MN_RC_L Medium 0.65 100 25 0.78 0.05 0.70 LN_IMN_RC_R Medium 0.65 323 81 0.86 0.03 0.81 LN_Inguinofem_L <th>Larynx Glottic</th> <th>Medium</th> <th>0.65</th> <th>438</th> <th>110</th> <th>0.81</th> <th>0.04</th> <th>0.74</th>	Larynx Glottic	Medium	0.65	438	110	0.81	0.04	0.74
Larynz,SG Medium 0.65 413 104 0.78 0.04 0.71 LN_Ax_L1_L Large 0.8 437 110 0.81 0.06 0.71 LN_Ax_L1_R Large 0.8 437 110 0.81 0.06 0.71 LN_Ax_L2_L Medium 0.65 203 51 0.79 0.04 0.72 LN_Ax_L2_L3_L Medium 0.65 437 110 0.82 0.06 0.72 LN_Ax_L2_R Medium 0.65 203 51 0.74 0.07 0.62 LN_Ax_L3_R Medium 0.65 203 51 0.74 0.07 0.62 LN_MAR_L Medium 0.65 100 25 0.78 0.05 0.70 LN_IMN_RC_R Medium 0.65 100 25 0.78 0.05 0.77 LN_Inguinofem_L Large 0.8 310 78 0.85 0.05 0.77 LN_Neck_II-V_L								
LN_AX_L1_L Large 0.8 437 110 0.81 0.06 0.71 LN_AX_L1_R Large 0.8 437 110 0.81 0.06 0.71 LN_AX_L2_L Medium 0.65 203 51 0.79 0.04 0.72 LN_AX_L2_L3_L Medium 0.65 437 110 0.82 0.06 0.72 LN_AX_L2_R Medium 0.65 203 51 0.79 0.04 0.72 LN_AX_L3_R Medium 0.65 203 51 0.74 0.07 0.62 LN_AX_L3_R Medium 0.65 203 51 0.74 0.07 0.62 LN_MAX_L3_R Medium 0.65 100 25 0.78 0.05 0.70 LN_IMN_RC_R Medium 0.65 323 81 0.86 0.03 0.81 LN_Inguinofem_R Large 0.8 310 78 0.85 0.05 0.77 LN_Neck_II-V_L </td <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td>								
LI_A_L_1_R Large 0.8 437 110 0.81 0.06 0.71 LN_AX_L2_L Medium 0.65 203 51 0.79 0.04 0.72 LN_AX_L2_L3_L Medium 0.65 437 110 0.82 0.06 0.72 LN_AX_L2_R Medium 0.65 203 51 0.79 0.04 0.72 LN_AX_L3_R Medium 0.65 203 51 0.74 0.07 0.62 LN_AX_L3_R Medium 0.65 203 51 0.74 0.07 0.62 LN_MN_RC_L Medium 0.65 100 25 0.78 0.05 0.70 LN_IMN_RC_R Medium 0.65 100 25 0.78 0.05 0.77 LN_Inguinofem_L Large 0.8 310 78 0.85 0.05 0.77 LN_Neck_II-V_L Medium 0.65 323 81 0.86 0.03 0.81 LN_Neck_LI-V								
LN_Ax_L2_L Medium 0.65 203 51 0.79 0.04 0.72 LN_Ax_L2_L3_L Medium 0.65 437 110 0.82 0.06 0.72 LN_Ax_L2_R Medium 0.65 437 110 0.82 0.06 0.72 LN_Ax_L2_R Medium 0.65 203 51 0.74 0.07 0.62 LN_Ax_L3_R Medium 0.65 203 51 0.74 0.07 0.62 LN_MN_RC_L Medium 0.65 100 25 0.78 0.05 0.70 LN_IMN_RC_R Medium 0.65 100 25 0.78 0.05 0.70 LN_Inguinofem_L Large 0.8 310 78 0.85 0.05 0.77 LN_Neck_II-V_L Medium 0.65 323 81 0.86 0.03 0.81 LN_Neck_V_L Medium 0.65 267 67 0.8 0.08 0.67 LN_Neck_V_		-						
LN_Ax_L2_L3_L Medium 0.65 437 110 0.82 0.06 0.72 LN_Ax_L2_L3_R Medium 0.65 437 110 0.82 0.06 0.72 LN_Ax_L3_L Medium 0.65 203 51 0.79 0.04 0.72 LN_Ax_L3_L Medium 0.65 203 51 0.74 0.07 0.62 LN_AX_L3_R Medium 0.65 203 51 0.74 0.07 0.62 LN_MR_C_L Medium 0.65 100 25 0.78 0.05 0.70 LN_Inguinofem_L Large 0.8 310 78 0.85 0.05 0.77 LN_Neck_II-V_L Medium 0.65 323 81 0.86 0.03 0.81 LN_Neck_II-V_L Medium 0.65 323 81 0.86 0.03 0.81 LN_Neck_V_L Medium 0.65 267 67 0.8 0.08 0.67 LN_Ne		-						
LN_Ax_L2_L3_R Medium 0.65 437 110 0.82 0.06 0.72 LN_Ax_L3_R Medium 0.65 203 51 0.79 0.04 0.72 LN_Ax_L3_L Medium 0.65 203 51 0.74 0.07 0.62 LN_Ax_L3_R Medium 0.65 203 51 0.74 0.07 0.62 LN_IMN_RC_L Medium 0.65 100 25 0.78 0.05 0.70 LN_Inguinofem_L Large 0.8 310 78 0.85 0.05 0.77 LN_loguinofem_R Large 0.8 310 78 0.85 0.03 0.81 LN_Neck_II-V_L Medium 0.65 323 81 0.86 0.03 0.81 LN_Neck_II-V_R Medium 0.65 267 67 0.8 0.08 0.67 LN_Neck_V_R Medium 0.65 200 51 0.8 0.04 0.73 LNPel								
LN_Ax_L2_R Medium 0.65 203 51 0.79 0.04 0.72 LN_Ax_L3_L Medium 0.65 203 51 0.74 0.07 0.62 LN_Ax_L3_R Medium 0.65 203 51 0.74 0.07 0.62 LN_IMN_RC_L Medium 0.65 100 25 0.78 0.05 0.70 LN_INN_RC_R Medium 0.65 100 25 0.78 0.05 0.70 LN_Inguinofem_L Large 0.8 310 78 0.85 0.05 0.77 LN_Inguinofem_R Large 0.8 310 78 0.86 0.03 0.81 LN_Neck_II-V_L Medium 0.65 323 81 0.86 0.03 0.81 LN_Neck_V_L Medium 0.65 267 67 0.8 0.08 0.67 LN_Neck_V_R Medium 0.65 200 50 0.89 0.04 0.82 LN_Pelvics								
LN_Ax_L3_L Medium 0.65 203 51 0.74 0.07 0.62 LN_Ax_L3_R Medium 0.65 203 51 0.74 0.07 0.62 LN_IMN_RC_L Medium 0.65 100 25 0.78 0.05 0.70 LN_Inguinofem_L Large 0.8 310 78 0.85 0.05 0.77 LN_Inguinofem_R Large 0.8 310 78 0.85 0.05 0.77 LN_Inguinofem_R Large 0.8 310 78 0.85 0.05 0.77 LN_Induck_II-V_L Medium 0.65 323 81 0.86 0.03 0.81 LN_Neck_II-V_L Medium 0.65 267 67 0.8 0.08 0.67 LN_Neck_V_R Medium 0.65 200 50 0.89 0.04 0.82 LN_Perios_NRG Large 0.8 149 38 0.88 0.02 0.85 LN								
LN_Ax_L3_R Medium 0.65 203 51 0.74 0.07 0.62 LN_IMN_RC_L Medium 0.65 100 25 0.78 0.05 0.70 LN_IMN_RC_R Medium 0.65 100 25 0.78 0.05 0.70 LN_Inguinofem_L Large 0.8 310 78 0.85 0.05 0.77 LN_Inguinofem_R Large 0.8 310 78 0.85 0.05 0.77 LN_Neck_II-V_L Medium 0.65 323 81 0.86 0.03 0.81 LN_Neck_II-V_R Medium 0.65 323 81 0.86 0.03 0.81 LN_Neck_U-R Medium 0.65 267 67 0.8 0.08 0.67 LN_Paraaortic Large 0.8 149 38 0.88 0.02 0.85 LN_Sclav_RC_R Medium 0.65 200 51 0.8 0.04 0.73 Lobe								
Ln_IMN_RC_L Medium 0.65 100 25 0.78 0.05 0.70 LN_IMN_RC_R Medium 0.65 100 25 0.78 0.05 0.70 LN_Inguinofem_L Large 0.8 310 78 0.85 0.05 0.77 LN_Inguinofem_R Large 0.8 310 78 0.85 0.05 0.77 LN_Neck_II-V_L Medium 0.65 323 81 0.86 0.03 0.81 LN_Neck_II-V_R Medium 0.65 323 81 0.86 0.03 0.81 LN_Neck_U_R Medium 0.65 267 67 0.8 0.08 0.67 LN_PekcV_R Medium 0.65 267 67 0.8 0.04 0.82 LN_Paraaortic Large 0.8 149 38 0.88 0.02 0.85 LN_Sclav_RC_R Medium 0.65 200 51 0.8 0.04 0.73 Lobe_								
LN_IMN_RC_R Medium 0.65 100 25 0.78 0.05 0.70 LN_Inguinofem_L Large 0.8 310 78 0.85 0.05 0.77 LN_Inguinofem_R Large 0.8 310 78 0.85 0.05 0.77 LN_Neck_II-V_L Medium 0.65 323 81 0.86 0.03 0.81 LN_Neck_II-V_R Medium 0.65 323 81 0.86 0.03 0.81 LN_Neck_VL Medium 0.65 267 67 0.8 0.08 0.67 LN_Neck_VR Medium 0.65 267 67 0.8 0.08 0.67 LN_Pelvics_NRG Large 0.8 200 50 0.89 0.04 0.82 LN_Sclav_RC_L Medium 0.65 200 51 0.8 0.04 0.73 Lobe_Temporal_L Large 0.8 174 44 0.88 0.03 0.83 Mac								
LN_norm Large 0.8 310 78 0.85 0.05 0.77 LN_Inguinofem_R Large 0.8 310 78 0.85 0.05 0.77 LN_Neck_II-V_L Medium 0.65 323 81 0.86 0.03 0.81 LN_Neck_II-V_R Medium 0.65 323 81 0.86 0.03 0.81 LN_Neck_V_L Medium 0.65 267 67 0.8 0.08 0.67 LN_Neck_V_R Medium 0.65 267 67 0.8 0.04 0.82 LN_Pelvics_NRG Large 0.8 200 50 0.89 0.04 0.82 LN_Pelvics_NRG Large 0.8 149 38 0.88 0.02 0.85 LN_Sclav_RC_L Medium 0.65 200 51 0.8 0.04 0.73 Lobe_Temporal_L Large 0.8 174 44 0.88 0.03 0.83 Macula_L								
LN_lnguinofem_R Large 0.8 310 78 0.85 0.05 0.77 LN_Neck_II-V_L Medium 0.65 323 81 0.86 0.03 0.81 LN_Neck_II-V_R Medium 0.65 323 81 0.86 0.03 0.81 LN_Neck_V_R Medium 0.65 267 67 0.8 0.08 0.67 LN_Neck_V_R Medium 0.65 267 67 0.8 0.08 0.67 LN_Pekc_V_R Medium 0.65 267 67 0.8 0.08 0.67 LN_Perivics_NRG Large 0.8 200 50 0.89 0.04 0.82 LN_Sclav_RC_R Medium 0.65 200 51 0.8 0.04 0.73 Lobe_Temporal_L Large 0.8 174 44 0.88 0.03 0.83 Macula_L Small 0.5 120 31 0.64 0.1 0.48 Macula_R<		Large	0.8	310	78	0.85	0.05	
LN_Neck_II-V_R Medium 0.65 323 81 0.86 0.03 0.81 LN_Neck_V_L Medium 0.65 267 67 0.8 0.08 0.67 LN_Neck_V_R Medium 0.65 267 67 0.8 0.08 0.67 LN_Neck_V_R Medium 0.65 267 67 0.8 0.08 0.67 LN_Paraaortic Large 0.8 200 50 0.89 0.04 0.82 LN_Pelvics_NRG Large 0.8 149 38 0.88 0.02 0.85 LN_Sclav_RC_L Medium 0.65 200 51 0.8 0.04 0.73 Lobe_Temporal_L Large 0.8 174 44 0.88 0.03 0.83 Macula_L Small 0.5 120 31 0.64 0.1 0.48 Macula_R Small 0.55 91 23 0.74 0.1 0.58 Nipple_L		-	0.8	310	78	0.85	0.05	0.77
LN_Neck_V_L Medium 0.65 267 67 0.8 0.08 0.67 LN_Neck_V_R Medium 0.65 267 67 0.8 0.08 0.67 LN_Paraaortic Large 0.8 200 50 0.89 0.04 0.82 LN_Paraaortic Large 0.8 149 38 0.88 0.02 0.85 LN_Sclav_RC_L Medium 0.65 200 51 0.8 0.04 0.73 LN_Sclav_RC_R Medium 0.65 200 51 0.8 0.04 0.73 Lobe_Temporal_L Large 0.8 174 44 0.88 0.03 0.83 Macula_L Small 0.5 120 31 0.64 0.1 0.48 Macula_R Small 0.5 120 31 0.64 0.1 0.58 Nipple_L Medium 0.65 91 23 0.74 0.1 0.58 Pancreas Me	LN_Neck_II-V_L	Medium	0.65	323	81	0.86	0.03	0.81
LN_Neck_V_R Medium 0.65 267 67 0.8 0.08 0.67 LN_Paraaortic Large 0.8 200 50 0.89 0.04 0.82 LN_Pelvics_NRG Large 0.8 149 38 0.88 0.02 0.85 LN_Sclav_RC_L Medium 0.65 200 51 0.8 0.04 0.73 LN_Sclav_RC_R Medium 0.65 200 51 0.8 0.04 0.73 Lobe_Temporal_L Large 0.8 174 444 0.88 0.03 0.83 Lobe_Temporal_R Large 0.8 174 44 0.88 0.03 0.83 Macula_L Small 0.5 120 31 0.64 0.1 0.48 Macula_R Small 0.55 91 23 0.74 0.1 0.58 Nipple_R Medium 0.65 91 23 0.74 0.1 0.58 Pancreas <	LN_Neck_II-V_R	Medium	0.65	323	81	0.86	0.03	0.81
LN_ParaaorticLarge0.8200500.890.040.82LN_Pelvics_NRGLarge0.8149380.880.020.85LN_Sclav_RC_LMedium0.65200510.80.040.73LN_Sclav_RC_RMedium0.65200510.80.040.73Lobe_Temporal_LLarge0.8174440.880.030.83Lobe_Temporal_RLarge0.81744440.880.030.83Macula_LSmall0.5120310.640.10.48Macula_RSmall0.5120310.640.10.58Nipple_LMedium0.6591230.740.10.58PancreasMedium0.65706450.920.10.76PericardiumLarge0.8160410.940.020.91Rectum_FMedium0.65252640.910.020.88Retina_LSmall0.5N/A*N/A*N/A*N/A*	LN_Neck_V_L	Medium	0.65	267	67	0.8	0.08	0.67
Ln_Pelvics_NRG Large 0.8 149 38 0.88 0.02 0.85 LN_Sclav_RC_L Medium 0.65 200 51 0.8 0.04 0.73 LN_Sclav_RC_R Medium 0.65 200 51 0.8 0.04 0.73 LN_Sclav_RC_R Medium 0.65 200 51 0.8 0.04 0.73 Lobe_Temporal_L Large 0.8 174 444 0.88 0.03 0.83 Lobe_Temporal_R Large 0.8 174 444 0.88 0.03 0.83 Macula_L Small 0.5 120 31 0.64 0.1 0.48 Macula_R Small 0.5 120 31 0.64 0.1 0.48 Nipple_L Medium 0.65 91 23 0.74 0.1 0.58 Pancreas Medium 0.65 706 45 0.92 0.1 0.76 Pericardium	LN_Neck_V_R	Medium	0.65	267	67	0.8	0.08	0.67
LN_Sclav_RC_L Medium 0.65 200 51 0.8 0.04 0.73 LN_Sclav_RC_R Medium 0.65 200 51 0.8 0.04 0.73 LN_Sclav_RC_R Medium 0.65 200 51 0.8 0.04 0.73 Lobe_Temporal_L Large 0.8 174 44 0.88 0.03 0.83 Lobe_Temporal_R Large 0.8 174 44 0.88 0.03 0.83 Lobe_Temporal_R Large 0.8 174 44 0.88 0.03 0.83 Macula_L Small 0.5 120 31 0.64 0.1 0.48 Macula_R Small 0.5 91 23 0.74 0.1 0.58 Nipple_R Medium 0.65 91 23 0.74 0.1 0.58 Pancreas Medium 0.65 706 45 0.92 0.1 0.76 Pericardium <	LN_Paraaortic	Large	0.8	200	50	0.89	0.04	0.82
LN_Sclav_RC_R Medium 0.65 200 51 0.8 0.04 0.73 Lobe_Temporal_L Large 0.8 174 44 0.88 0.03 0.83 Lobe_Temporal_R Large 0.8 174 44 0.88 0.03 0.83 Lobe_Temporal_R Large 0.8 174 44 0.88 0.03 0.83 Macula_L Small 0.5 120 31 0.64 0.1 0.48 Macula_R Small 0.5 120 31 0.64 0.1 0.48 Nipple_L Medium 0.65 91 23 0.74 0.1 0.58 Nipple_R Medium 0.65 91 23 0.74 0.1 0.58 Pancreas Medium 0.65 706 45 0.92 0.1 0.76 Pericardium Large 0.8 160 41 0.94 0.02 0.88 Retina_L Small <td>LN_Pelvics_NRG</td> <td>Large</td> <td>0.8</td> <td>149</td> <td>38</td> <td>0.88</td> <td>0.02</td> <td>0.85</td>	LN_Pelvics_NRG	Large	0.8	149	38	0.88	0.02	0.85
Lobe_Temporal_L Large 0.8 174 44 0.88 0.03 0.83 Lobe_Temporal_R Large 0.8 174 44 0.88 0.03 0.83 Lobe_Temporal_R Large 0.8 174 44 0.88 0.03 0.83 Macula_L Small 0.5 120 31 0.64 0.1 0.48 Macula_R Small 0.5 120 31 0.64 0.1 0.48 Nipple_L Medium 0.65 91 23 0.74 0.1 0.58 Nipple_R Medium 0.65 91 23 0.74 0.1 0.58 Pancreas Medium 0.65 706 45 0.92 0.1 0.76 Pericardium Large 0.8 160 41 0.94 0.02 0.91 Rectum_F Medium 0.65 252 64 0.91 0.02 0.88 Retina_R Small	LN_Sclav_RC_L	Medium	0.65	200	51	0.8	0.04	0.73
Lobe_Temporal_R Large 0.8 174 44 0.88 0.03 0.83 Macula_L Small 0.5 120 31 0.64 0.1 0.48 Macula_R Small 0.5 91 23 0.74 0.1 0.58 Nipple_R Medium 0.65 91 23 0.74 0.1 0.58 Pancreas Medium 0.65 706 45 0.92 0.1 0.76 Pericardium Large 0.8 160 41 0.94 0.02 0.91 Rectum_F Medium 0.65 252 64 0.91 0.02 0.88 Retina_R Small 0.5 <td>LN_Sclav_RC_R</td> <td>Medium</td> <td>0.65</td> <td>200</td> <td>51</td> <td>0.8</td> <td>0.04</td> <td>0.73</td>	LN_Sclav_RC_R	Medium	0.65	200	51	0.8	0.04	0.73
Macula_L Small 0.5 120 31 0.64 0.1 0.48 Macula_R Small 0.5 120 31 0.64 0.1 0.48 Nipple_L Medium 0.65 91 23 0.74 0.1 0.58 Nipple_R Medium 0.65 91 23 0.74 0.1 0.58 Pancreas Medium 0.65 91 23 0.74 0.1 0.58 Pancreas Medium 0.65 91 23 0.74 0.1 0.58 Pericardium Large 0.65 706 45 0.92 0.1 0.76 Pericardium Large 0.8 160 41 0.94 0.02 0.91 Rectum_F Medium 0.65 252 64 0.91 0.02 0.88 Retina_L Small 0.5 N/A* N/A* N/A* N/A*	Lobe_Temporal_L	Large	0.8	174	44	0.88	0.03	0.83
Macula_R Small 0.5 120 31 0.64 0.1 0.48 Nipple_L Medium 0.65 91 23 0.74 0.1 0.58 Nipple_R Medium 0.65 91 23 0.74 0.1 0.58 Pancreas Medium 0.65 91 23 0.74 0.1 0.58 Pancreas Medium 0.65 706 45 0.92 0.1 0.76 Pericardium Large 0.8 160 41 0.94 0.02 0.91 Rectum_F Medium 0.65 252 64 0.91 0.02 0.88 Retina_L Small 0.5 N/A* N/A* N/A* N/A*	Lobe_Temporal_R	Large	0.8	174	44	0.88	0.03	0.83
Nipple_L Medium 0.65 91 23 0.74 0.1 0.58 Nipple_R Medium 0.65 91 23 0.74 0.1 0.58 Pancreas Medium 0.65 91 23 0.74 0.1 0.58 Pancreas Medium 0.65 706 45 0.92 0.1 0.76 Pericardium Large 0.8 160 41 0.94 0.02 0.91 Rectum_F Medium 0.65 252 64 0.91 0.02 0.88 Retina_L Small 0.5 N/A* N/A* N/A* N/A*	Macula_L	Small	0.5	120	31	0.64	0.1	0.48
Nipple_R Medium 0.65 91 23 0.74 0.1 0.58 Pancreas Medium 0.65 706 45 0.92 0.1 0.76 Pericardium Large 0.8 160 41 0.94 0.02 0.91 Rectum_F Medium 0.65 252 64 0.91 0.02 0.88 Retina_L Small 0.5 N/A* N/A* N/A* N/A* Retina_R Small 0.5 N/A* N/A* N/A* N/A*	Macula_R	Small	0.5	120	31	0.64	0.1	0.48
Pancreas Medium 0.65 706 45 0.92 0.1 0.76 Pericardium Large 0.8 160 41 0.94 0.02 0.91 Rectum_F Medium 0.65 252 64 0.91 0.02 0.88 Retina_L Small 0.5 N/A* N/A* N/A* N/A* Retina_R Small 0.5 N/A* N/A* N/A* N/A*	Nipple_L	Medium	0.65	91	23	0.74	0.1	0.58
Pericardium Large 0.8 160 41 0.94 0.02 0.91 Rectum_F Medium 0.65 252 64 0.91 0.02 0.88 Retina_L Small 0.5 N/A* N/A* N/A* N/A* N/A* Retina_R Small 0.5 N/A* N/A* N/A* N/A*	Nipple_R	Medium	0.65	91	23	0.74	0.1	0.58
Rectum_F Medium 0.65 252 64 0.91 0.02 0.88 Retina_L Small 0.5 N/A* N/A* N/A* N/A* N/A* Retina_R Small 0.5 N/A* N/A* N/A* N/A*	Pancreas	Medium	0.65	706	45	0.92	0.1	0.76
Retina_L Small 0.5 N/A* N/A* N/A* N/A* N/A* Retina_R Small 0.5 N/A* N/A* N/A* N/A* N/A*	Pericardium	Large	0.8	160	41	0.94	0.02	0.91
Retina_R Small 0.5 N/A* N/A* N/A* N/A*	Rectum_F	Medium	0.65	252	64	0.91	0.02	0.88
	Retina_L	Small	0.5	N/A*	N/A*	N/A*	N/A*	N/A*
Rib_L Large 0.8 N/A* N/A* N/A* N/A*	Retina_R	Small	0.5	N/A*	N/A*	N/A*	N/A*	N/A*
	Rib_L	Large	0.8	N/A*	N/A*	N/A*	N/A*	N/A*

Dih D	Lorgo	0.8	N/A*	N/A*	N/A*	N/A*	N/A*
Rib_R	Large						
Spleen	Large	0.8	160	41	0.92	0.08	0.79
UteroCervix	Medium	0.65	143	36	0.81	0.11	0.63
V_Venacava_I	Medium	0.65	399	100	0.81	0.09	0.66
VB	Large	0.8	1051	64	0.99	0.01	0.98
VB_C1	Medium	0.65	196	14	0.90	0.25	0.48
VB_C2	Medium	0.65	202	13	0.99	0.01	0.97
VB_C3	Medium	0.65	214	14	0.97	0.03	0.92
VB_C4	Medium	0.65	234	15	0.90	0.24	0.51
VB_C5	Medium	0.65	328	20	0.89	0.21	0.55
VB_C6	Medium	0.65	520	33	0.87	0.24	0.49
VB_C7	Medium	0.65	651	36	0.98	0.01	0.95
VB_L1	Large	0.8	731	49	0.97	0.09	0.82
VB_L2	Large	0.8	650	44	0.99	0.03	0.93
VB_L3	Large	0.8	582	44	0.99	0.03	0.93
VB_L4	Large	0.8	569	44	0.99	0.01	0.97
VB_L5	Large	0.8	550	43	0.99	0.02	0.95
VB_T01	Medium	0.65	663	37	0.98	0.01	0.96
VB_T02	Medium	0.65	682	38	0.96	0.12	0.76
VB_T03	Medium	0.65	700	38	0.96	0.14	0.73
VB_T04	Medium	0.65	695	39	0.95	0.16	0.68
VB_T05	Medium	0.65	687	39	0.97	0.08	0.84
VB_T06	Medium	0.65	682	36	0.97	0.06	0.86
VB_T07	Medium	0.65	686	39	0.94	0.16	0.69
VB_T08	Medium	0.65	729	45	0.96	0.14	0.73
VB_T09	Medium	0.65	778	49	0.98	0.04	0.92
VB_T10	Medium	0.65	803	49	0.97	0.09	0.82
 VB_T11	Medium	0.65	814	48	0.98	0.07	0.86
 VB_T12	Medium	0.65	795	50	0.97	0.13	0.75
—							

*N/A: Structures are generated based on a post-processing/boolean operation from previously released structure models (Eye, Rib) rather than generated from a CNN model. Quantitative and Qualitative testing for these structures is still performed in the external clinical testing below to validate appropriate contour generation and clinical acceptability of these derived structure models.

Additional external clinical testing was performed in order to validate the accuracy of the models on image sets acquired that were unique to the training datasets. Publically available CT datasets from The Cancer Imaging Archive (TCIA archive) were used and both AutoContour and manually added ground truth contours following the

same structure guidelines used for structure model training were added to the image sets.

	Table	5: CT External Clinical Dataset References
Model Group	Data Source ID	Data Citation
CT Pelvis	TCIA - Pelvic-Ref	Afua A. Yorke, Gary C. McDonald, David Solis Jr., Thomas Guerrero. (2019) Pelvic Reference Data. The Cancer Imaging Archive. DOI: <u>10.7937/TCIA.2019.woskq500</u>
CT Head and Neck	TCIA - Head-Neck-PET-CT	Martin Vallières, Emily Kay-Rivest, Léo Jean Perrin, Xavier Liem, Christophe Furstoss, Nader Khaouam, Phuc Félix Nguyen-Tan, Chang-Shu Wang, Khalil Sultanem. (2017). Data from Head-Neck-PET-CT. The Cancer Imaging Archive. doi: <u>10.7937/K9/TCIA.2017.80je5q00</u>
CT Abdomen	TCIA - Pancreas-CT-CB	Hong, J., Reyngold, M., Crane, C., Cuaron, J., Hajj, C., Mann, J., Zinovoy, M., Yorke, E., LoCastro, E., Apte, A. P., & Mageras, G. (2021). Breath-hold CT and cone-beam CT images with expert manual organ-at-risk segmentations from radiation treatments of locally advanced pancreatic cancer [Data set]. The Cancer Imaging Archive. https://doi.org/10.7937/TCIA.ESHQ-4D90
CT Thorax:	TCIA - NSCLC	Aerts, H. J. W. L., Wee, L., Rios Velazquez, E., Leijenaar, R. T. H., Parmar, C., Grossmann, P., Carvalho, S., Bussink, J., Monshouwer, R., Haibe-Kains, B., Rietveld, D., Hoebers, F., Rietbergen, M. M., Leemans, C. R., Dekker, A., Quackenbush, J., Gillies, R. J., Lambin, P. (2019). Data From NSCLC-Radiomics [Data set]. The Cancer Imaging Archive. https://doi.org/10.7937/K9/TCIA.2015.PF0M9REI
CT Thorax	TCIA - LCTSC	Yang, J., Sharp, G., Veeraraghavan, H., Van Elmpt, W., Dekker, A., Lustberg, T., & Gooding, M. (2017). Data from Lung CT Segmentation Challenge (Version 3) [Data set]. The Cancer Imaging Archive. https://doi.org/10.7937/K9/TCIA.2017.3R3FVZ08

DSC values were calculated between ground truth contour data and AutoContour structures and rated on the same DSC passing criteria used for the Training DSC validation. All structures passed the minimum DSC criteria for small, medium and large structures with an mean DSC of 0.79+/-0.11, 0.83+/-0.12, and 0.90+/-0.09 respectively: Additionally, the qualitative clinical appropriateness of AutoContour structures generated on these scans was graded by clinical experts. Autocontour structures were graded on a scale from 1 to 5 where 5 refers to contour requiring no additional edits, and 1 refers to a score in which full manual re-contour of the structure would be required. An average score >= 3 was used to determine whether a structure model would ultimately be beneficial clinically. An average rating of 4.5 was found across all CT structure models demonstrating that only minor edits would be required in order to make the structure models acceptable for clinical use.

Table 6: CT External Reviewer Results for AutoContour Model RADAC V3										
CT Structure	Size	Pass Criteria	# Testings Sets	Average DSC	Average DSC Std. Dev	Lower Bound 95% Confidence Interval	External Reviewer Average Rating (1-5)			
A_Pulmonary	Medium	0.65	20	0.93	0.02	0.89	4.6			
Bladder_F	Large	0.8	20	0.87	0.22	0.52	4.8			
Bone_Pelvic	Large	0.8	41	0.93	0.04	0.86	4.4			
Bone_Skull	Large	0.8	23	0.98	0.01	0.97	4.5			
Bone_Sternum	Medium	0.65	20	0.92	0.03	0.88	4.8			
Bowel	Medium	0.65	46	0.86	0.07	0.75	4.3			
Bowel_Large	Medium	0.65	46	0.81	0.07	0.70	4.1			
Bowel_Small	Medium	0.65	46	0.77	0.07	0.66	4.0			
BuccalMucosa	Medium	0.65	23	0.68	0.10	0.51	4.2			
Cavity_Oral_Ext	Medium	0.65	23	0.97	0.01	0.95	4.9			
Chestwall_L	Large	0.8	20	0.88	0.10	0.71	4.4			
Chestwall_R	Large	0.8	20	0.90	0.05	0.82	4.4			
Chestwall_OAR	Large	0.8	20	0.94	0.03	0.89	4.8			
Chestwall_RC_L	Large	0.8	20	0.91	0.05	0.83	5.0			
Chestwall_RC_R	Large	0.8	20	0.91	0.03	0.85	5.0			
Colon_Sigmoid	Medium	0.65	40	0.68	0.20	0.35	4.2			
Cornea_L	Small	0.5	23	0.80	0.04	0.73	4.8			
Cornea_R	Small	0.5	23	0.79	0.06	0.69	4.8			
Duodenum	Medium	0.65	25	0.72	0.15	0.48	4.8			
Femur_Head_L	Medium	0.65	41	0.94	0.06	0.85	4.9			
Femur_Head_R	Medium	0.65	41	0.93	0.07	0.82	4.8			
Gallbladder	Medium	0.65	21	0.86	0.05	0.78	4.8			
Genitals_F	Large	0.8	20	0.89	0.04	0.83	4.5			

r							
Genitals_M	Large	0.8	21	0.96	0.03	0.92	4.7
Hippocampus_L	Medium	0.65	23	0.85	0.08	0.71	4.6
Hippocampus_R	Medium	0.65	23	0.87	0.08	0.74	4.6
Larynx_Glottic	Medium	0.65	23	0.83	0.10	0.66	5.0
Larynx_NRG	Medium	0.65	23	0.93	0.05	0.85	4.7
Larynx_SG	Medium	0.65	23	0.81	0.10	0.65	4.8
LN_Ax_L1_L	Large	0.8	20	0.86	0.10	0.70	4.4
LN_Ax_L1_R	Large	0.8	20	0.87	0.08	0.74	4.3
LN_Ax_L2_L	Medium	0.65	20	0.74	0.05	0.65	4.1
LN_Ax_L2_L3_L	Medium	0.65	20	0.86	0.07	0.75	4.3
LN_Ax_L2_L3_R	Medium	0.65	20	0.87	0.04	0.80	4.5
LN_Ax_L2_R	Medium	0.65	20	0.75	0.05	0.67	4.0
LN_Ax_L3_L	Medium	0.65	20	0.73	0.09	0.57	4.2
LN_Ax_L3_R	Medium	0.65	20	0.75	0.08	0.61	4.2
LN_IMN_RC_L	Medium	0.65	20	0.84	0.09	0.68	3.9
LN_IMN_RC_R	Medium	0.65	19	0.81	0.14	0.58	3.9
LN_Inguinofem_L	Large	0.8	40	0.91	0.09	0.77	4.1
LN_Inguinofem_R	Large	0.8	38	0.90	0.08	0.77	4.1
LN_Neck_II-V_L	Medium	0.65	23	0.98	0.01	0.97	4.5
LN_Neck_II-V_R	Medium	0.65	23	0.98	0.01	0.96	4.5
LN_Neck_V_L	Medium	0.65	23	0.79	0.04	0.72	4.6
LN_Neck_V_R	Medium	0.65	23	0.76	0.07	0.65	4.4
LN_Paraaortic	Large	0.8	23	0.88	0.05	0.79	4.6
LN_Pelvics_NRG	Medium	0.65	41	0.79	0.18	0.48	4.4
LN_Sclav_RC_L	Medium	0.65	20	0.72	0.16	0.45	4.3
LN_Sclav_RC_R	Medium	0.65	20	0.70	0.14	0.47	4.2

Lobe_Temporal_L	Large	0.8	23	0.87	0.05	0.79	4.6
Lobe_Temporal_R	Large	0.8	23	0.88	0.05	0.81	4.6
Macula_L	Small	0.5	23	0.70	0.26	0.27	5.0
Macula_R	Small	0.5	23	0.72	0.23	0.34	5.0
Nipple_L	Medium	0.65	35	0.78	0.11	0.60	5.0
Nipple_R	Medium	0.65	38	0.78	0.17	0.50	5.0
Pancreas	Medium	0.65	25	0.71	0.15	0.47	4.1
Pericardium	Large	0.8	20	0.99	0.00	0.98	4.7
Rectum_F	Medium	0.65	20	0.90	0.09	0.76	4.2
Retina_L	Small	0.5	23	0.88	0.03	0.83	5.0
Retina_R	Small	0.5	23	0.88	0.03	0.84	5.0
Rib_L	Large	0.8	20	0.85	0.04	0.78	4.9
Rib_R	Large	0.8	20	0.85	0.04	0.79	4.8
Spleen	Large	0.8	24	0.95	0.07	0.84	4.9
Stomach (Update)	Large	0.8	25	0.90	0.06	0.80	4.9
CaudaEquina (Update)	Medium	0.65	42	0.88	0.06	0.77	4.5
UteroCervix	Medium	0.65	14	0.74	0.23	0.36	3.8
V_Venacava_I	Medium	0.65	41	0.81	0.08	0.68	4.8
VB	Large	0.8	63	0.96	0.03	0.91	4.4
VB_C1	Medium	0.65	23	0.95	0.02	0.92	4.6
VB_C2	Medium	0.65	24	0.96	0.02	0.92	4.7
VB_C3	Medium	0.65	26	0.93	0.10	0.77	4.5
VB_C4	Medium	0.65	33	0.87	0.20	0.54	4.6
VB_C5	Medium	0.65	40	0.87	0.16	0.61	4.6
VB_C6	Medium	0.65	42	0.90	0.07	0.78	4.6
VB_C7	Medium	0.65	42	0.93	0.04	0.87	4.8

VB_L1	Large	0.8	32	0.88	0.22	0.53	4.7
VB_L2	Large	0.8	25	0.89	0.21	0.55	4.7
VB_L3	Large	0.8	20	0.92	0.19	0.60	4.9
VB_L4	Large	0.8	20	0.93	0.21	0.58	4.8
VB_L5	Large	0.8	20	0.92	0.20	0.60	4.4
VB_T01	Medium	0.65	43	0.94	0.13	0.72	4.7
VB_T02	Medium	0.65	43	0.95	0.05	0.86	4.5
	Medium	0.65	43	0.93	0.12	0.73	4.9
 VB T04	Medium	0.65	43	0.93	0.14	0.69	4.6
 VB T05	Medium	0.65	44	0.90	0.20	0.58	4.6
VB T06	Medium	0.65	43	0.88	0.22	0.51	4.7
VB_T07	Medium	0.65	37	0.83	0.29	0.36	4.7
VB T08	Medium	0.65	27	0.79	0.33	0.25	4.7
VB T09	Medium	0.65	23	0.79	0.32	0.26	4.7
VB_T10	Medium	0.65	23	0.81	0.31	0.29	4.7
VB_T11	Medium	0.65	24	0.82	0.31	0.30	4.7
VB_T12	Medium	0.65	31	0.86	0.24	0.46	4.8

The MR training data set used for initial testing of the Brain models (Cerebellum, Hypothalamus, Hypo_True, OpticTract_L, OpticTract_R, Pituitary, Eye_L, Eye_R) had an average of 274 training image sets and 92 testing image sets and were acquired from the Cancer Imaging Archive GLIS-RT dataset. These data sets consisted primarily of glioblastoma and astrocytoma patients. Images were acquired on either a GE Signa HDxT (3T) or Siemens Skyra (3T) scanner and had an average slice thickness of 1mm, In-plane resolution between 0.5-1.0 mm, and acquisition parameters of TR=2.3-8.9ms, TE=3.0-3.2s.

The MR training data used for initial testing of the MR Pelvis models (Prostate, GInd_Prostate, and SeminalVes) was taken from the Cancer Imaging Archive Prostate-MRI-US-Biopsy dataset which consisted patients with a of suspicion of prostate cancer due to elevated PSA and/or suspicious imaging findings. The images used were T2-Axial image sets acquired on a 3T Siemens Skyra scanner. The majority of pulse

sequences used are 3D T2:SPC, with TR/TE 2200/203, Matrix/FOV 256 × 205/14 × 14 cm, and 1.5mm slice spacing.

	Table 7: MR Initial Testing Dataset References						
Model Group	Data Source ID	Data Citation					
MR Brain	MR - Renown	Shusharina, N., & Bortfeld, T. (2021). <i>Glioma Image Segmentation for</i> <i>Radiotherapy: RT targets, barriers to cancer spread, and organs at risk</i> [Data set]. The Cancer Imaging Archive. <u>https://doi.org/10.7937/TCIA.T905-ZQ20</u>					
MR Pelvis	Gold Atlas Pelvis	Natarajan, S., Priester, A., Margolis, D., Huang, J., & Marks, L. (2020). Prostate MRI and Ultrasound With Pathology and Coordinates of Tracked Biopsy (Prostate-MRI-US-Biopsy) [Data set]. The Cancer Imaging Archive. DOI: <u>10.7937/TCIA.2020.A61IOC1A</u>					

Datasets used for testing were removed from the training dataset pool before model training began, and used exclusively for testing.

Ground truthing of each test data set was generated manually using consensus (NRG/RTOG) guidelines as appropriate by three clinically experienced experts consisting of 2 radiation therapy physicists and 1 radiation dosimetrist. For MR Structure models, a mean training DSC of 0.87+/-0.07 was found for medium models and 0.74+/-0.07 for small models.

Table 8: MR Trai	Table 8: MR Training Data Results for AutoContour Model RADAC V3										
MR Models	Size	Pass Criteria	DSC (Avg)	DSC Std Dev (Avg)	Lower Bound 95% Confidence Interval						
Cerebellum	Medium	0.65	0.93	0.01	0.91						
GInd_Prostate	Medium	0.65	0.87	0.04	0.80						
Hypothalamus	Small	0.50	0.79	0.04	0.72						
Hypo_True	Small	0.50	0.71	0.06	0.61						
OpticTract_L	Small	0.50	0.72	0.08	0.59						
OpticTract_R	Small	0.50	0.72	0.08	0.59						
Pituitary	Small	0.50	0.75	0.11	0.57						
Prostate	Medium	0.65	0.89	0.03	0.84						
SeminalVes	Medium	0.65	0.74	0.13	0.53						
Eye_L	Medium	0.65	0.90	0.10	0.74						
Eye_R	Medium	0.65	0.90	0.10	0.74						

Additional external clinical testing was performed in order to validate the accuracy of the models on image sets acquired that were unique to the training datasets.

	Table 9: MR External Clinical Dataset References							
Model Group	Data Source ID	Data Citation						
MR Brain	MR - Renown	N/A						
MR Pelvis	Gold Atlas Pelvis	Nyholm, Tufve, Stina Svensson, Sebastian Andersson, Joakim Jonsson, Maja Sohlin, Christian Gustafsson, Elisabeth Kjellén, et al. 2018. "MR and CT Data with Multi Observer Delineations of Organs in the Pelvic Area - Part of the Gold Atlas Project." <i>Medical Physics</i> 12 (10): 3218–21. doi:10.1002/mp.12748.						

For the Brain models, datasets acquired via data-use agreement from a clinical partner were acquired containing 20 MR T1 Ax post (BRAVO) image scans acquired with a GE MR750w scanner. Images had an average slice thickness of 1.6mm, In-plane resolution between 0.94 mm, and acquisition parameters of TR=5.98ms, TE=96.8s. Data for testing of the MR Pelvis structure models were acquired from a publicly available Gold Atlas Data set which contained 19 images of patients with prostate or rectal cancer. Images were acquired on a GE DISCOVERY MR750w (Ax T2 FRFSE) with in-plane resolution of 0.9mm, Slice thickness of 2.5mm, TR=5.988s and TE=96.8ms.

DSC values were calculated between ground truth contour data and AutoContour structures and rated on the same DSC passing criteria as was used for the training DSC validation. All structures passed the minimum DSC criteria for small and medium structures with a mean DSC of 0.74+/-0.07, 0.87+/-0.07 respectively: Additionally, the qualitative clinical appropriateness of AutoContour structures generated on these scans was graded by clinical experts. Autocontour structures were graded on a scale from 1 to 5 where 5 refers to contour requiring no additional edits, and 1 refers to a score in which full manual re-contour of the structure would be required. An average score >= 3 was used to determine whether a structure model would ultimately be beneficial clinically. An average rating of 4.4 was found across all MR structure models demonstrating that only minor edits would be required in order to make the structure models acceptable for clinical use.

Table 10:	Table 10: MR External Reviewer Results for AutoContour Model RADAC V3											
MR Models	Size	Pass Criteria	# External Test Data Sets	Average DSC	Average DSC Std. Dev	Lower Bound 95% Confidence Interval	External Reviewer Average Rating (1-5)					
Cerebellum	Medium	0.65	20	0.93	0.01	0.91	4					
GInd_Prostate	Medium	0.65	18	0.87	0.04	0.80	4.3					
Hypothalamus	Small	0.50	20	0.79	0.04	0.72	4.2					
Hypo_True	Small	0.50	19	0.71	0.06	0.61	4.3					
OpticTract_L	Small	0.50	20	0.72	0.08	0.59	4.4					
OpticTract_R	Small	0.50	19	0.72	0.08	0.59	4.5					
Pituitary	Small	0.50	19	0.75	0.11	0.57	4.2					
Prostate	Medium	0.65	19	0.89	0.03	0.84	4.2					
SeminalVes	Medium	0.65	19	0.74	0.13	0.53	4.2					
Eye_L	Medium	0.65	19	0.90	0.10	0.74	4.9					
Eye_R	Medium	0.65	20	0.90	0.10	0.74	4.9					

Validation testing of the AutoContour application demonstrated that the software meets user needs and intended uses of the application.

Mechanical and Acoustic Testing Not Applicable (Standalone Software)

Not Applicable (Standalone Software)

Animal Study

No animal studies were conducted using the Subject Device, AutoContour.

Clinical Studies

No clinical studies were conducted using the Subject Device, AutoContour

5.9. Conclusion

AutoContour Model RADAC V3 is deemed substantially equivalent to the Predicate Device, AutoContour Model RADAC V2 (K220598). Verification and Validation testing and the Risk Management Report demonstrate that AutoContour Model RADAC V3 is as safe and effective as the Predicate Device. The technological characteristics table demonstrates the similarity between AutoContour Model RADAC V3 and the Predicate Device and does not raise any questions on the safety and effectiveness of the Subject Device.