

Nanjing Jusha Display Technology Co., Ltd Donny Lee Certification Engineer 8A, Block 1. Nanjing International Service Outsourcing Mansion, No. 301 Hanzhongmen Street NANJING, JIANGSU 210036 CHINA

May 18, 2023

Re: K230723

Trade/Device Name: C510G/JUSHA-C510G/C510/JUSHA- C510 LCD Monitor Regulation Number: 21 CFR 892.2050 Regulation Name: Medical Image Management And Processing System Regulatory Class: Class II Product Code: PGY Dated: March 16, 2023 Received: March 16, 2023

Dear Donny Lee:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

essica damb

Jessica Lamb, Assistant Director Imaging Software Team DHT8B: Division of Radiological Imaging Devices and Electronic Products OHT8: Office of Radiological Health Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number *(if known)* K230723

Device Name C510G/JUSHA-C510G/C510/JUSHA-C510

Indications for Use (Describe)

C510G/JUSHA-C510G/C510/JUSHA-C510 LCD Monitor is intended to be used in displaying and viewing digital images, including standard and multi-frame digital mammography, for view, analysis, and diagnosis by trained medical practitioners. It is specially designed for breast tomosynthesis applications.

Type of Use (Select one or both, as applicable)	
Rescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov*

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

(K230723)

In accordance with 21 CFR 807.92 the following summary of information is provided:

Date:	March 16, 2023
Submitter:	Nanjing Jusha Display Technology Co., Ltd
	Add: 301, 8F Block A, No.1, Nanjing International Service Outsourcing Mansion, Hanzhongmen Street, Nanjing, 210036 China
Contact Person:	Donny Lee
	Certification Manager
	Nanjing Jusha Display Technology Co., Ltd
	Tel: +86-25- 83305050
	Fax: +86-25- 58783273
	lidongdong@jusha.com.cn
Device Trade Name:	C510G/JUSHA-C510G/C510/JUSHA-C510 LCD Monitor
Common/Usual Name:	5MP Color LCD Monitor
Classification Name:	Medical image management and processing system 21CFR 892.2050
Product Code:	PGY
Predicate Device(s):	JUSHA-M550G; K190848
Device Description:	C510G/JUSHA-C510G/JUSHA-C510/C510 LCD Monitor is the
1	display system with the high resolution (2560*2048), high luminance
	(1150 cd/m2), and 16-bit (65536 level), built-in DICOM standard
	LUT. In particular, C510G has ambient brightness adapting, real-time DICOM automatic calibration, full-screen brightness equalization and presence induction system, with these this display can automatic adjustment according to different requirements in order to achieve the best results.
	The product is consisted of the following components:
	- 21.3" Color-TFT LCD Panel
	- DMX0704AR1/Main Board/REV: 1.2
	- C510G LCD Monitor software

	- Power Adapter	
	- Data Cable.	
	- Data Cable.	
	The LCD Monitor is designed, tested, and will be manufactured in	
	accordance with both mandatory and voluntary standards:	
	1. IEC 60601-1:2012, EN 60601-1:2013, ANSI/AAMI	
	ES60601-1:2005+A1:2012+C1:2009+A2:2010,	
	CAN/CSA C22.2 NO.60601-1:14, Medical equipment	
	medical electrical equipment - Part 1: General	
	requirements for basic safety and essential performance.	
	2. IEC 60601-1-2 Edition 4:2014, EN 60601-1-2:2015,	
	Medical electrical equipment - Part 1-2: General	
	requirements for basic safety and essential performance -	
	Collateral standard: Electromagnetic disturbances -	
	Requirements and tests.	
Intended Use:	C510G/JUSHA-C510G/JUSHA-C510/C510 LCD Monitor is intended	
	to be used in displaying and viewing digital images, including	
	standard and multi-frame digital mammography, for view, analysis,	
	and diagnosis by trained medical practitioners. It is specially designed	
	for breast tomosynthesis applications.	
Technology:	C510G/JUSHA-C510G/JUSHA-C510/C510 LCD Monitor is the	
	display system with the high resolution (2560*2048), high luminance	
	(1150 cd/m2), and 16-bit grayscale (65536 level), built-in DICOM	
	standard LUT. In particular, C510G has ambient brightness adapting,	
	real-time DICOM automatic calibration, full-screen brightness	
	equalization and presence induction system, with these this display	
	can automatic adjustment according to different requirements in order	
	to achieve the best results.	
Determination of Substantial	Summary of Non-Clinical Tests:	
Equivalence:	The LCD Marites could be desired as following	
	The LCD Monitor complies with voluntary standards as following:	
	1 IEC 60601-1:2012, EN 60601-1:2013, ANSI/AAMI ES60601-	
	1:2005+A1:2012+C1:2009+A2:2010, CAN/CSA C22.2	
	NO.60601-1:14, Medical equipment medical electrical equipment	
	- Part 1: General requirements for basic safety and essential	
	performance.	
	2 IEC 60601-1-2 Edition 4:2014, EN 60601-1-2:2015, Medical	
	electrical equipment - Part 1-2: General requirements for basic	
	safety and essential performance - Collateral standard:	
	Electromagnetic disturbances - Requirements and tests	

	C510G is substantially equivalent to JUSHA-M550G. C510G employs the maximum resolution values same as that of JUSHA-M550G.
	The following quality assurance measures were applied to the development of the system:
	• Risk Analysis
	Requirements Reviews
	• Design Reviews
	• Raw materials verification
	• Testing on unit level (Module verification)
	• Integration testing (System verification)
	• Final acceptance testing (Validation)
	• Performance testing (Verification)
	• Safety testing (Verification)
	Summary of Clinical Tests:
	The subject of this premarket submission, LCD Monitor, did not
	require clinical studies to support substantial equivalence.
	The proposed device is Substantially Equivalent (SE) to the predicate
	device which is US legally market device. Therefore, the subject
	device is determined as safe and effectiveness.
Conclusion:	Nanjing Jusha Display Technology Co., Ltd Considers the C510G
	LCD Monitor to be as safe, as effective, and performance is
	substantially equivalent to the predicate device JUSHA-M550G.

Product Comparison

This comparison identifies the similarities and differences of the proposed C510G LCD Monitor device to the legally marketed predicate JUSHA-M550G LCD Monitor device to which substantial equivalency is claimed.

Attributes	Predicate Device	Proposed Device	Discussion of Differences
Product	JUSHA-M550G LCD Monitor	C510G LCD Monitor	
510(k) Number	K190848	K230723	
	Display Performa	nce/Specifications	
Screen	21.3inches, Mono-TFT LCD	21.3inches, Color-TFT LCD	Same
technology	Panel	Panel	
Viewing angle (H, V)	Horizontal 178°,	Horizontal 178°,	Same
	Vertical 178°	Vertical 178°	
Resolution	2560 x 2048/2048x 2560	2560 x2048/2048x 2560	Same
Display area	422.4(H) x 377.9 (V) mm	422.4(H) x337.9 (V) mm	The two monitors have different front frame style while the M550G employ a Slim Bezel design so that the display area is different.
Contrast Ratio	1700:1	2000:1	C510G has better contrast ratio for the panel, while the difference does not make difference to the diagnosis.
DICOM calibrated luminance	1000cd/m2	500cd/m2	The two models have different panels with different brightness but they have the same resolution for 5MP for Mammo image diagnosis.
Pixel Pitch	0.165x0.165 mm	0.165x0.165 mm	Same
Backlight	LED	LED	Same.

Attributes	Predicate Device	Proposed Device	Discussion of Differences
Product	JUSHA-M550G LCD Monitor	C510G LCD Monitor	
510(k) Number	K190848	K230723	
DICOM LUT	16-bit: 65536	16-bit: 65536	Same
Luminance	Built in calibration sensor	Built in calibration sensor	Same
calibration	provided	provided	
	Video S	Signal Input	
Input signals	DVI standard 1.0,	DVI standard 1.0,	Same
	DisplayPort 1.2a	DisplayPort 1.2a	
Input	DVI-D x 1,	DVI-D x 1,	Same
terminational	DisplayPort x 1	DisplayPort x 1	
Output signals	DisplayPort 1.2a	DisplayPort 1.2a	Same
Output	DisplayPort x 1	DisplayPort x 1	Same
Terminational			
Display controller	Off the shelf	Off the shelf	Same
	Power Relat	ed Specification	
Power Requirement	DC 24V	DC 24V	Same
Power Consumption/Sa ve Mode	55W/less than 0.5W	55W/less than 0.5W	Same
Power	DVI DMPM	DVI DMPM	Same
Management			
	DisplayPort 1.2a	DisplayPort 1.2a	
Miscellaneous Fea	atures/Specifications		
USB Ports/standard	1 upstream (endpoint),	1 upstream (endpoint),	Same
1 516/ Stundard	2 downstream/ Rev. 2.0	2 downstream/ Rev. 2.0	

Attributes	Predicate Device	Proposed Device	Discussion of Differences
Product	JUSHA-M550G LCD	C510G LCD Monitor	
	Monitor		
510(k) Number	K190848	K230723	
Dimensions w/o	Without stand:	Without stand:	Different
stand			Dimensions
	363mm x475mm x66mm	399mm x491mm x74mm	design due to the
(W x H x D)	With stand:	With stand:	different
			structure size.
	363mm x522mm x238mm	399mm x530mm x238mm	
Indication for	JUSHA-M550G/JUSHA-	C510G/JUSHA-	Same
use	M550/M550G/M550 LCD	C510G/C510/JUSHA-C510	
	Monitor is intended to be	LCD Monitor is intended to	
	used in displaying and	be used in displaying and	
	viewing digital images,	viewing digital images,	
	including standard and	including standard and	
	multi-frame digital	multi-frame digital	
	mammography, for view,	mammography, for view,	
	analysis, and diagnosis by	analysis, and diagnosis by	
	trained medical practitioners.	trained medical practitioners.	
	It is specially designed for	It is specially designed for	
	breast tomosynthesis	breast tomosynthesis	
	applications.	applications.	

CD C510G LCD Monitor K230723 K230723 2, EN 1 IEC 60601-1:2012, EN SI/AAMI 60601-1:2013, ANSI/AAMI
2, EN 1 IEC 60601-1:2012, EN Same
ES60601-C1:2009+1:2005+A1:2012+C1:2009+A C22.2A2:2010, CAN/CSA C22.2edicalNO.60601-1:14, Medicalelectricalequipment medical electricalequipmentPart 1: Generalsic safetyrequirements for basic safetyand essentialperformance.dition2 IEC 60601-1-2 Edition1-2:2015,4:2014, EN 60601-1-2:2015,wedipmentPart 1-2: Generalsic safetyrequirements for basic safetyand essential performance -Collateral standard:sturbancesElectromagnetic disturbancestests- Requirements and tests

PERFORMANCE DATA:

The following performance data were provided in support of the substantial equivalence determination.

Bench testing:

Bench testing was conducted to demonstrate the C510G meets all performance standards as follows:

- Measurement of the angular dependency of luminance response in horizontal, vertical and diagonal directions
- Measurement of the luminance non-uniformity characteristics of the display screen as specified in TGI18 guideline.
- Measurement of the chromaticity non-uniformity characteristics of the display screen as specified in TG18 guideline.
- Measurement of small-spot contrast ratio.
- Measurement of temporal response

• Performance data on luminance stability

Electrical safety and electromagnetic compatibility (EMC)

Electrical safety and EMC testing were conducted on the C510G The device complies with the IEC 60601-1 standard for safety and the IEC 60601-1-2 standard for EMC. Animal and clinical study

The subject of this premarket submission, C510G, does not require animal or clinical studies to support substantial equivalence.

CONCLUSIONS

C510G Medical Display is substantially equivalent to the predicate device with respect to technical characteristics, performance, application and intended use. The non-clinical data support the safety of the device. The device should perform as intended in the specified use conditions. Nanjing Jusha Display Technology Co., Ltd considers the C510G Medical Display does not raise any new issues of safety or effectiveness.