

April 11, 2023

Ray Co., Ltd. % Changhwan Lee RA Manager 1F~3F, 4F(Part), 5F, 265, Daeji-Ro, Suji-gu Yongin-si, Gyeonggi-do 16882 SOUTH KOREA

Re: K230753

Trade/Device Name: RCT800

Regulation Number: 21 CFR 892.1750

Regulation Name: Computed Tomography X-Ray System

Regulatory Class: Class II

Product Code: OAS Dated: March 15, 2023 Received: March 17, 2023

Dear Changhwan Lee:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part

801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Lu Jiang, Ph.D.

Assistant Director
Diagnostic X-Ray Systems Team
DHT8B: Division of Imaging Devices

and Electronic Products

OHT8: Office of Radiological Health Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

Form Approved: OMB No. 0910-0120 Expiration Date: 06/30/2023

Expiration Date: 06/30/2023 See PRA Statement below.

510(k) Number (if known) K230753	
Device Name RCT800	
Indications for Use (Describe) RCT800 is CBCT and panoramic x-ray imaging system with co. Which is intended to radiographic examination of the dento-ma adult and pediatric patients. And a model scan is included as a Cephalometric image is also includes wrist to obtain carpus impreatment. The device is to be operated and used by dentists or other legal	axillofacial, sinus, TMJ, Airway for diagnostic support for n option. lages for growth and maturity assessment for orthodontic
The device is to be operated and used by definists of other legal	ny quantited fleatur care professionalis.
Type of Use (Select one or both, as applicable)	
Prescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)
CONTINUE ON A SEPARA	ATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

K230753

1. 510(k) Summary

The 510(k) summary of safety and effectiveness information is being submitted in accordance with the requirements of 21 CFR Part 807.92.

2. Date

March 6, 2023

3. Administrative Information

Applicant Ray Co., Ltd.		Ray Co., Ltd.	
Address		1F~3F, 4F(Part), 5F, 265, Daeji-ro, Suji-gu, Yongin-si, 16882, Korea	
Name		Ray Co., Ltd.	
Manufacturer	Address	1F~3F, 4F(Part), 5F, 265, Daeji-ro, Suji-gu, Yongin-si, 16882, Korea	
	Tel	+82-31-605-1000	
	Fax	+82-2-6280-5534	
Contact Person	Name	Changhwan Lee	
Contact Person	Email	ch0406.lee@raymedical.co.kr	

4. Device Information

Trade/Proprietary Name		RCT800	
Common Name		Dental Panoramic/Tomography and Cephalometric X-ray System	
	Device	X-ray tomography, computed, dental	
Classification Name Regulation Number Class Product Code	_	21 CFR 892.1750	
	Class	2	
	Product Code	OAS	
	Review Panel	Radiology	

5. Predicate device

Parameter	Predicate Device-1	Predicate Device-2
Device Name	RCT800	RCT700
Manufacturer	Ray Co., Ltd	Ray Co., Ltd
510(K) Number	K192737	K213226
Classification name	Computed tomography x-ray system	Computed tomography x-ray system
Regulation number	892.1750	892.1750
Primary product code	OAS	OAS

6. Device Description

The system's purpose is RCT800, a 3D computed tomography scanner for scanning hard tissues like bones and teeth. By rotating the c-arm, which is embedded with an all-in-one x-ray tube and a detector on each end, CBCT images of the dental maxillofacial area can be attained by recombining data from the same level that is scanned from different angles. Additionally, the system includes a panoramic image scanning function for attaining images of the whole teeth, a cephalometric scanning option for attaining a cephalic image, and a Model Scan option for attaining a dental model CBCT image.

7. Indication for use

RCT800 is CBCT and panoramic x-ray imaging system with cephalometric. Which is intended to radiographic examination of the dento-maxillofacial, sinus, TMJ, Airway for diagnostic support for adult and pediatric patients. And a model scan is

Allway for diagnostic support for adult and pediatric patients. And a model scan is included as an option.

Cephalometric image is also includes wrist to obtain carpus images for growth and maturity assessment for orthodontic treatment.

The device is to be operated and used by dentists or other legally qualified heath care professionals.

8. Patient population

This device is intended to acquire diagnostic x-ray images of adult and pediatric individuals/patients without restriction on ethnic group, gender, weight, health status, or condition.

We recommend that patients who undergo X-ray diagnostic radiation exposure be over 5 years old.

9. Comparison with predicate device

The following table provides the summary of the technological characteristics of RCT800 compared to the predicate device.

Parameter	Proposed Device	Predicate Device-1	Predicate Device-2
Manufacturer	Ray Co., Ltd.	Ray Co., Ltd.	Ray Co., Ltd.
Device name	RCT800	RCT800	RCT700
510(K) Number	(Special 510K)	K192737 (Traditional 510K)	K213226 (Traditional 510K)
Common Name	Dental panoramic/tomography and cephalometric x-ray system	Dental panoramic/tomography and cephalometric x-ray system	Dental panoramic/tomography and cephalometric x-ray system
Indications for use	RCT800 is CBCT and panoramic x-ray imaging system with cephalometric. Which is intended to radiographic examination of the dento-maxillofacial, sinus, TMJ, Airway for diagnostic support for adult and pediatric patients. And a model scan is included as an option. Cephalometric image is also includes wrist to obtain carpus images for growth and maturity assessment for orthodontic treatment. The device is to be operated and used by dentists or other legally qualified health care professionals.	RCT800 is CBCT and panoramic x-ray imaging system with cephalometric. Which is intended to radiographic examination of the dento-maxillofacial, sinus, TMJ, Airway for diagnostic support for adult and pediatric patients. And a model scan is included as an option. Cephalometric image is also includes wrist to obtain carpus images for growth and maturity assessment for orthodontic treatment. The device is to be operated and used by dentists or other legally qualified health care professionals.	RCT700 is CBCT and panoramic x-ray imaging system with cephalometric. Which is intended to radiographic examination of the dento-maxillofacial, sinus, TMJ, Airway and ENT structure for diagnostic support for adult and pediatric patients. And a model scan is included as an option. Cephalometric image is also includes wrist to obtain carpus images for growth and maturity assessment for orthodontic treatment. The device is to be operated and used by dentists or other legally qualified health care professionals.
Mode of Operation	Same as predicate device #1	Continuous operation with intermittent, stated permissible loading	Same as predicate device #1
3D technology	Same as predicate device #1	CBCT Cone beam Computed Tomography	Same as predicate device #1
reconstruction algorithm	Same as predicate device #1	FBP(Filtered Back Projection)	Same as predicate device #1

Performance Specification		Same as predicate device #1	1) CBCT Computed tomography - Patient - Dental Model Scan(Optional) 2) Panoramic 3) Cephalometric(optional) - One shot type - Scan type	1) CBCT Computed tomography - Patient 2) Panoramic 3) Cephalometric(optional) - One shot type - Scan type
Functional Option		Same as predicate device #1	Base CT+PANO Option(CEPH) CT + PANO + SCAN CEPH CT + PANO + One shot(One shot, Standard Type) CT + PANO + One shot(One shot, Large Type).	Same as predicate device #1
	СТ	FXDD-1724RA	FXDD-0606CA	Same as predicate device #1
		Same as predicate device #1	FXDD-1012CHA	Jupi0606X
	PANO	FXDD-1724RA	FXDD-0606CA	Same as predicate device #1
Detector Type	PANO	Same as predicate device #1	FXDD-1012CHA	Jupi0606X
	Ceph (Scan)	Same as predicate device #1	XID-C24DC	Same as predicate device #1
	Conh (One shot)	Same as predicate device #2	1717SCC	FXRD-1717VA
	Ceph (One shot)	Same as predicate device #2	PaxScan 2530C	FXDD-1012CA
Exposure switch Type		Same as predicate device #1	"Deadman" Button type	Same as predicate device #1
Main Components		Same as predicate device #1	Ceph Apparatus	Same as predicate device #1

Same as predicate device #1	Vertical Carriage	Same as predicate device #1
Same as predicate device #1	Rotator	Same as predicate device #1
Same as predicate device #1	X-RAY Generator	Same as predicate device #1
Same as predicate device #1	X-ray tube	Same as predicate device #1
Same as predicate device #1	High Frequency Generator	Same as predicate device #1
Same as predicate device #1	Column	Same as predicate device #1
Same as predicate device #1	Touch monitor (panel)	Same as predicate device #1
Detector	Detector	Detector
-CT	- CT	-CT
FXDD-1724RA	FXDD-0606CA	FXDD-0606CA
FXDD-1012CHA	FXDD-1012CHA	Jupi0606X
-PANO	- PANO	-PANO
FXDD-1724RA	FXDD-0606CA	FXDD-0606CA
FXDD-1012CHA	FXDD-1012CHA	Jupi0606X
-Ceph	- Ceph	-Ceph
XID-C24DC(Scan)	XID-C24DC(Scan)	XID-C24DC(Scan)
, ,	` '	` ,
FXRD-1717VA(Oneshot, Large Size)	1717SCC(One shot, Large Size)	FXRD-1717VA(Oneshot, Large Size)
FXDD-1012CA(One shot, Standard	PaxScan 2530C(One shot, Standard	FXDD-1012CA(One shot, Standard
Size)	Size)	Size)
Same as predicate device #1	Chinrest	Chinrest
Same as predicate device #1	Head rest	Head rest
Same as predicate device #1	Automatic Collimator	Automatic Collimator
Same as predicate device #1	Exposure switch	Exposure switch
Same as predicate device #1	Emergency stop switch	Emergency stop switch
Same as predicate device #1	Console PC set	Console PC set
 -	-	

Automatic	Collimator	Same as predicate device #1	CT exams Panoramic exams Cephalometric exams	Same as predicate device #1
Display Ty	ре	Same as predicate device #1	TFT LCD type(Normally black) *1280x800 pixel	Same as predicate device #1
Class		Same as predicate device #1	Class I with type B applied parts according to IEC 60601-1	Same as predicate device #1
Focal size		Same as predicate device #2	Patient 0.5 Model scan 0.04 (Optional)	0.5
Field of Vie	ew(CT)	FXDD-1724RA : Max.180x160 mm FXDD-1012CHA : Max. 200x200 mm	FXDD-0606CA : Max.160x100 mm FXDD-1012CHA : Max. 200x200 mm	Max. 160x100 mm
X-ray Volta	age	Same as predicate device #1	60~100kVp	Same as predicate device #1
X-ray Curre	ent	Same as predicate device #2	4~17mA	1~17mA
Total Filtra	tion	Same as predicate device #1	Min. 2.8 mm Al equivalent	Same as predicate device #1
	СТ	FXDD-1724RA: 95μm	FXDD-0606CA: 119μm	Same as predicate device #1
	CI	Same as predicate device #1	FXDD-1012CHA: 124μm	Jupi0606X : 100μm
	PANO	FXDD-1724RA: 95μm	FXDD-0606CA: 119μm	Same as predicate device #1
Detector Pixel size	FANO	Same as predicate device #1	FXDD-1012CHA: 124μm	Jupi0606X : 100μm
	Ceph (Scan)	Same as predicate device #1	XID-C24DC: 100μm	Same as predicate device #1
	Ceph(One shot)	Same as predicate device #2	1717SCC: 127μm	FXRD-1717VA : 140μm
	Cepti(Offe shot)	Same as predicate device #2	PaxScan 2530C: 139μm	FXDD-1012CA : 124μm
	СТ	Same as predicate device #2	1.44(Patient) 1.91(Model Scan)	1.44
Magnifica	PANO	Same as predicate device #2	1.3	1.31
tion	Ceph (Scan)	Same as predicate device #1	1.11	Same as predicate device #1
	Ceph(One shot)	Same as predicate device #1	Large Type: 1.13	Same as predicate device #1

		Same as predicate device #1	Standard Type: 1.12	Same as predicate device #1
		Same as predicate device #1		Same as predicate device #1
		CT : below 20sec	CT : below 20sec(Patient) CT : below 180sec(Model Scan)	CT : below 14sec
Scan time		Same as predicate device #1	Pano : below 14sec	Same as predicate device #1
		Same as predicate device #2	Ceph[Scan size] : below 19sec	Ceph[Scan size] : below 20sec
		Same as predicate device #1	Ceph[One shot size]: below 2sec	Same as predicate device #1
Format comp	patible	Same as predicate device #1	DICOM 3.0 Format compatible	Same as predicate device #1
Image Viewin	ng Software	RayScan (Cleared under K192737)	RayScan (Cleared under K192737)	RayScan (Cleared under K182614)
Image acquis	sition	Same as predicate device #1	Giga-Ethernet Network	Same as predicate device #1
Total Height		Same as predicate device #1	Max 2,296mm	Same as predicate device #1
			1) Computed Tomography(CT) + Panoramic(PANO)=189kg(416.6lb) ± 10%	1) Computed Tomography(CT) + Panoramic(PANO)=185kg(407.9lb) ± 10%
			2) Computed Tomography(CT) + Panoramic(PANO) + Ceph (Scan type)= 219kg (482.8lb) ± 10%	2) Computed Tomography(CT) + Panoramic(PANO) + Ceph (Scan type)= 212.5kg (468.5lb) ± 10%
Weight		Same as predicate device #1	3) Computed Tomography(CT) + Panoramic(PANO) + Ceph (One shot type, installed in Standard size)= 217kg (478.4lb) ± 10%	3) Computed Tomography(CT) + Panoramic(PANO) + Ceph (One shot type, installed in Standard size)= 211kg (465.2lb) ± 10%
			4) Computed Tomography(CT) + Panoramic(PANO) + Ceph (One shot type, installed in Large size) 212kg (467.3lb) ± 10%	4) Computed Tomography(CT) + Panoramic(PANO) + Ceph (One shot type, installed in Large size) 211kg (465.2lb) ± 10%
Type of instal	llation	Same as predicate device #1	Wall or floor mount	Same as predicate device #1

Patient position	Same as predicate device #1	Standing / Wheelchair	Same as predicate device #1
Applicable Standards	Same as predicate device #1	IEC 60601-1 IEC 60601-1-3 IEC 60601-2-63 IEC 60601-1-2	Same as predicate device #1

The product is principally just the same as in the previous 510(k) #K192737 and #K213226.

The complete of differences of the subject device to the predicate device #K192737 is as follows

- The minimum X-ray current of the tube has been changed from 4mA to 1mA
- The Magnification of Panorama has been changed from 1.3 to 1.31
- Detector (using for CT, Pano, One-shot Ceph), Field of view.

 However, X-ray current, magnification and one-shot ceph detector was identified in #K213226.

The 510(k) for the existing detector used in our equipment is provided below.

Modality	Detector Model	Cleared	510(k) No.
2.7	FXDD-1724RA	No PMA	K222219
СТ	FXDD-1012CHA	No PMA	K182805
Pano	FXDD-1724RA	No PMA	K222219
	FXDD-1012CHA	No PMA	K182805
Scan Ceph	XID-C24DC	No PMA	K181452
One shot Ceph	FXRD-1717VA	No PMA	K213226
One shot Ceph	FXDD-1012CA	No PMA	K213226

10. Safety and Effectiveness Information

The RCT800 system described in this 510(k) is similar to the predicate device in terms of indications for use, materials, safety characteristics and X-ray source.

The following information further substantiates the substantial equivalence between the subject device and predicate device. The fundamental technological characteristics of the subject and predicate devices are similar. The imaging modes are similar; PANO, CEPH (Optional), CBCT. All viewing software programs have been cleared with previous 510k submissions; RAYSCAN (K192737). The sponsor tested the subject device in a laboratory and provided a non-clinical performance report. The same test protocol was used to test the performance of the subject and predicate devices for comparison. The sponsor certifies that adequate design and development controls (according to 21 CFR 820.30) were in place for manufacturing the subject device.

The differences are as follows.

- The minimum X-ray current of the tube has been changed from 4mA to 1mA
- The Magnification of Panorama has been changed from 1.3 to 1.31
- Detector (using for CT, Pano and One-shot Ceph), Field of view.

 However, X-ray current, magnification and one-shot ceph detector was identified in #K213226.

Electrical, mechanical and environmental safety testing according to the standard of IEC 60601-1:2005/AMD1:2012 (3.1 Edition), IEC 60601-1-3:2008/AMD1:2013 (Second Edition), IEC 60601-1-6:2010 (Third Edition), and IEC 60601-2-63:2012 (First Edition) were performed.

EMC testing was conducted in accordance with the standard IEC 60601-1-2:2014 (Edition 4.0).

The software of RCT800 saves patient and image data and offers an inquiry function. In addition, it supports the image generate function intended to obtain images using the RCT800 equipment and various sensors for diagnosis. That has been validated according to the FDA "Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices" and "Content of Premarket Submissions for Management of Cybersecurity in Medical Devices" to assure substantial equivalence. The software for this device was considered a "moderate" level of concern since a failure or latent flaw in the software would not directly result in serious injury or death to the patient or operator. As a result, we identified the level of concern associated with a new device and provided documentation consistent with that level. Based on our risk analysis of software, the difference does not affect its safety and effectiveness.

Bench testing was conducted according to the FDA Guidance "Format for Traditional and Abbreviated 510(k)s," Section 18, "Performance Testing – Bench." Bench testing is used to assess whether the parameters required to describe functionalities related to imaging properties of the dental X-ray device and patient dosage satisfy the designated tolerance.

Performance (Imaging performance) testing was conducted according to the standards of IEC 61223-3-4 and IEC 61223-3-7. All test results were satisfactory.

Non-clinical considerations were conducted in accordance with the FDA Guidance

"Guidance for the Submissions of 510(k)'s for Solid State X-ray Imaging Devices." Because the subject device uses the same detector as the predicate device, there are no significant differences between the two devices as a result of non-clinical testing.

Clinical considerations were conducted according to the FDA Guidance "Format for Traditional and Abbreviated 510(k)s," Section 20. Clinical images were provided, and while these images were not necessary to establish substantial equivalence based on the modifications to the device, they provide further evidence, in addition to the laboratory performance data, to show that the complete system works as intended.

The features of RCT800 were clinically tested and approved by two licensed practitioners/clinicians. Clinical imaging samples were collected from new detectors on the proposed device at the two offices where the predicate device was installed for the clinical test images. These images were gathered from all detectors installed with RCT800 using protocols with random patient age, gender, and size. A licensed practitioner reviewed the sample clinical images and deemed them to be of acceptable quality for the intended use.

11. Conclusions

In accordance with the Federal Food, Drug and Cosmetic Act, 21 CFR Part 807 and based on the information provided in this premarket notification. Ray Co., Ltd. concludes that the newly RCT800 is safe, effective and substantially equivalent to the predicate device as described herein.