

April 20, 2023

GE Healthcare Japan Corporation % He Haibo Regulatory Affairs Leader 7-127, 4-chome, Asahigaoka Hino, Tokyo 191-8503 JAPAN

Re: K230807

Trade/Device Name: Deep Learning Image Reconstruction

Regulation Number: 21 CFR 892.1750

Regulation Name: Computed tomography x-ray system

Regulatory Class: Class II

Product Code: JAK Dated: March 23, 2023 Received: March 23, 2023

#### Dear He Haibo:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <a href="https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm">https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm</a> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for

K230807 - He Haibo Page 2

devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <a href="https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products">https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products</a>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <a href="https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems">https://www.fda.gov/medical-device-problems</a>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<a href="https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance">https://www.fda.gov/training-and-continuing-education/cdrh-learn</a>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<a href="https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</a>) for more information or contact DICE by email (<a href="DICE@fda.hhs.gov">DICE@fda.hhs.gov</a>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Lu Jiang, Ph.D. Assistant Director

Diagnostic X-Ray Systems Team

Lu Jiang

DHT8B: Division of Radiological Imaging Devices

and Electronic Products

OHT8: Office of Radiological Health Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

# DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

# **Indications for Use**

510(k) Number (if known)

Form Approved: OMB No. 0910-0120
Expiration Date: 06/30/2023

Expiration Date: 06/30/2023 See PRA Statement below.

K230807			
Device Name Deep Learning Image Reconstruction			
dications for Use (Describe) The Deep Learning Image Reconstruction software is a deep learning based reconstruction method intended to produce ross-sectional images of the head and whole body by computer reconstruction of X-ray transmission data taken at ifferent angles and planes, including Axial, Helical (Volumetric), and Cardiac acquisitions, for all ages.			
Deep Learning Image Reconstruction software can be used for head, whole body, cardiac, and vascular CT applications.			
Type of Use (Select one or both, as applicable)			
Prescription Use (Part 21 CFR 801 Subpart D) Over-The-Counter Use (21 CFR 801 Subpart C)			
CONTINUE ON A SEPARATE PAGE IF NEEDED.			

This section applies only to requirements of the Paperwork Reduction Act of 1995.

#### \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



# K230807 510(k) SUMMARY OF SAFETY AND EFFECTIVNESS

This 510(k) summary of Safety and Effectiveness information is submitted in accordance with the requirement of 21 CFR Part 807.92:

**Date:** April 20, 2023

**Submitter:** GE Healthcare Japan Corporation

7-127, Asahigaoka, 4-chome

Hino-shi, Tokyo, 191-8503, Japan

Primary Contact: He Haibo

Regulatory Affairs Leader Phone: 86-010-5708-3413 Email: haibo.he1@ge.com

**Secondary Contacts:** Helen Peng

Senior Regulatory Affairs Director

Phone: 1-262-424-8222

Email: hong.peng@med.ge.com

Laura Turner

Regulatory Affairs Manager Phone: 1-262-200-1044 Email: laura.turner@ge.com

**Subject Device Name:** Deep Learning Image Reconstruction

**Device Classification** Class II

**Regulation Number/** 

**Product Code:** 

21 CFR 892.1750 Computed tomography x-ray system / JAK

#### **Predicate Device Information**

**Device Name:** Deep Learning Image Reconstruction

Manufacturer: GE Medical Systems, LLC

**510(k) Number:** K213999 cleared on February 18, 2022

**Regulation Number/** 

**Product Code:** 

21 CFR 892.1750 Computed tomography x-ray system / JAK

# 510(k) Premarket Notification Submission - DLIR



#### **Reference Devices Information**

**Device Name:** ASiR-V

Manufacturer: GE Medical Systems, LLC

**510(k) Number:** K133640 cleared on March 25, 2014

**Regulation Number/** 

21 CFR 892.1750 Computed tomography x-ray system / JAK

**Product Code:** 

**Device Name:** Deep Learning Image Reconstruction **Manufacturer:** GE Healthcare Japan Corporation

**510(k) Number:** K212067 cleared on September 17, 2021

**Regulation Number/** 

21 CFR 892.1750 Computed tomography x-ray system / JAK

**Product Code:** 

#### **Device Description**

Deep Learning Image Reconstruction is an image reconstruction method that uses a dedicated Deep Neural Network (DNN) that has been designed and trained specifically to generate CT Images to give an image appearance, as shown on axial NPS plots, similar to traditional FBP images while maintaining the performance of ASiR-V in the following areas: image noise (pixel standard deviation), low contrast detectability, high-contrast spatial resolution, and streak artifact suppression.

The images produced are branded as "TrueFidelity™ CT Images". Reconstruction times with Deep Learning Image Reconstruction support a normal throughput for routine CT.

The deep learning technology is integrated into the scanner's existing raw data-based image reconstruction chain to produce DICOM compatible "TrueFidelity<sup>TM</sup> CT Images".

The system allows user selection of three strengths of Deep Learning Image Recon: Low, Medium or High. The strength selection will vary with individual users' preferences and experience for the specific clinical need.

The DLIR algorithm was modified on the Revolution CT/Apex platform for improved reconstruction speed and image quality and cleared in February 2022 with K213999. The same modified DLIR is now being ported to Revolution EVO (K131576) /Revolution Maxima (K192686), Revolution Ascend (K203169, K213938) and Discovery CT750 HD family CT systems including Discovery CT750 HD, Revolution Frontier and Revolution Discovery CT (K120833).

Compared to the predicate device, the intended use and indications for use of Deep Learning Image Reconstruction are identical.

# 510(k) Premarket Notification Submission - DLIR



#### **Intended Use**

The Deep Learning Image Reconstruction software is intended for head, whole body, cardiac, and vascular CT scans.

#### **Indications for Use**

The Deep Learning Image Reconstruction software is a deep-learning based reconstruction method intended to produce cross-sectional images of the head and whole body by computer reconstruction of X-ray transmission data taken at different angles and planes, including Axial, Helical (Volumetric), and Cardiac acquisitions, for all ages.

Deep Learning Image Reconstruction software can be used for head, whole body, cardiac, and vascular CT applications.

#### **Comparisons**

The GEHC Deep Learning Image Reconstruction (DLIR) option for Revolution EVO/Revolution Maxima, Revolution Ascend and Discovery CT750 HD family CT systems is substantially equivalent to the predicate device K213999 for Revolution CT/Apex platform. The fundamental technology, i.e., the DLIR algorithm, remains unchanged from the predicate. The table below summarizes the substantive feature/technological similarities and differences between the predicate device and the proposed device.

Specification/ Attribute	Predicate Device  Deep Learning Image Reconstruction (K213999)	Proposed Device Deep Learning Image Reconstruction
Technology	DLIR uses a dedicated Deep Neural Network (DNN) which is trained on the CT scanner and therefore models the propagation of noise through the system to identify and remove the noise	Same
System statistics - Noise modeling of the data collection imaging chain (photon noise and electronic noise)	Characterization of the photon statistics as it propagates through the preprocessing and calibration imaging chain	Same
System statistics – Noise characteristics of the reconstructed images	DLIR uses a trained DNN which models the scanned object using information obtained from extensive phantom and clinical data to identify the noise characteristics and remove it	Same



Specification/ Attribute	Predicate Device  Deep Learning Image Reconstruction (K213999)	Proposed Device  Deep Learning Image  Reconstruction
Clinical Workflow	Select recon type and strength (Low, Medium, High).	Same
Reference protocols/dose	Using the same reference protocols provided on the <b>Revolution CT/Apex platform</b> systems for ASiR-V	Using the same reference protocols provided on the Revolution EVO, Revolution Maxima, Revolution Ascend and Discovery CT750 HD family systems for ASiR-V
Deployment Environment	On CT Console	On CT Console On GE's Edison Platform.

The subject device Deep Learning Image Reconstruction does not introduce any additional risks/hazards, warnings, or limitations.

## **Determination of Substantial Equivalence**

## **Summary of Non-Clinical Testing**

Deep Learning Image Reconstruction has successfully completed the design control testing per our quality system. No additional hazards were identified, and no unexpected test results were observed. Deep Learning Image Reconstruction was designed under the Quality System Regulations of 21CFR 820 and ISO 13485. GEHC believes that the extensive bench testing is sufficient for FDA's substantial equivalence determination.

The following quality assurance measures have been applied to the development of the system:

- Requirement Definition
- Risk Analysis and Control
- Technical Design Reviews
- Formal Design Reviews
- Software Development Lifecycle
  - o Code Review
  - Software Unit Implementation
  - o Software Integrations and Integration Testing
- System Testing
  - Safety Testing (Verification)
  - Image Performance Testing (Verification)
  - Simulating Use Testing (Validation)
- Software Release

# 510(k) Premarket Notification Submission - DLIR



The testing and results did not raise different questions of safety and effectiveness than associated with predicate device. We consider the proposed device is substantially equivalent to the predicate device, DLIR for Revolution CT/Apex platform systems.

The substantial equivalence is also based on the software documentation for a "Moderate" level of concern.

Engineering bench testing was also performed to support substantial equivalence and the product performance claims. The evaluation and analysis used the same test methodologies and acceptance criteria with the identical raw datasets obtained on GEHC's Revolution EVO, Revolution Ascend, Revolution Frontier and Discovery CT750 HD systems and then applying the Deep Learning Image Reconstruction or ASiR-V reconstruction (hence the dose (CTDIvol) is identical for both). The resultant images were then compared for:

- Low Contrast Detectability (LCD)
- Image Noise (pixel standard deviation)
- High-Contrast Spatial Resolution (MTF)
- Streak Artifact Suppression
- Spatial Resolution, longitudinal (FWHM slice sensitivity profile)
- Noise Power Spectrum (NPS) and Standard Deviation of noise
- CT Number Uniformity
- CT Number Accuracy
- Contrast to Noise (CNR) ratio
- Artifact analysis metal objects, unintended motion, truncation
- Pediatric Phantom IQ Performance Evaluation
- Low Dose Lung Cancer Screening Protocol IQ Performance Evaluation

#### **Substantial Equivalence**

Deep Learning Image Reconstruction for Revolution EVO/Revolution Maxima, Revolution Ascend and Discovery CT750 HD family was developed under GEHC Healthcare's quality system. The changes associated with Deep Learning Image Reconstruction (DLIR) do not change the Intended Use from the predicate, and represent equivalent technology characteristics, with no impact on control mechanism, operating principle, and energy type. Design verification, along with bench testing provided in this submission demonstrates that the Deep Learning Image Reconstruction (DLIR) met all of its design requirement and performance criteria. GEHC's quality system's design, verification, and risk management processes did not identify any additional hazards, unexpected results, or adverse effects stemming from the changes to the predicate.

GE HealthCare believes that Deep Learning Image Reconstruction software is substantially equivalent to the legally marketed predicate device and hence is safe and effective for its intended use.