

June 1, 2023

Candela Corporation Danielle Gibboney Sr. Regulatory Affairs Specialist 251 Locke Drive Marlborough, Massachusetts 01752

Re: K230990

Trade/Device Name: Candela Vbeam Family of Pulsed Dye Lasers (Vbeam Prima, Vbeam Perfecta)
Regulation Number: 21 CFR 878.4810
Regulation Name: Laser Surgical Instrument For Use In General And Plastic Surgery And In Dermatology
Regulatory Class: Class II
Product Code: GEX
Dated: April 5, 2023
Received: April 6, 2023

Dear Danielle Gibboney:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,



Jianting Wang Acting Assistant Director DHT4A: Division of General Surgery Devices OHT4: Office of Surgical and Infection Control Devices Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K230990

Device Name

Candela Vbeam Family of Pulsed Dye Lasers (Vbeam Prima, Vbeam Perfecta)

Indications for Use (Describe)

Vbeam Prima is indicated for the following:

595 nm

• General Surgery:

• Photocoagulation of benign cutaneous vascular lesions and benign cutaneous lesions.

• Dermatology/Plastic Surgery:

o For treatment of benign cutaneous vascular lesions, such as facial and leg telangiectasia, rosacea, port wine stains, hemangiomas, angioma, spider angioma, Poikiloderma of Civatte, and benign cutaneous lesions, such as warts, scars, striae and Psoriasis and the treatment of wrinkles.

- Treatment of Benign Epidermal Pigmented Lesions.
- Treatment of Inflammatory Acne Vulgaris.
- Gynecology:

o Photocoagulation of benign cutaneous lesions and benign vascular lesions in gynecology.

• Podiatry:

- o Treatment of benign cutaneous lesions, such as warts.
- Pediatric Population

o Treatment of cutaneous capillary malformations, also known as port wine stains (PWS), and infantile hemangiomas (IH) / congenital hemangiomas

1064 nm

The Vbeam Prima laser system is intended for the coagulation and hemostasis of benign vascular lesions such as, but not limited to, port wine stains, hemangiomas, warts, telangiectasia, rosacea, Venus lakes, leg veins, spider veins, and poikiloderma of Civatte and treatment of benign cutaneous lesions such as, but not limited to lentigos (age spots), solar lentigos (sun spots), café-au-lait macules, seborrheic keratoses, nevi, chloasma, verrucae, skin tags, and keratoses. The laser is also indicated for the treatment of wrinkles such as, but not limited to, peri-ocular and peri-oral wrinkles.

Type of Use	(Select one	or both,	as applicable)	
-------------	-------------	----------	----------------	--

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

510(k) Number (if known)

K230990

Device Name

Candela Vbeam Family of Pulsed Dye Lasers (Vbeam Prima, Vbeam Perfecta)

Indications for Use (Describe)

Vbeam Perfecta is indicated for the following:

595 nm

• General Surgery:

- Photocoagulation of benign cutaneous vascular lesions and benign cutaneous lesions.
- Dermatology/Plastic Surgery:

o For treatment of benign cutaneous vascular lesions, such as facial and leg telangiectasia, rosacea, port wine stains, hemangiomas, angioma, spider angioma, Poikiloderma of Civatte, and benign cutaneous lesions, such as warts, scars, striae and Psoriasis and the treatment of wrinkles.

- Treatment of Benign Epidermal Pigmented Lesions.
- Treatment of Inflammatory Acne Vulgaris.
- Gynecology:

o Photocoagulation of benign cutaneous lesions and benign vascular lesions in gynecology.

- Podiatry:
- o Treatment of benign cutaneous lesions, such as warts.
- Pediatric Population

Treatment of cutaneous capillary malformations, also known as port wine stains (PWS), and infantile hemangiomas (IH) / congenital hemangiomas

Type of Use (Select one or both, as applicable)	
Prescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number." This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov*

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



510(k) SUMMARY K230990 Candela Vbeam Family of Pulsed Dye Lasers (Vbeam Prima, Vbeam Perfecta)

This summary of 510(k) submitted in accordance with the requirements of 21 CFR 807.92.

1. DATE PREPARED

April 5, 2023

2. APPLICANT NAME

Candela Corporation 251 Locke Drive Marlborough MA 01752 USA

3. OFFICIAL CORRESPONDENT

Danielle Gibboney Sr. Regulatory Affairs Specialist Candela Corporation 251 Locke Drive Marlborough MA 01752 USA Phone: 617-904-3820 Email: danielleg@candelamedical.com

4. PRODUCT INFORMATION

Name of Device: Candela Vbeam Family of Pulsed Dye Lasers (Vbeam Prima, Vbeam Perfecta)

Common/Usual Name: Powered Laser Surgical Instrument

Classification Name: Laser surgical instrument for use in general and plastic surgery and in dermatology (per 21 CFR Part 878.4810)

Device Classification: Class II (per 21 CFR Part 878.4810)

Product Code: GEX

5. LEGALLY MARKETED PREDICATE DEVICE FOR CLAIMED EQUIVALENCE:

Predicate Device: Vbeam Prima Laser System (K183452) Predicate Device: Candela Family of Pulse Laser Systems (K050673)

6. DEVICE DESCRIPTION:

The Candela Vbeam Family of Pulsed Dye Lasers includes both the Vbeam Prima and Vbeam Perfecta Laser Systems.

The Vbeam Prima Laser System has been previously cleared for both port wine stains, and hemangiomas under K183452.

The Vbeam Perfecta (Candela Family of Pulse Laser Systems) has been previously cleared for both port wine stains, and hemangiomas under K050673.

This 510(k) Premarket Notification is to expand the indications for use for the Candela Vbeam Family of Pulsed Dye Lasers for the Vbeam Prima Laser System and Vbeam Perfecta to include **the pediatric population for treatment of cutaneous capillary malformations, also known as port wine stains (PWS), and infantile hemangiomas (IH) / congenital hemangiomas** for the 595 nm wavelength. There is no new technology being introduced than what has been previously cleared Vbeam Prima Laser System under its predicate K183452 and Vbeam Perfecta under its predicate K050673.

7. INTENDED USE AND INDICATIONS FOR USE:

Vbeam Prima:

Vbeam Prima is indicated for the following: <u>595 nm</u>

- General Surgery:
- Photocoagulation of benign cutaneous vascular lesions and benign cutaneous lesions.
- Dermatology/Plastic Surgery:
 - For treatment of benign cutaneous vascular lesions, such as facial and leg telangiectasia, rosacea, port wine stains, hemangiomas, angioma, spider angioma, Poikiloderma of Civatte, and benign cutaneous lesions, such as warts, scars, striae and Psoriasis and the treatment of wrinkles.
- Treatment of Benign Epidermal Pigmented Lesions.
- Treatment of Inflammatory Acne Vulgaris.
- Gynecology:
 - Photocoagulation of benign cutaneous lesions and benign vascular lesions in gynecology.
- Podiatry:
 - Treatment of benign cutaneous lesions, such as warts.
- Pediatric Population
 - Treatment of cutaneous capillary malformations, also known as port wine stains (PWS), and infantile hemangiomas (IH) / congenital hemangiomas

<u>1064 nm</u>

The Vbeam Prima laser system is intended for the coagulation and hemostasis of benign vascular lesions such as, but not limited to, port wine stains, hemangiomas, warts, telangiectasia, rosacea, Venus lakes, leg veins, spider veins, and poikiloderma of Civatte and treatment of benign cutaneous lesions such as, but not limited to lentigos (age spots), solar lentigos (sun spots), café-au-lait macules, seborrheic keratoses, nevi, chloasma, verrucae, skin tags, and keratoses. The laser is also indicated for the treatment of wrinkles such as, but not limited to, peri-ocular and peri-oral wrinkles.

Vbeam Perfecta:

Vbeam Perfecta is indicated for the following:

<u>595 nm</u>

- General Surgery:
- Photocoagulation of benign cutaneous vascular lesions and benign cutaneous lesions.
- Dermatology/Plastic Surgery:
 - For treatment of benign cutaneous vascular lesions, such as facial and leg telangiectasia, rosacea, port wine stains, hemangiomas, angioma, spider angioma, Poikiloderma of Civatte, and benign cutaneous lesions, such as warts, scars, striae and Psoriasis and the treatment of wrinkles.
- Treatment of Benign Epidermal Pigmented Lesions.
- Treatment of Inflammatory Acne Vulgaris.
- Gynecology:
 - Photocoagulation of benign cutaneous lesions and benign vascular lesions in gynecology.
- Podiatry:
 - Treatment of benign cutaneous lesions, such as warts.
- Pediatric Population
 - Treatment of cutaneous capillary malformations, also known as port wine stains (PWS), and infantile hemangiomas (IH) / congenital hemangiomas

8. TECHNOLOGICAL COMPARISON:

The subject device Candela Vbeam Family of Pulsed Dye Lasers (Vbeam Prima, Vbeam Perfecta) is substantially equivalent and identical in the design, function, and intended use to the Vbeam Prima Laser System under K183452 and Vbeam Perfecta under K050673. The difference between the subject Candela Vbeam Family of Pulsed Dye Lasers (Vbeam Prima, Vbeam Perfecta) and its predicates is the additional indications for **Pediatric Population: Treatment of cutaneous** capillary malformations, also known as port wine stains (PWS), and infantile hemangiomas (IH) / congenital hemangiomas that this Premarket Notification is proposing. The expanded indications between the subject device and its predicates does not raise any new concerns of safety or effectiveness of the device. Thus, based on the information presented in this Premarket Notification, Candela Vbeam Family of Pulsed Dye Lasers (Vbeam Prima, Vbeam Perfecta) is substantially equivalent to its predicates Vbeam Prima Laser System under K183452 and Vbeam Perfecta under K050673. Please refer to specification comparison tables in Table 1 and Table 2 for comparisons between intended use/indications for use, and technological & biological characteristic comparison below.

Name of Device: 510(k) Product Code Device Class	Candela Vbeam Family of Pulsed Dye Lasers (Vbeam Prima) Proposed <u>Not Assigned</u> <u>GEX</u>	Vbeam Prima Laser System Predicate <u>K183452</u> <u>GEX</u> <u>Class II</u>	Vbeam Perfecta Laser System Predicate <u>K050673</u> <u>GEX</u> <u>Class II</u>
Intended use / Indications:	Class IIVbeam Prima is indicated for the following:595 nm: General Surgery: 	 Vbeam Prima is indicated for the following: <u>595 nm:</u> General Surgery: Photocoagulation of benign cutaneous vascular lesions and benign cutaneous lesions. Dermatology/Plastic Surgery: For treatment of benign cutaneous vascular lesions, such as facial and leg telangiectasia, rosacea, port wine stains, hemangiomas, angioma, spider angioma, Poikiloderma of Civatte, and benign cutaneous lesions, such as warts, scars, striae and Psoriasis and the treatment of wrinkles. 	 595 nm: General Surgery: Photocoagulation of benign cutaneous vascular lesions and benign cutaneous lesions. Dermatology/Plastic Surgery: For treatment of benign cutaneous vascular lesions, such as facial and leg telangiectasia, rosacea, port wine stains, hemangiomas, angioma, spider angioma, Poikiloderma of Civatte, and benign cutaneous lesions, such as warts, scars, striae and Psoriasis and the treatment of wrinkles. Treatment of Benign Epidermal Pigmented Lesions. Treatment of Inflammatory Acne Vulgaris.

Table 1: Intended/Indication for use comparison table.

Name of Device: 510(k) Product Code Device Class	Candela Vbeam Family of Pulsed Dye Lasers (Vbeam Prima) Proposed <u>Not Assigned</u> <u>GEX</u> <u>Class II</u>	Vbeam Prima Laser System Predicate <u>K183452</u> <u>GEX</u> <u>Class II</u>	Vbeam Perfecta Laser System Predicate <u>K050673</u> <u>GEX</u> <u>Class II</u>
	 Treatment of Benign Epidermal Pigmented Lesions. Treatment of Inflammatory Acne Vulgaris. 	 Treatment of Benign Epidermal Pigmented Lesions. Treatment of Inflammatory Acne Vulgaris. 	Gynecology: Photocoagulation of benign cutaneous lesions and benign vascular lesions in gynecology.
	Gynecology: ○ Photocoagulation of benign cutaneous lesions and benign vascular lesions in gynecology.	Gynecology: • Photocoagulation of benign cutaneous lesions and benign vascular lesions in gynecology.	 Treatment of benign cutaneous lesions, such as warts.
	Podiatry: o Treatment of benign cutaneous lesions, such as warts.	Podiatry: o Treatment of benign cutaneous lesions, such	
	 Treatment of Treatment of cutaneous capillary malformations, also known as port wine stains (PWS), and infantile hemangiomas (IH) / congenital hemangiomas 1064 nm The Vbeam Prima laser system is intended for the coagulation 	as warts. <u>1064 nm</u> The Vbeam Prima laser system is intended for the coagulation and hemostasis of benign vascular lesions such as, but not limited to, port wine stains, hemangiomas, warts, telangiectasia, rosacea, Venus lakes, leg veins, spider veins, and poikiloderma of Civatte and treatment of benign cutaneous	

Name of Device: 510(k) Product Code Device Class	Candela Vbeam Family of Pulsed Dye Lasers (Vbeam Prima) Proposed <u>Not Assigned</u> <u>GEX</u> <u>Class II</u>	Vbeam Prima Laser System Predicate <u>K183452</u> <u>GEX</u> <u>Class II</u>	Vbeam Perfecta Laser System Predicate <u>K050673</u> <u>GEX</u> <u>Class II</u>
	and hemostasis of benign vascular lesions such as, but not limited to, port wine stains, hemangiomas, warts, telangiectasia, rosacea, Venus lakes, leg veins, spider veins, and poikiloderma of Civatte and treatment of benign cutaneous lesions such as, but not limited to lentigos (age spots), solar lentigos (sun spots), café-au-lait macules, seborrheic keratoses, nevi, chloasma, verrucae, skin tags, and keratoses. The laser is also indicated for the treatment of wrinkles such as, but not limited to, peri-ocular and peri- oral wrinkles.	lesions such as, but not limited to lentigos (age spots), solar lentigos (sun spots), café-au-lait macules, seborrheic keratoses, nevi, chloasma, verrucae, skin tags, and keratoses. The laser is also indicated for the treatment of wrinkles such as, but not limited to, peri-ocular and peri-oral wrinkles.	
	Vbeam Perfecta is indicated for the following: <u>595 nm:</u> General Surgery:		

Name of Device: 510(k) Product Code Device Class	Candela Vbeam Family of Pulsed Dye Lasers (Vbeam Prima) Proposed <u>Not Assigned</u> <u>GEX</u> <u>Class II</u>	Vbeam Prima Laser System Predicate <u>K183452</u> <u>GEX</u> <u>Class II</u>	Vbeam Perfecta Laser System Predicate <u>K050673</u> <u>GEX</u> <u>Class II</u>
	 lesions. Dermatology/Plastic Surgery: For treatment of benign cutaneous vascular lesions, such as facial and leg telangiectasia, rosacea, port wine stains, hemangiomas, angioma, spider angioma, Poikiloderma of Civatte, and benign cutaneous lesions, such as warts, scars, striae and Psoriasis and the treatment of wrinkles. Treatment of Benign Epidermal Pigmented Lesions. Treatment of Inflammatory Acne Vulgaris. Gynecology: Photocoagulation of benign cutaneous lesions and benign vascular lesions in gynecology. 		

Name of Device: 510(k) Product Code Device Class	Candela Vbeam Family of Pulsed Dye Lasers (Vbeam Prima) Proposed <u>Not Assigned</u> <u>GEX</u> <u>Class II</u>	Vbeam Prima Laser System Predicate <u>K183452</u> <u>GEX</u> <u>Class II</u>	Vbeam Perfecta Laser System Predicate <u>K050673</u> <u>GEX</u> <u>Class II</u>
	 Podiatry: Treatment of benign cutaneous lesions, such as warts. Pediatric Population: Treatment of cutaneous capillary malformations, also known as port wine stains (PWS), and infantile hemangiomas (IH) / congenital 		
Similarities/Differences	Identical to Predicate device, but with expanded indications identified in BOLD.	Identical to the subject device, but without the expanded indications.	Identical to the subject device, but without the expanded indications.

General Specifications	Candela Vbeam Family of Pulsed Dye Lasers (Vbeam Prima, Vbeam Perfecta) Proposed <u>Not Assigned</u>		Vbeam Prima Laser System Predicate <u>K183452</u>		Vbeam Perfecta Predicate <u>K050673</u>
Technical Characteristics	Identical		Identical		Identical
Wavelength	595 nm IDENTICAL	1064 nm (Vbeam Prima only) IDENTICAL	595 nm IDENTICAL	1064 nm IDENTICAL	595 nm IDENTICAL
Laser Type	Flash lamp-excited, pulsed dye laser IDENTICAL	Flashlamp-excited, Nd:YAG laser (Vbeam Prima Only) IDENTICAL	Flash lamp- excited, pulsed dye laser IDENTICAL	Flashlamp-excited, Nd:YAG laser IDENTICAL	Flash lamp-excited, pulsed dye laser IDENTICAL
Pulse Energy	Up to 12 J (Vbeam Prima only) IDENTICAL Up to 8 J (Vbeam Perfecta only) IDENTICAL	Up to 45 J (Vbeam Prima only) IDENTICAL	Up to 12 J IDENTICAL	Up to 45 J IDENTICAL	Up to 8 J IDENTICAL
Pulse Width	0.45-40 ms IDENTICAL	Up to 60 ms (Vbeam Prima Only) IDENTICAL	0.45-40 ms IDENTICAL	Up to 60 ms IDENTICAL	0.45-40 ms IDENTICAL
Laser Pulse Repetition Rate	Up to 1.5 Hz IDENTICAL	Up to 10 Hz (Vbeam Prima only) IDENTICAL	Up to 1.5 Hz IDENTICAL	Up to 10 Hz IDENTICAL	Up to 1.5 Hz IDENTICAL

Table 2: Technological & Biological specification comparison

Page **9** of **16**

General Specifications	Candela Vbeam Family of Pulsed Dye Lasers (Vbeam Prima, Vbeam Perfecta) Proposed <u>Not Assigned</u>	Vbeam Prima Laser System Predicate <u>K183452</u>	Vbeam Perfecta Predicate <u>K050673</u>
Spot Size (mm)	DCD (Vbeam Prima only): 3-15mm HP 3x10mm HP 1.5 Zoom HP <u>CC (Vbeam Prima only):</u> 3-15mm Zoom HP IDENTICAL <u>Vbeam Perfecta Only:</u> 3, 5, 7, 10, 12 millimeters and 3x10 (elliptical)	DCD: 3-15mm HP 3x10mm Zoom HP 1.5 ZHP <u>CC:</u> 3-15mm Zoom HP IDENTICAL	3, 5, 7, 10, 12 millimeters and 3x10 (elliptical) IDENTICAL
Electrical Power	200-240 VAC~; 4600 VA; 50/60 Hz ; single phase (Vbeam Prima only) IDENTICAL 20 - 230 V~ , 50/60 Hz, single phase, 4,000 VA or 17.4 A at 230 V~ (Vbeam Perfecta only)	200-240 VAC~; 4600 VA; 50/60 Hz ; single phase IDENTICAL	220 - 230 V~ , 50/60 Hz, single phase, 4,000 VA or 17.4 A at 230 V~ IDENTICAL
Physical Dimensions /Weight (Console)	 280 lbs. (Vbeam Prima only) IDENTICAL 290 lbs. (Vbeam Perfecta only) IDENTICAL 	280 lbs. IDENTICAL	290 lbs. IDENTICAL
Patient Contacting Material	Distance Gauges: ULTEM 1000-1000	Distance Gauges: ULTEM 1000- 1000 IDENTICAL	Distance Gauges: ULTEM 1000-1000 IDENTICAL

9. PERFORMANCE DATA:

Performance Testing: Bench:

The performance testing of the subject Candela Vbeam of Pulsed Dye Lasers (Vbeam Prima, Vbeam Perfecta) is based on the established testing previous cleared under Vbeam Prima Laser System under its predicate K183452 and Vbeam Perfecta under its predicate K050673. There are no changes in the design therefore the subject Candela Vbeam of Pulsed Dye Lasers is based on the established performance testing of the device's predicates.

Performance Testing-Clinical

A systematic literature search using PubMed, Embase and Cochrane databases, and a supplemental search in Clinicaltrials.gov, was conducted to identify peerreviewed articles in which the Vbeam PDL System with the 595-nm wavelength was used to treat cutaneous capillary malformations, also known as port wine stains (PWS), and infantile hemangiomas (IH) / congenital hemangiomas in the pediatric population. A total of 33 articles were identified that reported on randomized controlled, prospective, open label, evaluator-blinded clinical trials, or retrospective evaluator-blinded studies that treated at least 10 individuals in each study using the Vbeam 595-nm PDL. The studies were conducted globally in the USA, UK, Europe, and Asia, which included China, Japan, Taiwan, India, and the Middle East. The 33 studies included a total of 7,725 patients of which 5,692 were pediatric (neonates, infants, children, and adolescents) Fitzpatrick Skin Type (FST) I-VI patients treated for hemangioma (4,782 IH, 910 congenital) and 1,354 were pediatric and 679 were pediatric and adult patients (up to 79 years) with Fitzpatrick Skin Type (FST) I- VI treated for PWS / cutaneous capillary malformation.

In conclusion, the data from these published reports support the intended indications for use of the Vbeam 595-nm PDL to treat cutaneous capillary malformations, also known as port wine stains (PWS), and infantile hemangiomas (IH) / congenital hemangiomas in the pediatric population. The articles are identified in Table 3 below.

<u>No.</u>	<u>Citation</u>	Indication(s):
1	Zhang W, Li F, Yang Y, Xue L, Cao M, Wang L. Hemangioma treatment with pulsed dye laser- distinct parameters used between neonatal and non-neonatal patients. J Cosmet Laser Ther. 2016 Nov;18(7):389- 392. doi: 10.1080/14764172.2016.1197402.	Infantile Hemangiomas (pediatric)
2	Yang B, Li L, Zhang LX, Sun YJ, Ma L. Clinical Characteristics and Treatment Options of Infantile Vascular Anomalies. Medicine (Baltimore). 2015 Oct;94(40):e1717. doi:	Vascular Anomalies in children including Infantile Hemangioma & Port Wine Stain

Table 3. Identification of Clinical Articles in Candela Vbeam Pulse Dyed Lasers literature search analysis

<u>No.</u>	Citation	Indication(s):
	10.1097/MD.0000000000001717.	
3	Sadeghinia A, Moghaddas S, Tavakolpour S,	Port Wine Stain (pediatric)
	Teimourpour A, Danespazhooh M, Mahmoudi H.	
	I reatment of port wine stains with 595-nm pulsed	
	dye laser in 27 pediatric patients: A prospective	
	study in the Iranian population of Cosmet Laser Ther	
	2019;21(7-8):373-377. doi:	
	10.1080/14764172.2019.1661489.	
4	Therapy of Pulsed Dye Laser With Intense Pulsed	Port Wine Stain (pediatric and
	Light in Port-Wine Stain Treatment: A Prospective	adult)
	Side-by-Side Comparison. Dermatol Surg. 2021	
	10.1097/DSS.000000000003114.	
5	Fallahi M, Hallaji Z, Tavakolpour S, Niknam S,	Port Wine Stain (pediatric and
	Salehi Farid A, Nili A, Teimourpour A,	adult)
	Daneshpazhooh M, Rahmati J, Haddady Abianeh	
	S, Mahmoudi H. Evaluating the efficacy and safety	
	of topical sirolimus 0.2% cream as adjuvant	
	therapy with pulsed dye laser for the treatment of	
	port wine stain: A randomized, double-blind,	
	placebo-controlled trial.	
	doi: 10.1111/jocd.13867.	
6	Yu W, Ma G, Qiu Y, Chen H, Jin Y, Yang X, Chang	Port Wine Stain (pediatric)
	L, Wang T, Hu X, Li W, Lin X. Prospective	
	comparison treatment of 595-nm pulsed-dye	
	lasers for virgin port-wine stain. Br J Dermatol.	
	2015 Mar;172(3):684- 91. doi: 10.1111/bjd.13356.	
7	Bernstein EF. High-energy 595 nm pulsed dye laser	Port Wine Stain (pediatric and
	port-wine stains. Dermatol Surg. 2006	adult)
	Jan;32(1):26-33. doi: 10.1111/1524-	
8	4/25.2006.32003.	Dout Mine Otain (nediatric and
0	Yu W, Wang T, Zhu J, Qiu Y, Chen H, Jin Y, Yang	Port wine Stain (pediatric and
	ream does not influence officacy and pain	adult)
	reduction during pulsed due laser treatment of	
	nort-wine stain: a	
	prospective side-by-side comparison.	
	Lasers Med Sci. 2018 Apr;33(3):573-579.	
9	001. 10. 1007/S10103-017-2415-3.	Dort Wing Stain (nadiatria and
	Wang T Chang L Chen V Ma G Lin Y Shorter	adult)
	Treatment Intervals of Fast Asians with Port-Wine	addity
	Stain with Pulsed Dve Laser Are Safe and	
	Effective-A Prospective Side-by-Side	
	Comparison, Photomed Laser Surg. 2018	
	Jan;36(1):37-43. doi: 10.1089/pho.2017.4315.	

<u>No.</u>	Citation	Indication(s):
10	Zhu J, Yu W, Wang T, Chen Y, Lyu D, Chang L, Ma	Port Wine Stain (pediatric)
	G, Lin X. Less is more: similar efficacy in three	
	sessions and seven sessions of pulsed dye laser	
	treatment in infantile port-wine stain patients. Lasers	
	2018 Nov;33(8):1707-1715. doi: 10.1007/s10103- 018-2525-6.	
11	Brauer JA, Farhadian JA, Bernstein LJ, Bae YS, Geronemus RG. Pulsed Dye Laser at Subpurpuric Settings for the Treatment of Pulsed Dye Laser- Induced Ecchymoses in Patients With Port-Wine Stains. Dermatol Surg. 2018 Feb;44(2):220-226. doi: 10.1097/DSS.000000000001255.	Port Wine Stain (pediatric and adult)
12	Yu W, Zhu J, Wang L, Qiu Y, Chen Y, Yang X,	Virgin Port Wine Stain
	Chang L, Ma G, Lin X. Double Pass 595 nm Pulsed	(pediatric and adult)
	Dye Laser Does Not Enhance the Efficacy of Port	
	Wine Stains Compared with Single Pass: A	
	Comparison with Histological Examination. Photomed Laser Surg. 2018 Jun;36(6):305-312. doi: 10.1089/pho.2017.4392.	
13	Yu W, Ying H, Chen Y, Qiu Y, Chen H, Jin Y, Yang	Port Wine Stain (pediatric and
	X, Wang T, Ma G, Lin X. In Vivo Investigation of the	adult)
	Safety and Efficacy of Pulsed Dye Laser with Two	
	Spot Sizes in Port-Wine Stain Treatment: A	
	Prospective Side-by-Side Comparison. Photomed Laser Surg. 2017 Sep;35(9):465-	
	471. doi: 10.1089/pho.2016.4186.	
14	Faurschou A, Togsverd-Bo K, Zachariae C,	Port Wine Stain (pediatric and
	Haedersdal M. Pulsed dye laser vs. intense pulsed	adult)
	light for port-wine stains: a randomized side-by- side	
	trial with blinded response evaluation. Br J	
	Dermatol. 2009	
	Feb;160(2):359-64. doi: 10.1111/J.1365-	
15	Z133.2008.08993.X. Yung A. Sheehan-Dare R. A comparative study of a	Defractory Dart Mina Stain
10	595-nm with a 585-	(podiatric and adult)
	nm pulsed dye laser in refractory port wine stains.	(pediatric and addit)
	Br J Dermatol. 2005 Sep;153(3):601-6. doi:	
10	10.1111/j.1365-2133.2005.06707.x.	
10	Yu W, Ma G, Qiu Y, Chen H, Jin Y, Yang X, Hu X,	Port Wine Stain (pediatric and
	Wang T, Chang L, Zhou H, Li W, Lin X. Why do	adult)
	port-wine stains (PWS) on the lateral face respond	
	better to pulsed dye laser (PDL) than those	
	central face2 I Am Acad Dermatol 2016	
	Mar:74(3):527-35 doi: 10 1016/i ipad 2015 08 026	
17	Lil Kono T Groff WF Chan HH Kitazawa V	Port Wine Stain (nediatric and
-	Nozaki M Comparison study of a long-pulse	adult)
	pulsed dve laser and a long-pulse	
	pulsed alexandrite laser in the treatment of port	

<u>No.</u>	<u>Citation</u>	Indication(s):
	wine stains. J Cosmet Laser Ther. 2008 Mar;10(1):12-5. doi: 10.1080/14764170701817023.	
18	Tomson N, Lim SP, Abdullah A, Lanigan SW. The treatment of port- wine stains with the pulsed-dye laser at 2-week and 6-week intervals: a comparative study. Br J Dermatol. 2006 Apr;154(4):676-9. doi: 10.1111/j.1365-2133.2005.07113.x.	Port Wine Stain (pediatric and adult)
19	Asahina A, Watanabe T, Kishi A, Hattori N, Shirai A, Kagami S, Watanabe R, Le Pavoux A, Maekawa T, Tamaki K, Ohara K. Evaluation of the treatment of port-wine stains with the 595-nm long pulsed dye laser: a large prospective study in adult Japanese patients. J Am Acad Dermatol. 2006 Mar;54(3):487-93. doi: 10.1016/j.jaad.2005.11.1034.	Port Wine Stain (pediatric and adult)
20	Liu S, Yang C, Yang S, Wang Z, Luo D, Zhang X. Topical application of 5-aminolevulinic acid followed by 595-nm pulsed dye laser irradiation for the treatment of recalcitrant port-wine stains: a primary study. J Cosmet Laser Ther. 2012 Aug;14(4):189-92. doi: 10.3109/14764172.2012.699677.	Recalcitrant Port Wine Stain (pediatric and adult)
21	Woo WK, Jasim ZF, Handley JM. Evaluating the efficacy of treatment of resistant port-wine stains with variable-pulse 595-nm pulsed dye and 532-nm Nd:YAG lasers. Dermatol Surg. 2004 Feb;30(2 Pt 1):158-62; discussion 162. doi: 10.1046/j.1076-0512.2003.30055.x.	Port Wine Stain (pediatric and adult)
22	Peters MA, van Drooge AM, Wolkerstorfer A, van Gemert MJ, van der Veen JP, Bos JD, Beek JF. Double pass 595 nm pulsed dye laser at a 6 minute interval for the treatment of port-wine stains is not more effective than single pass. Lasers Surg Med. 2012 Mar;44(3):199-204. doi: 10.1002/lsm.22011.	Port Wine Stain (pediatric and adult)
23	Tremaine AM, Armstrong J, Huang YC, Elkeeb L, Ortiz A, Harris R, Choi B, Kelly KM. Enhanced port- wine stain lightening achieved with combined treatment of selective photothermolysis and imiquimod. J Am Acad Dermatol. 2012 Apr;66(4):634-41. doi: 10.1016/j.jaad.2011.11.958.	Port Wine Stain (pediatric and adult)
24	Frohm Nilsson M, Passian S, Wiegleb Edstrom D. Comparison of two dye lasers in the treatment of port-wine stains. Clin Exp Dermatol. 2010 Mar;35(2):126-30. doi: 10.1111/j.1365-	Port Wine Stain (pediatric and adult)

<u>No.</u>	Citation	Indication(s):
	2230.2009.03399.x.	
25	Thajudheen CP, Jyothy K, Priyadarshini A. Treatment of port-wine stains with flash lamp	Port Wine Stain (pediatric and adult)
	pumped pulsed dye laser on Indian skin: a six year study. J Cutan Aesthet Surg. 2014	
26	Jan;7(1):32-6. doi: 10.4103/0974-2077.129973.	
20	exploration of optimal time and safety of 595-nm	(pediatric)
	pulsed dye laser for the treatment of early	
	2022	
27	May;35(5):e15406. doi: 10.1111/dth.15406.	Homongiama (padiatria)
	M, Yamaki T, Soejima K, Nozaki M. Comparison	riemangionia (pediatric)
	study of a traditional pulsed dye laser versus a	
	childhood hemangiomas. Lasers Surg Med. 2006	
28	Feb;38(2):112-5. doi: 10.1002/lsm.20257.	Infantila Llamangiama
20	Haenssle HA, Zutt M, von Fintel H, Kühnle I, Schön	(pediatric)
	MP, Thoms KM. Nd:YAG and pulsed dye laser	
	analysis of 271 treated haemangiomas: a retrospective	
	children. J Eur Acad Dermatol Venereol. 2017	
29	Aug;31(8):1372-1379. doi: 10.1111/jdv.14074.	Vascular Anomalies including
	Zhang ZY, Chang SC, Huo YP. The application of	Infantile Hemangioma
	595-nm pulsed dye laser for vascular anomalies in	(pediatric) & Port Wine Stain
	Cosmet Laser Ther. 2019;21(3):171-178. doi:	
20	10.1080/14764172.2018.1502450.	
50	pass laser therapy for	Port Wine Stain (pediatric)
	recalcitrant facial port wine stains. Lasers Med Sci. 2022 Apr:37(3):1643-1650. doi: 10.1007/s10103-	
	021-03414-x.	
31	Li D, Chen B, Zhang H, Yuan Y, Fan W, Ying Z. Retrospective study of the treatment of port-wine	Port Wine Stain (pediatric and adult)
	stains with 595-nm pulsed dye laser in 261	
	Chinese patients. Lasers Med Sci. 2020	
	03013-2.	
32	Jeon H, Bernstein LJ, Belkin DA, Ghalili S,	Port Wine Stain (pediatric)
	Port-Wine Stains in Infancy Without the Need for	
	General Anesthesia. JAMA Dermatol. 2019 Apr	

<u>No.</u>	Citation	Indication(s):
	1;155(4):435-441. doi: 10.1001/jamadermatol.2018.5249.	
33	Hartmann F, Lockmann A, Himpel O, Kühnle I, Hensen J, Schön MP, Thoms KM. Combination therapy of oral propranolol and combined Nd:YAG/pulsed dye laser therapy in infantile hemangiomas: a retrospective analysis of 48 treated hemangiomas in 30 children. J Dtsch Dermatol Ges. 2020 Sep;18(9):984-993. doi: 10.1111/ddg.14184.	Infantile Hemangioma (pediatric)

10. SUBSTANTIAL EQUIVALENCE COMPARISON

When comparing the subject Candela Vbeam of Pulsed Dye Lasers (Vbeam Prima, Vbeam Perfecta) is identical to the previously cleared Vbeam Prima Laser System under K183452 and Vbeam Perfecta under K050673. The additional indications for use do not raise any new issues of safety and effectiveness. There are no technological changes between the subject and predicate devices. The subject Candela Vbeam of Pulsed Dye Lasers (Vbeam Prima, Vbeam Perfecta) is substantially equivalent, in terms of technological characteristics, performance, and intended use to the predicate devices Vbeam Prima Laser System under K183452 and Vbeam Perfecta under K050673 as they are identical.