

May 19, 2023

Xenta Biomedical Science Co., Ltd.
Huang Ling
Product Manager
Building C5, 9th Floor, Rm 901, No. 11 Kaiyuan Avenue
Huangpu District
Guangzhou, Guangdong 510535
China

Re: K231137

Trade/Device Name: Xenta Drug Screen Cup, Xenta Drug Screen Dipcard

Regulation Number: 21 CFR 862.3100 Regulation Name: Amphetamine test system

Regulatory Class: Class II

Product Code: DKZ, DIS, JXM, DJR, DJG, LCM

Dated: April 21, 2023 Received: April 21, 2023

### Dear Huang Ling:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <a href="https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm">https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm</a> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801 and Part 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <a href="https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products">https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products</a>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <a href="https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems">https://www.fda.gov/medical-device-problems</a>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<a href="https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</a>) and CDRH Learn (<a href="https://www.fda.gov/training-and-continuing-education/cdrh-learn">https://www.fda.gov/training-and-continuing-education/cdrh-learn</a>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<a href="https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</a>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Paula V. Caposino -S

Paula Caposino, Ph.D.
Acting Deputy Director
Division of Chemistry
and Toxicology Devices
OHT7: Office of In Vitro Diagnostics
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

# DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

### Indications for Use

Form Approved: OMB No. 0910-0120

Expiration Date: 06/30/2023 See PRA Statement below.

510(k) Number (if known) K231137		<u>'</u>	
Device Name Xenta Drug Screen Cup Xenta Drug Screen Dipcard			
Indications for Use (Describe)  Xenta Drug Screen Cup and Xenta	-		
qualitatively detect the presence of o	drugs and drug metabolites	in human urine at the following c	cut-off concentrations:
Test	Calibrator		Cut-off level
Barbiturates (BAR)	Secobarbital		300 ng/mL
Benzodiazepines (BZO)	Oxazepam		300 ng/mL
Amphetamine (AMP)	D-Amphetam	ine	1000 ng/mL
Methadone (MTD)	Methadone		300 ng/mL
Oxycodone (OXY)	Oxycodone		100 ng/mL
Phencyclidine (PCP)	Phencyclidine		25 ng/mL
The tests contain two formats: 1) Te multiple drug tests in any combinati diagnostics use. They are intended f	on of the drug analytes list	•	2 2
The assays provide only a prelimina preferred confirmatory method. Clir test result, particularly when preliminations of the present the stress of the stres	nical consideration and pro-	fessional judgment should be appl	
Type of Use (Select one or both, as app	licable)		
Prescription Use (Par	t 21 CFR 801 Subpart D)	Over-The-Counter Use (21 C	FR 801 Subpart C)
	CONTINUE ON A SEPAR	ATE PAGE IF NEEDED.	

# This section applies only to requirements of the Paperwork Reduction Act of 1995. \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

## Section 5 - 510(k) Summary

Date of Summary Preparation: 5/17/2023

510(k) Number: K231137

1. Submitter's Identifications

Submitter: Xenta Biomedical Science Co., Ltd.

Address: Building C5, 9th Floor, Rm 901, No. 11 Kaiyuan Avenue, Huangpu District, Guangzhou,

510535, P.R. China

Contact Person: Huang Ling

Contact Email Address: 1433969171@qq.com

Telephone: +86-20-31707187

Fax: +86-20-31707187

2. Correspondent's Identifications

Correspondent's Name: Xenta Biomedical Science Co., Ltd.

Address: Building C5, 9th Floor, Rm 901, No. 11 Kaiyuan Avenue, Huangpu District, Guangzhou,

510535, P.R. China

Contact Person: Huang Ling

Contact Email Address: 1433969171@qq.com

Telephone: 86-20-31707187

Fax: 86-20-31707187

#### 3. Name of the Device

Proprietary names:

Xenta Drug Screen Cup Xenta Drug Screen Dipcard

Recommended classification regulation:

21 CFR 862.3150 Barbiturate test system 21 CFR 862.3170 Benzodiazepine test system 21 CFR 862.3100 Amphetamine test system 21 CFR 862.3620 Methadone test system 21 CFR 862.3650 Opiate test system

Unclassified Enzyme Immunoassay Phencyclidine

Device class: Class II Panel: Toxicology

Product code: DIS,JXM,DKZ,DJR,DJG,LCM

#### 4. The Predicate Devices

K122809 Advin Multi-Drug Screen Test Dip Card

Advin Multi-Drug Screen Test Cup Advin Multi-Drug Screen Test Cassette

### 5. Device Description

Xenta Drug Screen Cup and Xenta Drug Screen Dipcard are competitive binding, lateral flow immunochromatographic assays for the qualitative detection of Barbiturate, Benzodiazepine, Amphetamine, Methadone, Oxycodone, Phencyclidine at or above the cut-off levels as indicated. The tests are performed without the use of an instrument.

The test cup and test dipcard formats use identical test strips made with the same chemical formulation and manufacturing procedures.

#### 6. Indications for Use

Xenta Drug Screen Cup and Xenta Drug Screen Dipcard are lateral flow chromatographic immunoassays designed to qualitatively detect the presence of drugs and drug metabolites in human urine at the following cut-off concentrations:

Test	Calibrator	Cut-off level
Barbiturates (BAR)	Secobarbital	300 ng/mL
Benzodiazepines (BZO)	Oxazepam	300 ng/mL
Amphetamine (AMP)	D-Amphetamine	1000 ng/mL
Methadone (MTD)	Methadone	300 ng/mL
Oxycodone (OXY)	Oxycodone	100 ng/mL
Phencyclidine (PCP)	Phencyclidine	25 ng/mL

The tests contain two formats:1) Test Cup and 2) Test Dipcard. The tests may be configured as single drug tests or multiple drug tests in any combination of the drug analytes listed in the table above. The tests are intended for in vitro diagnostics use. They are intended for prescription use.

The assays provide only a preliminary analytical test result. Gas Chromatography/Mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are indicated.

### 7. Comparison to Predicate Devices:

A summary comparison of features of the Xenta Drug Screen Cup and Xenta Drug Screen Dipcard and the predicate devices is provided in the following Table:

Item	Device	Predicate (K122809)
Indication for use	Qualitative detection of	Same (but the number of

	drugs-of-abuse in urine	drugs detected different)		
Intended Users	Prescription Use	Over the Counter (OTC) Use		
		and Prescription Use		
Specimen	Urine	Same		
Cutoff	Barbiturates:300 ng/mL	Barbiturate:300 ng/mL		
	Benzodiazepines:300 ng/mL	Benzodiazepine:300 ng/mL		
	Amphetamine:1000 ng/mL	Amphetamine:500 ng/mL		
	Methadone:300 ng/mL	Methadone:300 ng/mL		
	Oxycodone:100 ng/mL	Oxycodone:100 ng/mL		
	Phencyclidine: 25 ng/mL	Phencyclidine: 25 ng/mL		
Read time	5 minutes	Same		
Results	Qualitative	Same		
Methodology	Competitive binding, Lateral flow	Same		
	immunochromatographic assay based			
	on the principle of antigen antibody			
	immunochemistry			
Configuration	Dipcard and Cup	Cassette,Dip Card and Cup		

### 8. Performance Data:

### **8.1** Cross-reactivity with structurally similar compounds

To test the cross reactivity of the test, test Dipcard and test Cup was used to test with drug metabolites and drug structurally similar compounds in urine. All the components were added to drug-free normal human urine. Each sample was tested in 5 replicates using Test Cup and Test Dipcard. If any positive result was observed, the compounds were further diluted with known drug-free urine specimen sequentially to different concentrations and tested in quintuplicate, until the highest concentration that generates a negative result was identified.

The cross reacting substances with the lowest concentration that produced a positive result was identified and is listed in the table below.

Amphetamine (AMP)	Lowest Concentrat ion (ng/mL)	% Cross-react ivity	Methadone (MTD)	Lowest Concentrat ion (ng/mL)	% Cross-react ivity
d-Amphetamine	1000	100%	Methadone	300	100%
L-Amphetamine	50000	2%	(±)2-Ethy1-1,5-dimethy1-3, 3-diphenylpyrrolinium	50000	0.6%
MDA	2000	50%	Doxylamine	50000	0.6%
Phentermine	45000	2.2%	Oxycodone (OXY)		
Benzodiazepines (BZO)			Oxycodone	100	100%
Oxazepam	300	100%	Hydrocodone	3000	3.3%
Alpha-hydroxyalprazolam	1900	15.8%	Hydromorphone	75000	0.1%

Alprazolam	200	150%	Oxymorphone	1000	10%
Bromazepam	1000	30%	Ethymorphine	50000	0.2%
Clobazam	100	300%	Codeine	50000	0.2%
Clonazepam	800	37.5%	Barbiturates (BAR)		
Desalkylflurazepam	400	75%	Secobarbital	300	100%
Diazepam	200	150%	Amobarbital	500	60%
Flunitrazepam	350	85.7%	Aprobarbital	200	150%
Lorazepam	1200	25%	Butabarbital	75	400%
Lorazepam-glucuronide	5000	6%	Butalbital	1500	20%
Nitrazepam	100	300%	Cyclopentobarbital	600	50%
Nordiazepam	400	75%	Pentobarbital	700	42.9%
Norchlordiazepapoxide	500	60%	Phenobarbital	300	100%
Nordiazepoxide	400	75%	Phencyclidine (PCP)		
Temazepam	120	250%	Phencyclidine	25	100%
Triazolam	1000	30%	4-Hydroxy-PCP	15000	0.2%

#### 8.2 Interference

Clinical urine samples may contain substances that could potentially interfere with the test. The following compounds were added to drug-free urine or drug positive urine containing AMP, BAR, BZO, MTD, OXY, PCP with the concentration 50% below the cutoff and the concentration 50% above the cutoff, respectively. All potential interfering substances were added at a concentration of  $100\mu g/mL$  (All concentrations of the drugs were confirmed with GC/MS). The urine specimens were tested with two lots of the corresponding Single/Multi-drug Test Cup and Test Dipcard. None of the compounds listed below were shown to interfere.

Acetaminophen	Estrone-3-sulfate	d,l-Octopamine
Acetophenetidin	Ethyl-p-aminobenzoate	Oxalic acid
Amoxicillin	Erythromycin	Oxolinic acid
Ampicillin	Fenoprofen	Oxymetazoline
Aspirin	Flucloxacillin	Oxytetracycline
Atenolol	Fluoxetine	Papaverine
Atorvastatin	Furosemide	Penicillin-G
Azlocillin	Gentisic acid	Pentazocine
Benzilic acid	Hemoglobin	Perphenazine
Benzylpenicillin	Hydralazine	Phenelzine
Benzoic acid	Hydrochlorothiazide	Prednisolone
Bilirubin	Hydrocortisone	Prednisone
Benzydamine	o-Hydroxyhippuric acid	d,l-Propanolol
Caffeine	p-Hydroxytyramine	d-Pseudoephedrine
Carbamazepine	Ibuprofen	Quinacrine
Cephalexin	Indomethacin	Quinine
Chloralhydrate	Iproniazid	Quindine
Chloramphenicol	d,l-Isoproterenol	Ranitidine
Chlorothiazide	Isoxsuprine	Salicylic acid
Chlorpheniramine	Ketamine	Serotonin
d,l-Chlorpromazine	Ketoprofen	Sulfamethazine
Cholesterol	Labetalol	Sulindac
Clonidine	Lisinopril	Tetracycline

Cimetidine Loperamide Tetrahydrozoline

Citalopram Meperidine Thiamine Thioridazine Cortisone Meprobamate Creatinine Methoxyphenamine d, 1-Thyroxine Deoxycorticosterone Methylphenidate Tolbutamine Dexamethasone Nadolol Tolbutamide Dextromethorphan Nalidixic acid Trifluoperazine Diclofenac Naproxen Tryptamine Diflunisal Niacinamide Uric acid Digoxin Nicotine Verapamil Diphenhydramine Nifedipine Zomepirac

Ephedrine Norethindrone

β-Estradiol Noscapine

Acetone Acetylsalicylic acid Albumin
Ascorbic Acid Aspartame Ascorbic Acid
Atropine Benzocaine Benzoylecgonine

Chlorquine (±) Chlolrpheniramine Creatine

Dexbrompheniramine Dophenhydramine

Dophenhydramine (+/-)-Isoproterenol

1R,2S(+)-EphedrineEthanolGlucoseGuaiacol glyceryl etherLevorphanolLidocaineLysergic acidMethadoneMethanol

Methaqualone Morphine (1R,2S)-(-)-n-Methyl-ephe

drine

Nortriptyline

(+)-Naproxen

Nordiazepam (+/-)-Norephedrine

Pheniramine Phenothiazine L-Phenylephrine

B-Phenylethylamin Phencyclidine Procaine Propoxyphene Ranitidine Riboflavin

Salicylic acid Secobarbital Sodium Chloride

Theophyline Tyramine Uric acid

Vitamin(L-Ascorbic Acid) 4-Dimethylaminoantipyrine d-Amphetamine

### 8.3 Effect of urinary pH

The pH of an aliquot negative urine pool is adjusted to a pH range of 3 to 9 in 1 pH unit increments and spiked with each drug at 50% below and 50% above cutoff levels (All concentrations were confirmed with GC/MS). Each sample was tested by two lots of the corresponding Xenta Drug Screen Cup and Xenta Drug Screen Dipcard. The results demonstrate that varying ranges of pH do not interfere with the performance of the test.

#### 8.4 Effect of Urinary specific gravity

The specific gravity studies were conducted on different specific gravity including 1.002,1.010, 1.020, 1.030, 1.040 specimens with drug free urine containing AMP, BAR, BZO, MTD, OXY, PCP at 50% below and 50% above cutoff level (All concentrations were confirmed with GC/MS).

Each sample was tested by two lots of the corresponding Xenta Drug Screen Cup and Xenta Drug Screen Dipcard. The results demonstrate that varying ranges of urinary specific gravity do not affect the test result.

#### 8.5 Precision

Precision studies were performed using the single drug and multi-drug test formats. Drug free specimens were spiked with analytes at 0,  $\pm 75\%$  cutoff,  $\pm 50\%$  cutoff,  $\pm 25\%$  cutoff and  $\pm 100\%$  cutoff of drug. The concentrations of the target drugs were confirmed with GC/MS. In both the single drug test and multi-drug test precision studies each concentration of the urine specimen was divided into aliquots. Each aliquot was blindly labeled by a nonparticipant. Separate sets of blinded coded samples were assigned and randomized prior to testing. The study was conducted by 6 operators at 3 Point-of-Care sites. Two operators per location tested 3 aliquots at each concentration for each lot per day (3 runs/day) for 10 non-consecutive days using one device lot per location. One operator tested the test dipcard format and the second operator tested the test cup format. There were 1620 observations by 3 sites at 9 concentrations.

#### Single drug Test Cup:

_	Approximate		Number of			Res	ult			
Drug	concentration	% of cutoff	determinatio	L	Lot 1		Lot 2		Lot 3	
test	of sample		ns per lot	Positive	Negative	Positive	Negative	Positive	Negative	
	0ng/ml	Negative	60	0	60	0	60	0	60	
	250ng/ml	-75%cutoff	60	0	60	0	60	0	60	
	500ng/ml	-50%cutoff	60	0	60	0	60	0	60	
	750ng/ml	-25%cutoff	60	6	54	8	52	8	52	
AMP	1000ng/ml	cutoff	60	34	26	36	24	32	28	
	1250ng/ml	+25%cutoff	60	52	8	52	8	54	6	
	1500ng/ml	+50%cutoff	60	60	0	60	0	60	0	
	1750ng/ml	+75%cutoff	60	60	0	60	0	60	0	
	2000ng/ml	+100%cutoff	60	60	0	60	0	60	0	
	0ng/ml	Negative	60	0	60	0	60	0	60	
	75ng/ml	-75%cutoff	60	0	60	0	60	0	60	
	150ng/ml	-50%cutoff	60	0	60	0	60	0	60	
	225ng/ml	-25%cutoff	60	4	56	6	54	4	56	
BAR	300ng/ml	cutoff	60	38	22	38	22	36	24	
	375ng/ml	+25%cutoff	60	54	6	56	4	58	2	
	450ng/ml	+50%cutoff	60	60	0	60	0	60	0	
	525ng/ml	+75%cutoff	60	60	0	60	0	60	0	
	600ng/ml	+100%cutoff	60	60	0	60	0	60	0	
	0ng/ml	Negative	60	0	60	0	60	0	60	
BZO	75ng/ml	-75%cutoff	60	0	60	0	60	0	60	
	150ng/ml	-50%cutoff	60	0	60	0	60	0	60	

	225ng/ml	-25%cutoff	60	6	54	4	56	6	54
	300ng/ml	cutoff	60	38	22	36	24	34	26
	375ng/ml	+25%cutoff	60	56	4	54	6	52	8
	450ng/ml	+50%cutoff	60	60	0	60	0	60	0
	525ng/ml	+75%cutoff	60	60	0	60	0	60	0
	600ng/ml	+100%cutoff	60	60	0	60	0	60	0
	0ng/ml	Negative	60	0	60	0	60	0	60
	75ng/ml	-75%cutoff	60	0	60	0	60	0	60
	150ng/ml	-50%cutoff	60	0	60	0	60	0	60
	225ng/ml	-25%cutoff	60	6	54	4	56	4	56
MTD	300ng/ml	cutoff	60	34	26	38	22	36	24
	375ng/ml	+25%cutoff	60	54	6	58	2	56	4
	450ng/ml	+50%cutoff	60	60	0	60	0	60	0
	525ng/ml	+75%cutoff	60	60	0	60	0	60	0
	600ng/ml	+100%cutoff	60	60	0	60	0	60	0
	0ng/ml	Negative	60	0	60	0	60	0	60
	25ng/ml	-75%cutoff	60	0	60	0	60	0	60
	50ng/ml	-50%cutoff	60	0	60	0	60	0	60
	75ng/ml	-25%cutoff	60	6	54	4	56	8	52
OXY	100ng/ml	cutoff	60	36	24	38	22	34	26
	125ng/ml	+25%cutoff	60	56	4	56	4	58	2
	150ng/ml	+50%cutoff	60	60	0	60	0	60	0
	175ng/ml	+75%cutoff	60	60	0	60	0	60	0
	200ng/ml	+100%cutoff	60	60	0	60	0	60	0
	0ng/ml	Negative	60	0	60	0	60	0	60
	6.3ng/ml	-75%cutoff	60	0	60	0	60	0	60
	12.5ng/ml	-50%cutoff	60	0	60	0	60	0	60
	18.8ng/ml	-25%cutoff	60	4	56	4	56	6	54
PCP	25ng/ml	cutoff	60	36	24	38	22	34	26
	31.3ng/ml	+25%cutoff	60	54	6	56	4	56	4
	37.5ng/ml	+50%cutoff	60	60	0	60	0	60	0
	43.8ng/ml	+75%cutoff	60	60	0	60	0	60	0
	50ng/ml	+100%cutoff	60	60	0	60	0	60	0

# Multi-drug Test Cup:

Drug	Approximate		Number of			Re	sult		
test	concentration	% of cutoff	determinations	Lo	ot 1	Lo	ot 2	Lo	ot 3
	of sample		per lot	Positive	Negative	Positive	Negative	Positive	Negative
AMP	0ng/ml	Negative	60	0	60	0	60	0	60
AMP	250ng/ml	-75%cutoff	60	0	60	0	60	0	60

	500ng/ml	-50%cutoff	60	0	60	0	60	0	60
	750ng/ml	-25%cutoff	60	6	54	6	54	8	52
	1000ng/ml	cutoff	60	34	26	38	22	36	24
	1250ng/ml	+25%cutoff	60	54	6	58	2	56	4
	1500ng/ml	+50%cutoff	60	60	0	60	0	60	0
	1750ng/ml	+75%cutoff	60	60	0	60	0	60	0
	2000ng/ml	+100%cutoff	60	60	0	60	0	60	0
	0ng/ml	Negative	60	0	60	0	60	0	60
	75ng/ml	-75%cutoff	60	0	60	0	60	0	60
	150ng/ml	-50%cutoff	60	0	60	0	60	0	60
	225ng/ml	-25%cutoff	60	6	54	2	58	4	56
BAR	300ng/ml	cutoff	60	38	22	36	24	34	26
	375ng/ml	+25%cutoff	60	58	2	54	6	56	4
	450ng/ml	+50%cutoff	60	60	0	60	0	60	0
	525ng/ml	+75%cutoff	60	60	0	60	0	60	0
	600ng/ml	+100%cutoff	60	60	0	60	0	60	0
	0ng/ml	Negative	60	0	60	0	60	0	60
	75ng/ml	-75%cutoff	60	0	60	0	60	0	60
	150ng/ml	-50%cutoff	60	0	60	0	60	0	60
	225ng/ml	-25%cutoff	60	8	52	4	56	6	54
BZO	300ng/ml	cutoff	60	36	24	36	24	38	22
	375ng/ml	+25%cutoff	60	52	8	54	6	56	4
	450ng/ml	+50%cutoff	60	60	0	60	0	60	0
	525ng/ml	+75%cutoff	60	60	0	60	0	60	0
	600ng/ml	+100%cutoff	60	60	0	60	0	60	0
	0ng/ml	Negative	60	0	60	0	60	0	60
	75ng/ml	-75%cutoff	60	0	60	0	60	0	60
	150ng/ml	-50%cutoff	60	0	60	0	60	0	60
	225ng/ml	-25%cutoff	60	2	58	4	56	4	56
MTD	300ng/ml	cutoff	60	38	22	34	26	36	24
	375ng/ml	+25%cutoff	60	58	2	54	6	56	4
	450ng/ml	+50%cutoff	60	60	0	60	0	60	0
	525ng/ml	+75%cutoff	60	60	0	60	0	60	0
	600ng/ml	+100%cutoff	60	60	0	60	0	60	0
	0ng/ml	Negative	60	0	60	0	60	0	60
OXY	25ng/ml	-75%cutoff	60	0	60	0	60	0	60
OAI	50ng/ml	-50%cutoff	60	0	60	0	60	0	60
	75ng/ml	-25%cutoff	60	6	54	8	52	4	56

	100ng/ml	cutoff	60	34	26	36	24	38	22
	125ng/ml	+25%cutoff	60	54	6	58	2	56	4
	150ng/ml	+50%cutoff	60	60	0	60	0	60	0
	175ng/ml	+75%cutoff	60	60	0	60	0	60	0
	200ng/ml	+100%cutoff	60	60	0	60	0	60	0
	0ng/ml	Negative	60	0	60	0	60	0	60
	6.3ng/ml	-75%cutoff	60	0	60	0	60	0	60
	12.5ng/ml	-50%cutoff	60	0	60	0	60	0	60
	18.8ng/ml	-25%cutoff	60	4	56	4	56	6	54
PCP	25ng/ml	cutoff	60	38	22	36	24	34	26
	31.3ng/ml	+25%cutoff	60	56	4	54	6	52	8
	37.5ng/ml	+50%cutoff	60	60	0	60	0	60	0
	43.8ng/ml	+75%cutoff	60	60	0	60	0	60	0
	50ng/ml	+100%cutoff	60	60	0	60	0	60	0

## Single drug Test Dipcard:

D	Approximate		Number of			Re	sult		
Drug	concentration	% of cutoff	determinations	Lo	ot 1	Lo	ot 2	Lo	ot 3
test	of sample		per lot	Positive	Negative	Positive	Negative	Positive	Negative
	0ng/ml	Negative	60	0	60	0	60	0	60
	250ng/ml	-75%cutoff	60	0	60	0	60	0	60
	500ng/ml	-50%cutoff	60	0	60	0	60	0	60
	750ng/ml	-25%cutoff	60	8	52	6	54	10	50
AMP	1000ng/ml	cutoff	60	36	24	34	26	36	24
	1250ng/ml	+25%cutoff	60	50	10	52	8	54	6
	1500ng/ml	+50%cutoff	60	60	0	60	0	60	0
	1750ng/ml	+75%cutoff	60	60	0	60	0	60	0
	2000ng/ml	+100%cutoff	60	60	0	60	0	60	0
	0ng/ml	Negative	60	0	60	0	60	0	60
	75ng/ml	-75%cutoff	60	0	60	0	60	0	60
	150ng/ml	-50%cutoff	60	0	60	0	60	0	60
	225ng/ml	-25%cutoff	60	6	54	4	56	4	56
BAR	300ng/ml	cutoff	60	38	22	36	24	34	26
	375ng/ml	+25%cutoff	60	56	4	54	6	56	4
	450ng/ml	+50%cutoff	60	60	0	60	0	60	0
	525ng/ml	+75%cutoff	60	60	0	60	0	60	0
	600ng/ml	+100%cutoff	60	60	0	60	0	60	0
	0ng/ml	Negative	60	0	60	0	60	0	60
BZO	75ng/ml	-75%cutoff	60	0	60	0	60	0	60
	150ng/ml	-50%cutoff	60	0	60	0	60	0	60

	225ng/ml	-25%cutoff	60	8	52	6	54	4	56
	300ng/ml	cutoff	60	36	24	38	22	36	24
	375ng/ml	+25%cutoff	60	54	6	54	6	56	4
	450ng/ml	+50%cutoff	60	60	0	60	0	60	0
	525ng/ml	+75%cutoff	60	60	0	60	0	60	0
	600ng/ml	+100%cutoff	60	60	0	60	0	60	0
	0ng/ml	Negative	60	0	60	0	60	0	60
	75ng/ml	-75%cutoff	60	0	60	0	60	0	60
	150ng/ml	-50%cutoff	60	0	60	0	60	0	60
	225ng/ml	-25%cutoff	60	4	56	2	58	6	54
MTD	300ng/ml	cutoff	60	32	28	36	24	34	26
MTD	375ng/ml	+25%cutoff	60	56	4	54	6	58	2
	450ng/ml	+50%cutoff	60	60	0	60	0	60	0
	525ng/ml	+75%cutoff	60	60	0	60	0	60	0
	600ng/ml	+100%cutoff	60	60	0	60	0	60	0
	0ng/ml	Negative	60	0	60	0	60	0	60
	0ng/ml	Negative	60	0	60	0	60	0	60
	25ng/ml	-75%cutoff	60	0	60	0	60	0	60
	50ng/ml	-50%cutoff	60	0	60	0	60	0	60
	75ng/ml	-25%cutoff	60	2	58	6	54	4	56
OXY	100ng/ml	cutoff	60	34	26	36	24	34	26
	125ng/ml	+25%cutoff	60	54	6	56	4	58	2
	150ng/ml	+50%cutoff	60	60	0	60	0	60	0
	175ng/ml	+75%cutoff	60	60	0	60	0	60	0
	200ng/ml	+100%cutoff	60	60	0	60	0	60	0
	0ng/ml	Negative	60	0	60	0	60	0	60
	6.3ng/ml	-75%cutoff	60	0	60	0	60	0	60
	12.5ng/ml	-50%cutoff	60	0	60	0	60	0	60
	18.8ng/ml	-25%cutoff	60	6	54	6	54	4	56
PCP	25ng/ml	cutoff	60	38	22	36	24	38	22
	31.3ng/ml	+25%cutoff	60	56	4	56	4	58	2
	37.5ng/ml	+50%cutoff	60	60	0	60	0	60	0
	43.8ng/ml	+75%cutoff	60	60	0	60	0	60	0
	50ng/ml	+100%cutoff	60	60	0	60	0	60	0

## Multi-drug Test Dipcard:

Denza	Approximate		Number of	Number of Result					
Drug	concentration	% of cutoff	determinations	Lot 1		Lot 2		Lot 3	
		per lot	Positive	Negative	Positive	Negative	Positive	Negative	
AMP	0ng/ml	Negative	60	0	60	0	60	0	60
AMP	250ng/ml	-75%cutoff	60	0	60	0	60	0	60

	500ng/ml	-50%cutoff	60	0	60	0	60	0	60
	750ng/ml	-25%cutoff	60	10	50	8	52	6	54
	1000ng/ml	cutoff	60	36	24	32	28	34	26
	1250ng/ml	+25%cutoff	60	56	4	54	6	54	6
	1500ng/ml	+50%cutoff	60	60	0	60	0	60	0
	1750ng/ml	+75%cutoff	60	60	0	60	0	60	0
	2000ng/ml	+100%cutoff	60	60	0	60	0	60	0
	0ng/ml	Negative	60	0	60	0	60	0	60
	75ng/ml	-75%cutoff	60	0	60	0	60	0	60
	150ng/ml	-50%cutoff	60	0	60	0	60	0	60
	225ng/ml	-25%cutoff	60	4	56	6	54	8	52
BAR	300ng/ml	cutoff	60	36	24	34	26	38	22
	375ng/ml	+25%cutoff	60	56	4	54	6	54	6
	450ng/ml	+50%cutoff	60	60	0	60	0	60	0
	525ng/ml	+75%cutoff	60	60	0	60	0	60	0
	600ng/ml	+100%cutoff	60	60	0	60	0	60	0
	0ng/ml	Negative	60	0	60	0	60	0	60
	75ng/ml	-75%cutoff	60	0	60	0	60	0	60
	150ng/ml	-50%cutoff	60	0	60	0	60	0	60
	225ng/ml	-25%cutoff	60	6	54	4	56	6	54
BZO	300ng/ml	cutoff	60	38	22	34	26	36	24
	375ng/ml	+25%cutoff	60	56	4	54	6	56	4
	450ng/ml	+50%cutoff	60	60	0	60	0	60	0
	525ng/ml	+75%cutoff	60	60	0	60	0	60	0
	600ng/ml	+100%cutoff	60	60	0	60	0	60	0
	0ng/ml	Negative	60	0	60	0	60	0	60
	75ng/ml	-75%cutoff	60	0	60	0	60	0	60
	150ng/ml	-50%cutoff	60	0	60	0	60	0	60
	225ng/ml	-25%cutoff	60	6	54	4	56	2	58
MTD	300ng/ml	cutoff	60	36	24	34	26	36	24
	375ng/ml	+25%cutoff	60	56	4	58	2	56	4
	450ng/ml	+50%cutoff	60	60	0	60	0	60	0
	525ng/ml	+75%cutoff	60	60	0	60	0	60	0
	600ng/ml	+100%cutoff	60	60	0	60	0	60	0
	0ng/ml	Negative	60	0	60	0	60	0	60
	25ng/ml	-75%cutoff	60	0	60	0	60	0	60
OXY	50ng/ml	-50%cutoff	60	0	60	0	60	0	60
	75ng/ml	-25%cutoff	60	8	52	4	56	6	54
	100ng/ml	cutoff	60	38	22	36	24	34	26

	125ng/ml	+25%cutoff	60	56	4	58	2	54	6
	150ng/ml	+50%cutoff	60	60	0	60	0	60	0
	175ng/ml	+75%cutoff	60	60	0	60	0	60	0
	200ng/ml	+100%cutoff	60	60	0	60	0	60	0
	0ng/ml	Negative	60	0	60	0	60	0	60
	6.3ng/ml	-75%cutoff	60	0	60	0	60	0	60
	12.5ng/ml	-50%cutoff	60	0	60	0	60	0	60
	18.8ng/ml	-25%cutoff	60	2	58	4	56	6	54
PCP	25ng/ml	cutoff	60	36	24	38	22	34	26
	31.3ng/ml	+25%cutoff	60	58	2	56	4	54	6
	37.5ng/ml	+50%cutoff	60	60	0	60	0	60	0
	43.8ng/ml	+75%cutoff	60	60	0	60	0	60	0
	50ng/ml	+100%cutoff	60	60	0	60	0	60	0

### 8.6 Accuracy

80 clinical urine samples collected all from sample place several hospitals and drug relief reformatory. All clinical urine specimens for each drug were analyzed by GC/MS and by two lots of the corresponding Single/Multi-drug Test Cup and Test Dipcard. Samples were divided by concentration into five categories: drug free, less than half the cutoff, near cutoff negative, near cutoff positive, and high positive. All samples were blindly labeled by a nonparticipant. Separate sets of the blind coded were assigned. Samples were also randomized prior to testing. The study was conducted by 4 nurses at two Point-of-Care sites. The test dipcard format was performed at one site and the test cup format at the second site. Each operator only performed one test format and different nurses tested each format. Results were as follows:

### Single drug Test Cup:

Drug Test	Xenta Result	Drug free by GC/MS analysis	Less than half the cutoff concentration by GC/MS analysis	Near Cutoff Negative (Between 50% below the cutoff and the cutoff concentration)	Near Cutoff Positive (Between the cutoff and 50% above the cutoff concentration)	High Positive (greater than 50% above the cutoff concentration)	Total
AMP	+	0	0	0	6	34	80
AMP	_	33	2	5	0	0	80
BAR	+	0	0	0	6	34	80
BAK	-	33	0	7	0	0	80
D70	+	0	0	1	7	33	80
BZO	-	31	0	8	0	0	80
MTD	+	0	0	1	5	35	80
MTD	-	32	2	5	0	0	80
OVV	+	0	0	0	6	34	90
OXY	-	35	0	5	0	0	80
DCD.	+	0	0	1	5	35	90
PCP	ı	35	0	4	0	0	80

## Analysis of Discordant Results with Single drug Test Cup

Sin	gle drug Test Cup		GC/MS Analysis		
Drug Test	Cutoff(ng/mL)	Test Result	Drug Concentration (ng/mL)	Drug in Urine	
BZO	300	Positive	188	Oxazepam	
MTD	300	Positive	209	Methadone	
PCP	25	Positive	23	Phencyclidine	

## Multi-drug Test Cup:

Drug Test	Xenta Result	Drug free by GC/MS analysis	Less than half the cutoff concentration by GC/MS analysis	Near Cutoff Negative (Between 50% below the cutoff and the cutoff concentration)	Near Cutoff Positive (Between the cutoff and 50% above the cutoff concentration)	High Positive (greater than 50% above the cutoff concentration)	Total
AMP	+	0	0	0	6	34	80
Alvii	ı	33	2	5	0	0	80
BAR	+	0	0	0	6	34	80
DAK	-	33	0	7	0	0	80
BZO	+	0	0	1	7	33	80
	_	31	0	8	0	0	80
MTD	+	0	0	1	5	35	80
	_	32	2	5	0	0	80
OXY	+	0	0	0	6	34	80
	_	35	0	5	0	0	
PCP	+	0	0	1	5	35	80
	_	35	0	4	0	0	80

## Analysis of Discordant Results with Multi-drug Test Cup

M	Iulti-drug Test C	up	GC/MS Analysis		
Drug Test	Cutoff(ng/m L)	Test Result	Drug Concentration (ng/mL)	Drug in Urine	
BZO	300	Positive	188	Oxazepam	
MTD	300	Positive	209	Methadone	
PCP	25	Positive	23	Phencyclidine	

## Single drug Test Dipcard:

Drug Test	Co-Innovation Result	Drug free by GC/MS analysis	Less than half the cutoff concentration by GC/MS analysis	Near Cutoff Negative (Between 50% below the cutoff and the cutoff concentration)	Near Cutoff Positive (Between the cutoff and 50% above the cutoff concentration)	High Positive (greater than 50% above the cutoff concentration)	Total
AMP	+	0	0	0	6	34	80
AIVIP	-	33	2	5	0	0	80
BAR	+	0	0	0	6	34	80

	-	33	0	7	0	0	
BZO	+	0	0	1	7	33	80
	_	31	0	8	0	0	80
MTD	+	0	0	1	5	35	80
	_	32	2	5	0	0	80
OXY	+	0	0	0	6	34	80
	_	35	0	5	0	0	80
PCP	+	0	0	1	5	35	80
	_	35	0	4	0	0	30

### Analysis of Discordant Results with Single drug Test Dipcard

Sir	ngle drug Test Dipc	ard		GC/MS Analysis
Drug Test	Cutoff(ng/mL)	Test Result	Drug Concentration (ng/mL)	Drug in Urine
BZO	300	Positive	188	Oxazepam
MTD	300	Positive	209	Methadone
PCP	25	Positive	23	Phencyclidine

### Multi-drug Test Dipcard:

Drug Test	Co-Innovation Result	Drug free by GC/MS analysis	Less than half the cutoff concentration by GC/MS analysis	Near Cutoff Negative (Between 50% below the cutoff and the cutoff concentration)	Near Cutoff Positive (Between the cutoff and 50% above the cutoff concentration)	High Positive (greater than 50% above the cutoff concentration)	Total	
AMP	+	0	0	0	6	34	80	
	-	33	2	5	0	0		
BAR	+	0	0	0	6	34	80	
	-	33	0	7	0	0	80	
BZO	+	0	0	1	7	33	00	
	_	31	0	8	0	0	80	
MTD	+	0	0	1	5	35	00	
	_	32	2	5	0	0	80	
OXY	+	0	0	0	6	34	0.0	
	_	35	0	5	0	0	80	
PCP	+	0	0	1	5	35	00	
	_	35	0	4	0	0	80	

### Analysis of Discordant Results with Rapid Multi-drug Test Dipcard

M	ulti-drug Test Dipca	ard	GC/MS Analysis		
Drug Test	Cutoff(ng/mL)	Test Result	Drug Concentration (ng/mL)	Drug in Urine	
BZO	300	Positive	188	Oxazepam	
MTD	300	Positive	209	Methadone	
PCP	25	Positive	23	Phencyclidine	

#### Stability of the test line studies:

Stability of the test line studies were performed using the single drug and multi-drug test formats. Drug free specimens were spiked with analytes at different concentration containing  $\pm$ 0% cutoff and  $\pm$ 100% cutoff. Each sample contain all of eleven drugs and the concentrations of target drugs were confirmed with GC/MS. Spiked specimens and five drug free specimens were test.

The results show that the color T line of Xenta Drug Screen Cup and Xenta Drug Screen Dipcard are stable from 3 to 50 minutes.

We suggest to read the test results within 5 to 30 minutes.

#### 9. Conclusion:

The data collected in the performance and accuracy studies demonstrate that the Xenta Drug Screen Cup and Xenta Drug Screen Dipcard are substantially equivalent to the predicate device.

--- End of this section ---