

August 2, 2023

Simplivia Healthcare LTD. Shay Shaham VP QA/RA North Industrial Zone Kiryat Shmona, 1101801 Israel

Re: K231286

Trade/Device Name: Chemfort<sup>®</sup> Catheter Adaptor Regulation Number: 21 CFR 880.5440 Regulation Name: Intravascular Administration Set Regulatory Class: Class II Product Code: ONB Dated: May 3, 2023 Received: May 4, 2023

Dear Shay Shaham:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <a href="https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm">https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm</a> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part

801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <u>https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products</u>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Davil Walloscher

David Wolloscheck, Ph.D. Assistant Director DHT3C: Division of Drug Delivery and General Hospital Devices, and Human Factors OHT3: Office of GastroRenal, ObGyn, General Hospital and Urology Devices Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

# Indications for Use

510(k) Number *(if known)* K231286

Device Name Chemfort® Catheter Adaptor

#### Indications for Use (Describe)

for up to 7 days.

Chemfort® Catheter Adaptor is a single use, sterile Closed System Transfer Device (CSTD) that mechanically prohibits the release of drugs, including antineoplastic and hazardous drugs, in vapor, aerosol or liquid form during administration, thus minimizing exposure of individuals, healthcare personnel, and the environment to hazardous drugs. Chemfort® Catheter Adaptor prevents the introduction of microbial and airborne contaminants into the drug or fluid path

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

#### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

#### \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov* 

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



# K231286 510(K) SUMMARY

Preparation Date:	July 31, 2023
Device name:	Chemfort <sup>®</sup> Catheter Adaptor
Type of 510(k) submission:	Traditional
Date of Submission:	May 3, 2023
Applicant's name:	Simplivia Healthcare LTD. North Industrial Zone Kiryat Shmona, 1101801 Israel
Phone: Fax:	(972) 4 6908826 (972) 74 7652161
FDA Registration Number	9611423
Contact Person:	Shay Shaham VP QA / RA
FDA Product Code:	ONB
FDA Regulation Number:	21 CFR 880.5440
FDA Regulation Name:	Intravascular administration set
<b>Classification Panel:</b>	General Hospital
Common Name:	Closed Antineoplastic and Hazardous Drug Reconstitution and Transfer System
FDA Classification:	Class II
Predicate Device:	Tevadaptor <sup>®</sup> Catheter Adaptor (K180489)



### **Indications for Use**

Chemfort<sup>®</sup> Catheter Adaptor is a single use, sterile Closed System Transfer Device (CSTD) that mechanically prohibits the release of drugs, including antineoplastic and hazardous drugs, in vapor, aerosol or liquid form during administration, thus minimizing exposure of individuals, healthcare personnel, and the environment to hazardous drugs.

Chemfort<sup>®</sup> Catheter Adaptor prevents the introduction of microbial and airborne contaminants into the drug or fluid path for up to 7 days.

### **Device Description**

The Chemfort<sup>®</sup> Catheter Adaptor enables drug transfer to the catheter, thus allowing drug administration to the patient's urinary bladder. The use of elastomeric seals of the Chemfort<sup>®</sup> Catheter Adaptor prevents hazardous drug contamination of healthcare professionals, the patient, and the environment.

The Chemfort<sup>®</sup> Catheter Adaptor is an addition to the cleared Chemfort<sup>®</sup> system (K192866). The Catheter Adaptor provides closed system protection during the following procedures:

- a) Drug transfer from a standard luer lock syringe to the Catheter Adaptor through the Chemfort<sup>®</sup> Syringe Adaptor (K192866).
- b) Closed system drug administration to the urinary bladder, through a urinary catheter. The Catheter Adaptor fits a wide range of standard catheter sizes and converts an open catheter connection to a closed Chemfort<sup>®</sup> connection.

The Chemfort<sup>®</sup> Catheter Adaptor allows the healthcare professional to have the option for safe drug administration to the urinary catheter and safe disconnection of the Chemfort<sup>®</sup> Syringe Adaptor (K192866) from the patient's urinary catheter.

## Summary of Technological Characteristics:

The following table (**Table 2**) compares the Chemfort<sup>®</sup> Catheter Adaptor to the predicate device with respect to intended use, technological characteristics and principles of operation, providing detailed information regarding the basis for the determination of substantial equivalence.

	Proposed Device Chemfort <sup>®</sup> Catheter Adaptor	Predicate Device Tevadaptor® Catheter Adaptor (K180489)	Equivalence to predicate
Device Class	Class II	Class II	Same
Classification Panel	General Hospital	General Hospital	Same
Product Code	ONB	ONB	Same
Regulation	Intravascular	Intravascular	Same
Description	Administration Set	Administration Set	
<b>Regulation No.</b>	21 C.F.R. §880.5440	21 C.F.R. §880.5440	Same
Indications for use	Chemfort <sup>®</sup> Catheter Adaptor is a single use, sterile Closed System Transfer Device (CSTD) that mechanically prohibits the release of drugs, including antineoplastic and hazardous drugs, in vapor, aerosol or liquid form during administration, thus minimizing exposure of individuals, healthcare personnel, and the environment to hazardous drugs. Chemfort <sup>®</sup> Catheter Adaptor prevents the introduction of microbial and airborne contaminants into the drug or fluid path for up to 7 days.	Tevadaptor <sup>®</sup> is a Closed System Drug Transfer Device (CSTD) that mechanically prohibits the release of the drug in vapor, aerosol or liquid form during preparation and administration, and prevents the introduction of microbial and airborne contaminants into the drug or fluid path, allowing the system to minimize exposure of individuals, healthcare personnel, and the environment to hazardous drugs.	First part: Same meaning. Second part: Tevadaptor <sup>®</sup> was tested and proved to prevent contaminants from entering the drug or fluid path for up to 3 days, Chemfort <sup>®</sup> has been tested and approved for 7 days
Components	Part of Chemfort <sup>®</sup> , a multi- components system	Part of Tevadaptor <sup>®</sup> , a multi-components system	Same
	including Catheter Adaptor	including Catheter Adaptor	
Interaction with	The distal end connects to	The distal end connects to	Same
other devices	the urinary catheter.	the urinary catheter.	
	The proximal end connects	The proximal end connects	

Table 1.	Proposed	Device,	Refefance	<b>Device and</b>	Predicate	Device	Comparation
----------	----------	---------	-----------	-------------------	-----------	--------	-------------

	Proposed Device Chemfort <sup>®</sup> Catheter Adaptor	Predicate Device Tevadaptor <sup>®</sup> Catheter Adaptor (K180489)	Equivalence to predicate
	to Chemfort <sup>®</sup> Syringe	to Tevadaptor <sup>®</sup> Syringe	
	Adaptor / Syringe Adaptor	Adaptor / Syringe Adaptor	
	Lock.	Lock.	~
Re-use	Distal end: to maintain the	Distal end: to maintain the	Same
capability	closed system, the Catheter	closed system, the Catheter	
	Adaptor should not be	Adaptor should not be	
	uringry antheter	uringry astheter	
	Provimal end: The	Provimal end: The	
	Chemfort <sup>®</sup> port of the	Tevadaptor <sup>®</sup> port of the	
	Catheter Adaptor can be	Catheter Adaptor can be	
	connected and	connected and	
	disconnected from the	disconnected from the	
	Syringe Adaptor port up to	Syringe Adaptor port up to	
	10 times.	10 times.	
Principles of	Multi-component system,	Multi-component system,	Same
Operation	components are intended	components are intended	
	to be used as a system,	to be used as a system,	
	manually manipulated.	manually manipulated.	
Interaction with	No direct interaction-	No direct interaction-	Same
patient	device interaction with the	device interaction with the	
	patient is achieved through	patient is achieved through	
	through the urinery	through the urinery	
	catheter	catheter	
	cameter.	cameter.	
Interconnecting	Mechanical snap	Mechanical snap	Same
features	connections.	connections.	
Safety features	Vented cap	Vented cap	Same
-	• Septum to septum	• Septum to septum	
	contact	contact	
Target users	Nurses or other healthcare	Nurses or other healthcare	Same
	professionals.	professionals.	
Technology	All of the Chemfort <sup>®</sup>	All of the Tevadaptor <sup>®</sup>	Same
	devices ports are sealed	devices ports are sealed	
	with resealing Septum.	with resealing Septum.	
	When Syringe Adaptor	When Syringe Adaptor	
	and Chemfort <sup>®</sup> port are	and levadaptor port are	
	Joined, the two septums	Joined, the two septums	
	then nierced by a needle	then nierced by a needle	
	(from the Syringe Adaptor	(from the Syringe Adaptor	
	or Svringe Adaptor Lock).	or Syringe Adaptor Lock).	

	Proposed Device Chemfort <sup>®</sup> Catheter Adaptor	Predicate Device Tevadaptor <sup>®</sup> Catheter Adaptor (K180489)	Equivalence to predicate
	thus creating a secured	thus creating a secured	
Environment of use	Hospitals, compounding centers and clinics	Hospitals, compounding centers and clinics	Same
Sterilization method	Ethylene Oxide validated cycle SAL 10 <sup>-6</sup>	Ethylene Oxide validated cycle SAL 10 <sup>-6</sup>	Same
Biocompatibility	All Catheter Adaptor parts that are in contact with patient comply with the requirements of ISO 10993-1	All Catheter Adaptor parts that are in contact with patient comply with the requirements of ISO 10993-1	Same
Prescription use	Rx only	Rx only	Same

### **Performance Data**

Simplivia Healthcare conducted several performance tests to demonstrate that the Chemfort<sup>®</sup> Catheter Adaptor complies with the following standards and that it functions as intended.

- ISO 10993-1:2018, Biological Evaluation of Medical Devices. Part 1: Evaluation and testing within a risk management process.
- ISO 10993-4:2017, Biological Evaluation of Medical Devices. Part 4: Selection of tests for interactions with blood.
- ISO 10993-5:2009, Biological Evaluation of Medical Devices. Part 5: Tests for in vitro cytotoxicity.
- ISO 10993-7:2008/Amd 1:2019, Biological Evaluation of Medical Devices. Part 7: Ethylene oxide sterilization residuals.
- ISO 10993-10:2021, Biological Evaluation of Medical Devices. Part 10: Tests for irritation and skin sensitization.
- ISO 10993-11:2017, Biological Evaluation of Medical Devices. Part 11: Tests for systemic toxicity.
- ISO 10993-18:2020, Biological Evaluation of Medical Devices. Part 18: Chemical characterization of medical device materials within a risk management process.
- ISO 11135:2014 + Amd.1:2018, Sterilization of health-care products Ethylene oxide Requirements for the development, validation and routine control of a sterilization process for medical devices.
- ISO 11607-1:2019, Packaging for terminally sterilized medical devices Part 1: Requirements for materials, sterile barrier systems and packaging systems.

- ISO 14971:2019 –Medical devices Medical devices —Application of risk management to medical devices
- USP <85>, Bacterial Endotoxins Test.
- USP <161>, Transfusion and Infusion Assemblies and Similar Medical Devices.
- USP <788>, Particulate Matter in Injections

#### Substantial Equivalence

Simplivia Healthcare's Chemfort<sup>®</sup> Catheter Adaptor has similar indications for use, and similar technological characteristics and principles of operation as the predicate device, K180489. Performance data demonstrated that the Chemfort<sup>®</sup> Catheter Adaptor is as safe and as effective as its predicate and does not raise any new safety and effectiveness issues. Thus, Simplivia Healthcare's Chemfort<sup>®</sup> Catheter Adaptor is substantially equivalent to its predicate device, K180489.