



June 6, 2023

Nephron Nitrile, LLC.
Lou Kennedy
Chief Executive Officer
4777 12th Street Extension
West Columbia, South Carolina 29172

Re: K231349

Trade/Device Name: Nephron Nitrile Powder-Free Nitrile Examination Gloves (Tested For Use With
Chemotherapy Drugs and Fentanyl)
Regulation Number: 21 CFR 880.6250
Regulation Name: Non-Powdered Patient Examination Glove
Regulatory Class: Class I, reserved
Product Code: LZA, LZC, OPJ, QDO
Dated: May 9, 2023
Received: May 9, 2023

Dear Lou Kennedy:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,


Bifeng Qian -S

Bifeng Qian, M.D., Ph.D.
Assistant Director
DHT4B: Division of Infection Control
and Plastic Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K231349

Device Name

Nephron Nitrile Powder-Free Nitrile Examination Gloves (Tested For Use With Chemotherapy Drugs and Fentanyl)

Indications for Use (Describe)

Nephron Nitrile Powder-Free Nitrile Examination Gloves (Tested For Use With Chemotherapy Drugs and Fentanyl) is a disposable device intended for medical purpose that is worn on the examiner's hand to prevent contamination between patient and examiner.

In addition, these gloves were tested for use with Chemotherapy drugs and Fentanyl in accordance with ASTM D6978-05 (2019) Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs.

The following chemotherapy drugs and concentration had NO breakthrough detected up to 240 minutes:

| | | |
|--------------------------------|-----------------------------|---------------------------------|
| Bleomycin Sulfate (15.0 mg/ml) | Doxorubicin HCl (2.0 mg/ml) | Mechlorethamine HCl (1.0 mg/ml) |
| Busulfan (6.0 mg/ml) | Epirubicin HCl (2.0 mg/ml) | Melphalan (5.0 mg/ml) |
| Carboplatin (10.0 mg/ml) | Etoposide (20.0 mg/ml) | Methotrexate (25.0 mg/ml) |
| Cisplatin (1.0 mg/ml) | Fludarabine (25.0 mg/ml) | Mitomycin C (0.5 mg/ml) |
| Cyclophosphamide (20.0 mg/ml) | Fluorouracil (50.0 mg/ml) | Mitoxantrone HCl (2.0 mg/ml) |
| Cytarabine (100.0 mg/ml) | Gemcitabine (38.0 mg/ml) | Paclitaxel (6.0 mg/ml) |
| Dacarbazine (10.0 mg/ml) | Idarubicin HCl (1.0 mg/ml) | Rituximab (10.0 mg/ml) |
| Daunorubicin HCl (5.0 mg/ml) | Ifosfamide (50.0 mg/ml) | Trisenox (1.0 mg/ml) |
| Docetaxel (10.0 mg/ml) | Irinotecan (20.0 mg/ml) | Vincristine Sulfate (1.0 mg/ml) |

The following chemotherapy drugs have low permeation times:

Carmustine (3.3 mg/ml) : 33.8 minutes Thiotepea (10.0 mg/ml) : 128.1 minutes

Warning: Not for Use with: Carmustine, Thiotepea

The tested Opioid is:

Fentanyl Citrate Injection (100 mcg/2 mL). Permeation: no breakthrough up to 240 minutes

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(K) SUMMARY

K231349

AS REQUIRED BY: 21CFR§807.92

A. APPLICANT INFORMATION

| | |
|---------------------|--|
| 510(K) Owner's Name | Nephron Nitrile, LLC. |
| Address | 4777 12th Street Extension, West Columbia, SC 29172. |
| Phone | 844-937-3888 |
| Fax | 1-803-926-9853 |
| E-mail | lkennedy@nephronpharm.com |
| Contact Person | Lou Kennedy |
| Designation | Chief Executive Officer |
| Contact Number | 1-803-569-3110 |
| Contact Email | lkennedy@nephronpharm.com |
| | |
| Date Prepared | 31 May 2023 |

B. DEVICE IDENTIFICATION

| | |
|-----------------------------------|--|
| Name of the device | Nephron Nitrile Powder-Free Nitrile Examination Gloves (Tested For Use With Chemotherapy Drugs and Fentanyl) |
| Product proprietary or trade name | Nephron Nitrile |
| Common or usual name | Nitrile Examination Gloves (Tested for use with Chemotherapy drugs and Fentanyl) |
| 510(k) Number | K231349 |
| Classification name | Patient Examination Glove, Specialty |
| Device Classification | Class-1 |
| Product Code | LZA, LZC, OPJ, QDO |
| Regulation Number | 21 CFR 880.6250 |
| Review Panel | General Hospital |

C. PREDICATE DEVICE

| | |
|-------------------------|---|
| Predicate Device | Nephron Nitrile Powder-Free Nitrile Examination Gloves (Tested For Use With Chemotherapy Drugs) |
| 510(k) Number | K223559 |
| Regulatory Class | Class-1 |
| Product code | LZA, LZC, OPJ |
| Owner | Nephron Nitrile, LLC. |

510(K) SUMMARY

K231349

AS REQUIRED BY: 21CFR§807.92

D. DESCRIPTION OF THE DEVICE:

Nephron Nitrile Powder-Free Nitrile Examination Gloves (Tested For Use With Chemotherapy Drugs and Fentanyl) is a Class I patient examination gloves bearing the product codes LZA, LZC, OPJ, QDO (21 CFR 880.6250). They meet all the current specifications listed under the ASTM D6319-19, Standard Specification for Nitrile Examination Gloves for Medical Application and also complies with requirements for Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs as per ASTM D6978-05 (2019). They are made from Nitrile (NBR). These gloves are blue in color and are powder free. The product is non-sterile, fingertip textured, ambidextrous with beaded cuff and single use only.

Nephron Nitrile Powder-Free Nitrile Examination Gloves (Tested For Use With Chemotherapy Drugs and Fentanyl) with sizes Medium, Large, X-Large and XX-Large are included in the submission.

E. INDICATION FOR USE OF THE DEVICE:

Nephron Nitrile Powder-Free Nitrile Examination Gloves (Tested For Use With Chemotherapy Drugs and Fentanyl) is a disposable device intended for medical purpose that is worn on the examiner's hand to prevent contamination between patient and examiner. In addition, these gloves were tested for use with chemotherapy drugs and fentanyl in accordance with ASTM D6978-05 (2019) Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs.

| | | | |
|---|---|--|---|
| The following chemotherapy drugs and concentration had NO breakthrough detected up to 240 minutes: | | | The tested Opioid is: Fentanyl Citrate Injection (100 mcg/2 mL) Permeation: no breakthrough up to 240 minutes |
| Bleomycin Sulfate (15.0 mg/ml) Busulfan (6.0 mg/ml) Carboplatin (10.0 mg/ml) Cisplatin (1.0 mg/ml) Cyclophosphamide (20.0 mg/ml) Cytarabine (100.0 mg/ml) Dacarbazine (10.0 mg/ml) Daunorubicin HCl (5.0 mg/ml) Docetaxel (10.0 mg/ml) | Doxorubicin HCl (2.0 mg/ml) Epirubicin HCl (2.0 mg/ml) Etoposide (20.0 mg/ml) Fludarabine (25.0 mg/ml) Fluorouracil (50.0 mg/ml) Gemcitabine (38.0 mg/ml) Idarubicin HCl (1.0 mg/ml) Ifosfamide (50.0 mg/ml) Irinotecan (20.0 mg/ml) | Mechlorethamine HCl (1.0 mg/ml) Melphalan (5.0 mg/ml) Methotrexate (25.0 mg/ml) Mitomycin C (0.5 mg/ml) Mitoxantrone HCl (2.0 mg/ml) Paclitaxel (6.0 mg/ml) Rituximab (10.0 mg/ml) Trisenox (1.0 mg/ml) Vincristine Sulfate (1.0 mg/ml) | |
| The following chemotherapy drugs have low permeation times: Carmustine (3.3 mg/ml) : 33.8 minutes Thiotepa (10.0 mg/ml) : 128.1 minutes | | | |
| Warning: Not for Use with: Carmustine, Thiotepa | | | |

510(K) SUMMARY

K231349

AS REQUIRED BY: 21CFR§807.92

F. SUMMARY OF THE TECHNOLOGICAL CHARACTERISTICS OF THE DEVICE COMPARED TO THE PREDICATE DEVICE

| CHARACTERISTICS | STANDARDS | DEVICE PERFORMANCE | | Comparison |
|--------------------|-----------|--|--|---|
| | | PREDICATE | SUBJECT | |
| 510(K) Number | --- | K223559 | K231349 | |
| Name of device | --- | Nephron Nitrile Powder-Free Nitrile Examination Gloves (Tested For Use With Chemotherapy Drugs) | Nephron Nitrile Powder-Free Nitrile Examination Gloves (Tested For Use With Chemotherapy Drugs and Fentanyl) | Different: Except for the addition of Fentanyl for subject device, it is identical |
| Product Code | --- | LZA, LZC, OPJ | LZA, LZC, OPJ, QDO | Different: Added one product code |
| Indication for use | --- | Nephron Nitrile Powder-Free Nitrile Examination Gloves (Tested For Use With Chemotherapy Drugs) is a disposable device intended for medical purpose that is worn on the examiner's hand to prevent contamination between patient and examiner. In addition, these gloves were tested for use with chemotherapy drugs in accordance with ASTM D6978-05 (2019) Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs. | Nephron Nitrile Powder-Free Nitrile Examination Gloves (Tested For Use With Chemotherapy Drugs and Fentanyl) is a disposable device intended for medical purpose that is worn on the examiner's hand to prevent contamination between patient and examiner. In addition, these gloves were tested for use with chemotherapy drugs and fentanyl in accordance with ASTM D6978-05 (2019) Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs. | Different: Except for the addition of Fentanyl for subject device, it is identical Note: Adding the Fentanyl to the indication does not affect the intended use. |
| Regulation Number | --- | 21 CFR 880.6250 | 21 CFR 880.6250 | Identical |
| Material | --- | Nitrile | Nitrile | Identical |
| Color | --- | Blue | Blue | Identical |
| Size | --- | M, L, XL, XXL | M, L, XL, XXL | Identical |
| Single Use | --- | Single-use | Single-use | Identical |

510(K) SUMMARY

K231349

AS REQUIRED BY: 21CFR§807.92

| CHARACTERISTICS | STANDARDS | DEVICE PERFORMANCE | | Comparison |
|---|--------------------------|---|--|------------|
| | | PREDICATE | SUBJECT | |
| 510(K) Number | --- | K223559 | K231349 | |
| Sterile/non sterile | --- | Non Sterile | Non Sterile | Identical |
| Rx Only or OTC | --- | OTC | OTC | Identical |
| Dimensions - Length | ASTM D6319-19 | Minimum 230 mm (sizes M – XXL) | Minimum 230 mm (sizes M – XXL) | Identical |
| Dimensions - Width | ASTM D6319-19 | M: 95±10 mm L: 110±10 mm XL: 120±10 mm XXL: 130±10 mm | M: 95±10 mm L: 110±10 mm XL: 120±10 mm XXL: 130±10 mm | Identical |
| Physical Properties- Tensile Strength | ASTM D63192019 | Before aging 14MPa, min | Before aging 14MPa, min | Identical |
| | | After aging 14MPa, min | After aging 14MPa, min | Identical |
| Physical Properties- Ultimate Elongation | ASTM D63192019 | Before aging 500%, min | Before aging 500%, min | Identical |
| | | After aging 400%, min | After aging 400%, min | Identical |
| Thickness | ASTM D6319-19 | Palm: Minimum 0.05 mm Finger: Minimum 0.05 mm | Palm: Minimum 0.05 mm Finger: Minimum 0.05 mm | Identical |
| Powder Free Residue | ASTM D6319- 19 | ≤ 2 mg per glove | ≤ 2 mg per glove | Identical |
| Freedom from holes | ASTM D5151- 2019 | In accordance with ASTM D 5151-19, following ASTM D6319-19, G-I, AQL 2.5 | In accordance with ASTM D 5151-19, following ASTM D6319- 19, G-I, AQL 2.5 | Identical |
| Chemotherapy Drugs Tested with Minimum Breakthrough Detection Time | ASTM D6978- 05 (2019) | Bleomycin Sulfate 15 mg/ml (15,000 ppm) >240 Minutes | Bleomycin Sulfate 15 mg/ml (15,000 ppm) >240 Minutes | Identical |
| | | Busulfan 6 mg/ml (6,000 ppm) >240 Minutes | Busulfan 6 mg/ml (6,000 ppm) >240 Minutes | Identical |
| | | Carboplatin 10 mg/ml (10,000 ppm) >240 Minutes | Carboplatin 10 mg/ml (10,000 ppm) >240 Minutes | Identical |
| | | Carmustine 3.3 mg/ml (3,300 ppm) 33.8 Minutes | Carmustine 3.3 mg/ml (3,300 ppm) 33.8 Minutes | Identical |
| | | Cisplatin 1 mg/ml (1,000 ppm) >240 Minutes | Cisplatin 1 mg/ml (1,000 ppm) >240 Minutes | Identical |

510(K) SUMMARY

K231349

AS REQUIRED BY: 21CFR§807.92

| CHARACTERISTICS | STANDARDS | DEVICE PERFORMANCE | | Comparison |
|-----------------|-----------|--|--|------------|
| | | PREDICATE | SUBJECT | |
| 510(K) Number | --- | K223559 | K231349 | |
| | | Cyclophosphamide 20 mg/ml (20,000 ppm) >240 Minutes | Cyclophosphamide 20 mg/ml (20,000 ppm) >240 Minutes | Identical |
| | | Cytarabine 100 mg/ml (100,000 ppm) >240 Minutes | Cytarabine 100 mg/ml (100,000 ppm) >240 Minutes | Identical |
| | | Dacarbazine 10 mg/ml (10,000 ppm) >240 Minutes | Dacarbazine 10 mg/ml (10,000 ppm) >240 Minutes | Identical |
| | | Daunorubicin HCl 5 mg/ml (5,000 ppm) >240 Minutes | Daunorubicin HCl 5 mg/ml (5,000 ppm) >240 Minutes | Identical |
| | | Docetaxel 10 mg/ml (10,000 ppm) >240 Minutes | Docetaxel 10 mg/ml (10,000 ppm) >240 Minutes | Identical |
| | | Doxorubicin HCl 2 mg/ml (2,000 ppm) >240 Minutes | Doxorubicin HCl 2 mg/ml (2,000 ppm) >240 Minutes | Identical |
| | | Epirubicin HCl 2 mg/ml (2,000 ppm) >240 Minutes | Epirubicin HCl 2 mg/ml (2,000 ppm) >240 Minutes | Identical |
| | | Etoposide 20 mg/ml (20,000 ppm) >240 Minutes | Etoposide 20 mg/ml (20,000 ppm) >240 Minutes | Identical |
| | | Fludarabine 25 mg/ml (25,000 ppm) >240 Minutes | Fludarabine 25 mg/ml (25,000 ppm) >240 Minutes | Identical |
| | | Fluorouracil 50 mg/ml (50,000 ppm) >240 Minutes | Fluorouracil 50 mg/ml (50,000 ppm) >240 Minutes | Identical |
| | | Gemcitabine 38 mg/ml (38,000 ppm) >240 Minutes | Gemcitabine 38 mg/ml (38,000 ppm) >240 Minutes | Identical |
| | | Idarubicin HCl 1 mg/ml (1,000 ppm) >240 Minutes | Idarubicin HCl 1 mg/ml (1,000 ppm) >240 Minutes | Identical |
| | | Ifosfamide 50 mg/ml (50,000 ppm) >240 Minutes | Ifosfamide 50 mg/ml (50,000 ppm) >240 Minutes | Identical |
| | | Irinotecan 20 mg/ml (20,000 ppm) >240 Minutes | Irinotecan 20 mg/ml (20,000 ppm) >240 Minutes | Identical |
| | | Mechlorethamine HCl 1 mg/ml (1,000 ppm) >240 Minutes | Mechlorethamine HCl 1 mg/ml (1,000 ppm) >240 Minutes | Identical |

510(K) SUMMARY

K231349

AS REQUIRED BY: 21CFR§807.92

| CHARACTERISTICS | STANDARDS | DEVICE PERFORMANCE | | Comparison |
|---|--|---|---|------------|
| | | PREDICATE | SUBJECT | |
| 510(K) Number | --- | K223559 | K231349 | |
| | | Melphalan 5 mg/ml (5,000 ppm) >240 Minutes | Melphalan 5 mg/ml (5,000 ppm) >240 Minutes | Identical |
| | | Methotrexate 25 mg/ml (25,000 ppm) >240 Minutes | Methotrexate 25 mg/ml (25,000 ppm) >240 Minutes | Identical |
| | | Mitomycin C 0.5 mg/ml (500 ppm) >240 Minutes | Mitomycin C 0.5 mg/ml (500 ppm) >240 Minutes | Identical |
| | | Mitoxantrone HCl 2 mg/ml (2,000 ppm) >240 Minutes | Mitoxantrone HCl 2 mg/ml (2,000 ppm) >240 Minutes | Identical |
| | | Paclitaxel 6 mg/ml (6,000 ppm) >240 Minutes | Paclitaxel 6 mg/ml (6,000 ppm) >240 Minutes | Identical |
| | | Rituximab 10 mg/ml (10,000 ppm) >240 Minutes | Rituximab 10 mg/ml (10,000 ppm) >240 Minutes | Identical |
| | | Thiotepa 10 mg/ml (10,000 ppm) 128.1 Minutes | Thiotepa 10 mg/ml (10,000 ppm) 128.1 Minutes | Identical |
| | | Trisenox 1 mg/ml (1,000 ppm) >240 Minutes | Trisenox 1 mg/ml (1,000ppm) >240 Minutes | Identical |
| | | Vincristine Sulfate 1 mg/ml (1,000 ppm) >240 Minutes | Vincristine Sulfate 1 mg/ml (1,000 ppm) >240 Minutes | Identical |
| Opioid Drugs Tested with Minimum Breakthrough Detection Time | ASTM D6978- 05 (2019) | Not tested | Fentanyl Citrate Injection (100mcg/2mL) >240 Minutes | *Different |
| Biocompatibility | Primary Skin Irritation- ISO 10993-23: First Edition 2021-01 | Under the conditions of the study, the test article met the requirements of the test | Under the conditions of the study, the test article met the requirements of the test | Identical |
| | Dermal Sensitization- ISO 10993-10: Fourth Edition 2021-11 | Under the conditions of the study, the test article was not considered a sensitizer | Under the conditions of the study, the test article was not considered a sensitizer | Identical |

510(K) SUMMARY

K231349

AS REQUIRED BY: 21CFR§807.92

| CHARACTERISTICS | STANDARDS | DEVICE PERFORMANCE | | Comparison |
|----------------------|--|---|---|------------|
| | | PREDICATE | SUBJECT | |
| 510(K) Number | --- | K223559 | K231349 | |
| | In vitro cytotoxicity- ISO 10993-5: Third Edition 2009-06-01 | Under the Conditions of the study, the undiluted test article extract and 50% test article extract dilution did not meet the requirements of the test and the 25%, 12.5%, 6.25%, and 3.13% test article extract dilutions met the requirements of the test. | Under the conditions of the study, the undiluted test article extract and 50% test article extract dilution did not meet the requirements of the test and the 25%, 12.5%, 6.25%, and 3.13% test article extract dilutions met the requirements of the test. | Identical |
| | Acute Systemic Toxicity- ISO 10993-11: Third Edition 2017-09 | Under the conditions of the study, there was no mortality or evidence of systemic toxicity | Under the conditions of the study, there was no mortality or evidence of systemic toxicity | Identical |

* Identical except for the addition of Fentanyl, which is the subject of this submission.

There are no significant differences between the products and are identical in terms of intended use, materials, design and manufacturing methods. The devices meet the ASTM standard D6319-19 and D6978-05 (2019).

G. NON-CLINICAL TESTING SUMMARY PERFORMANCE DATA

BENCH TEST DATA

| TEST METHOD | PURPOSE | ACCEPTANCE CRITERIA | RESULT | | |
|--|--|--|--|---|---|
| ASTM D6319-19 Standard Specification for Nitrile Examination Gloves for Medical Application. | To determine the length of the gloves | Medium : 230 mm min Large : 230 mm min X-Large : 230 mm min XX-Large : 230 mm min | Medium : 235 mm Large : 237 mm X-Large : 250 mm XX-Large : 238 mm | | |
| ASTM D6319-19 Standard Specification for Nitrile Examination Gloves for Medical Application. | To determine the width of the gloves | Medium : 95+/-10 mm Large : 110+/-10 mm X-Large : 120+/-10 mm X-Large : 130+/-10 mm | Medium : 95 Large : 113 X-Large : 121 X-Large : 129 | | |
| ASTM D6319-19 Standard Specification for Nitrile Examination Gloves for Medical Application. | To determine the thickness of the gloves | Palm: 0.05 mm min for all sizes Finger: 0.05 mm min for all sizes | <u>Size</u> Medium Large X-Large XX-Large | <u>Palm</u> 0.077 mm 0.106 mm 0.089 mm 0.113 mm | <u>Finger</u> 0.111 mm 0.109 mm 0.115 mm 0.107 mm |

510(K) SUMMARY

K231349

AS REQUIRED BY: 21CFR§807.92

| TEST METHOD | PURPOSE | ACCEPTANCE CRITERIA | RESULT | | |
|--|--|---|---------------------------|--------------------------------------|-------------------------------------|
| ASTM D6319-19 Standard Specification for Nitrile Examination Gloves for Medical Application. | To determine the physical properties- Tensile strength | <u>Before Ageing</u> Tensile Strength 14MPa min for all sizes <u>After Ageing</u> Tensile Strength 14MPa min for all sizes | <u>Size</u> Medium | <u>Before Ageing</u> 34.0 MPa | <u>After Ageing</u> 37.3 MPa |
| | To determine the physical properties- Ultimate Elongation | <u>Before Ageing</u> Ultimate Elongation 500% min for all sizes <u>After Ageing</u> Ultimate Elongation 400% min for all sizes | <u>Size</u> Medium | <u>Before Ageing</u> 542% | <u>After Ageing</u> 503% |
| ASTM D5151-19 Standard Test Method for Detection of Holes in Medical Gloves | To determine the holes in the gloves | AQL 2.5 | Gloves Pass AQL 2.5 | | |
| ASTM D6124-06 (Reapproved 2017) Standard Test Method for Residual Powder on Medical Gloves | To determine the residual powder in the gloves | ≤ 2 mg/glove | Medium: 0.3516 mg/glove | | |

| TEST METHOD | PURPOSE | ACCEPTANCE CRITERIA | RESULT |
|---|--|--|---|
| ASTM D6978-05 (Reapproved 2019) Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs. | To determine the breakthrough detection time of chemotherapy drugs | Bleomycin Sulfate 15mg/ml (15,000 ppm) >240 Minutes | Bleomycin Sulfate 15 mg/ml (15,000 ppm) >240 Minutes |
| | | Busulfan 6 mg/ml (6,000 ppm) >240 Minutes | Busulfan 6 mg/ml (6,000 ppm) >240 Minutes |
| | | Carboplatin 10 mg/ml (10,000 ppm) >240 Minutes | Carboplatin 10 mg/ml (10,000 ppm) >240 Minutes |
| | | Cisplatin 1 mg/ml (1,000 ppm) >240 Minutes | Cisplatin 1 mg/ml (1,000 ppm) >240 Minutes |
| | | Cyclophosphamide 20 mg/ml (20,000 ppm) >240 Minutes | Cyclophosphamide 20 mg/ml (20,000 ppm) >240 Minutes |
| | | Cytarabine 100 mg/ml (100,000 ppm) >240 Minutes | Cytarabine 100 mg/ml (100,000 ppm) >240 Minutes |
| | | Dacarbazine 10 mg/ml (10,000 ppm) >240 Minutes | Dacarbazine 10 mg/ml (10,000 ppm) >240 Minutes |
| | | Daunorubicin HCl 5 mg/ml (5,000 ppm) >240 Minutes | Daunorubicin HCl 5 mg/ml (5,000 ppm) >240 Minutes |
| | | Docetaxel 10 mg/ml (10,000 ppm) >240 Minutes | Docetaxel 10 mg/ml (10,000 ppm) >240 Minutes |
| | | Doxorubicin HCl 2 mg/ml (2,000 ppm) >240 Minutes | Doxorubicin HCl 2 mg/ml (2,000 ppm) >240 Minutes |
| Epirubicin HCl 2 mg/ml (2,000 ppm) >240 Minutes | Epirubicin HCl 2 mg/ml | | |

510(K) SUMMARY

K231349

AS REQUIRED BY: 21CFR§807.92

| TEST METHOD | PURPOSE | ACCEPTANCE CRITERIA | RESULT |
|---|--|---|---|
| | | | (2,000 ppm) >240 Minutes |
| | | Etoposide 20 mg/ml (20,000 ppm) >240 Minutes | Etoposide 20 mg/ml (20,000 ppm) >240 Minutes |
| | | Fludarabine 25 mg/ml (25,000 ppm) >240 Minutes | Fludarabine 25 mg/ml (25,000 ppm) >240 Minutes |
| | | Fluorouracil 50 mg/ml (50,000 ppm) >240 Minutes | Fluorouracil 50 mg/ml (50,000 ppm) >240 Minutes |
| | | Gemcitabine 38 mg/ml (38,000 ppm) >240 Minutes | Gemcitabine 38 mg/ml (38,000 ppm) >240 Minutes |
| | | Idarubicin HCl 1 mg/ml (1,000 ppm) >240 Minutes | Idarubicin HCl 1 mg/ml (1,000 ppm) >240 Minutes |
| | | Ifosfamide 50 mg/ml (50,000 ppm) >240 Minutes | Ifosfamide 50 mg/ml (50,000 ppm) >240 Minutes |
| | | Irinotecan 20 mg/ml (20,000 ppm) >240 Minutes | Irinotecan 20 mg/ml (20,000 ppm) >240 Minutes |
| | | Mechlorethamine HCl 1 mg/ml (1,000 ppm) >240 Minutes | Mechlorethamine HCl 1 mg/ml (1,000 ppm) >240 Minutes |
| | | Melphalan 5 mg/ml (5,000 ppm) >240 Minutes | Melphalan 5 mg/ml (5,000 ppm) >240 Minutes |
| | | Methotrexate 25 mg/ml (25,000 ppm) >240 Minutes | Methotrexate 25 mg/ml (25,000 ppm) >240 Minutes |
| | | Mitomycin C 0.5 mg/ml (500 ppm) >240 Minutes | Mitomycin C 0.5 mg/ml (500 ppm) >240 Minutes |
| | | Mitoxantrone HCl 2 mg/ml (2,000 ppm) >240 Minutes | Mitoxantrone HCl 2 mg/ml (2,000 ppm) >240 Minutes |
| | | Paclitaxel 6 mg/ml (6,000 ppm) >240 Minutes | Paclitaxel 6 mg/ml (6,000 ppm) >240 Minutes |
| | | Rituximab 10 mg/ml (10,000 ppm) >240 Minutes | Rituximab 10 mg/ml (10,000 ppm) >240 Minutes |
| | | Trisenox 1 mg/ml (1,000 ppm) >240 Minutes | Trisenox 1 mg/ml (1,000 ppm) >240 Minutes |
| | | Vincristine Sulfate 1 mg/ml (1,000 ppm) >240 Minutes | Vincristine Sulfate 1 mg/ml (1,000 ppm) >240 Minutes |
| ASTM D6978-05 (Reapproved 2019) Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs. | To determine the breakthrough detection time of Opioid drugs | Fentanyl Citrate Injection (100mcg/2mL) >240 Minutes | Fentanyl Citrate Injection (100mcg/2mL) >240 Minutes |

510(K) SUMMARY

K231349

AS REQUIRED BY: 21CFR§807.92

BIOCOMPATIBILITY DATA

| TEST METHOD | PURPOSE | ACCEPTANCE CRITERIA | RESULT |
|--|---|--|--|
| ISO 10993-23 First edition 2021-01 Biological Evaluation of Medical Devices - Part 23, Tests for Irritation. | To evaluate the local dermal irritation of a test article extract following intracutaneous injection in rabbits. | Under the condition of study not an irritant | Under the conditions of the study, the test article met the requirements of the test |
| 10993-10 Fourth edition 2021-11 Biological Evaluation of Medical Devices - Part 10, Tests for Skin Sensitization. | To evaluate the test item, for the skin sensitization in Guinea pigs by maximization test. | Under the conditions of the study, not a sensitizer | Under the conditions of the study, the test article was not considered a sensitizer |
| ISO 10993-5 Third edition 2009-06-01 Biological Evaluation of Medical Devices - Part 5, Tests for In Vitro Cytotoxicity. | To determine the potential of a test article to cause cytotoxicity | Under the conditions of the study, non-cytotoxic | The undiluted test article extract and 50% test article extract dilution did not meet the requirements of the test and the 25%, 12.5%, 6.25%, and 3.13% test article extract dilutions met the requirements of the test. Cytotoxicity concern was addressed by acute systemic toxicity testing. |
| ISO 10993-11 Third edition 2017-09 Biological Evaluation of Medical Devices - Part 11, Tests for Systemic Toxicity. | To evaluate the acute systemic toxicity of a test article extract following injection in mice. | Under the conditions of study, the device extracts do not pose a systemic toxicity concern | Under the conditions of study, there was no mortality or evidence of systemic toxicity. |

The performance test data of the non-clinical tests meet following standards:

ASTM D6319-19 Standard Specification for Nitrile examination Gloves for Medical Application.

ASTM D5151-19 Standard Test Method for Detection of Holes in Medical Gloves.

ASTM D6124-06 (Reapproved 2017) Standard Test Method for Residual Powder on Medical Gloves.

ASTM D6978-05 (Reapproved 2019) Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs.

510(K) SUMMARY

K231349

AS REQUIRED BY: 21CFR§807.92

ISO 10993-23 First Edition 2021-01 Biological Evaluation of Medical Devices - Part 23, Tests for Irritation.

ISO 10993-10 Fourth Edition 2021-11 Biological Evaluation of Medical Devices - Part 10, Tests for Skin Sensitization.

ISO 10993-5 Third Edition 2009-06-01 Biological Evaluation of Medical Devices - Part 5, Tests for In Vitro Cytotoxicity.

ISO 10993-11 Third Edition 2017-09 Biological Evaluation of Medical Devices - Part 11, Tests for Systemic Toxicity.

H. CLINICAL TESTING SUMMARY

Not applicable - Clinical data is not needed for gloves.

I. CONCLUSION

The conclusions drawn from the non-clinical test demonstrate that the subject device in 510(K) K231349 submission, Nephron Nitrile Powder-Free Nitrile Examination Gloves (Tested For Use With Chemotherapy Drugs and Fentanyl) is as safe, as effective, and performs as well as or better than the legally marketed predicate device K223559.