

August 11, 2023

Shandong Intco Medical Products Co, Ltd. % Deze Wang
Official Correspondent
Intco Medical Industries, Inc
805 Barrington Ave.
Ontario, California 91764

Re: K231408

Trade/Device Name: Basic Synguard Nitrile Exam Gloves, Powder Free, Blue Colored, Non-Sterile

Regulation Number: 21 CFR 880.6250

Regulation Name: Non-Powdered Patient Examination Glove

Regulatory Class: Class I, reserved

Product Code: LZA Dated: May 11, 2023 Received: May 15, 2023

Dear Deze Wang:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

K231408 - Deze Wang Page 2

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance) and CDRH Learn (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,



BiFeng Qian, M.D., Ph.D
Assistant Director
DHT4B: Division of Infection Control
and Plastic Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

510(k) Number (if known)

Form Approved: OMB No. 0910-0120
Expiration Date: 06/30/2023

Expiration Date: 06/30/2023 See PRA Statement below.

K231408				
Device Name				
Basic Synguard Nitrile Exam Gloves, Powder Free, Blue Colored, Non-Sterile				
Indications for Use (Describe)				
A patient examination glove is a disposable device intended for medical purposes that is worn upon the examiner's hands or fingers to prevent contamination between patient and examiner.				
Type of Use (Select one or both, as applicable)				
☐ Prescription Use (Part 21 CFR 801 Subpart D) ☐ Over-The-Counter Use (21 CFR 801 Subpart C)				
CONTINUE ON A SEPARATE PAGE IF NEEDED				

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



510(K) Summary

K231408

(As requirement by 21 CFR 807.92)

K231408

Date: August 8, 2023

A. Applicant:

Name: Shandong Intco Medical Products Co., Ltd.

Address: Qiwang Road, Naoshan Industrial Park, Qingzhou, Shandong, China

Contact: Ms. Min Gao

Title: QA

Tel: +86-536-6136888
Email: gaomin@intco.com

Submission Correspondent: Primary contact: Max Li

Shandong Intco Medical Products Co., Ltd.

Qiwang Road, Naoshan Industrial Park, Qingzhou, Shandong, China

Tel: +86-536-6136888
Email: maxli@intco.com
Secondary contact: Min Gao

Qiwang Road, Naoshan Industrial Park, Qingzhou, Shandong, China

Tel: +86-536-6136888 Email: gaomin@intco.com

B. Device:

Trade Name: Basic Synguard Nitrile Exam Gloves, Powder Free, Blue Colored, Non-Sterile.

Common Name: Non-Powdered Patient Examination Glove

Regulatory Information

Classification Name: Non-Powdered Patient Examination Glove

Classification: Class I Product code: LZA

Regulation Number: 21 CFR 880.6250

Review Panel: General Hospital and Personal Use

C. Predicate device:

K221169

Name of the device: Nitrile Examination Gloves Powder Free-Black Owner's Name: JR Engineering & Medical Technologies (M) Sdn. Bhd.

Classification: Class I Product code: LZA

Regulation Number: 21 CFR 880.6250

D. Indications for use of the device:

The Basic Synguard Nitrile Exam Gloves, Powder Free, Blue Colored, Non-Sterile is a disposable device intended for medical purposes that is worn on the examiners' hand or fingers to prevent contamination between patient and examiner.

E. Device Description:

The subject device is a patient examination glove made of synthetic nitrile latex compound. It is Non-Sterile, Powder free and is Blue in color. The device is ambidextrous and can be worn on either the left or right hand. The device meets all the specifications in ASTM D6319-19, Standard specification for Nitrile Examination Gloves. Additionally, the gloves have been tested for biocompatibility per ISO 10993-10, ISO 10993-11.

F. Summary of Technological Characteristics

The Basic Synguard Nitrile Exam Gloves, Powder Free, Blue Colored, Non-Sterile are compared with the predicate device in terms of intended use, design, material, specifications, and performance.

The following table shows similarities and differences of use, design, and material between the proposed device and the predicate device.

Table 1 General Comparison of Proposed and Predicate Devices

Froduct Name K231408 K221169 Product Name Basic Synguard Nitrile Exam Gloves, Powder Free, Blue Colored, Non-Sterile JR MEDIC Nitrile Examination Gloves Powder Free- Black	-
Powder Free, Blue Colored, Non- Powder Free- Black	-
	-
Sterile	-
	-
Manufacturer(s) Shandong Intco Medical Products Co., JR Engineering & Medical	
Ltd. Technologies (M) SDN.BHD.	
Product Code LZA LZA	Same
Regulation 21 CFR 880.6250 21 CFR 880.6250	Same
Number	
Indications for The Basic Synguard Nitrile Exam Nitrile Examination Gloves Powder	Same
use Gloves, Powder Free, Blue Colored, free is a disposable device intended	
Non-Sterile is a disposable device for medical purpose that is worn on	
intended for medical purposes that is the examiner's hand to prevent	
worn on the examiners' hand or finger contamination between patient and	
to prevent contamination between examiner.	
patient and examiner.	
Powder free Yes Yes	Same
Design feature Ambidextrous. Ambidextrous.	Same
Material Nitrile Nitrile	Same
Color Blue Black	Different
Size XS, S, M, L, XL, XXL XS, S, M, L, XL	Similar
Sterile Non-Sterile Non-Sterile	Same
Use Singe use Single use	Same

Analysis: The proposed device has product size XS, S, M, L, XXL while the predicate device has product size XS, S, M, L, XL. The proposed device is Blue color, while the predicate device is black color. But safety and performance testing has been done to the proposed device and the results showed that the device meets the requirements of standard ASTM D6319-19. Therefore, this difference do not raise any new safety or performance questions.

Table 2 Comparison of Physical, Biocompatibility and Performance Testing

CHARACTERSTICS	STANDARDS	DEVICE PERFORMANCE		Remarks
		SUBJECT	PREDICATE	
510(k) Number		K231408	K221169	
Name of device		Basic Synguard Nitrile	Nitrile Examination	
		Exam Gloves, Powder	Gloves Powder Free-	
		Free, Blue Colored, Non-	Black	
		Sterile		
Dimensions	ASTMD6319- 19	Length Min 230 m	Length Min 230 mm	Same
		Width Min 95+/- 10 mm	Width Min 95+/- 10	
		(for medium size)	mm	
			(for medium size)	
Physical Properties	ASTMD6319- 19	Before Aging	Before Aging	Same
		Tensile Strength	Tensile Strength	
		min 14 MPa	min 14 MPa	
		Ultimate Elongation	Ultimate Elongation	
		Min 500%	Min 500%	
		After Aging	After Aging	
		Tensile Strength	Tensile Strength	
		min 14 MPa	min 14 MPa	
		Ultimate Elongation Min 400%	Ultimate Elongation Min 400%	
Thickness	ASTMD6319- 19	Palm min 0.05 mm	Palm min 0.05 mm	Same
		Finger min 0.05 mm	Finger min 0.05 mm	
Powder Free	ASTMD6319- 19	≤2 mg/glove	≤2 mg/glove	Similar
Water leak	ASTM D5151-19	3/315 AQL2.5	Meets ASTM D6319 requirements	Similar
	Primary Skin	Under the condition of	Under the condition	Same
	Irritation-ISO	study not an irritant	of study not an	
	10993-10:2010		irritant	
	Dermal	Under the conditions of	Under the conditions of	Same
Biocompatibility	Sensitization- ISO	the study not a sensitizer	the study not a sensitizer	
	10993- 10:2010			

In vitro cytotoxicity	Data Not available	Under the conditions	
ISO10993-5 :2009		of the study, cytotoxic.	
Acute Systemic Toxicity Test ISO	Under the conditions of this study, there was no evidence of systemic toxicity.		Same
10993- 11:2017	or systemic toxicity.	toxicity.	

Analysis: The proposed device has the Skin Sensitization test, Primary Skin Irritation test and the Acute Systemic Toxicity test. The result is the same as the PREDICATE device's.

G. Summary of Non-Clinical Testing

Biocompatibility

The following tests for the subject device were conducted to evaluate the biocompatibility of Basic Synguard Nitrile Exam Gloves, Powder Free, Blue Colored, Non-Sterile:

- ISO 10993-10: Skin Sensitization test
- ISO 10993-10: Intracutaneous Reactivity test
- ISO 10993-11: Acute Systemic Toxicity test

> Performance Testing

Physical performance testing of the proposed device was conducted as per ASTM D6319-19 *Standard Specification for Nitrile Examination Gloves*.

To summarize, the performance testing of the subject device was conducted to adequately demonstrate the effectiveness of the device in accordance with the relevant test methods cited below:

- ASTM D6319-19 Standard Specification for Nitrile Examination Gloves.
- · ASTM D6124-06 (Reapproved 2017) Standard Test Method for Residual Powder on Medical Gloves
- ASTM D5151-19 Standard Test Method for Detection of Holes in Medical Gloves

Test Method	Purpose	Acceptance Criteria	Results
Dimensions (length)	The purpose of the	Length	Pass
(width)	test is to evaluate	230 mm min	240 mm min
(thickness)	the physical	Width(mm)	Pass
	dimension of the	XS: 70±10	XS: average 78.3mm
	glove	S: 80±10	S: average 86.5mm
		M: 95±10	M: average 97.6mm
		L: 110±10	L: average 108.7mm
		XL: 120±10	XL: average 115.4mm
		XXL: 130±10	XXL: average 123.7mm
		Thickness(mm):	Pass
		Palm:_Minimum 0.05	Palm - 0.06mm min.

		Finger: Minimum 0.05	Finger - 0.082mm min
Physical properties	The purpose of the	Before Aging:	Pass
	test is to evaluate	Tensile Strength: 14 MPa, min.	Before Aging:
	the tensile strength	Elongation: 500%, min.	Tensile Strength: 22.6MPa, min.
	and ultimate	After Aging:	Elongation: 531%, min.
	elongation before	Tensile Strength: 14 MPa, min.	After Aging:
	and after aging	Elongation: 400%, min.	Tensile Strength: 21.9MPa, min.
			Elongation: 437%, min.
Freedom from holes	The purpose of the	No leakage at sampling level of	Pass
	test is to detect	G-1, AQL 2.5	No leakage, 312 of 315 passed
	holes in the gloves		
Residual Powder	The purpose of the	<2mg per glove	Pass
	test is to detect the		average 0.17 mg per glove
	powder residue in		
	the glove		

H. Clinical Test Conclusion

No clinical study is included in this submission.

I. Conclusion

The conclusion drawn from the nonclinical tests demonstrates that the subject device in 510(K) submission, the Basic Synguard Nitrile Exam Gloves, Powder Free, Blue Colored, Non-Sterile is as safe, as effective, and performs as well as or better than the legally marketed predicate device cleared under K221169.