

August 10, 2023

Abbott Medical Alyssa Timmers Senior Regulatory Affairs Specialist One St. Jude Medical Drive St. Paul, Minnesota 55117

Re: K231415

Trade/Device Name: EnSite[™] X EP System Regulation Number: 21 CFR 870.1425 Regulation Name: Programmable Diagnostic Computer Regulatory Class: Class II Product Code: DQK Dated: May 15, 2023 Received: May 16, 2023

Dear Alyssa Timmers:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Aneesh S. Deoras -S

Aneesh Deoras Assistant Director Division of Cardiac Electrophysiology, Diagnostics and Monitoring Devices Office of Cardiovascular Devices Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number *(if known)* K231415

Device Name EnSite[™] X EP System

Indications for Use (Describe)

EnSite[™] X EP System

The EnSite[™] X EP System is a suggested diagnostic tool in patients for whom electrophysiology studies have been indicated.

The EnSite[™] X EP System provides information about the electrical activity of the heart and displays catheter location during conventional electrophysiological (EP) procedures.

EnSiteTM X EP System Contact Force Software License:

When used with the TactiSysTM Quartz Equipment, the EnSiteTM X EP System Contact Force Module is intended to provide visualization of force information from compatible catheters.

EnSite[™] X EP System Surface Electrode Kit:

The EnSiteTM X EP Surface Electrode Kit is indicated for use with the EnSiteTM X EP System in accordance with the EnSiteTM X EP System indications for use. The EnSiteTM X EP System TactiFlexTM Ablation Catheter, Sensor EnabledTM Software Module is indicated for use with the EnSiteTM X EP System in accordance with the EnSiteTM X EP System indications for use.

EnSiteTM X EP System, TactiFlexTM Ablation Catheter, Sensor EnabledTM, Software Upgrade and EnSiteTM X EP System, TactiFlexTM Ablation Catheter, Sensor EnabledTM, Software License:

The EnSiteTM X EP System TactiFlexTM Ablation Catheter, Sensor EnabledTM Software Module is indicated for use with the EnSiteTM X EP System in accordance with the EnSiteTM X EP System indications for use.

Type of Use (Select one or both, as applicable)		
☑ Prescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)	

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov*

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



510(k) Information	
510(k) Number	K231415
510(k) Type	Traditional 510(k)
Date Prepared	15 May 2023
Submitter Information	
Manufacturer Name & Address	Abbott Medical One St. Jude Medical Drive, St. Paul, Minnesota, 55117, USA Manufacturer of the EnSite X EP System
	Abbott Medical Costa Rica Ltda. Edificio #44 Calle 0, Ave. 2 Zona Franca Coyol El Coyol Alajuela, Costa Rica 1897-4050 Manufacturer of the EnSite X EP System Surface Electrode Kit
Contact Person	Alyssa Timmers Senior Regulatory Affairs Specialist 651-756-3706 <u>alyssa.timmers@abbott.com</u>
EnSite™ X EP System	Device Information
Trade Name	EnSite™ X EP System
Common Name	Programmable Diagnostic Computer
Class	II
Classification Name	870.1425, computer, diagnostic, programmable
Product Code	DQK
Predicate Device	EnSite™ X EP System (K223094)
Device Description	The EnSite [™] X EP System is a catheter navigation and mapping system. A catheter navigation and mapping system is capable of displaying the 3-dimensional (3-D) position of conventional and Sensor Enabled [™] (SE) electrophysiology catheters, as well as displaying cardiac electrical activity as waveform traces and as three-dimensional (3D) isopotential and isochronal maps of the cardiac chamber.
	cardiac chamber. The system creates a model by collecting and labeling the anatomic locations within the chamber. A surface is created by moving a selected catheter to locations within a cardiac structure. As the catheter moves, points are collected at and between all electrodes on the catheter. A surface is wrapped around the outermost points.



Indications for Use	EnSite™ X EP Svstem
	The EnSite™ X EP System is a suggested diagnostic tool in patients for whom
	electrophysiology studies have been indicated.
	The EnSite [™] X EP System provides information about the electrical activity of the heart and displays catheter location during conventional electrophysiological (EP) procedures.
	En Site TM X ED System Contact Force Software License
	Ensite I''' X EP System Contact Force Software License
	Force Module is intended to provide visualization of force information from compatible catheters.
	EnSite™ X EP System Surface Electrode Kit
	The EnSite™ X EP Surface Electrode Kit is indicated for use with the EnSite™ X EP
	System in accordance with the EnSite™ X EP System indications for use.
	EnSite™ X EP System, TactiFlex™ Ablation Catheter, Sensor Enabled™, Software Upgrade and Software License
	The EnSite™ X EP System TactiFlex™ Ablation Catheter, Sensor Enabled™ Software Module is indicated for use with the EnSite™ X EP System in accordance with the EnSite™ X EP System indications for use.
Predicate Comparison	
Comparison	The subject device, EnSite [™] X EP System v3.0, and the predicate device, EnSite [™] X EP System v2.0.1 have the same intended use and indications for use. They use the same fundamental scientific technology to facilitate catheter position and orientation, as well as cardiac mapping and model creation. There were no changes to the hardware. The subject device software was revised to include:
	 A display of the distance between two AutoMarks, or an AutoMark and the projected ablation catheter distal electrode,
	 A duplicate selection method to allow for the display of map points based on the highest frequency,
	 An expansion of the data able to be imported and exported between the subject device and third-party systems,
	 A workflow enhancement that allows physicians to change navigation modes without restarting a study,
	The introduction of a new AutoMark metric,
	A minor update to the Force Direction Indicator,
	A minor update to the Sandpaper tool, and
	Fixes to minor known software issues.
	All risks associated with these modifications were mitigated to acceptable levels. No new questions of safety or effectiveness were raised.
Non-Clinical Testing Summary	Design verification activities were performed and met the respective acceptance criteria to ensure that the devices in scope of this submission are safe and effective.
	Testing The EnSite™ X EP System v3.0 was developed and tested in accordance with the following industry guidance documents and standards:
	 Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices
	 IEC 62304:2015-06 Edition 1.1, Medical Device Software - Software Life Cycle Processes
	 ISO 14971:2019 Medical Devices – Application of Risk Management to Medical Devices



Non-Clinical Testing Summary (Continued)	 ANSI AAMI IEC 62366-1:2015, Medical devices – Part 1: Application of usability engineering to medical devices
	Types of Testing Performed
	Software Verification at unit, software and system level
	Performance testing of updated feature functionality
	 Preclinical Validation Testing to confirm the system could meet user requirements and its intended use after modifications
	 Human Factors Evaluations to confirm the user interface of the subject device can be used as intended by the defined user groups
Statement of Equivalence	The subject and predicate devices have the same intended use, and the same indications for use. The devices operate using the same fundamental scientific technology to facilitate catheter position and orientation, as well as cardiac mapping and model creation. The testing completed and submitted in this Traditional 510(k) provides objective evidence the subject device is at least as safe and effective as the predicate device.