

August 7, 2023

STERIS Corporation Gregory Land Senior Regulatory Affairs Specialist 5960 Heisley Rd Mentor, Ohio 44060

#### Re: K231490

Trade/Device Name: Celerity 20 HP Biological Indicator; VERIFY V24 Self-Contained Biological Indicator Regulation Number: 21 CFR 880.2800 Regulation Name: Sterilization Process Indicator Regulatory Class: Class II Product Code: FRC Dated: May 22, 2023 Received: May 23, 2023

Dear Gregory Land:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <a href="https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm">https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm</a> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <a href="https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products">https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products</a>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely, Christopher K. Dugard -S

for Clarence W. Murray, III, Ph.D. Assistant Director DHT4B: Division of Infection Control and Plastic Surgery Devices OHT4: Office of Surgical and Infection Control Devices Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

### Indications for Use

Form Approved: OMB No. 0910-0120 Expiration Date: 06/30/2023 See PRA Statement below.

510(k) Number (*if known*) K231490

#### Device Name

VERIFY V24 Self-Contained Biological Indicator

Indications for Use (Describe)

The VERIFY® V24 Self-Contained Biological Indicator is intended for routine monitoring of the following sterilizer cycles:

• Lumen, Non Lumen, Flexible, Fast, Fast Non Lumen and Specialty Cycles of the V-PRO® 1, 1 Plus, maX, maX 2, 60 and s2 Low Temperature Sterilization Systems

• STERRAD® 100S Sterilizer (Default Cycle)

• STERRAD® 200 Sterilizer (Default Cycle)

• Standard and Advanced Cycles of the STERRAD® NX Sterilizer to include sterilizers with ALLClear® Technology

• Express, Flex Scope and Standard Cycles of the STERRAD® 100NX Sterilizer to include sterilizers with ALLClear® Technology.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

#### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

#### \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov* 

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

### Indications for Use

510(k) Number (*if known*) K231490

Device Name Celerity 20 HP Biological Indicator

The Celerity 20 HP Biological indicator is intended for the routine monitoring of the following sterilizer cycles:

• Lumen, Non Lumen, Fast Non Lumen, Fast, Flexible and Specialty Cycles of the V-PRO 1, 1 Plus, maX, maX 2, 60 and s2 Low Temperature Sterilization Systems

• STERRAD® 100S Sterilizer (Default Cycle)

• Standard and Advanced Cycles of the STERRAD® NX Sterilizer with or without ALLClear® Technology

• Standard, Flex Scope, Express, and DUO Cycles of the STERRAD® 100NX Sterilizer with or without ALLClear® Technology

When used in conjunction with the Celerity HP Incubator or the Celerity Incubator, the Celerity HP BI provides a fluorescent result within 20 minutes.

Type of Use (Select one or both, as applicable)	
---	--

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

#### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

#### \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov* 

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



# 510(K) Summary For K231490 Celerity 20 HP Biological Indicator

#### **Sponsor Facility**

STERIS Corporation 5960 Heisley Road Mentor, OH 44060 Phone: (440) 354-2600 Fax: (440)357-9198

#### **Manufacturing Facility**

 STERIS Corporation

 9325 Pinecone Drive

 Mentor, OH 44060

 Phone: (440) 392-7800

 Fax: (440) 392-7896

#### Contact

Gregory Land Lead Regulatory Affairs Specialist Phone: (440) 352-7424 Fax: (440) 357-9198 Email: greg land@steris.com

#### **Submission Date:**

May 22, 2023

#### 1 Device Name

Trade Name:CELERITY 20 HP Biological IndicatorCommon/usual Name:Biological IndicatorDevice Classification:Class IIClassification Name:Indicator, Biological Sterilization Process [21 CFR 880.2800(a), FRC]

### 2 Predicate Device

Proprietary Name Common/usual Name Classification Name: 510(k) Submitter/Holder 510(k) Number: Celerity 20 HP Biological Indicator Biological indicator Indicator, Biological Sterilization Process STERIS Corporation K183294

## 3 Description of Device

The product is intended to monitor the vapor phased hydrogen peroxide sterilization cycles described in the indication for use. It produces an optical change in signal that is detected by the STERIS proprietary reader, Celerity 20 HP Incubator, within 20 minutes to confirm the viability of the biological indicator at the end of a sterilization process. The product consists of *Geobacillus stearothermophilus* spores and a defined nutrient media in a plastic vial. A reporter enzyme, which is produced by the organism, reacts with a fluorogenic substrate within the defined nutrient media to produce a fluorescent moiety.

### 4 Intended Use/Indications for Use

The Celerity 20 HP Biological indicator is intended for the routine monitoring of the following sterilizer cycles:

- Lumen, Non Lumen, Fast Non Lumen, Fast, Flexible and Specialty Cycles of the V-PRO 1, 1 Plus, maX, maX 2, 60 and s2 Low Temperature Sterilization Systems
- STERRAD<sup>®</sup> 100S Sterilizer (Default Cycle)
- Standard and Advanced Cycles of the STERRAD<sup>®</sup> NX Sterilizer with or without ALLClear<sup>®</sup> Technology
- Standard, Flex Scope, Express, and DUO Cycles of the STERRAD<sup>®</sup> 100NX Sterilizer with or without ALLClear<sup>®</sup> Technology

When used in conjunction with the Celerity HP Incubator or the Celerity Incubator, the Celerity HP BI provides a fluorescent result within 20 minutes.

### 5 <u>Summary of Technical Characteristics</u>

5.1

Feature	CELRITY 20 HP Biological Indicator (proposed)	Celerity 20 HP Biological Indicator (K183294)	Comparison
Intended Use	The Celerity 20 HP Biological Indicator is intended for	The Celerity 20 HP Biological Indicator is intended for	The proposed and predicate devices are identical. The

## STERIS TRADITIONAL 510(k) PREMARKET NOTIFICATION Celerity 20 HP Biological Indicator

Feature	CELRITY 20 HP Biological Indicator (proposed)	Celerity 20 HP Biological Indicator (K183294)	Comparison
	<ul> <li>routine monitoring of the following sterilizer cycles:</li> <li>Lumen, Non Lumen, Fast Non Lumen, Fast, Flexible and Specialty Cycles of the V-PRO 1, 1 Plus, maX, maX 2, 60 and s2 Low Temperature Sterilization Systems.</li> <li>STERRAD® 100S Sterilizer (Default Cycle)</li> <li>Standard and Advanced Cycles of the STERRAD® NX Sterilizer with or without ALLClear® Technology</li> <li>Standard, Flex Scope, Express and DUO Cycles of the STERRAD® 100NX Sterilizer with or without ALLClear Technology</li> <li>When used in conjunction with the Celerity HP Incubator or the Celerity Incubator, the Celerity 20 HP BI provides a fluorescent result within 20 minutes.</li> </ul>	<ul> <li>routine monitoring of the following sterilizer cycles:</li> <li>Lumen, Non Lumen, Fast Non Lumen, Fast and Flexible Cycles of the V-PRO: 1, 1 Plus, maX, maX2 and s2 Low Temperature Sterilization Systems.</li> <li>STERRAD® 100S Sterilizer (Default Cycle)</li> <li>Standard and Advanced Cycles of the STERRAD® NX Sterilizer with or without ALLClear</li> <li>Standard, Flex Scope, Express and DUO Cycles of the STERRAD® 100NX Sterilizer with or without ALLClear</li> <li>When used in conjunction with the Celerity HP Incubator, the Celerity 20 HP BI provides a fluorescent result within 20 minutes.</li> </ul>	Specialty Cycle is a new cycle in the V- PRO maX 2 Low Temperature Sterilizer, which has been submitted in this premarket notification.
Indicator organism	Geobacillus stearothermophilus	Geobacillus stearothermophilus	Same
Mechanism of action	An enzyme, which is produced by the organism, reacts with a fluorogenic substrate within the defined nutrient media to produce a fluorescent moiety.	An enzyme, which is produced by the organism, reacts with a fluorogenic substrate within the defined nutrient media to produce a fluorescent moiety.	Same
Accessories	Automated incubator / reader	Automated incubator / reader	Same
Viable spore population	1.0 – 4.0 x 10 <sup>6</sup> spore/BI	1.0 – 4.0 x 10 <sup>6</sup> spore/BI	Same
Resistance	Resistance @ 9.1 mg/L H <sub>2</sub> O <sub>2</sub> : • <u>D-value <math>\geq</math> 6 sec • <u>Survival Time</u> <math>\geq</math> 4 sec • <u>Kill Time</u> <math>\leq</math> 7 min</u>	Resistance @ 9.1 mg/L H <sub>2</sub> O <sub>2</sub> : • <u>D-value <math>\geq</math> 6 sec • <u>Survival Time</u> <math>\geq</math> 4 sec • <u>Kill Time</u> <math>\leq</math> 7 min</u>	Same

#### **STERIS TRADITIONAL 510(k) PREMARKET NOTIFICATION** Celerity 20 HP Biological Indicator

Feature	CELRITY 20 HP Biological Indicator (proposed)	Celerity 20 HP Biological Indicator (K183294)	Comparison
Culture Conditions	55- 59°C, media included in BI, 20-minute incubation time.	55- 59°C, media included in BI, 20-minute incubation time.	Same
Primary Packaging	Direct inoculum on plastic vial, cap with recovery media.	Direct inoculum on plastic vial, cap with recovery media.	Same
Process indicator	VERIFY V-PRO Chemical Indicator (K140515); magenta to yellow color change.	VERIFY V-PRO Chemical Indicator (K140515); magenta to yellow color change.	Same
Shelf-life	10 months	10 months	Same

## 6 <u>Summary of Non-clinical Tests</u>

Performance testing to demonstrate substantial equivalence to the predicate has been completed and is summarized in Table 5-3 below.

Table 5-3. Summary	of Non-clinical Testing
--------------------	-------------------------

Test	Acceptance Criteria	Conclusion
Verification of Celerity HP BI	Celerity BIs exposed to the Specialty	Pass
Performance after V-PRO Specialty	Cycle demonstrate essentially	
Cycle Extended Aeration	equivalent or not significantly	
	different resistance as compared to	
	SCBIs exposed to the 136L V-PRO	
	Sterilizer Fast Non Lumen Cycle	
Final Process Qualification of the	The V-PRO maX 2 Sterilizer	Pass
VPRO maX2 Sterilizer Specialty	Specialty Cycle final process	
Cycle	qualification was successful.	

## 7 <u>Conclusion</u>

Based on the intended use, technological characteristics and the non-clinical tests performed, the subject device is as safe, as effective, and performs at least as safely and effectively as the legally marketed predicate device, the Celerity 20 HP Biological Indicator, K183294 Class II [21 CFR 880.2800(a)], product code FRC.



# 510(K) Summary For K231490 VERIFY V24 Self-Contained Biological Indicator

#### **Sponsor Facility**

STERIS Corporation 5960 Heisley Road Mentor, OH 44060 Phone: (440) 354-2600 Fax: (440)357-9198

#### **Manufacturing Facility**

 STERIS Corporation

 9325 Pinecone Drive

 Mentor, OH 44060

 Phone: (440) 392-7800

 Fax: (440) 392-7896

#### Contact

Gregory Land Lead Regulatory Affairs Specialist Phone: (440) 352-7424 Fax: (440) 357-9198 Email: greg\_land@steris.com

#### **Submission Date:**

May 22, 2023

#### **STERIS TRADITIONAL 510(k) PREMARKET NOTIFICATION VERIFY V24 Self-Contained Biological Indicator**

#### 8 Device Name

Trade Name:	VERIFY V24 Self-Contained Biological Indicator
Common/usual Name:	Biological Indicator
Device Classification:	Class II
Classification Name:	Indicator, Biological Sterilization Process [21 CFR 880.2800(a), FRC]

### 9 Predicate Device

Proprietary Name:VERIFY V24 Self-Contained Biological IndicatorCommon/usual NameBiological indicatorClassification Name:Indicator, Biological Sterilization Process510(k) Submitter/Holder:STERIS Corporation510(k) Number:K183300

### 10 Description of Device

The VERIFY V24 Self-Contained Biological Indicator (SCBI) is used by healthcare facilities to monitor the V-PRO<sup>®</sup> Low Temperature Sterilization Systems. It is designed to accompany medical devices placed in the sterilizer.

The user places the VERIFY V24 SCBI into the V-PRO<sup>®</sup> Low Temperature Sterilization System and performs a sterilization cycle. After cycle completion, the SCBI can either be immediately activated or it can be held at room temperature for a maximum of 72 hours (3 days) prior to activation.

The SCBI is activated by sealing the vial and rupturing the medium ampoule using the STERIS VERIFY SCBI HP activator. The activator automatically seals the SCBI vial and releases the growth media.

The activated SCBI is incubated at 55-60°C for  $\geq 24$  hours. The SCBI indicates a pass if the media remains orange and non-turbid. The SCBI indicates a failure of sterilization if the media changes from orange to yellow and/or if the media is turbid.

### 11 Intended Use/Indications for Use

The VERIFY<sup>®</sup> V24 Self-Contained Biological Indicator is intended for routine monitoring of the following sterilizer cycles:

- Lumen, Non Lumen, Flexible, Fast, Fast Non Lumen and Specialty Cycles of the V-PRO<sup>®</sup> 1, 1 Plus, maX, maX 2, 60 and s2 Low Temperature Sterilization Systems
- STERRAD<sup>®</sup> 100S Sterilizer (Default Cycle)
- STERRAD<sup>®</sup> 200 Sterilizer (Default Cycle)
- Standard and Advanced Cycles of the STERRAD<sup>®</sup> NX Sterilizer to include sterilizers with ALLClear<sup>®</sup> Technology
- Express, Flex Scope and Standard Cycles of the STERRAD<sup>®</sup> 100NX Sterilizer to include sterilizers with ALLClear<sup>®</sup> Technology

### 12 <u>Summary of Technical Characteristics</u>

Feature	VERIFY V24 SCBI	VERIFY V24 SCBI Predicate	Comparison
reature	(proposed)	(K183300)	-
Intended Use	<ul> <li>The VERIFY® V24 Self- Contained Biological Indicator is intended for routine monitoring of the following sterilizer cycles:</li> <li>Lumen, Non Lumen, Flexible, Fast, Fast Non Lumen and Specialty Cycles of the V-PRO<sup>®</sup> 1, 1 Plus, maX, maX 2, 60 and s2 Sterilizers</li> <li>STERRAD<sup>®</sup> 100S Sterilizer (Default Cycle)</li> <li>STERRAD<sup>®</sup> 200 Sterilizer (Default Cycle)</li> <li>STERRAD<sup>®</sup> 200 Sterilizer (Default Cycle)</li> <li>STERRAD<sup>®</sup> NX Sterilizer to include sterilizers with ALLClear® Technology</li> <li>Express, Flex Scope and Standard Cycles of the STERRAD<sup>®</sup> 100NX Sterilizer to include sterilizers with ALLClear® technology.</li> </ul>	<ul> <li>The VERIFY® V24 Self- Contained Biological Indicator is intended for routine monitoring of the following sterilizer cycles:</li> <li>Lumen, Non Lumen, Flexible, Fast and Fast Non Lumen Cycles of the V- PRO® 1, 1 Plus, maX, maX 2, 60 and s2 Sterilizers</li> <li>STERRAD® 100S Sterilizer (Default Cycle)</li> <li>STERRAD® 200 Sterilizer (Default Cycle)</li> <li>Standard and Advanced Cycles of the STERRAD® NX Sterilizer to include sterilizers with ALLClear® Technology</li> <li>Express, Flex Scope and Standard Cycles of the STERRAD® 100NX Sterilizer to include sterilizers with ALLClear® technology.</li> </ul>	The proposed and predicate devices are identical. The Specialty Cycle is a new cycle in the V-PRO maX 2 Low Temperature Sterilizer, which has been submitted in this premarket notification.
organism	stearothermophilus	stearothermophilus	Same
Mechanism of action	Visual detection of growth based on media color change in the presence of surviving indicator organisms.	Visual detection of growth based on media color change in the presence of surviving indicator organisms.	Same
Accessories	VERIFY Incubator and VERIFY SCBI HP Activator (optional)	VERIFY Incubator and VERIFY SCBI HP Activator (optional)	Same
Viable spore population	2.0 – 3.4 x 10 <sup>6</sup> spore/BI	2.0 – 3.4 x 10 <sup>6</sup> spore/BI	Same
Resistance characteristics	Resistance @ 2.7 mg/L H <sub>2</sub> O <sub>2</sub> : • <u>D-value</u> $4.0 - 8.0$ sec • <u>Survival Time</u> $4 - 30$ sec • <u>Kill Time</u> $\leq 16$ min	Resistance @ 2.7 mg/L H <sub>2</sub> O <sub>2</sub> :         • <u>D-value</u> 4.0 - 8.0 sec         • <u>Survival Time</u> 4 - 30 sec         • <u>Kill Time</u> ≤ 16 min	Same
Culture Conditions	55- 60°C, media included in SCBI, 24 hour incubation time.	55- 60°C, media included in SCBI, 24 hour incubation time.	Same
Primary Packaging	Direct inoculum on plastic vial, glass ampoule with recovery media.	Direct inoculum on plastic vial, glass ampoule with recovery media.	Same
Process indicator	VERIFY V24 Self-Contained Biological Indicator Vail	VERIFY V24 Self-Contained Biological Indicator Vail	Same

### 12.1 <u>VERIFY V24 Self-Contained Biological Indicator</u>

Feature	VERIFY V24 SCBI (proposed)	VERIFY V24 SCBI Predicate (K183300)	Comparison
	Label; magenta to	Label; magenta to	
	orange/yellow color change.	orange/yellow color change.	
Shelf-life	15 Months	15 Months	Same

## 13 <u>Summary of Non-clinical Tests</u>

Performance testing to demonstrate substantial equivalence to the predicate has been completed and is summarized in Table 5-3 below.

Test	Acceptance Criteria	Conclusion
Verification of V24 SCBI	SCBIs exposed to the Specialty Cycle	Pass
Performance after V-PRO Specialty	demonstrate essentially equivalent or	
Cycle Extended Aeration	not significantly different resistance as compared to SCBIs exposed to the	
	136L V-PRO Sterilizer Fast Non	
	Lumen Cycle	
Final Process Qualification of the	The V-PRO maX 2 Sterilizer	Pass
VPRO maX2 Sterilizer Specialty	Specialty Cycle final process	
Cycle	qualification was successful.	

Table 5-3. Summary of Non-clinical Testing

## 14 Conclusion

Based on the intended use, technological characteristics and the non-clinical tests performed, the subject device is as safe, as effective, and performs at least as safely and effectively as the legally marketed predicate device, the VERIFY V24 Self-Contained Biological Indicator, K183300 Class II [21 CFR 880.2800(a)], product code FRC.