



December 6, 2023

NeuroOne Medical Technologies Corp.
% John Doucet, Ph.D.
Vice President, Regulatory Affairs - Neuromodulation and Evolving Technologies
MCRA, LLC
803 7th Street, NW, 3rd Floor
Washington, District of Columbia 20001

Re: K231675

Trade/Device Name: OneRF Ablation System
Regulation Number: 21 CFR 882.4400
Regulation Name: Radiofrequency Lesion Generator
Regulatory Class: Class II
Product Code: GXD, GXI
Dated: November 6, 2023
Received: November 6, 2023

Dear Dr. John Doucet:

We have reviewed your section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (the Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database available at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Additional information about changes that may require a new premarket notification are provided in the FDA guidance documents entitled "Deciding When to Submit a 510(k) for a Change to an Existing Device"

(<https://www.fda.gov/media/99812/download>) and "Deciding When to Submit a 510(k) for a Software Change to an Existing Device" (<https://www.fda.gov/media/99785/download>).

Your device is also subject to, among other requirements, the Quality System (QS) regulation (21 CFR Part 820), which includes, but is not limited to, 21 CFR 820.30, Design controls; 21 CFR 820.90, Nonconforming product; and 21 CFR 820.100, Corrective and preventive action. Please note that regardless of whether a change requires premarket review, the QS regulation requires device manufacturers to review and approve changes to device design and production (21 CFR 820.30 and 21 CFR 820.70) and document changes and approvals in the device master record (21 CFR 820.181).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR Part 803) for devices or postmarketing safety reporting (21 CFR Part 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR Part 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR Parts 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Adam D. Pierce -S
Digitally signed by
Adam D. Pierce -S
Date: 2023.12.06
15:41:45 -05'00'

Adam D. Pierce, Ph.D.

Assistant Director

DHT5A: Division of Neurosurgical,
Neurointerventional

and Neurodiagnostic Devices

OHT5: Office of Neurological

and Physical Medicine Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K231675

Device Name

NeuroOne OneRF Radiofrequency Ablation System

Indications for Use (Describe)

The NeuroOne OneRF Radiofrequency Ablation System is indicated for creation of radiofrequency lesions in nervous tissue for functional neurosurgical procedures.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary
K231675

This summary of 510(k) safety and effectiveness information is submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

510(k) Number:	K231675
Date Prepared:	December 6, 2023
Applicant:	NeuroOne Medical Technologies Corp. 7599 Anagram Drive Eden Prairie, MN 55344 Phone: (952) 426-1383 E-mail: debk@nmtcl.com
Contact Person:	Debra Kridner Regulatory Affairs Consultant 7599 Anagram Drive Eden Prairie, MN 55344

SUBJECT DEVICE

Trade/Device Name:	OneRF™ Ablation System
Device Regulation Number:	21 CFR§882.4400 21 CFR§882.4725
Device/Regulation Name:	Radiofrequency Lesion Generator, Radiofrequency Lesion Probe
Product Code:	GXD GXI
Device Class/ Regulation Classification:	Class II

DEVICE DESCRIPTION (For the Device Subject to this 510(k) Premarket Notification)

OneRF™ Ablation System

The OneRF™ Ablation System components consist of the:

1. **Radiofrequency (RF) Generator** and Accessories
 - a. Generator Interface Cable (GIC), Cart and Foot Pedal (optional)
2. **sEEG-RF Probe** (with Universal Cable Assembly) and Ablation Accessories
 - a. Temperature Accessory, Spacer Tubes, Stylet and RF Connector Box

The OneRF™ Ablation System uses radiofrequency ablation to create lesion (s) in an area of nervous tissue that the neurosurgeon has identified for ablation. This is accomplished by diagnostically locating the area to ablate using the implanted Stereoelectroencephalography (sEEG). Once the area of ablation is located, the generator and ablation accessories are taken to the patient with implanted sEEG(s). The sEEG(s) now function as an sEEG-RF Probe for RF ablation.

To perform the ablation the Temperature Accessory (TA) is inserted into the sEEG-RF Probe to the predetermined electrode contact location using the specified spacer tube. A stylet may be used to ensure the lumen is patent prior to inserting the TA. The universal cable assembly remains attached to the sEEG-RF Probe while its cables are disconnected from the diagnostic head box. One (monopolar) or two (bipolar) cables from the universal cable assembly are inserted into the radiofrequency connector box (RFCB) depending on the contact area(s) to be ablated. The TA is also connected to the RFCB. The RFCB is attached to the Generator Interface Cable which is attached to the Generator. Temperature and time are used to create lesion (s).

PREDICATE DEVICE

NeuroOne has chosen a Radiofrequency Lesion Generator as our primary predicate device; specifically, the Cosman G4 Radiofrequency Generator that received marketing authorization after FDA review of K082051 and was cleared under 21 CFR 882.4400 (product code GXD). This RF Generator was chosen as the primary predicate because it has the same intended use, namely lesioning nervous tissue for neurosurgical procedures and similar technological characteristics as the subject device generator.

A second predicate device Radiofrequency Lesion Probe was chosen as RF lesions can only be created with a RF Generator and RF Probe. This device is the Diros OWL RF Probe which received marketing authorization after FDA review of K010202 and was cleared under 21 CFR 882.4725 (product code GXI). To create an RF Ablation System this RF Probe was chosen as an additional predicate because it has the same intended use as NeuroOne's sEEG-RF Probe; namely, to be connected to an RF Generator for lesioning nervous tissue for neurosurgical procedures and similar technological characteristics as the subject device probe.

The predicate devices were chosen as the Intended Use is the same. The Indications for Use, Fundamental Scientific Technology, and Principles of Operation are either the same or if different did not raise new questions of safety and/or effectiveness when compared to the subject device.

INDICATIONS FOR USE

The NeuroOne OneRF Radiofrequency Ablation System is indicated for creation of radiofrequency lesions in nervous tissue for functional neurosurgical procedures.

COMPARISON OF SUBJECT DEVICE (SYSTEM) TO PRIMARY PREDICATE DEVICE (RF GENERATOR) and PREDICATE (RF PROBE)

DEVICE CLASSIFICATION COMPARISON – RF GENERATOR AND RF PROBE

Device Classification Comparison - Generator			
	Subject Device	Primary Predicate Device K082051	Comparison Same/Different
Trade/Device Name:	NeuroOne OneRF™ Ablation System Generator	Cosman RF Generator	NA
Device Regulation Number	21 CFR§882.4400	21 CFR§882.4400	Same
Device / Regulation Name	Radiofrequency lesion generator	Radiofrequency lesion generator	Same
Device Regulation Identification	A radiofrequency lesion generator is a device used to produce lesions in the nervous system or other tissue by the direct application of radiofrequency currents to selected sites.	A radiofrequency lesion generator is a device used to produce lesions in the nervous system or other tissue by the direct application of radiofrequency currents to selected sites.	Same
Product Codes	GXD	GXD	Same
Device Class / Regulation Classification	Class II	Class II	Same

Device Classification Comparison - Probe			
	Subject Device	Predicate Device K010202	Comparison Same/Different
Trade/Device Name:	NeuroOne OneRF™ Ablation System sEEG-RF Probe	Diros RF Probe	NA
Device Regulation Number	21 CFR§882.4725	21 CFR§882.4725	Same
Device / Regulation Name:	Radiofrequency lesion probe	Radiofrequency lesion probe	Same
Device Regulation Identification	A radiofrequency lesion probe is a device connected to a radiofrequency (RF) lesion generator to deliver the RF energy to the site within the nervous system where a lesion is desired.	A radiofrequency lesion probe is a device connected to a radiofrequency (RF) lesion generator to deliver the RF energy to the site within the nervous system where a lesion is desired.	Same
Product Codes	GXI	GXI	Same
Device Class / Regulation Classification	Class II	Class II	Same

INTENDED USE/INDICATIONS FOR USE COMPARISON TABLE – RF GENERATOR/RF PROBE (ABLATION SYSTEM)

Intended Use/Indications for Use Comparison RF Generator and RF Probe (Ablation System)			
	Subject Device (RF Ablation System) K231675	Primary Predicate RF Generator K082051 and Predicate RF Probe K010202	Comparison Same/Different
Intended Use	RF Ablation System consists of a generator and probe used for lesioning nervous tissue for neurosurgical procedures.	RF Ablation System consists of a generator and probe used for lesioning nervous tissue for neurosurgical procedures.	Same
Indications for Use	The NeuroOne OneRF™ Radiofrequency Ablation System is indicated for creation of radiofrequency lesions in nervous tissue for functional neurosurgical procedures.	<p><u>RF Generator</u></p> <p>The Cosman G4 Radiofrequency Generator is indicated for use in procedures to create radiofrequency lesions for the treatment of pain, or for lesioning nerve tissue for functional neurosurgical procedures.</p> <p><i>The Cosman G4 Radiofrequency Generator is used with separately approved Cosman Radiofrequency Probes</i></p> <p><u>RF Probe</u></p> <ol style="list-style-type: none"> 1. Lesioning nerve tissue for functional neurosurgical procedures such as thalamotomies, pallidotomies, tractotomies, and myelotomies; or 2. radiofrequency heat lesion procedures for the relief of pain 	<p>Different</p> <p>Removed reference to relief/treatment of pain.</p> <p>The subject device indications for use does not separately define a specific probe or generator as they are defined as part of the NeuroOne system.</p>

TECHNOLOGICAL AND PERFORMANCE CHARACTERISTICS COMPARISON

Technological and Performance Characteristics Comparison – RF Generator			
	Subject Device	Primary Predicate Device K082051	Comparison Same/Different
Device name	NeuroOne OneRF™ Generator	Cosman RF Generator	N/A
RF Output Frequency	465.1 kHz	480 kHz	Different
Maximum output power (Hardware)	50 watts	50 watts	Same
Maximum output power (Software)	2 watts	50 watts	Different
Power delivery modes	Continuous and pulsed	Continuous and pulsed	Same
AC Power Compatibility	AC Line 100-240V	AC Line 100-240V	Same
Lesion Creation	Yes	Yes	Same
Temperature Monitoring	Yes	Yes	Same
Temperature Range (°C)	20 to 90	0 to 110	Different
Temperature Shut Off (°C)	90	>100	Different
Auto shutdown for temperature exceeding safe levels	Yes	Yes	Same
Electrical safety/EMC compliant	IEC 60601-1 and IEC 60601-1-2 compliant	IEC 60601-1 and IEC 60601-1-2 compliant	Same
RF energy delivery modes:	<ul style="list-style-type: none"> • Manual control • Temperature control 	<ul style="list-style-type: none"> • Manual control • Temperature control 	Same
RF energy delivery channel types	Monopolar Bipolar	Monopolar Bipolar	Same
User Touchscreen Interface	Yes	Yes	Same
USB Port updating software and downloading logs	Yes	Yes	Same

Technological and Performance Characteristics Comparison - RF Probe			
	Subject Device	Predicate Device K010202	Comparison Same/Different
Device name	NeuroOne OneRF™ sEEG-RF Probe	Diros RF Probe (Tasker Intracranial Lesion Electrode)	N/A
RF Probe Configuration	Single open lumen probe (closed at tip)	Single closed lumen probe	Same
Temperature Measurement Capabilities	Yes	Yes	Same
Location of Temperature Sensor	On Temperature Accessory (inserted into lumen of probe)	Embedded in probe tip	Different
Patient- Contact Materials	Polyimide - electrode Platinum - contact	Stainless steel Insulation – (medical grade, abrasion resistant)	Different
Compatibility with RF Generator	Yes, NeuroOne RF Generator (using adapters)	Yes, Diros and Cosman RF Generators (using legally marketed adapters)	Same
Key Dimensions:			
Electrode Active Length	2 mm	2 mm	Same
Diameter	0.8 mm	0.8 mm	Same
Number of Contacts	4 – 15 on each probe enabled for Ablation	1 on each probe	Different
Creation of Lesions in Nervous Tissue	Yes	Yes	Same
Comparative Lesion Size Testing in Ex-vivo Tissue	Yes	Yes	Same

System Characteristics	
System Characteristic	Identification
User	Neurosurgeons familiar with RF lesion techniques
Anatomical Site of Use	Nervous tissue
Access Method	Intracranial
Energy Type	Radiofrequency
Procedure Type	Ablation
Principle of Operation	Operator controlled; RF delivered from RF generator to compatible sEEG - RF probes to create lesions in nervous tissue
Mechanism of Action	Cellular necrosis through thermal coagulation
System Feedback Mechanism	Temperature controlled
Ability to Make Multiple Lesions	Yes

SUMMARY OF PERFORMANCE TESTING AND STANDARDS

The following performance data were provided in support of the substantial equivalence determination for the subject NeuroOne OneRF™ Ablation System. Performance testing of the subject device (system) was conducted to demonstrate that the device (system) meets its performance specifications. Results of design verification and validation activities did not raise any new or different questions of safety or effectiveness. The risk management process was used throughout the non-clinical verification activities in accordance with ISO 14971. The following tests/analysis were conducted.

Non-Clinical Performance Tests		
Test	Overview Summary	Results and Conclusions
Lesion Size Testing	This testing was performed to provide a measurement of the size (width and length) of the lesion for each active electrode contact: temperature, time, configuration (monopolar, bipolar), mode (temperature control, manual).	Lesion sizes were determined based on time and temperature. Lesion size is comparable to predicate.
Dimensional Verification and RFCB Cable Durability	This testing was performed to evaluate the dimensional characteristics and to demonstrate compatibility between components. In addition, the flexural durability of the RFCB cables was tested.	Pass – The test results indicate that the sEEG-RF Probe / Temperature Accessory / Stylet / Spacer Tube and Radio Frequency Connector Box designs meet the dimensional and cable durability requirements.
Mechanical Performance	This testing was performed to verify specifications related to the mechanical interaction between the sEEG-RF Probe and Accessories	Pass – The test results indicate that the sEEG-RF Probe / Temperature Accessory / Stylet / Spacer Tube and Radio Frequency Connector Box designs meet the mechanical performance requirements
Mechanical Integrity	This testing was performed to evaluate the mechanical durability of the sEEG-RF Probe Accessories.	Pass - The test results indicate that the sEEG-RF Probe Accessories (Temperature Accessory and Radio Frequency Connector Box) designs meet the mechanical integrity requirements
Ablation System Performance	The purpose of this testing is to evaluate specifications related to energy delivery and temperature accuracy of the OneRF Ablation System including durability after use.	Pass - The test results indicate that the sEEG-RF Probe / Temperature Accessory / Stylet / Spacer Tube and Radio Frequency Connector Box designs meet the system performance requirements.
Generator System	This testing was performed to verify specifications related to the Generator and UI Software.	Pass – The test results indicate that the Generator and UI Software designs meet the system performance requirements.
Electrical Safety	Product shall meet the applicable requirements of Electromagnetic Compatibility and Electrical Safety standards - IEC 60601-1, -2, -6 and IEC 60601-2-2	Pass – Met applicable requirements
Temperature Accessory Kit and RFCB Package Integrity	The packaged device and labeling shall withstand the conditions of packaging, shelf life, and distribution testing to ISO 11607-1, ISTA 3A, ASTM D4169, ASTM F1980-16, ASTM 2096, ASTM F88 without loss of function, sterility, or legibility.	Pass - The test results indicate that the sEEG-RF Probe Accessories (Temperature Accessory / Stylet / Spacer Tube and Radio Frequency Connector Box) packaging designs meet the integrity requirements (i.e., seal strength, bubble leak, label inspection, and no damage that impacts device sterility).

Sterilization	The sterilization process shall be validated to demonstrate a minimum of SAL of 10^{-6} for the product using Ethylene Oxide per ISO 11135	Pass – All criteria passed and the sterilization cycle was validated.
Usability – Summative Validation	This testing was performed in accordance with FDA guidance, “ <i>Applying Human Factors and Usability Engineering to Medical Devices</i> ”, February 3, 2016	Pass – The NeuroOne OneRF™ Ablation System has been found to be safe and effective for the intended users, uses, and use environments.
Software	Software was presented in accordance with FDA Guidance “ <i>Content of Premarket Submissions for Software Contained in Medical Devices</i> ” (issued on: May 11, 2005). In addition, IEC 62304 was followed, as applicable	Software analysis addressed applicable requirements
Cybersecurity	Cybersecurity analysis was performed on the OneRF™ Ablation System in accordance with “ <i>Content of Premarket Submissions for Management of Cybersecurity in Medical Devices</i> ” – Final Guidance October 2, 2014 and “ <i>Cybersecurity in Medical Devices: Quality System Considerations and Content of Premarket Submissions Draft Guidance for Industry and Food and Drug Administration Staff</i> ” – Draft Guidance April 7, 2022.	Cybersecurity analysis addressed applicable requirements

Biocompatibility	
Biocompatibility Overview	Conclusions
<p>Components tested: sEEG Electrode (sEEG-RF Probe), Anchor Bolt, Cap with prolonged (>24 hours to 30 days) contact with tissue/bone.</p> <p>Components tested: Strain Relief, Lock Band, Stylet Assembly, Electrode Tail, and Electrode Connector with limited (<24 hours) contact with intact skin.</p>	Passed – Reference (K211367/K222404)
OneRF Ablation System Sterile Components: Temperature Accessory/Spacer Tubes/Stylet and Radio Frequency Connector Box	No testing - there is No direct or indirect patient contact

BASIS FOR DETERMINATION OF SUBSTANTIAL EQUIVALENCE

The subject device and predicate devices have the same Intended Use, Device Regulation Number/Name/Identification and Product Codes. The differences in the Indications for Use and Fundamental Scientific Technology between the subject device, OneRF™ Ablation System (RF Generator and RF Probe) and the predicate RF Generator (Cosman) and RF Probe (Diros) do not raise new questions regarding safety and effectiveness when compared. Conclusions drawn from the nonclinical testing demonstrate the device is as safe, as effective, and performs as well as the legally marketed device predicates, per 21 CFR 807.92(b)(3). The OneRF™ Ablation System is substantially equivalent to the predicate devices.