

August 11, 2023

Prevest Denpro Limited % Angela Blackwell Senior Consultant Blackwell Device Consulting P.O. Box 718 Gresham, Oregon 97030-0172

Re: K231696

Trade/Device Name: Fusion Bond 5, Fusion Bond 7, Fusion Bond DC, Renew MDP, Renew Universal Regulation Number: 21 CFR 872.3200 Regulation Name: Resin Tooth Bonding Agent Regulatory Class: Class II Product Code: KLE Dated: May 22, 2023 Received: June 12, 2023

Dear Angela Blackwell:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely, Michael E. Adjodha -S

Michael E. Adjodha, M.ChE. Assistant Director DHT1B: Division of Dental and ENT Devices OHT1: Office of Ophthalmic, Anesthesia, Respiratory, ENT and Dental Devices Office of Product Evaluation and Quality Center for Devices and Radiological Health

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

510(k) Number *(if known)* K231696

Device Name

Fusion Bond 5, Fusion Bond 7, Fusion Bond DC, Renew MDP, Renew Universal

Indications for Use (Describe)

Fusion Bond 5 is indicated for direct composite or compomer restorations, adhesive cementation, and composite repair.

Fusion Bond 7 is indicated for bonding of composites to tooth structure, core build up, and adhesive cementation of crown & bridges, including inlays and onlays.

Fusion Bond DC is indicated for direct light cure composite or compomer restorations, core build up, and adhesive cementation of crown & bridges, including inlays and onlays.

Renew MDP is indicated for bonding of dual cure, light cure or self cure composite or compomer restorations to tooth structure, treatment of hypersensitive teeth, and intraoral repairs of fractured restorations.

Renew Universal is indicated for direct bonding of light-cured composites, and compomers to tooth structure, bonding of dual-cured core build up composites to tooth structure as long as these materials are light-cured, intraoral repair of composite, metal-based and zirconia /alumina-based restorations, intraoral repair of ceramic restorations in combination with a silane coupling agent, treatment of hypersensitive teeth, and cavity sealing as a pretreatment for indirect restorations.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

K231696

Fusion Bond 5, Fusion Bond 7, Fusion Bond DC, Renew MDP, Renew Universal 510(k) Summary August 9, 2023

Name and Address: Prevest Denpro Limited Unit II, Export Promotion Industrial Park Bari Brahmana, Jammu 181133 India Contact Person: Atul Modi Email: prevestindia@gmail.com Telephone: (941) 919 4280

Name of device: Fusion Bond 5, Fusion Bond 7, Fusion Bond DC, Renew MDP, Renew Universal Trade Name: Fusion Bond 5, Fusion Bond 7, Fusion Bond DC, Renew MDP, Renew Universal Common Name: dental bonding agent Classification Name: Resin tooth bonding agents CFR: 21 CFR 872.3200 Primary Product Code: KLE Regulatory Class: II

Submission Contact:

Angela Blackwell Blackwell Device Consulting P.O. Box 718 Gresham, OR 97030-0172 (704)450-9934 angela@blackwelldevice.com

Device Description:

Prevest Denpro Tooth Bonding Agents attach different types of restorations to teeth.

Fusion Bond 5 is a light curing single component bonding agent. It primes and bonds in one step to attach restorations of composites, copomers, amalgams, and porcelains to tooth structures. Tooth structures should have a 37% phosphoric acid etchant applied before application of Fusion Bond 5.

Fusion Bond 7 is a light curing single component bonding agent. It etches, primes, and bonds in one step to attach restorations of composites, copomers, amalgams, and porcelains to tooth structures.

Fusion Bond DC is a dual cure bonding agent for bonding direct and indirect restorations. Enamel and dentin surfaces should be conditioned with 37% phosphoric acid etchant before application of Fusion Bond DC. It is a two component adhesive with bonding agent and activator.

Renew MDP is a light curing single component bonding agent. It combines etching, priming and bonding in one bottle. It is an ethanol-water based dental adhesive containing 10-Methacryloyloxydececyl Dihydrogen Phosphate (10-MDP) a functional monomer which helps in bonding to dentin and cut and un-cut enamel. Renew MDP works with light-cured, self-cured and dual-cured composite materials. Renew Universal is a light curing single component bonding agent. It combines etching, priming and bonding. Renew Universal contains a combination of functional monomers such as 10-MDP and 4-META. Renew Universal works with light-cured and dual-cured composite materials.

Device Name	Indications
Fusion Bond 5	Fusion Bond 5 is indicated for direct composite or
	compomer restorations, adhesive cementation,
	and composite repair.
Fusion Bond 7	Fusion Bond 7 is indicated for bonding of
	composites to tooth structure, core build up, and
	adhesive cementation of crown & bridges,
	including inlays and onlays.
Fusion Bond DC	Fusion Bond DC is indicated for direct light cure
	composite or compomer restorations, core build
	up, and adhesive cementation of crown &
	bridges, including inlays and onlays.
Renew MDP	Renew MDP is indicated for bonding of dual cure,
	light cure or self cure composite or compomer
	restorations to tooth structure, treatment of
	hypersensitive teeth, and intraoral repairs of
	fractured restorations.
Renew Universal	Renew Universal is indicated for direct bonding of
	light-cured composites, and compomers to tooth
	structure, bonding of dual-cured core build up
	composites to tooth structure as long as these
	materials are light-cured, intraoral repair of
	composite, metal-based and zirconia /alumina-
	based restorations, intraoral repair of ceramic
	restorations in combination with a silane coupling
	agent, treatment of hypersensitive teeth, and
	cavity sealing as a pretreatment for indirect
	restorations.

Indications for Use:

Testing Summary:

Fusion Bond 5, Fusion Bond 7, Fusion Bond DC, Renew MDP and Renew Universal were tested for appearance, shear bond strength and curing time according to protocols based on ISO 29022:2013.

All test results met the criteria in standard.

Shelf life for Fusion Bond 5 and Fusion Bond 7 is 3 years. Shelf life for Fusion Bond DC, Renew MDP and Renew Universal is 2 years. All shelf life determinations use the same testing protocols as the characterization testing which are based on ISO 29022:2013. The predicate and reference devices use the same ISO standard for their testing. Their shelf lives are 3 years or are not given.

Primary Predicate Device: Prolink and Prolink SE K110403 from SIlmet

Additional Predicate Devices: Clearfil Photo Bond K943165 and Clearfil SE Bond 2 K131432 from Kuraray

Reference Devices Used for Ingredients:

Fusion Bond 5 - Adhese Universal (K133318) XP Bond Dual Cure Universal Total Etch Adhesive (K070538) Cal LC (K212457)

Fusion Bond 7 - Adper TM Prompt TM (K060684) Adhese Universal (K133318) Adhese Universal DC-(K210804) XP Bond Dual Cure Universal Total Etch Adhesive (K070538) Cal LC (K212457) Fusion Bond DC - XP Bond Dual Cure Universal Total Etch Adhesive (K070538) Cal LC (K212457) Renew MDP - Adhese Universal (K133318) Cal LC (K212457)

Renew Universal - Adhese Universal (K133318) Gluma Comfort Bond (K992985) Cal LC (K212457)

Substantial Equivalence:

The bonding agents have similar ingredients to the predicate and reference devices, the same indications for use, and similar physical parameter testing.

Resin Tooth Bonding Agents from Prevest Denpro

	Fusion Bond 5	Prolink K110403 from Silmet	
		Predicate Device	
Product	KLE	KLE	
Code			
Indications	Fusion Bond 5 is indicated for	Direct Composite or Compomer	
for Use	direct composite or compomer	restorations	
	restorations, adhesive	Adhesive cementation	
	cementation, and composite	Composite repair	
	repair.		
Mechanism	Bond restorations to teeth.	Bond restorations to teeth.	
of Action			
Applicable	ISO 29022 - Dentistry – Adhesive	ISO 29022 - Dentistry – Adhesive –	
Standards	 Notched-edge shear bond 	Notched-edge shear bond	
	strength test	strength test	
Compositio	 Bisphenol A Glycidyl 	Bisphenol diglycidyl	
n	Methacrylate	methacrylate	
	 Ethoxylated Bisphenol A 	 Ethoxylated Bisphenol A 	
	Dimethacrylate	Dimethacrylate	
	 Urethane Dimethacrylate 	Urethane Dimethacrylate	
	 Triethylene Glycol 	Triethylene Glycol	
	Dimethacrylate	Dimethacrylate	
	 2-Hydroxethylmethyl 	 2-Hydroxethylmethyl 	
	acrylate	acrylate	
	Camphorquinone	 camphoroquinone 	
	 Ethyl-4 Di methyl amino 	 Ethyl-4 Di methyl amino 	
	benzoate	benzoate	
	 Tertiary butanol 	Ethanol	

	 Butylated hydroxy toluene (BHT) 	Acetone
	2-dimethyl amino ethyl Methyacrylate	
Shear Bond	Complies with ISO 29022	Complies with ISO 29022
Strength	10 MPA	14 MPa
Curing	Complies with ISO 29022	Complies with ISO 29022
Time	20-30 sec	20 sec
Shelf Life	3 years	Not identified in 510k summary,
		package information says 3 years

	Fusion Bond 7	Prolink SE K110403 from Silmet	
		Predicate Device	
Product Code	KLE	KLE	
Indications for Use	Fusion Bond 7 is indicated for bonding of composites to tooth structure, core build up, and adhesive cementation of crown & bridges, including	Bonding of composites to tooth structure Core Build up Adhesive cementation of crown & bridges, including inlays and	
		onlays	
Mechanism of Action	Bond restorations to teeth.	Bond restorations to teeth.	
Composition	 Bisphenol A Glycidyl Methacrylate Ethoxylated Bisphenol A Dimethacrylate Urethane Dimethacrylate Triethylene Glycol Dimethacrylate Hydroxethylmethyl Phosphate Camphorquinone Ethyl-4 Di methyl amino benzoate Tertiary butanol Butylated hydroxy toluene (BHT) 2-dimethyl amino ethyl Methyacrylate 	 Bisphenol diglycidyl methacrylate Ethoxylated Bisphenol A Dimethacrylate Urethane Dimethacrylate Triethylene Glycol Dimethacrylate 2-Hydroxethylmethyl acrylate camphoroquinone Ethyl-4 Di methyl amino benzoate Ethanol Acetone 	
Applicable	ISO 29022 - Dentistry –	ISO 29022 - Dentistry – Adhesive –	
Standards	Adhesive – Notched-edge	Notched-edge shear bond	
	shear bond strength test	strength test	

Shear Bond	Complies with ISO 29022	Complies with ISO 29022
Strength	25 MPA	29MPa
Curing Time	Complies with ISO 29022	Complies with ISO 29022
	20-30 sec	20 sec
Shelf Life	3 years	Not identified in 510k summary,
		package information says 3 years

	Fusion Bond DC	Clearfil Photo Bond K943165	Prolink K110403 from
		from Kuraray	Silmet
		Additional Predicate Device	Predicate Device
Product	KLE	KLE	KLE
Code			
Indications	Fusion Bond DC is	No 510k database file	Direct Composite or
for Use	indicated for direct light		Compomer restorations
	cure composite or		Adhesive cementation
	compomer restorations,		Composite repair
	core build up, and		
	adhesive cementation of		
	crown & bridges,		
	including inlays and		
	oniays.	Devel veste vetieves to to the	Developerations to
of Action	Bond restorations to	Bond restorations to teeth.	Bond restorations to
OF ACTION Commonsistin	Leeln.		teetn.
Compositio	Fusion Bona DC Bonaing	Universal Liquid:	Bisphenol dialusidud
n	Agent	- NI NI disetta di s	algiyciayi
	- Dischargel A	• N-N dietnyi-p	methacrylate
	Bisphenol A	loludiene	Ethoxylated
	Glycidyi	Sodium benzene	Bisphenol A
		suifinate	Dimethacrylate
	Ethoxylated Bisphonol A	Ethyl Alcohol	Oretnane Dimotheonylate
	Dimothacrylate		
			Inethylene Glycol Dimothacrylata
	Oreunane Dimethacrylate		
			• 2-
	Glycol		
	Dimethacrylate		
	 Diffective via te Di Hydroxy Etbyl- 	Catalyst Liquid:	 Ethyl_4 Di methyl
	P-Toluene		amino benzoate
	Camphorquinon	Bisphenol A Glucidul	Fthanol
	e	Methacrylate	
	Ethyl-4 Di methyl	• 10-	
	amino benzoate	Methacryloyloxydecy	
	Tertiary butanol	l di hydrogen	
	,	phosphate	

	 Butylated hydroxy toluene (BHT) Fusion Bond DC Activator Bisphenol A Glycidyl Methacrylate Ethoxylated Bisphenol A Dimethacrylate Urethane Dimethacrylate Urethane Bimethacrylate Triethylene Glycol Benzyl Peroxide Tertiary butanol Butylated hydroxy toluene (BHT) 	 2-Hydroxethylmethyl acrylate Hydrophobic Dimethacrylate Benzyl Peroxide DL-Camphorquinone 	
Applicable	ISO 29022 - Dentistry –	ISO 29022 - Dentistry –	ISO 29022 - Dentistry –
Stanuarus	edge shear bond strength test	shear bond strength test	shear bond strength test
Shear Bond	Complies with ISO 29022	Complies with ISO 29022	Complies with ISO 29022
	LJ IVIPa	LD.9 IVIPa	14 IVIPd Complies with ISO 20022
	20-30 sers	10 sec	20 sec
Shelf Life	2 vears	Unknown	Not identified in 510k
	_ ,		summary, package information says 3 years

	Renew MDP	Prolink K110403 from	Clearfil SE Bond 2 K131432
		Silmet	from Kuraray
		Additional Predicate	Predicate Device
		Device	
Product	KLE	KLE	KLE
Code			
Indications	Renew MDP is indicated	Direct Composite or	[1] Direct restorations
for Use	for bonding of dual cure,	Compomer restorations	using light-cured
	light cure or self cure	Adhesive cementation	composite resin
	composite or compomer	Composite repair	[2] Cavity sealing as a
	restorations to tooth	I III	pretreatment for indirect
			restorations

	structure, treatment of hypersensitive teeth, and intraoral repairs of fractured restorations.		 [3] Treatment of exposed root surfaces [4] Treatment of hypersensitive teeth [5] Intraoral repairs of fractured restorations [6] Post cementation using a dual- or self-cured composite resin [7] Core build-ups using a light-, dual- or self-cured core material [8] Cementing inlays, onlays, crowns, bridges and veneers using a composite resin cement
Mechanism of Action	Bond restorations to teeth.	Bond restorations to teeth.	Bond restorations to teeth.
Composition	 Bisphenol A Glycidyl Methacrylate Urethane Dimethacrylate Triethylene Glycol Dimethacrylate 2- Hydroxethylmethyl acrylate Ethanol DL Camphorquinone Ethyl-4 Di methyl amino benzoate Butylated hydroxy toluene (BHT) 2-dimethyl amino ethyl Methyacrylate DM water 10- Methacryloxydecyl Dihydogen Phosphate 	 Bisphenol diglycidyl methacrylate Ethoxylated Bisphenol A Dimethacrylate Urethane Dimethacrylate Triethylene Glycol Dimethacrylate 2- Hydroxethylmethyl acrylate camphoroquinone Ethyl-4 Di methyl amino benzoate Ethanol Acetone 	 Bisphenol A Glycidyl Methacrylate 10- Methacryloxydecyl Dihydogen Phosphate 2- Hydroxethylmethyl acrylate Hydrophobic aliphatic dimethacrylate Dimethacrylate monomer dl- Camphorquinone Microfillers Initiators Accelerators
Applicable Standards	ISO 29022 - Dentistry – Adhesive – Notched-edge shear bond strength test	ISO 29022 - Dentistry – Adhesive – Notched-edge shear bond strength test	ISO 29022 - Dentistry – Adhesive – Notched-edge shear bond strength test

Shear Bond	Complies with ISO 29022	Complies with ISO 29022	Complies with ISO 29022
Strength	21 MPa	14 Mpa	12.45 MPa or 25-30 MPa
Curing Time	Complies with ISO 29022	Complies with ISO 29022	Complies with ISO 29022
	20-30 secs	20 sec	10 sec
Shelf Life	2 years	Not identified in 510k	Not identified in 510k
		summary, package	summary, package
		information says 3 years	information says 2 years

	Renew Universal	Prolink K110403 from	Clearfil SE Bond 2 K131432
		Silmet	from Kuraray
		Additional Predicate	Predicate Device
		Device	
Product	KLE	KLE	KLE
Code			
Indications	Renew Universal is	Direct Composite or	[1] Direct restorations
for Use	indicated for direct	Compomer restorations	using light-cured
	bonding of light-cured	Adhesive cementation	composite resin
	composites, and	Composite repair	[2] Cavity sealing as a
	compomers to tooth	1 1	pretreatment for indirect
	structure, bonding of dual-		restorations
	cured core build up		[3] Treatment of exposed
	composites to tooth		root surfaces
	structure as long as these		[4] I reatment of
	materials are light-cured,		s introductive teeth
	intraoral repair of		[5] Intraoral repairs of
	composite, metal-based		[6] Post cementation using
	and zirconia /alumina-		a dual- or self-cured
	based restorations,		composite resin
	intraoral repair of ceramic		[7] Core build-ups using a
	restorations in		light-, dual- or self-cured
	combination with a silane		core material
	coupling agent, treatment		[8] Cementing inlays,
	of hypersensitive teeth.		onlays, crowns, bridges
	and cavity sealing as a		and veneers using a
	pretreatment for indirect		composite resin cement
	restorations.		
Mechanism	Bond restorations to teeth.	Bond restorations to teeth.	Bond restorations to teeth.
of Action			
Composition	Bisphenol A	Bisphenol	Bisphenol A
	Glycidyl	diglycidyl	Glycidyl
	Methacrylate	methacrylate	Methacrylate
	Urethane	Ethoxylated	• 10-
	Dimethacrylate	Bisphenol A	Methacryloxydecyl
	Triethylene Glycol	Dimethacrylate	Dihydogen
	Dimethacrylate	Urethane	Phosphate
	,	Dimethacrylate	

	 2- Hydroxethylmethyl acrylate Ethanol DL Camphorquinone Ethyl-4 Di methyl amino benzoate Butylated hydroxy toluene (BHT) 2-dimethyl amino ethyl Methyacrylate DM water 4-2- Methacryloxyethyl Trimelltic 	 Triethylene Glycol Dimethacrylate 2- Hydroxethylmethyl acrylate camphoroquinone Ethyl-4 Di methyl amino benzoate Ethanol Acetone 	 2- Hydroxethylmethyl acrylate Hydrophobic aliphatic dimethacrylate Dimethacrylate monomer dl- Camphorquinone Microfillers Initiators Accelerators
Applicable	ISO 29022 - Dentistry –	ISO 29022 - Dentistry –	ISO 29022 - Dentistry –
Standards	Adhesive – Notched-edge	Adhesive – Notched-edge	Adhesive – Notched-edge
	shear bond strength test	shear bond strength test	shear bond strength test
Shear Bond	Complies with ISO 29022	Complies with ISO 29022	Complies with ISO 29022
Strength	21 MPa	14 MPa	12.45 MPa or 30 MPa
Curing Time	Complies with ISO 29022	Complies with ISO 29022	Complies with ISO 29022
	20-30 secs	20 sec	10 sec
Shelf Life	2 years	Not identified in 510k	Not identified in 510k
		summary, package	summary, package
		information says 3 years	information says 2 years

Conclusion: Prevest Denpro tooth bonding agents are substantially equivalent to the predicate devices Prolink and Prolink SE K110403 from Silmet. They have the same indications, similar testing, and very similar ingredients. Both the subject devices and the predicate device have physical parameters which meet requirements of the relevant ISO standards. Shelf life testing is similar to the shelf life testing of predicate or reference device. Reference devices are included to cover any ingredients, or indications not covered by the predicate devices. Any differences in ingredients are minor and do not change the substantial equivalence.