



June 4, 2024

Owens & Minor (O&M) Halyard, Inc.
Anureet Singh
Regulatory Affairs Manager
9120 Lockwood Blvd
Mechanicsville, Virginia 23116

Re: K232820

Trade/Device Name: HALYARD* FLUIDSHIELD* 1 Fog-Free Surgical Mask with So Soft Lining
(28806)
Regulation Number: 21 CFR 878.4040
Regulation Name: Surgical Apparel
Regulatory Class: Class II
Product Code: FXX
Dated: October 30, 2023
Received: October 30, 2023

Dear Anureet Singh:

We have reviewed your section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (the Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database available at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Additional information about changes that may require a new premarket notification are provided in the FDA guidance documents entitled "Deciding When to Submit a 510(k) for a Change to an Existing Device" (<https://www.fda.gov/media/99812/download>) and "Deciding When to Submit a 510(k) for a Software Change to an Existing Device" (<https://www.fda.gov/media/99785/download>).

Your device is also subject to, among other requirements, the Quality System (QS) regulation (21 CFR Part 820), which includes, but is not limited to, 21 CFR 820.30, Design controls; 21 CFR 820.90, Nonconforming product; and 21 CFR 820.100, Corrective and preventive action. Please note that regardless of whether a change requires premarket review, the QS regulation requires device manufacturers to review and approve changes to device design and production (21 CFR 820.30 and 21 CFR 820.70) and document changes and approvals in the device master record (21 CFR 820.181).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR Part 803) for devices or postmarketing safety reporting (21 CFR Part 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR Part 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR Parts 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,


Allan Guan -S

For Bifeng Qian, M.D., Ph.D.
Assistant Director

DHT4B: Division of Infection Control
and Plastic and Reconstructive Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K232820

Device Name

HALYARD* FLUIDSHIELD* 1 Fog-Free Surgical Mask with So Soft Lining (28806)

Indications for Use (Describe)

HALYARD* FLUIDSHIELD* 1 Fog-Free Surgical Mask with So Soft Lining (28806) is intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids, and particulate material. These face masks are intended for use in infection control practices to reduce the potential exposure of the wearer to blood and body fluids. The HALYARD* FLUIDSHIELD* 1 Fog-Free Surgical Mask with So Soft Lining (28806) is a single use, disposable device, provided non-sterile.

Type of Use (Select one or both, as applicable)

☐ Prescription Use (Part 21 CFR 801 Subpart D)

☒ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary – K232820

Date of Summary	June 04, 2024
510(k) Submitter	O&M Halyard, Inc. 9120 Lockwood Boulevard Mechanicsville, VA 23116 Phone: 804-723-7000/800-488-8850 Fax: 804-723-7100
Primary Contact for this 510(k) Submission	Anureet Singh Regulatory Affairs Manager
Marketed Common Name	Surgical Mask
Device Submission Trade Name and Description	HALYARD* FLUIDSHIELD* 1 Fog-Free Surgical Mask with So Soft Lining (28806)
Device Product Code, Class and Classification Name	FXX II Surgical Mask (21 CFR 878.4040)
Predicate Device	KC100 Face Masks (K110455)
Subject Device Description	The subject device is an ASTM F2100 Level 1 (green in color) Surgical Mask. The ASTM F2100 Level 1 Surgical Mask family of products are a three-layer mask, constructed of well-known non-woven materials used in facial protection. The Surgical Masks are provided with ties and a malleable nosepiece. The HALYARD* FLUIDSHIELD* 1 Fog-Free Surgical Mask with So Soft Lining (28806) is provided with fog-free foam. These masks are single use, disposable devices, provided non-sterile.
Indications for Use	HALYARD* FLUIDSHIELD* 1 Fog-Free Surgical Mask with So Soft Lining (28806) is intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids, and particulate material. These face masks are intended for use in infection control practices to reduce the potential exposure of the wearer to blood and body fluids. The HALYARD* FLUIDSHIELD* 1 Fog-Free Surgical Mask with So Soft Lining (28806) is a single use, disposable device, provided non-sterile.

Technological Characteristics Comparison Table			
	Subject Device HALYARD* FLUIDSHIELD* 1 Fog-Free Surgical Mask with So Soft Lining (28806)	Predicate Device KC100 Face Masks (K110455)	Comparison
FDA Product Code	FXX	FXX	Same
FDA Classification	Class II	Class II	Same

Regulation Number	21 CFR 878.4040	21 CFR 878.4040	Same
Common Name	Surgical Mask	Surgical Mask	Same
Device Trade Name	HALYARD* FLUIDSHIELD* 1 Fog-Free Surgical Mask with So Soft Lining (28806)	KC100 Surgical Mask with Fog-Free Strip	Different
Indications for Use	HALYARD* FLUIDSHIELD* 1 Fog-Free Surgical Mask with So Soft Lining (28806) is intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids, and particulate material. These face masks are intended for use in infection control practices to reduce the potential exposure of the wearer to blood and body fluids. The HALYARD* FLUIDSHIELD* 1 Fog-Free Surgical Mask with So Soft Lining (28806) is a single use, disposable device, provided non-sterile.	The Kimberly-Clark KC100 Face Mask(s) is intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids, and particulate material. These face masks are intended for use in infection control practices to reduce the potential exposure of the wearer to blood and body fluids. The Kimberly-Clark KC100 face mask(s) is a single use, disposable device(s), provided non-sterile.	Same The device trade names are different

Design Attributes			
	Subject Device HALYARD* FLUIDSHIELD* 1 Fog-Free Surgical Mask with So Soft Lining (28806)	Predicate Device KC100 Face Masks (K110455)	Comparison
Outer Layer	Spunbond polypropylene	Spunbond polypropylene	Same
Middle Layer	Meltblown polypropylene	Meltblown polypropylene	Same
Inner Layer	Polyethylene/polyester	Polyethylene/polyester	Same
Bindings	Spunbond polypropylene	Spunbond polypropylene or polyester spunlace	Same
Ties	Spunbond polypropylene	Spunbond polypropylene or polyester spunlace	Same
Foam	Polyester foam/polyester spunlace laminate, blue	Polyester foam/polyester spunlace laminate, blue	Same
Nose wire	Aluminum	Aluminum or polyethylene coated steel	Same
Number of Layers	3	3	Same
Color	Green	Yellow, blue, green, white, lavender	Similar
Style	Pleated with ties	Pleated with ties	Same
ASTM Level	1	1	Same

EN 14683 (Type II)	Yes	No	Different
Single Use	Yes	Yes	Same
Sterility	Non-sterile	Non-sterile	Same

Summary of Non-Clinical Performance Testing:

Performance testing of the HALYARD* FLUIDSHIELD* 1 Fog-Free Surgical Mask with So Soft Lining (28806) was evaluated, and the results showed that acceptance criteria were met.

Purpose	Test	Acceptance Criteria	Result
Face Mask Performance	ASTM F2100	ASTM F2100 Level 1	Pass
<i>Bacterial Filtration Efficiency</i>	<i>ASTM F2101</i>	$\geq 95\%$	<i>Pass</i>
<i>Particulate Filtration Efficiency</i>	<i>ASTM F2299</i>	$\geq 95\%$	<i>Pass</i>
<i>Differential Pressure</i>	<i>EN 14683</i>	$< 5.0 \text{ mmH}_2\text{O}/\text{cm}^2$	<i>Pass</i>
<i>Fluid Resistance</i>	<i>ASTM F1862</i>	<i>80 mmHg</i>	<i>Pass</i>
<i>Flammability</i>	<i>16 CFR Part 1610</i>	<i>Class 1</i>	<i>Pass</i>
Face Mask Performance	EN 14683	EN 14683 Type II	Pass
<i>Bacterial Filtration Efficiency</i>	<i>EN 14683</i>	$\geq 98\%$	<i>Pass</i>
<i>Differential Pressure</i>	<i>EN 14683</i>	$< 40 \text{ Pa}/\text{cm}^2$	<i>Pass</i>
<i>Microbial Cleanliness</i>	<i>ISO 11737-1</i>	$\leq 30 \text{ cfu/g}$	<i>Pass</i>
Biocompatibility	ISO 10993		Pass
<i>Cytotoxicity</i>	<i>ISO 10993-5</i>	<i>Non-cytotoxic</i>	<i>Pass</i>
<i>Sensitization</i>	<i>ISO 10993-10</i>	<i>Non-sensitizing</i>	<i>Pass</i>
<i>Irritation</i>	<i>ISO 10993-10</i>	<i>Non-irritant</i>	<i>Pass</i>

Summary of Clinical Performance Testing:

No clinical testing required.

Conclusion:

The conclusions drawn from the nonclinical tests demonstrate that the HALYARD* FLUIDSHIELD* 1 Fog-Free Surgical Mask with So Soft Lining (28806) is as safe, as effective and performs as well as or better than the legally marketed KC100 Face Masks (K110455).