



April 2, 2024

Canon Medical Systems Corporation  
% Orlando Tadeo Jr  
Sr. Manager, Regulatory Affairs  
Canon Medical Systems, USA  
2441 Michelle Drive  
TUSTIN, CA 92780

Re: K232835

Trade/Device Name: Aquilion ONE (TSX-308A/3) V1.4 with PIQE Reconstruction System  
Regulation Number: 21 CFR 892.1750  
Regulation Name: Computed Tomography X-Ray System  
Regulatory Class: Class II  
Product Code: JAK  
Dated: February 27, 2024  
Received: February 28, 2024

Dear Orlando Tadeo Jr:

We have reviewed your section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (the Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database available at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Additional information about changes that may require a new premarket notification are provided in the FDA guidance documents entitled "Deciding When to Submit a 510(k) for a Change to an Existing Device" (<https://www.fda.gov/media/99812/download>) and "Deciding When to Submit a 510(k) for a Software Change to an Existing Device" (<https://www.fda.gov/media/99785/download>).

Your device is also subject to, among other requirements, the Quality System (QS) regulation (21 CFR Part 820), which includes, but is not limited to, 21 CFR 820.30, Design controls; 21 CFR 820.90, Nonconforming product; and 21 CFR 820.100, Corrective and preventive action. Please note that regardless of whether a change requires premarket review, the QS regulation requires device manufacturers to review and approve changes to device design and production (21 CFR 820.30 and 21 CFR 820.70) and document changes and approvals in the device master record (21 CFR 820.181).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR Part 803) for devices or postmarketing safety reporting (21 CFR Part 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR Part 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR Parts 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

The image shows a stylized signature of 'Lu Jiang' in a cursive font, overlaid on a large, light blue 'FDA' logo.

Lu Jiang, Ph.D.  
Assistant Director  
Diagnostic X-Ray Systems Team  
DHT8B: Division of Radiological Imaging  
Devices and Electronic Products  
OHT8: Office of Radiological Health  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K232835

Device Name

Aquilion ONE (TSX-308A/3) V1.4 with PIQE Reconstruction System

Indications for Use (Describe)

This device is indicated to acquire and display cross sectional volumes of the whole body, to include the head, with the capability to image whole organs in a single rotation. Whole organs include but are not limited to brain, heart, pancreas, etc. The Aquilion ONE has the capability to provide volume sets of the entire organ. These volume sets can be used to perform specialized studies, using indicated software/hardware, of the whole organ by a trained and qualified physician.

FIRST is an iterative reconstruction algorithm intended to reduce exposure dose and improve high contrast spatial resolution for abdomen, pelvis, chest, cardiac, extremities and head applications.

AiCE is a noise reduction algorithm that improves image quality and reduces image noise by employing Deep Convolutional Neural Network methods for abdomen, pelvis, lung, cardiac, brain, extremities, head, and inner ear applications.

The spectral imaging system allows the system to acquire two nearly simultaneous CT images of an anatomical location using distinct tube voltages and/or tube currents by rapid KV switching. The X-ray dose will be the sum of the dose at each respective tube voltage and current in a rotation.

Information regarding the material composition of various organs, tissues, and contrast materials may be gained from the differences in X-ray attenuation between these distinct energies. When used by a qualified physician, a potential application is to determine the course of treatment.

PIQE is a Deep Learning Reconstruction method designed to enhance spatial resolution. By incorporating noise reduction into the Deep Convolutional Neural Network (DCNN), it is possible to achieve both spatial resolution improvement and noise reduction for cardiac, abdomen, and pelvis applications, in comparison to FBP and hybrid iterative reconstruction..

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

**CONTINUE ON A SEPARATE PAGE IF NEEDED.**

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRAStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

**510(k) SUMMARY**

- 1. SUBMITTER'S NAME:**  
Fumiaki Teshima  
Senior Manager, Quality Assurance Department  
Canon Medical Systems Corporation  
1385 Shimoishigami  
Otagawa-Shi, Tochigi-ken, Japan 324-8550
- 2. ESTABLISHMENT REGISTRATION:**  
9614698
- 3. OFFICIAL CORRESPONDENT/CONTACT PERSON:**  
Orlando Tadeo, Jr.  
Sr. Manager, Regulatory Affairs  
Canon Medical Systems USA, Inc  
2441 Michelle Drive  
Tustin, CA 92780  
(714) 669-7459
- 4. DATE PREPARED:**  
September 13, 2023
- 5. TRADE NAME(S):**  
Aquilion ONE (TSX-308A/3) V1.4 with PIQE Reconstruction System
- 6. COMMON NAME:**  
Computed Tomography X-ray System
- 7. DEVICE CLASSIFICATION:**
  - a) Classification Name: Computed Tomography X-ray system
  - b) Regulation Number: 21 CFR §892.1750
  - c) Regulatory Class: Class II
- 8. PRODUCT CODE:**  
JAK
- 9. PERFORMANCE STANDARD:**  
This device conforms to applicable Performance Standards for Ionizing Radiation Emitting Products [21 CFR, Subchapter J, Part 1020]

**10. PREDICATE DEVICE:**

Product	Marketed by	Regulation Number	Regulation Name	Product Code	510(k) Number	Clearance Date
Aquilion ONE (TSX-306A/3) V10.12 with Spectral Imaging System	Canon Medical Systems, USA	21 CFR §892.1750	Computed Tomography X-ray System	JAK: System, X-ray, Tomography, Computed	K213504	June 16, 2022

**11. REASON FOR SUBMISSION:**

New medical device

**12. DEVICE DESCRIPTION:**

**Aquilion ONE (TSX-308A/3) V1.4 with PIQE Reconstruction System** is a whole body multi-slice helical CT scanner, consisting of a gantry, couch and a console used for data processing and display. This device captures cross sectional volume data sets used to perform specialized studies, using indicated software/hardware, by a trained and qualified physician. This system is based upon the technology and materials of previously marketed Canon CT systems.

**Aquilion ONE (TSX-308A/3) V1.4 with PIQE Reconstruction System** is equipped with PIQE, a deep learning reconstruction technology designed to fully utilize the maximum resolution of the detector, intended to improve spatial resolution. Original image data is available to end users when PIQE images are used for diagnosis.

**13. INDICATIONS FOR USE:**

This device is indicated to acquire and display cross sectional volumes of the whole body, to include the head, with the capability to image whole organs in a single rotation. Whole organs include but are not limited to brain, heart, pancreas, etc. The Aquilion ONE has the capability to provide volume sets of the entire organ. These volume sets can be used to perform specialized studies, using indicated software/hardware, of the whole organ by a trained and qualified physician.

FIRST is an iterative reconstruction algorithm intended to reduce exposure dose and improve high contrast spatial resolution for abdomen, pelvis, chest, cardiac, extremities and head applications.

AiCE is a noise reduction algorithm that improves image quality and reduces image noise by employing Deep Convolutional Neural Network methods for abdomen, pelvis, lung, cardiac, extremities, head, and inner ear applications.

The spectral imaging system allows the system to acquire two nearly simultaneous CT images of an anatomical location using distinct tube voltages and/or tube currents by rapid KV switching. The X-ray dose will be the sum of the dose at each respective tube voltage and current in a rotation. Information regarding the material composition of various organs, tissues, and contrast materials may be gained from the differences in X-ray attenuation between these distinct energies. When used by a qualified physician, a potential application is to determine the course of treatment.

PIQE is a Deep Learning Reconstruction method designed to enhance spatial resolution. By incorporating noise reduction into the Deep Convolutional Neural Network (DCNN), it is possible to achieve both spatial resolution improvement and noise reduction for cardiac, abdomen, and pelvis applications, in comparison to FBP and hybrid iterative reconstruction.

**14. SUBSTANTIAL EQUIVALENCE:**

The **Aquilion ONE (TSX-308A/3) V1.4 with PIQE Reconstruction System** is substantially equivalent to Aquilion ONE (TSX-306A/3) V10.12 with Spectral Imaging System, which received premarket clearance under K213504, and is marketed by Canon Medical Systems USA. The intended use of the **Aquilion ONE** is the same as that of the predicate device. A comparison of the technological characteristics between the subject and the predicate device is included below.

	<b>Subject Device</b>	<b>Predicate Device</b>
<b>Device Name, Model Number</b>	<b>Aquilion ONE (TSX-308A/3) V1.4 with PIQE Reconstruction System</b>	<b>Aquilion ONE (TSX-306A/3) V10.12 with Spectral Imaging System</b>
<b>510(k) Number</b>	<b>This submission</b>	<b>K213504</b>
<b>PIQE Reconstruction System</b>	Available	Available
<ul style="list-style-type: none"> <li>▪ Scan Regions</li> <li>▪ Scan Type</li>   <li>▪ Intensity</li> <li>▪ Reconstruction image matrix</li> <li>▪ Image thickness</li> <li>▪ Reconstruction Interval</li> </ul>	<ul style="list-style-type: none"> <li>▪ Cardiac, <i>Abdomen and pelvis (Body)</i></li> <li>▪ Volume scan, Dynamic volume scan, <i>Helical scan</i></li> <li>▪ L1 / L2 / L3</li> <li>▪ 512 x 512, <i>1024 x 1024</i></li>   <li>▪ 0.5 and 1.0 mm, <i>2.0 mm (Helical only)</i></li> <li>▪ 0.25 (With option installed), 0.5, <i>1.0, 2.0 mm (Helical only)</i></li> </ul>	<ul style="list-style-type: none"> <li>▪ Cardiac</li> <li>▪ Volume scan, Dynamic volume scan</li>   <li>▪ MILD / STANDARD / STRONG</li> <li>▪ 512 x 512</li>   <li>▪ 0.5 and 1.0 mm</li> <li>▪ 0.25 (With option installed) and 0.5 mm</li> </ul>
<b>Noise reduction processing</b>	<ul style="list-style-type: none"> <li>▪ Adaptive Iterative Dose Reduction 3D (AIDR 3D)</li> <li>▪ AIDR 3D Enhanced</li> <li>▪ Advanced Intelligent Clear-IQ Engine (AiCE)</li> </ul>	<ul style="list-style-type: none"> <li>▪ Adaptive Iterative Dose Reduction 3D (AIDR 3D)</li> <li>▪ AIDR 3D Enhanced</li> <li>▪ Advanced Intelligent Clear-IQ Engine (AiCE)</li> </ul>
<b>Scan modes</b>	Conventional scan (Axial Scan) Volume, Dynamic volume, <i>Dynamic scan</i> Helical scan	Conventional scan (S&S, S&V) Volume, Dynamic volume scan Helical scan
<b>Positioning Scan</b>	Single, Dual, 3D Landmark Scan	Single, Dual, S-Helical (3D Landmark Scan)
<b>Anatomical Landmark Detection</b>	<i>Available</i>	Not Available
<b>Scan (Rotation) time (Some may require additional options)</b>	Half scan: (0.15, 0.18, 0.23) 0.24, 0.275, 0.3, 0.32, 0.35, 0.375, 0.4, 0.45, 0.5, 0.6, 0.75, 1.0, 1.5, 2.0, 3.0 s	Half scan: (0.18, 0.23) 0.275, 0.3, 0.32, 0.35, 0.375, 0.4, 0.45, 0.5, 0.6, 0.75, 1.0, 1.5, 2.0, 3.0 s
<b>Gantry opening diameter (aperture)</b>	<i>800 mm in diameter</i>	780 mm in diameter
<b>Wedge filter types</b>	Three (3) types: Small, Large, SilverBeam Filter (Dose reduction wedge)	Three (3) types: Medium, Large, SilverBeam Filter (Dose reduction wedge)
<b>Gantry internal cameras</b>	<i>Available</i>	Not Available
<b>Operating panel (Gantry monitor)</b>	<i>Operating panel</i> <i>Two operating panels on the front</i>	∅Station
<b>Couch Speed - Horizontal</b>	0.8 - 450 mm/s	0.8 - 300 mm/s
<b>Image display matrix</b>	<i>2428 x 1230 (max.)</i>	1024 x 1024 (max.)

**15. SAFETY:**

The device is designed and manufactured under the Quality System Regulations as outlined in 21 CFR § 820 and ISO 13485 Standards. This device is in conformance with the applicable parts of the following standards IEC60601-1, IEC60601-1-2, IEC60601-1-3, IEC60601-1-6, IEC60601-1-9, IEC60601-2-28, IEC60601-2-44, IEC60825-1, IEC62304, IEC62366-1, NEMA XR-25, NEMA XR-26 and NEMA XR-29. Additionally, this device complies with all applicable requirements of the radiation safety performance standards, as outlined in 21 CFR §1010 and §1020.

This device conforms to applicable Performance Standards for Ionizing Radiation Emitting Products [21 CFR, Subchapter J, Part 1020]

**16. TESTING**

Risk analysis and verification/validation activities conducted through bench testing demonstrate that the established specifications for the device have been met.

***Performance Testing - Bench***

Image Quality Evaluations

CT image quality assessments were performed, utilizing phantoms, to evaluate the general Image Quality, SilverBeam, AiCE, Spectral Reconstruction and PIQE performance of TSX-308A (Aquilion ONE) system relative to the predicate device, TSX-306A (Aquilion Prism) with regard to Contrast-to-Noise Ratios, CT Number Accuracy, Uniformity, Slice Sensitivity Profile, Modulation Transfer Function, Standard Deviation of Noise, Noise Power Spectra, Low Contrast Detectability. It was concluded that the subject device demonstrated equivalent or improved performance, compared to the predicate device, as demonstrated by the results of the above testing.

SilverBeam Dose Reduction

A phantom study was conducted which compared dose in Head/Body modes between normal scan mode and in DR-mode (with SilverBeam Filter) and it was determined that DR-mode resulted in dose reduction.

Low Contrast Detectability

A phantom study was conducted to measure low contrast detectability at 0.3%/2 mm and 0.3%/3 mm using AIDR 3D and AiCE. In this study the low contrast module of the CATPHAN 600 phantom was scanned and an image read was conducted for the 2mm and 3mm objects in order to calculate the minimum mA at which the low contrast object can be identified. Study results support the following claims of low contrast detectability:

- 2 mm at 0.3%, 15.3 mGy CTDIvol, using AIDR3D
- 2 mm at 0.3%, 14.7mGy CTDIvol, using AiCE
- 3 mm at 0.3%, 5.7mGy CTDIvol, using AiCE

***Performance Testing – Clinical Images***

Representative body, cardiac, chest, head, and extremity diagnostic images, reviewed by American Board-Certified Radiologists, were obtained using the subject device and it was confirmed that the reconstructed images using the subject device were of diagnostic quality.

A summary of the risk analysis and verification/validation testing conducted through bench testing is included in this submission which demonstrates that the requirements for the system have been met.

Software Documentation for a Basic Documentation Level, per the FDA guidance document, “Content of Premarket Submissions for Device Software Functions” issued on June 14, 2023, is included in this submission. This documentation includes justification for the Basic Documentation Level determination as well as testing which demonstrates that the verification and validation requirements have been met.

Cybersecurity documentation, per the FDA cybersecurity premarket guidance document “Cybersecurity in Medical Devices: Quality System Considerations and Content of Premarket Submissions” issued on September 27, 2023, is also included as part of this submission.

**17. CONCLUSION**

The **Aquilion ONE (TSX-308A/3) V1.4 with PIQE Reconstruction System** performs in a manner similar to and is intended for the same use as the predicate device, as indicated in product labeling. Based upon this information, conformance to standards, successful completion of software validation, application of risk management and design controls and the performance data presented in this submission it is concluded that the subject device has demonstrated substantial equivalence to the predicate device and is safe and effective for its intended use.