



November 21, 2023

Belega Inc.
% Nana Nagashita
Regulatory Consultant
Ken Block Consulting Co., Ltd.
1-9-9 Tsukiji
Chuo-ku, Tokyo 104-0045
Japan

Re: K233010

Trade/Device Name: Beagank 4T Plus
Regulation Number: 21 CFR 882.5890
Regulation Name: Transcutaneous Electrical Nerve Stimulator For Pain Relief
Regulatory Class: Class II
Product Code: NFO
Dated: September 15, 2023
Received: September 22, 2023

Dear Nana Nagashita:

We have reviewed your section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (the Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database available at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Additional information about changes that may require a new premarket notification are provided in the FDA guidance documents entitled "Deciding When to Submit a 510(k) for a Change to an Existing Device" (<https://www.fda.gov/media/99812/download>) and "Deciding When to Submit a 510(k) for a Software Change to an Existing Device" (<https://www.fda.gov/media/99785/download>).

Your device is also subject to, among other requirements, the Quality System (QS) regulation (21 CFR Part 820), which includes, but is not limited to, 21 CFR 820.30, Design controls; 21 CFR 820.90, Nonconforming product; and 21 CFR 820.100, Corrective and preventive action. Please note that regardless of whether a change requires premarket review, the QS regulation requires device manufacturers to review and approve changes to device design and production (21 CFR 820.30 and 21 CFR 820.70) and document changes and approvals in the device master record (21 CFR 820.181).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR Part 803) for devices or postmarketing safety reporting (21 CFR Part 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR Part 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR Parts 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Heather L. Dean -S

Heather Dean, PhD
Assistant Director, Acute Injury Devices Team
DHT5B: Division of Neuromodulation
and Rehabilitation Devices

OHT5: Office of Neurological
and Physical Medicine Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Indications for Use

510(k) Number (if known)

K233010

Device Name

BEAGANK 4T PLUS

Indications for Use (Describe)

The BEAGANK 4T PLUS is a handheld portable device for over-the-counter aesthetic use including facial and neck stimulation or body skin stimulation.

Type of Use (Select one or both, as applicable)

☐ Prescription Use (Part 21 CFR 801 Subpart D)

☒ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) SUMMARY

Owner/Submitter: Belega Co., Ltd.
1-4-20 Sonezaki Shinchu, Kita-ku
Osaka-shi, Osaka, Japan 530-0002

Contact Person: Ken Block Consulting Co., Ltd.
Ms. Nana Nagashita
Regulatory consultant

1-9-9 Tsukiji, Chuo-ku
Tokyo, Japan 104-0045
TEL: +81-3-6260-6688
FAX: +81-3-6260-6689
nagashita@kenblockconsulting.com

Date Prepared: September 20, 2023

Submission Type: Traditional 510(k)

Proposed Device

Manufacturer:	BELEGA Co., Ltd.
Trade Name:	BEAGANK 4T PLUS
Common Name:	Aesthetic TENS Device
Classification Name:	Transcutaneous Electrical Nerve Stimulator for Pain Relief
Product Code:	NFO
Regulatory Standard	21 CFR 882.5890

Predicate Device:

510(k) Number:	K191951
Manufacturer:	Avazzia, Inc.
Trade Name:	Avazzia OTC TENS for Aesthetics, model BEST-AV1: EZZI-LIFT Device
Common Name:	Aesthetic TENS Device
Classification Name:	Transcutaneous Electrical Nerve Stimulator for Pain Relief
Product Code	NFO
Regulatory Standard:	21 CFR 882.5890

Reference Device:

510(k) Number:	K011935
Manufacturer:	Salton Inc.
Trade Name:	Rejuvenique [®] system, Model RJV-10
Common Name:	Aesthetic TENS Device
Classification Name:	Transcutaneous Electrical Nerve Stimulator for Pain Relief
Product Code:	NFO, GYB (Media, Conductive)
Regulatory Standard:	21 CFR 882.5890, 21 CFR 882.1275 (Media, Conductive)

Reference Device:

510(k) Number:	K181008
Manufacturer:	Carol Cole Company dba NuFACE
Trade Name:	NuFACE Trinity
Common Name:	Aesthetic TENS Device
Classification Name:	Transcutaneous Electrical Nerve Stimulator for Pain Relief
Product Code:	NFO
Regulatory Standard:	21 CFR 882.5890

Device Description: BEAGANK 4T PLUS is a cosmetic device which is intended for use on the face and body skin. The device contains a rechargeable lithium ion battery and is provided with a charger. The device has four electrode metal pads on the head that deliver electrical impulses to the face and body skin and provide four different operational modes.

The device operates by applying Transcutaneous Electrical Nerve Stimulation to the face and body skin. The electrodes that are positioned on the skin must be covered with a wet cotton pad. Purified water or commercially available cosmetic formulation liquid can be used as a conductive media to generate the electric pulses. Two attachments are provided with the BEAGANK 4T PLUS:

- Scalp attachment
- Small attachment

Statement of Intended Use: The BEAGANK 4T PLUS is a handheld portable device for over-the-counter aesthetic use including facial and neck stimulation or body skin stimulation.

Summary of Technological Characteristics: The major differences between the proposed and predicate devices are summarized in the tables below:

Table 1: Basic characteristics				
	New Device	Predicate Device	Reference Device	Reference Device
Trade Name	BEAGANK 4T PLUS	EZZI-LIFT Device	Rejuvenique	NuFACE Trinity
510(k)	K233010	K191951	K011935	K181008
Indication for Use	The BEAGANK 4T PLUS is a handheld portable device for over-the-counter aesthetic use including facial and neck stimulation or body skin stimulation.	The Avazzia OTC TENS for aesthetics, model BEST-AV1™; EZZI-LIFT™ Device is indicated for over-the-counter aesthetic use including facial and neck stimulation or body skin stimulation.	Rejuvenique system is indicated for cosmetic use.	The NuFACE Trinity is intended for facial and neck stimulation and is indicated for over-the-counter cosmetic use.
Device Class	Class II	Class II	Class II	Class II
Operating Power	Internal rechargeable Lithium-ion battery	Two 1.5 V AA batteries	9 V battery	Internal rechargeable Lithium-ion battery
Method of Line Current Isolation	Type BF	Unknown	Unknown	Type BF

Patient Leakage Current	Normal	DC: <1µA AC: <1µA Tested per IEC 60601-1	Unknown	Unknown	N/A – Battery operated
	Single fault	DC: <1µA AC: <1µA Tested per IEC 60601-1	Unknown	Unknown	N/A – Battery operated
Ave. DC current through electrodes when device is on but no pulse are being applied (µA)		< 1µA	Unknown	Unknown	Unknown
Output Channels	Method of Channel Isolation	N/A – 1 output channel	Unknown	Electrode group selected by relays	N/A – 1 output channel
Regulated Current or Regulated Voltage		Voltage	Unknown	Voltage	Both
Automatic Overload Trip		Not required due to circuit design	Unknown	No	Not required due to circuit design
Automatic No-Load Trip		Not required due to circuit design	Unknown	No	Yes
Automatic Shut Off		Yes	Yes	Yes	Yes
Patient Override Control		Yes	Unknown	Yes	Yes
Weight		136g	153g (5.4 ounces)	80 g	9 oz. without charging base
Number of output modes		4 modes: maximum 5 minutes each	4 modes: 60 minutes	1	Unknown
Low Battery Indicator		Yes	Yes	Yes	Yes
Sterility		Non-Sterile	Non-Sterile	Non-Sterile	Non-Sterile
Electrical safety		- IEC 60601-1 Medical electrical equipment – Part 1 General requirements for basic safety and essential performance - IEC 60601-1-2 Medical electrical equipment –	- IEC 60601-1 Medical electrical equipment – Part 1 General requirements for safety - IEC 60601-1-2 Medical electrical equipment – Part 1-2: General requirements	Unknown	- IEC 60601-1 - IEC 60601-1-2

	Part 1-2: General requirements for basic safety and essential performance – Collateral Standard: Electromagnet ic disturbances – Requirements and tests	for basic safety and essential performance – Collateral Standard: Electromagnet ic disturbances – Requirements and tests - IEC 60601-2- 10 medical electrical equipment – Part 2:10; Particular requirements for the safety of nerve and muscle stimulators		
Housing Materials and Construction	ABS	PCBs inside plastic case housing	ABS Plastic, snap latch assembly	Thermoplastic
Accessories	Small and Scalp attachment	Built-in, Y, Brush, Pencil	Unknown	Unknown

Table 2: Technological Characteristics

Output Specifications: Mode 1				
	New Device	Predicate Device	Reference Device	Reference Device
Trade Name	BEAGANK 4T PLUS	EZZI-LIFT Device	Rejuvenique	NuFACE Trinity
Indications for Use	The BEAGANK 4T PLUS is a handheld portable device for over-the-counter aesthetic use including facial and neck stimulation or body skin stimulation.	The Avazzia OTC TENS for aesthetics, model BEST-AV1™; EZZI-LIFT™ Device is indicated for over-the-counter aesthetic use including facial and neck stimulation or body skin stimulation.	Rejuvenique system is indicated for cosmetic use.	The NuFACE Trinity is intended for facial and neck stimulation and is indicated for over-the-counter cosmetic use.

Waveform Shape	Rectangle, biphasic asymmetric	Positive square wave followed by a damped sinusoidal waveform of variable duration depending on damping and body loading	Pulsed biphasic, Rectangular (+phase), Spike (-phase)	Pulsed biphasic modulated square wave
Max output voltage	(+/-15%)	(+/-20%)	(+/-10%)	28VDC
@ 500Ω	127mV (0.12V)	-42V	18.8V	not publicly available
@ 2kΩ	515mV (0.51V)	-122V	24.8V	not publicly available
@ 10kΩ	1.83V	-348V	28.0V	not publicly available
Max output current	(+/-15%)	(+/-20%)	(+/-10%)	400μA (0.4mA)
@ 500Ω	0.19mA	363μA(0.363 mA)	37.6mA	not publicly available
@ 2kΩ	0.18mA	117μA (0.117 mA)	12.4mA	not publicly available
@ 10kΩ	0.16mA	38μA (0.038 mA)	2.8mA	not publicly available
Duration of primary pulse at 500Ω	Positive: 130μs (0.13msec) Negative: 130μs (0.13msec)	0.5msec	not publicly available	60msec
Pulse Duration (Cycle) at 500Ω	265μs (0.265msec)	1.1msec	300msec (+phase) 124.7msec (-phase, exponential)	60msec
Frequency	3.80kHz	15-121Hz	8Hz fixed	8.3Hz
Net Charge per pulse at 500Ω	0.025μC	4μC	0μC	not publicly available
Max phase charge	0.059μC	10μC	11.3μC	not publicly available
Max current density at 500Ω	- Standard: 0.15mA/cm ² - Scalp: 0.53mA/cm ² - Small:	Built-in, Y, Brush: 800μA/cm ² (0.80mA /cm ²)	46.4mA/cm ²	not publicly available

	1.36mA/cm ²	Pencil: 19,000µA/cm ² (19mA/cm ²)		
	With a wet cotton pad (Purified water) - Standard:0.15 mA/cm ² - Scalp: 0.53mA/cm ² Small: 1.36mA/cm ²			
Max power density at 500Ω (Smallest electrode)	- 8.51µW/cm ²	Built-in, Y, Brush: Average 500µW/cm ²	2.31mW/cm ² (2,310µW/cm ²)	not publicly available
	- With a wet cotton pad 11.39µW/cm ²	Pencil: Average 3,500µW/cm ²		
Max. Average Current (mA)	0.105	not publicly available	not publicly available	not publicly available
Burst mode				
Pulse per burst	N/A	N/A	160	20
Burst per second	N/A	N/A	1/240 (per electrode group)	8.3
Burst duration	N/A	N/A	20sec	2.4
Duty cycle	N/A	N/A	1/12	20.2sec
On Time (seconds)	N/A	N/A	20sec/electrode group	not publicly available
Off Time (seconds)	N/A	N/A	not publicly available	not publicly available
Output Specifications: Mode 2				
	New Device	Predicate Device	Reference Device	Reference Device
Trade Name	BEAGANK 4T PLUS	EZZI-LIFT Device	Rejuvenique	NuFACE Trinity
Indications for Use	The BEAGANK 4T PLUS is a handheld portable device for over-the-counter aesthetic use including facial and	The Avazzia OTC TENS for aesthetics, model BEST-AV1™; EZZI-LIFT™ Device is indicated for over-	Rejuvenique system is indicated for cosmetic use.	The NuFACE Trinity is intended for facial and neck stimulation and is indicated for over-

	neck stimulation or body skin stimulation.	the-counter aesthetic use including facial and neck stimulation or body skin stimulation.		the-counter cosmetic use.
Waveform Shape	Rectangle, biphasic asymmetric	Positive square wave followed by a damped sinusoidal waveform of variable duration depending on damping and body loading	Pulsed biphasic, Rectangular (+phase), Spike (-phase)	Pulsed biphasic modulated square wave
Max output voltage @ 500Ω @ 2kΩ @ 10kΩ	(+/-15%) 110mV (0.11V) 507mV (0.50V) 1.87V	(+/-20%) -42V -122V -348V	(+/-10%) 18.8V 24.8V 28.0V	28VDC not publicly available not publicly available
Max output current @ 500Ω @ 2kΩ @ 10kΩ	(+/-15%) 0.45mA 0.39mA 0.24mA	(+/-20%) 363μA (0.363mA) 117μA (0.117mA) 38μA (0.038mA)	(+/-10%) 37.6mA 12.4mA 2.8mA	400μA (0.4mA) not publicly available not publicly available
Duration of primary pulse at 500Ω	Positive: 130μs (0.13msec) Negative: 130μs (0.13msec)	0.5msec	not publicly available	60msec
Pulse Duration (Cycle) at 500Ω	265μs (0.265msec)	1.1msec	300msec (+phase) 124.7msec (-phase, exponential)	60msec
Frequency	3.80kHz	15-121Hz	8Hz fixed	8.3Hz
Net Charge per pulse at 500Ω	0.025μC	4μC	0μC	not publicly available

Max phase charge	0.025 μ C	10 μ C	11.3 μ C	not publicly available
Max current density at 500 Ω	- Standard: 0.36mA/ cm² - Scalp: 1.25mA/cm² - Small: 3.21mA/cm²	Built-in, Y, Brush: 800 μ A/cm ² (0.8mA/cm ²) Pencil: 19,000 μ A/cm ² (19mA/cm ²)	46.4mA/cm ²	not publicly available
	With a wet cotton pad (Purified water) - Standard: 0.10mA/cm² - Scalp: 0.36 mA/cm² - Small: 0.92 mA/cm²			
Max power density at 500 Ω (Smallest electrode)	- 17.46 μ W/cm ²	Built-in, Y, Brush: Average 500 μ W/cm ² Pencil: Average 3,500 μ W/cm ²	2.31mW/cm ² (2,310 μ W/cm ²)	not publicly available
	- With a wet cotton pad 33.87 μ W/cm ²			
Max. Average Current (mA)	0.105	not publicly available	not publicly available	not publicly available
Burst mode				
Pulse per burst	N/A	N/A	160	20
Burst per second	N/A	N/A	1/240 (per electrode group)	8.3
Burst duration	N/A	N/A	20sec	2.4
Duty cycle	N/A	N/A	1/12	20.2sec
On Time (seconds)	N/A	N/A	20sec/electrode group	not publicly available
Off Time (seconds)	N/A	N/A	not publicly available	not publicly available
Output Specifications: Mode 3				
	New Device	Predicate Device	Reference Device	Reference Device

Trade Name	BEAGANK 4T PLUS	EZZI-LIFT Device	Rejuvenique	NuFACE Trinity
Indications for Use	The BEAGANK 4T PLUS is a handheld portable device for over-the-counter aesthetic use including facial and neck stimulation or body skin stimulation.	The Avazzia OTC TENS for aesthetics, model BEST-AV1™. EZZI-LIFT™ Device is indicated for over-the-counter aesthetic use including facial and neck stimulation or body skin stimulation.	Rejuvenique system is indicated for cosmetic use.	The NuFACE Trinity is intended for facial and neck stimulation and is indicated for over-the-counter cosmetic use.
Waveform Shape	Rectangle, biphasic symmetrical	Positive square wave followed by a damped sinusoidal waveform of variable duration depending on damping and body loading	Pulsed biphasic, Rectangular (+phase), Spike (-phase)	Pulsed biphasic modulated square wave
Max output voltage	(+/-15%)	(+/-20%)	(+/-10%)	28VDC
@ 500Ω	16.6V	-42	18.8V	not publicly available
@ 2kΩ	21.0V	-122	24.8V	not publicly available
@ 10kΩ	22.5V	-348	28.0V	not publicly available
Max output current	(+/-15%)	(+/-20%)	(+/-10%)	400μA
@ 500Ω	30.8mA	363μA(0.363 mA)	37.6mA	not publicly available
@ 2kΩ	8.52mA	117μA (0.117 mA)	12.4mA	not publicly available
@ 10kΩ	2.26mA	38μA (0.038 mA)	2.8mA	not publicly available
Duration of primary pulse at 500Ω	Positive: 80μs (0.08msec) Negative: 80μs (0.08msec)	0.5msec	not publicly available	60msec
Pulse Duration (Cycle) at 500Ω	274μs (0.247msec)	1.1msec	300msec (+phase) 124.7msec (-phase, exponential)	60msec

Frequency	3.65kHz	15-121Hz	8Hz fixed	8.3Hz
Net Charge per pulse at 500Ω	-1.3μC	4μC	0μC	not publicly available
Max phase charge	2.2μC	10μC	11.3μC	not publicly available
Max current density at 500Ω	- Standard: 24.4mA/ cm² - Scalp: 85.6mA/cm² - Small: 220 mA/cm²	Built-in, Y, Brush: 800μA/cm ² (0.8mA/cm ²) Pencil: 19,000μA/cm ² (19mA/cm ²)	46.4mA/cm ²	not publicly available
	With a wet cotton pad (Purified water) - Standard: 4.26mA/ cm² - Scalp: 28.21mA/cm² - Small: 38.43mA/cm²			
Max power density at 500Ω (Smallest electrode)	- 106638.40μW/cm ²	Built-in, Y, Brush: Average 500μW/cm ² Pencil: Average 3,500μW/cm ²	2.31mW/ cm ² (2,310μW/cm ²)	not publicly available
	- With a wet cotton pad 1124.0μW/cm ²			
Max. Average Current (mA)	0	not publicly available	not publicly available	not publicly available
Burst mode				
Pulse per burst	N/A	N/A	160	20
Burst per second	N/A	N/A	1/240 (per electrode group)	8.3
Burst duration	N/A	N/A	20sec	2.4
Duty cycle	N/A	N/A	1/12	20.2sec

On Time (seconds)	N/A	N/A	20sec/electrode group	not publicly available
Off Time (seconds)	N/A	N/A	not publicly available	not publicly available
Output Specifications: Mode 4				
	New Device	Predicate Device	Reference Device	Reference Device
Trade Name	BEAGANK 4T PLUS	EZZI-LIFT Device	Rejuvenique	NuFACE Trinity
Indications for Use	The BEAGANK 4T PLUS is a handheld portable device for over-the-counter aesthetic use including facial and neck stimulation or body skin stimulation.	The Avazzia OTC TENS for aesthetics, model BEST-AV1™; EZZI-LIFT™ Device is indicated for over-the-counter aesthetic use including facial and neck stimulation or body skin stimulation.	Rejuvenique system is indicated for cosmetic use.	The NuFACE Trinity is intended for facial and neck stimulation and is indicated for over-the-counter cosmetic use.
Waveform Shape	Rectangle, Biphasic symmetric	Positive square wave followed by a damped sinusoidal waveform of variable duration depending on damping and body loading	Pulsed biphasic, Rectangular (+phase), Spike (-phase)	Pulsed biphasic modulated square wave
Max output voltage	(+/-15%)	(+/-20%)	(+/-10%)	28VDC
@ 500Ω	15.3V	-42V	18.8V	not publicly available
@ 2kΩ	21.2V	-122V	24.8V	not publicly available
@ 10kΩ	22.9V	-348V	28.0V	not publicly available
Max output current	(+/-15%)	(+/-20%)	(+/-10%)	400μA
@ 500Ω	31.4mA	363μA(0.363 mA)	37.6mA	not publicly available
@ 2kΩ	8.60mA	117μA (0.117 mA)	12.4mA	not publicly available
@ 10kΩ	2.15mA	38μA (0.038 mA)	2.8mA	not publicly available

Duration of primary pulse at 500Ω	Positive: 92μs (0.092msec) Negative: 92μs (0.092msec)	0.5msec	not publicly available	60msec
Pulse Duration (Cycle) at 500Ω	610μs (0.61msec)	1.1msec	300msec (+phase) 124.7msec (-phase, exponential)	60msec
Frequency	1.64kHz	15 to 121 Hz	8 Hz fixed	8.3Hz
Net Charge per pulse at 500Ω	0μC (biphasic symmetrical)	4μC	0μC	not publicly available
Max phase charge	0.29μC	10μC	11.3μC	not publicly available
Max current density at 500Ω	- Standard: 25.0mA/cm² - Head: 87.2mA/cm² - Small: 224.3mA/cm²	Built-in, Y, Brush: 800μA/cm ² (0.8mA/cm ²) Pencil: 19,000μA/cm ² (19mA/cm ²)	46.4 mA/cm ²	not publicly available
	With a wet cotton pad (Purified water) - Standard: 3.13mA/cm² - Scalp: 10.97mA/cm² - Small: 28.21mA/cm²			
Max power density at 500Ω (Smallest electrode)	- 51775.54μW/cm²	Built-in, Y, Brush: Average 500μW/cm ² Pencil: Average 3,500μW/cm ²	2.31mW/cm ² (2,310μW/cm ²)	not publicly available
	- With a wet cotton pad 556.96μW/cm²			
Max. Average Current (mA)	0	not publicly available	not publicly available	not publicly available
Burst mode				
Pulse per burst	N/A	N/A	160	20

Burst per second	N/A	N/A	1/240 (per electrode group)	8.3
Burst duration	N/A	N/A	20sec	2.4
Duty cycle	N/A	N/A	1/12	20.2sec
On Time (seconds)	N/A	N/A	20sec/electrode group	not publicly available
Off Time (seconds)	N/A	N/A	not publicly available	not publicly available

Clinical Test Data: No new clinical studies have been submitted as part of this Premarket Notification.

Summary of Non-Clinical / Test Data: Tests were performed on the device which demonstrated that the device is safe and effective, performs comparably to and is substantially equivalent to the predicate device. Testing confirmed that the device complies with safety standards for Electrical Safety and Electromagnetic Compatibility testing, specifically IEC standards 60601-1 and 60601-1-2.

Testing confirmed that the device electrodes complies with safety standards for Biocompatibility testing, specifically ISO standards 10993-5, 10993-10, and 10993-23.

Conclusion: Company considers the BEAGANK 4T PLUS device to be substantially equivalent to the predicate device and reference devices listed above. This conclusion is based on the similarities in primary intended use, principles of operation, functional design, and established medical use.