



December 23, 2025

Sim&Cure  
Colette Maurin  
Senior Director, Regulatory Affairs and Quality Assurance  
95 rue Pierre Flourens, Batiment H  
Montpellier, 34090  
France

Re: K253122  
Trade/Device Name: Sim&Size  
Regulation Number: 21 CFR 892.2050  
Regulation Name: Medical Image Management And Processing System  
Regulatory Class: Class II  
Product Code: PZO  
Dated: September 25, 2025  
Received: November 28, 2025

Dear Colette Maurin:

We have reviewed your section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (the Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database available at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Additional information about changes that may require a new premarket notification are provided in the FDA guidance documents entitled "Deciding When to Submit a 510(k) for a Change to an Existing Device"

(<https://www.fda.gov/media/99812/download>) and "Deciding When to Submit a 510(k) for a Software Change to an Existing Device" (<https://www.fda.gov/media/99785/download>).

Your device is also subject to, among other requirements, the Quality System (QS) regulation (21 CFR Part 820), which includes, but is not limited to, 21 CFR 820.30, Design controls; 21 CFR 820.90, Nonconforming product; and 21 CFR 820.100, Corrective and preventive action. Please note that regardless of whether a change requires premarket review, the QS regulation requires device manufacturers to review and approve changes to device design and production (21 CFR 820.30 and 21 CFR 820.70) and document changes and approvals in the device master record (21 CFR 820.181).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR Part 803) for devices or postmarketing safety reporting (21 CFR Part 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR Part 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR Parts 1000-1050.

All medical devices, including Class I and unclassified devices and combination product device constituent parts are required to be in compliance with the final Unique Device Identification System rule ("UDI Rule"). The UDI Rule requires, among other things, that a device bear a unique device identifier (UDI) on its label and package (21 CFR 801.20(a)) unless an exception or alternative applies (21 CFR 801.20(b)) and that the dates on the device label be formatted in accordance with 21 CFR 801.18. The UDI Rule (21 CFR 830.300(a) and 830.320(b)) also requires that certain information be submitted to the Global Unique Device Identification Database (GUDID) (21 CFR Part 830 Subpart E). For additional information on these requirements, please see the UDI System webpage at <https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/unique-device-identification-system-udi-system>.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

**SARA S. THOMPSON -S**

Sara S. Thompson, D.V.M.

Assistant Director

DHT5A: Division of Neurosurgical,

Neurointerventional, and

Neurodiagnostic Devices

OHT5: Office of Neurological and

Physical Medicine Devices

Office of Product Evaluation and Quality

Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K253122

Device Name

Sim&Size

### Indications for Use (Describe)

Sim&Size enables visualization of cerebral blood vessels for preoperational planning and sizing for neurovascular interventions and surgery. Sim&Size also allows for the ability to computationally model the placement of neurointerventional devices.

General functionalities are provided such as:

- Segmentation of neurovascular structures
- Automatic centerline detection
- Visualization of X-ray based images for 2D review and 3D reconstruction
- Placing and sizing tools
- Reporting tools

Information provided by the software is not intended in any way to eliminate, replace or substitute for, in whole or in part, the healthcare provider's judgment and analysis of the patient's condition.

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

#### **\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) Summary

### 1. Contact Details

**Applicant Name:** Sim&Cure  
**Applicant Address :** 95 rue Pierre Flourens, Bâtiment H  
34090 Montpellier  
FRANCE

**Applicant Contact Telephone:** +33953438809  
**Applicant Contact:** Mrs. Lucile AZEMA  
**Applicant Contact Email:** [l.azema@sim-and-cure.com](mailto:l.azema@sim-and-cure.com)

**Correspondent Name:** Sim&Cure  
**Correspondent Address:** 95 rue Pierre Flourens, Bâtiment H  
34090 Montpellier  
FRANCE

**Correspondent Contact Telephone:** +33953438809  
**Correspondent Contact:** Mrs. Colette MAURIN  
**Correspondent Contact Email:** [c.maurin@sim-and-cure.com](mailto:c.maurin@sim-and-cure.com)  
**Date prepared:** November 28, 2025

### 2. Device

Device Trade Name: Sim&Size

Common Name: Medical image management and processing system

Classification Name: Software for Visualization of Vascular Anatomy and Intravascular Devices

Regulation Number: 892.2050

Product Code: PZO

### 3. Legally Marketed Predicate Device

Predicate	Predicate Trade Name	Product Code
K242124	Sim&Size	PZO

## 4 Device Description Summary

Sim&Size is a Software as a Medical Device (SaMD) for simulating neurovascular implantable medical devices. The product enables visualization of cerebral blood vessels for preoperational planning for neurovascular interventions and surgery. It uses an image of the patient produced by 3D rotational angiography. It offers clinicians the possibility of simulating neurovascular implantable medical devices in the artery or in the aneurysm to be treated through endovascular surgery and provides support in the treatment for the sizing and positioning of implantable medical devices.

Each type of implant device is simulated in a simulation module of Sim&Size:

- FSize, a module that allows pre-operationally planning Flow-Diverter (FD) devices.
- ISize, a module that allows pre-operationally planning Intracascular (ID) devices.
- SSize, a module that allows pre-operationally planning Stent (ST) devices.
- FSize, a module that allows pre-operationally planning First and filling coils (FC) devices.

## 5 Intended Use / Indications for Use

Sim&Size enables visualization of cerebral blood vessels for preoperational planning and sizing for neurovascular interventions and surgery.

Sim&Size also allows for the ability to computationally model the placement of neurointerventional devices.

General functionalities are provided such as:

- Segmentation of neurovascular structures
- Automatic centerline detection
- Visualization of X-Ray based images for 2D review and 3D reconstruction
- Placing and sizing tools
- Reporting tools

Information provided by the software is not intended in any way to eliminate, replace or substitute for, in whole or in part, the healthcare provider's judgment and analysis of the patient's condition.

## 6 Comparison of Intended Use and Technological Characteristics

Characteristics and features of the Sim&Size software have been compared to the predicate device, Sim&Size (K242124).

Characteristics	Sim&Size K242124 (Predicate Device)	Sim&Size K253122 (Subject of this submission)	Comments
Intended Use / Indications for Use	<p>Sim&amp;Size enables visualization of cerebral blood vessels for preoperational planning and sizing for neurovascular interventions and surgery.</p> <p>Sim&amp;Size also allows for the ability to computationally model the placement of neurointerventional devices.</p> <p>General functionalities are provided such as:</p> <ul style="list-style-type: none"> <li>• Segmentation of neurovascular structures</li> <li>• Automatic centerline detection</li> <li>• Visualization of X-ray based images for 2D review and 3D reconstruction</li> <li>• Placing and sizing tools</li> <li>• Reporting tools</li> </ul> <p>Information provided by the software is not intended in any way to eliminate, replace or substitute for, in whole or in part, the healthcare provider's judgment and analysis of the patient's condition.</p>	<p>Sim&amp;Size enables visualization of cerebral blood vessels for preoperational planning and sizing for neurovascular interventions and surgery.</p> <p>Sim&amp;Size also allows for the ability to computationally model the placement of neurointerventional devices. General functionalities are provided such as:</p> <ul style="list-style-type: none"> <li>• Segmentation of neurovascular structures</li> <li>• Automatic centerline detection</li> <li>• Visualization of X-Ray based images for 2D review and 3D reconstruction</li> <li>• Placing and sizing tools</li> <li>• Reporting tools</li> </ul> <p>Information provided by the software is not intended in any way to eliminate, replace or substitute for, in whole or in part, the healthcare provider's judgment and analysis of the patient's condition.</p>	Same
Patient contact	No	No	Same
Human intervention for image interpretation	Yes	Yes	Same
Computer OS Compatibility	MS Windows Mac OS	MS Windows Mac OS	Same
Patient data management	<p>Import: manual through keyboard/mouse, and automatic import with an image file, study creation list.</p> <p>Export</p> <p>Deletion</p> <p>Anonymization</p> <p>Search</p>	<p>Import: manual through keyboard/mouse, and automatic import with an image file, study creation list.</p> <p>Export</p> <p>Deletion</p> <p>Anonymization</p> <p>Search</p>	Same
Data interchange	Local files, transfer through physical media (e.g. USB memory stick) and PACS connectivity (query/retrieve)	Local files, transfer through physical media (e.g. USB memory stick) and PACS connectivity (query/retrieve)	Same
DICOM Support	Compatible with DICOM image data from 3D rotational angiography stations.	Compatible with DICOM image data from 3D rotational angiography stations.	Same
Image Processing	Segmentation by user.	Segmentation by user.	Same
Image display	Orthogonal, color volume rendering, active presets, and 3D view of assemblies of devices.	Orthogonal, color volume rendering, active presets, and 3D view of assemblies of devices.	Same

3D assessment	<p>Assessment based on the 3D model of the simulated implant in the cerebrovascular:</p> <ul style="list-style-type: none"> <li>- Apposition indication of the implant along the arterial wall/aneurysm wall,</li> <li>- Volume embolization ratio indication of coil (FCsize module),</li> <li>- Metal surface coverage indication of the flow diverter (FSize module),</li> <li>- Pull-down maneuver of flow diverter deployment in fusiform aneurysms (FSize module),</li> <li>- Compression indication of the implant along the aneurysm wall (IDsize),</li> <li>- Reconstructed structure measurements,</li> <li>- Simulated implant dimensions measurements.</li> </ul>	<p>Assessment based on the 3D model of the simulated implant in the cerebrovascular:</p> <ul style="list-style-type: none"> <li>- Apposition indication of the implant along the arterial wall/aneurysm wall,</li> <li>- <b>Relative sizing, which represents the oversizing/undersizing of the flow diverter or braided stent along the arterial wall (FSize and STsize module),</b></li> <li>- Volume embolization ratio indication of the coil (FCsize module),</li> <li>- Metal surface coverage indication of the flow diverter (FSize module),</li> <li>- <b>Pore density of the flow diverter (FSize module),</b></li> <li>- Pull-down maneuver of flow diverter deployment in fusiform aneurysms (FSize module),</li> <li>- Compression indication of the implant along the aneurysm wall (IDsize),</li> <li>- <b>Spruce Index, which represents the heterogeneity of the device's deformation (IDsize module),</b></li> <li>- Reconstructed structure measurements,</li> <li>- Simulated implant length and diameter measurements.</li> </ul>	Similar
Implantable Medical Device Database	<p>In the FSize module:</p> <ul style="list-style-type: none"> <li>- Medtronic Pipeline Flex Embolization Device (PED – P100018/S15);</li> <li>- Medtronic Pipeline Flex Embolization Device with Shield Technology (PED2 – P100018/S026);</li> <li>- Medtronic Pipeline Vantage Embolization Device with Shield Technology (PED3 – P100018/S034);</li> <li>- Stryker Surpass Evolve Flow Diverter System (P170024/S003);</li> <li>- Stryker Surpass Elite Flow Diverter (P170024/S012);</li> <li>- MicroVention Flow Re-Direction Endoluminal Device System (FRED – P180027);</li> <li>- MicroVention Flow Re-Direction Endoluminal Device X System (FRED X – P180027/S002).</li> </ul> <p>In the IDsize module:</p> <ul style="list-style-type: none"> <li>- MicroVention Woven EndoBridge Aneurysm Embolization System (WEB – P170032).</li> </ul>	<p>In the FSize module:</p> <ul style="list-style-type: none"> <li>- Medtronic Pipeline Flex Embolization Device (PED – P100018/S15);</li> <li>- Medtronic Pipeline Flex Embolization Device with Shield Technology (PED2 – P100018/S026);</li> <li>- Medtronic Pipeline Vantage Embolization Device with Shield Technology (PED3 – P100018/S034);</li> <li>- Stryker Surpass Evolve Flow Diverter System (P170024/S003);</li> <li>- Stryker Surpass Elite Flow Diverter (P170024/S012);</li> <li>- MicroVention Flow Re-Direction Endoluminal Device System (FRED – P180027);</li> <li>- MicroVention Flow Re-Direction Endoluminal Device X System (FRED X – P180027/S002).</li> </ul> <p>In the IDsize module:</p> <ul style="list-style-type: none"> <li>- MicroVention Woven EndoBridge Aneurysm Embolization System (WEB – P170032).</li> </ul>	Similar



	<p>In the STsize module:</p> <ul style="list-style-type: none"> <li>- Stryker Neuroform Atlas Stent System (P180031/S001);</li> <li>- MicroVention Low-Profile Visualized Intraluminal Support and LVIS Jr (LVIS and LVIS Jr – P170013);</li> <li>- MicroVention Low-Profile Visualized Intraluminal Support EVO (LVIS EVO) (P170013/S004).</li> </ul> <p>In the FCsize module:</p> <ul style="list-style-type: none"> <li>- Medtronic Axium Detachable Coil and Axium Prime Detachable Coil (K233420);</li> <li>- MicroVention HydroCoil Embolic System (HES – K161367);</li> <li>- Stryker Target and Target XXL Detachable Coils (K161429);</li> <li>- Stryker Target Tetra Detachable Coils (K222533);</li> <li>- Balt Optima Coil System (K223386).</li> </ul>	<p>In the STsize module:</p> <ul style="list-style-type: none"> <li>- Stryker Neuroform Atlas Stent System (P180031/S001);</li> <li>- MicroVention Low-Profile Visualized Intraluminal Support and LVIS Jr (LVIS and LVIS Jr – P170013);</li> <li>- MicroVention Low-Profile Visualized Intraluminal Support EVO (LVIS EVO) (P170013/S004).</li> </ul> <p>In the FCsize module:</p> <ul style="list-style-type: none"> <li>- Medtronic Axium Detachable Coil and Axium Prime Detachable Coil (K233420);</li> <li>- MicroVention HydroCoil Embolic System (HES – K161367);</li> <li>- Stryker Target and Target XXL Detachable Coils (K161429);</li> <li>- Stryker Target Tetra Detachable Coils (K222533);</li> <li>- Balt Optima Coil System (K223386).</li> <li>- <b>Phenox-Wallaby Avenir Coil System (K173711)</b></li> </ul>	
Results output	Simulation report in PDF and DCM format. Arterial reconstruction with the deployed device in DCM and VTP format.	Simulation report in PDF and DCM format. Arterial reconstruction with the deployed device in DCM, VTP, and <b>DICOM</b> format.	Similar

## 7. Performance Data

Software verification and validation testing were conducted, and documentation was provided as recommended by the FDA's Guidance for Industry and FDA Staff, "*Content of Premarket Submissions for Device Software Functions*".

Non-clinical bench performance testing was conducted, and documentation was provided as recommended by the FDA's Guidance for Industry and FDA Staff, "*Recommended Content and Format of Non-Clinical Bench Performance Testing Information in Premarket Submissions*".

The testing program comprised the following elements:

- Verification testing, which compares the predictive behavior of the implantable medical device with its theoretical behavior.
- Experimental bench testing, which compares the device's geometrical parameters with the simulation model. This is done by performing optical acquisition of implantable medical device samples in both unconstrained and constrained configurations.
- Retrospective *in vivo* testing, which compares the *in vitro* retrospective cases with the virtual device simulation.

## 8 Conclusion

The subject and predicate devices are substantially equivalent in technical characteristics, general function, and application. The results of the verification and validation tests demonstrate that the Sim&Size device performs as intended in the specified use conditions. New features added in the subject device do not raise new questions of safety and effectiveness.