K030067

2.0 General Information

JAN 1 7 2003

2.1 510(k) Summary

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

The assigned 510(k) number is:

2.1.1 Submitter Name, Address, Contact

Ortho-Clinical Diagnostics, Inc. 100 Indigo Creek Drive Rochester, New York 14626-5101 (585) 453-3154

Contact Person:

Sarah Parsons

2.1.2 Preparation Date

Date 510(k) Summary Prepared:

2.1.3 Device Name

Trade or Proprietary Name: Vitros Immunodiagnostic Products HBsAg

Controls

Common Name:

HBsAg Controls

Classification Name: 21CFR 862.1660 Quality Control Material (Assayed

and Unassayed).

2.1.4 Predicate Device

The *Vitros* Immunodiagnostic Products HBsAg Controls are substantially equivalent to Ortho-Clinical Diagnostics, Inc. *Vitros* Immunodiagnostic Products HBsAg Controls (K011250).

2.1.5 Device Description

The *Vitros* Immunodiagnostic System uses luminescence as the signal in the qualitative detection of HBsAg in human plasma and serum. Coated microwells are used as the solid phase separation system.

The system is compromised of three main elements:

- 1. The *Vitros* Immunodiagnostic Products range of products, in this case *Vitros* Immunodiagnostic Products HBsAg Reagent Pack and *Vitros* Immunodiagnostic Products Calibrator, which are combined by the *Vitros* Immunodiagnostic System to perform a *Vitros* assay. The *Vitros* Immunodiagnostic Products HBsAg Reagent Pack and Calibrator have been approved for sale (PMA P000044).
- 2. The *Vitros* Immunodiagnostic System- instrumentation, which provides automated use of the immunoassay kits. The *Vitros* Immunodiagnostic System was cleared be market by a separate 510(k) pre-market notification (K962919).
- 3. Common reagents used by the *Vitros* System in each assay. The *Vitros* Immunodiagnostic Products Signal Reagent and the *Vitros* Immunodiagnostic Products Universal Wash Reagent were cleared as part of the *Vitros* Immunodiagnostic Products Total T3 510(k) premarket notification (K964310).

The *Vitros* System and common reagents are dedicated specifically only for use with the *Vitros* Immunodiagnostic Products range of immunoassay products.

2.1.6 Device Intended Use

The *Vitros* Immunodiagnostic Products HBsAg Controls are intended for use in monitoring the performance of the *Vitros* ECi Immunodiagnostic System when used for the *in vitro* qualitative detection of Hepatitis B Surface Antigen (HBsAg) in human serum and plasma (EDTA, heparin or citrate). The performance of the *Vitros* Immunodiagnostic Products HBsAg Controls has not been established with any other HBsAg assays.

2.1.7 Comparison to Predicate Device

The *Vitros* Immunodiagnostic Products HBsAg Controls are substantially equivalent to Ortho-Clinical Diagnostics Inc. *Vitros* Immunodiagnostic Products HBsAg Controls (K011250)

Table 1 lists the characteristics of the *Vitros* Immunodiagnostic Products HBsAg Controls (new device) and the *Vitros* Immunodiagnostic Products HBsAg Controls (predicate device).

Table 1. Comparison of New Device and Predicate Device

Device	Vitros Vitros		
Characteristic	Immunodiagnostic	Immunodiagnostic	
	Products HBsAg	Products HBsAg	
	Controls	Controls	
	(New device)	(Predicate device)	
Intended use	For use in monitoring the	For use in monitoring the	
	performance of the Vitros	performance of the Vitros	
	ECi Immunodiagnostic	ECi Immunodiagnostic	
	System when used for the	System when used for the	
	in vitro qualitative	<i>in vitro</i> qualitative	
	detection of Hepatitis B	detection of Hepatitis B	
	Surface Antigen (HBsAg)	Surface Antigen (HBsAg)	
	in human serum and	in human serum. The	
	plasma (EDTA, citrate	performance of the <i>Vitros</i>	
	and heparin). The	Immunodiagnostic	
	performance of the Vitros	Products HBsAg Controls	
	Immunodiagnostic	has not been established	
	Products HBsAg Controls	with any other HBsAg	
	has not been established	assays.	
	with any other HBsAg		
	assays.		
Matrix of controls	Human serum with added	Human serum with added	
	constituents of human	constituents of human	
	origin and antimicrobial	origin and antimicrobial	
	agents	agents	
Control level	Positive and negative	Positive and negative	
Expected values	Each control has a quoted	Each control has a quoted	
	mean value derived form	mean value derived form	
	a minimum of 10 assays	a minimum of 10 assays	
	and a standard deviation	and a standard deviation	
	anticipated for single	anticipated for single	
	determinations of each	determinations of each	
	control in a number of	control in a number of	
	different laboratories	different laboratories	
	using different reagent	using different reagent	
	lots. Values are lot	lots. Values are lot	
	specific.	specific.	

2.1.8 Summary of Assessment of Performance Data

The new device is the same as the predicate device with the addition of data that supports a modification in the Intended Use. Additional data demonstrate the use of the controls to assess the performance of the *Vitros* Immunodiagnostic Products HBsAg assay when determining the qualitative

detection of HBsAg in human plasma (EDTA, citrate and heparin) in addition to serum. Included in this 510(k) submission is the technical report that describes the assessment of multiple samples collected as either whole serum or in the presence of EDTA citrate or heparin. Each sample was tested either unspiked (negative) or spiked with known positive HBsAg plasma close to the weak positive detection level of the assay (positive). The results showed that all samples (serum, EDTA, citrate or heparin) behave similarly in the assay supporting that the controls can monitor the assay performance regardless of the tested sample matrix (serum or plasma).

DEPARTMENT OF HEALTH & HUMAN SERVICES



Food and Drug Administration 2098 Gaither Road Rockville MD 20850

JAN 1 7 2003

Ms. Sarah Parsons Associate, Regulatory Affairs Ortho-Clinical Diagnostics, Inc. Regulatory Affairs MC00882 100 Indigo Creek Drive Rochester, NY 14626-5101

Re:

k030067

Trade/Device Name: Vitros Immunodiagnostic Products HBsAg Controls

Regulation Number: 21 CFR 862.1660

Regulation Name: Quality Control Material (Assayed and Unassayed).

Regulatory Class: Class I Product Code: MJX Dated: January 6, 2003 Received: January 7, 2003

Dear Ms. Parsons:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

Page 2 -

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (301) 594-3084. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/dsma/dsmamain.html.

Sincerely yours,

Steven I. Gutman, M.D., M.B.A.

Director

Office of In Vitro Diagnostic Device

Steven Butman

Evaluation and Safety

Center for Devices and

Radiological Health

Enclosure

SUMMARY INFORMATION 3.0

3.1	Statemen	4 ~ £	T-4	I III
1	Statemen	1 ()1	intenaec	i i ise

Page 1 of 1

510(k) Number (if known): <u>KD 30067</u>

Device Name:

Vitros Immunodiagnostic Products HBsAg Controls

Indications for Use:

For use in monitoring the performance of the Vitros ECi Immunodiagnostic System when used for the *in vitro* qualitative detection of Hepatitis B Surface Antigen (HBsAg) in human serum and plasma (EDTA, citrate and heparin). The

performance of the Vitros Immunodiagnostic Products HBsAg Controls has not been established with any other HBsAg assays.

For in vitro diagnostic use.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

(Division Sign-Off)

Division of Clinical Laboratory Devices

510(k) Number <u>K030067</u>

Prescription Use

OR

Over-The-Counter Use _