



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

MAR 1 2 2004

Mr. Asher Kassel  
Director of Regulatory Affairs  
Mennen Medical, Ltd.  
4 Hayarden Street, Yavne 81228  
P.O. Box 102  
Rehovot 76100  
ISRAEL

Re: K030722

Trade/Device Name: Envoy Patient Monitor  
Regulation Number: 21 CFR 870.1025  
Regulation Name: Physiological Patient Monitor (with arrhythmia detection or alarms)  
Regulatory Class: III  
Product Code: MHX, BZK  
Dated: January 14, 2004  
Received: February 2, 2004

Dear Mr. Kassel:

This letter corrects our substantially equivalent letter of February 26, 2004, regarding the Envoy Patient Monitor. Our letter identified the regulation name as Continuous Ventilator. This is in error; the correct regulation name is Physiological Patient Monitor (with arrhythmia detection or alarms) as indicated above.

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting

your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050. This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,



Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

## Indications for Use

510(k) Number: K030722

Device Name: Envoy Patient Monitor

### Indications for Use:

The Envoy is intended for use as a multiparameter monitoring system.

The Envoy can monitor ECG/heart rate, two invasive blood pressure channels, two temperature channels, pulse oximetry, respiration, non-invasive blood pressure, etCO<sub>2</sub> and Spirometry. This effectively allows the Envoy to monitor a wide-range of adult, pediatric and neonatal patient conditions, in many different areas of the hospital.

The Envoy also includes a spirometer intended for continuous monitoring of adult and pediatric patients with tidal volumes greater than 100mL. Functions include display of multiparameter waveforms, vital signs, and status messages.

Prescription Use  AND/OR Over-The-Counter Use \_\_\_\_\_  
(Part 21 CFR 801 Subpart D) (21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER  
PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

FEB 26 2004

Food and Drug Administration

Page 1 of 6

ENVOY Patient Monitor – 510(k) for new Spirometry module



Mennen Medical Ltd.,  
 4 Hayarden Street, Yavne 81228  
 PO Box 102, Rehovot 76100  
 Israel

Tel.: +972-8-9323333

Fax: +972-8-9328510

Date: 24 February, 2004

Topic: **510(k) Safety and Effectiveness Summary as per 21 CFR Section 807.92(c)**  
**Envoy Patient Monitor - new Spirometry module:**

**Establishment Name, Registration Number and Address:**

Name: Mennen Medical Ltd.  
 Registration Number: 9611022  
 Operator Number: 9011766  
 Address: 4 Ha-yarden Street, Yavne, 81228, Israel  
 Postal Address: PO Box 102,  
 Rehovot, 76100, Israel

Tel: +972-8-9323333

Fax: +972-8-9328510

Contact person: Asher Kassel, Director of Regulatory Affairs

To: Food and Drug Administration  
 Center for Devices and Radiological Health  
 Document Mail Center (HFZ-401)  
 9200 Corporate Boulevard  
 Rockville MD, 20850

Attn.: Document Control Clerk  
 From: Asher Kassel, Director of Regulatory Affairs

**Product Name:**

Proprietary: ENVOY

Common: Physiological Patient Monitor

Mennen Medical Part Number: 550-010-000 (full system)

554-000-010 (CPU only)

New Envoy Spirometry module P/N: 551-137-000

ENVOY Patient Monitor – 510(k) for new Spirometry module

**FDA Classification of Envoy Patient Monitor:**

Classification Name: Arrhythmia Detector and Alarm  
Classification Number: 21 CFR 870.1025  
Classification: Class III  
Product Code: 74 DSI

**FDA Classification of new Spirometry module:**

Classification Name: Monitoring spirometer  
Classification Number: 21 CFR 868.1850  
Classification: Class II  
Product Code: BZK

**Performance Standards:**

None promulgated

**Voluntary Standards:**

**\*IEC 60601-1:**

General Requirement for Safety for Medical Electrical Systems - part 1, (1988);  
Amendment 1 – 1991-11;  
Amendment 2 – 1995-03

**\*IEC 60601-1-2 (2001):**

General Requirement for Safety  
Collateral Standard: Electromagnetic compatibility - Requirements and tests.

**Predicate Device:**

METEOR Respiratory Mechanics Monitor(s) - K011784.

**Device Description - Envoy Patient Monitor:**

The Envoy is a multiparameter physiological patient monitor consisting of a main processing unit, a mountable color monitor, and a module rack housing the various Mennen Medical plug-in *vital signs* modules. The modules monitor the patient's vital signs. Up to six internal modules can be plugged into a module rack. The Envoy can accommodate two module racks. The vital sign data derived from the modules by the Envoy are presented on the monitor as waveform and numeric displays.

The Envoy vital signs modules acquire vital signs data from the patient, and display their waveforms and alarms indications on the Envoy display unit. Vital signs and waveform information are displayed simultaneously on the Envoy Display Unit. Up to 8 traces can be displayed at any one time.

The Envoy is a reusable, software driven, patient monitor, intended for use as part of a physiological monitoring system in a hospital environment. As such, it is not a life supporting, nor life-sustaining device; nor is it implantable and therefore sterility is not a consideration.

ENVOY Patient Monitor – 510(k) for new Spirometry module

The Envoy is not a kit and does not contain any drug or biological products. The **Spirometry** module of the Envoy patient monitor is not sold as a stand-alone spirometry device, but as part of a multiparameter physiological patient monitoring system (ENVOY). In Chapter 1 of the Envoy Operating Manual, the following **Prescription Notice** appears: “Federal United States law restricts the sale and use of this instrument to qualified medical personnel only”.

### **Intended Use and Functional Description and of the New Envoy Spirometry module:**

#### **Intended Use of the Spirometry module:**

The Spirometry module is used to provide an objective measurement of lung function. The Spirometry module is intended for use in the hospital clinical environment only (e.g. in the RICU). The module is used for the continuous monitoring of mechanically ventilated **adult** and **pediatric** patients; or for **adult** and **pediatric** patients able to breathe spontaneously. In both cases, breathing is via a tracheal tube or mask. The module is to be used only for patients with tidal volumes greater than 100 ml.

#### **Functional Description:**

The spirometry module can measure the following:

- a) Expiratory Vital Capacity (EVC): The maximum volume of gas which can be expired from the lungs during a relaxed expiration from a position of full inspiration.
- b) Inspiratory Vital Capacity (IVC): The maximum volume of gas which can be inspired into the lungs during a relaxed inspiration from a position of full expiration. The expiratory phase is the one more commonly used to measure obstruction and restriction within the lungs.

### **Intended Use of the Envoy Patient Monitor:**

The Envoy is intended for use as a multiparameter monitoring system.

The Envoy can monitor ECG/heart rate, two invasive blood pressure channels, two temperature channels, pulse oximetry, respiration, non-invasive blood pressure, etCO<sub>2</sub> and Spirometry. This effectively allows the Envoy to monitor a wide-range of adult, pediatric and neonatal patient conditions, in many different areas of the hospital.

The Envoy also includes a spirometry module intended for continuous monitoring of adult and pediatric patients with tidal volumes greater than 100mL. Functions include display of multiparameter waveforms, vital signs, and status messages.

ENVOY Patient Monitor – 510(k) for new Spirometry module

**Summary of the technological characteristics of the new Envoy Spirometry module (incorporating OEM technology for a Respiratory Mechanics System):**

The following tables summarize data on the Mennen Medical new Envoy Spirometry module:

<b>Envoy Spirometry Module</b>	
<b>Part Number:</b>	551-137-000
<b>Monitored Parameters:</b>	*Air Flow signal *Volume signal *Airway pressure
<b>Module size:</b>	<b>Single slot</b> Height: 10.0cm (4.0 in) Width: 4.0 cm (1.6 in) Depth: 14.0 cm (5.5 in)

<b>Monitored parameters</b>	<b>Envoy Spirometry module</b>
Vt_e	Yes
PIP	Yes
PEEP	Yes
MAP	Yes
Plat	Yes
PF_i	Yes
PF_e	Yes
RR	Yes
Ve	Yes
I:E	Yes
COMP	Yes
RES	Yes
RSBI	Yes

ENVOY Patient Monitor - 510(k) for new Spirometry module

<b>Alarm Indications</b>	<b>Envoy Spirometry module</b>
Clinical Alarms for derived parameters	Visual & Sound Alarms

<b>Specifications</b>	<b>Envoy Spirometry module</b>
Degree of protection against electrical shock	Type BF applied part
Sampling	100 sps

<b>OEM Airway Adapter/ Flow sensor</b>	<b>Envoy Spirometry module</b>
Flow Range:	2 - 180 lpm (33 - 3000 ml/s)
Accuracy:	± 5 % reading, or 0.5 lpm
Dead Space:	6.9 ml
Connections:	Airway - 15 mm ID/22 mmOD patient end by 15 mm ID ventilator end (ISO 5356-1) Proprietary Smart connector. Tri-Tubing - triple 0.055" inch ID lumen
Length:	2.5 inches (adapter) - 6 feet (tubing)
Weight:	6.3 grams (minus tubing)
Material:	Sensor - Polycarbonate (Makrolon). Tubing - Medical grade Polyvinyl Chloride

ENVOY Patient Monitor – 510(k) for new Spirometry module

**Conclusion of comparison of technological characteristics:**

We consider the Envoy Spirometry module to be substantially equivalent to the METEOR 200 Respiratory Mechanics Monitor; we submit that any differences between the two monitors do not raise any new issues of safety and effectiveness.

**Testing**

The Envoy Spirometry module has been subject to extensive safety and performance testing. Final testing for the system included various performance tests designed to ensure that the device meets all functional requirements and performance specifications. Safety testing and EMC testing were performed by an independent testing laboratory to ensure that the device complies to applicable industry and safety standards.

Signature: *Asher Kassel*

Asher Kassel  
Director of Regulatory Affairs  
Mennen Medical Ltd.

Tel: +972-8-9323311 (direct)  
Fax: +972-8-9328510  
E-mail: asher@mimi.co.il