

AUG 27 2004

## 510(k) Summary

USFilter/Ionpure, Inc. in accordance with the requirements of 21 CFR§807.92 is submitting this summary of 510(k)-safety and effectiveness information.

1 **Submitter's Identification**

USFilter/Ionpure, Inc.  
725 Wooten Road  
Colorado Springs, CO 80915  
719-570-9600

**Date Summary Prepared**

26 March 2003 , (Rev.4 – 02 August 2004)

2 **Official Correspondent/Contact**

Douglas L. Rufenacht  
725 Wooten Road  
Colorado Springs, CO 80915  
719-550-2055`

3 **Name of the Device**

USFilter Bicarbonate Mixing and Dispensing (Bicarb) System

4 **Substantial Equivalence**

MAR COR Services, Inc.  
MAR COR Services Bicarb System  
K003560

5 **Device Description**

The consistent mixing and delivery of bicarbonate solution for use in hemodialysis

6 **Intended Use**

The USFilter Bicarbonate Mixing and Dispensing (Bicarb) System is intended to be used in a Hemodialysis facility to allow safe and effective delivery of concentrates as well as central mixing and delivery of sodium bicarbonate solution necessary for a hemodialysis treatment.

7 **Comparison to Substantial Equivalence Devices**

The USFilter Bicarbonate Mixing and Dispensing (Bicarb) System unit contains the same fundamental scientific technology as the substantial equivalent device. Please reference the following Chart of Comparison:

Specifications and Technical Information	USFilter BICARB50SYS & BICARB100SYS	MAR COR MCB 210-50 & MCB 210-100
Indications for Use:	The USFilter Bicarbonate Mixing and Dispensing (Bicarb) System is intended to be used in a Hemodialysis facility to allow safe and effective delivery of concentrates as well as central mixing and delivery of sodium bicarbonate solution necessary for a hemodialysis treatment.	The Bicarb Mixing and Distribution System (Bicarb System) designed and manufactured by Mar Cor Services is, safe, reliable and easy to use. This system will reduce staff time by providing consistent bicarb mixing, and through system cleaning and disinfection.
<b>Utility Requirements</b>		
Water Feed:	RO or DI water produced to AAMI standards for Hemodialysis	RO or DI water produced to AAMI standards for Hemodialysis
Electrical requirements:	115 VAC, Single Phase, 20 AMP	115 VAC, Single Phase, 20 AMP
Drain:	1" Minimum	1" Minimum
Physical Dimensions	50                      100	50                      100
Length:	70"                      83"	85"                      85"
Depth:	27"                      35"	32"                      32"
Height:	60"                      60"	57"                      61"
Operating Weight:	1800 LBS              2500 LBS	1300 LBS              2050 LBS
<b>Dimensional/Operating Room</b>		
Width:	6' 10"	6' 10"
Depth:	3' 3"	3' 3"
Height:	6' 0"	6' 0"
<b>Installation Piping</b>		
Feedwater:	½" schedule 80 PVC	¾" schedule 80 PVC
Loop Distribution:	¾" schedule 80 PVC	¾" schedule 80 PVC
Loop Return:	¾" schedule 80 PVC	¾" schedule 80 PVC
Drain:	1" schedule 80 PVC	1" schedule 80 PVC
<b>System Floats</b>		
Mix Tank Float:	The mix tank is equipped with a high level float that will close the fill valve if the water in the tank exceeds 100 gallons (50 gallons for 50 gallon systems). The fill valve controller will not reactivate unless the water level is below the float level.	The mix tank is equipped with a high level float that will close the fill valve if the water in the tank exceeds 100 gallons (50 gallons for 50 gallon systems). The fill valve controller will not reactivate unless the water level is below the float level.
Distribution Tank Float:	The distribution tank is equipped with a low-level float that lights a warning light if the bicarbonate in the tank drops below 25 gallons. If the distribution tank is allowed to go empty an audible alarm will sound.	The distribution tank is equipped with a low-level float that lights a warning light if the bicarbonate in the tank drops below 25 gallons. If the distribution tank is allowed to go empty an audible alarm will sound.
<b>Pumps</b>		
Mix Pump:	¾ hp Totally Enclosed Fan Cooled (TEFC) motor	¾ hp Totally Enclosed Fan Cooled (TEFC) motor
Distribution Pump:	1/9 hp TEFC motor	1/8 hp TEFC motor
Tank Material:	Polyethylene	Polyethylene

8 **Discussion of Non-Clinical Test Performed for Determination of Substantial Equivalence is as follows:**

The USFilter Bicarbonate Mixing and Dispensing (Bicarb) System is an existing device actively marketed by USFilter with a long performance record. This submission seeks to obtain approval so as to begin marketing this device as part of a water system for use in hemodialysis.

9 **Discussion of Clinical Test Performed:**

N/A

10 **Biocompatibility Testing:**

USFilter certifies that all materials utilized in the design, manufacture, and operation of the USFilter Bicarbonate Mixing and Dispensing (Bicarb) System are utilized in the predicate device as well as used in other registered devices and systems currently sold by USFilter and registered with the Food and Drug Administration. The tables in appendix D provide the detailed break down of these materials as well as a documented review in accordance with Appendix C of FDA document G95-1. It is the contention of USFilter that as per the results of the review documented in appendix D no further testing is required.

11 **Conclusions:**

As was true for the MAR COR Services Bicarb System, the USFilter Bicarbonate Mixing and Dispensing (Bicarb) System is intended to provide consistent bicarbonate mixing and distribution for Hemodialysis applications.



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

**AUG 27 2004**

Mr. Douglas Rufenacht  
Validation Manager and Regulatory Compliance  
USFilter  
Process Water Systems  
High Purity Systems  
725 Wooten Road  
COLORADO SPRINGS CO 80915

Re: K031502  
Trade/Device Name: USFilter Bicarbonate Mixing and Dispensing (Bicarb) System  
Regulation Number: 21 CFR §876.5820  
Regulation Name: Hemodialysis system and accessories  
Regulatory Class: II  
Product Code: 78 FIN  
Dated: August 20, 2004  
Received: August 23, 2004

Dear Mr. Rufenacht:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at one of the following numbers, based on the regulation number at the top of the letter:

8xx.1xxx	(301) 594-4591
876.2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4616
884.2xxx, 3xxx, 4xxx, 5xxx, 6xxx	(301) 594-4616
892.2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4654
Other	(301) 594-4692

Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97) you may obtain. Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,



Nancy C. Brogdon  
Director, Division of Reproductive,  
Abdominal and Radiological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

## Indication for Use

The USFilter Bicarbonate Mixing and Dispensing (Bicarb) System is intended to be used in a Hemodialysis facility to allow safe and effective delivery of concentrates as well as central mixing and delivery of sodium bicarbonate solution necessary for a hemodialysis treatment.

*Prescription Use* \_\_\_\_\_ ✓

*David A. Reynolds*  
\_\_\_\_\_  
(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices  
510(k) Number K031502