

SEP 23 2003

Eye-deas, LLC

EyesOPen  
Manual Tonometer

510(k) Submission

### 510(k) Summary

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of 21 CFR 807.92.

The assigned 510(k) number is: K032769

Submitted by: Steven E. Feldon, M.D.  
President  
Eye-deas, LLC  
39 Sandringham Road  
Rochester, New York 14610

Date Prepared: July 21, 2003

### Device Name

Proprietary Name:	EyesOPen Tonometer
Common Name:	Ocular Tonometer
Classification Name:	Tonometer, Manual

### Substantial Equivalence

The Eye-deas EyesOPen tonometer is similar in indications, optical design and IOP calculation methodology to the Perkins tonometer.

The Eye-deas EyesOPen tonometer is similar in indications, size, user interface, microprocessor operation and solid-state force transducer design to the Tono-Pen 3 tonometer.

### Device Description

The Eye-deas EyesOPen tonometer is a precision electronic tonometer, which measures intraocular pressure (IOP). The body of the instrument is specially designed to fit comfortably in the user's hand, facilitating fast and accurate measurements.

The EyesOPen tonometer contains a solid-state force transducer and a CMOS image transducer that together produces electrical signals from which intraocular pressure is derived. When sufficient valid readings are obtained an IOP reading is displayed on a liquid crystal display.

Under each of the four digits are bars. The bars correspond to the goodness of fit to a linear regression of pressure against applanated diameter. The regression coefficient is 95% or better if the first bar is

illuminated, 90%-95% if the second bar is illuminated, 80%-90% if the third bar is illuminated, or less than 80% if the fourth bar is illuminated.

**Intended Use**

The Eye-deas EyesOPen tonometer is used to measure the intraocular pressure (IOP) during routine eye examinations or when an increased IOP is suspected.

**Substantial Equivalence Comparison – Perkins Tonometer**

	<b>EyesOPen tonometer</b>	<b>Perkins tonometer</b>
<b>Indication</b>	<b>Intraocular pressure (IOP) measurement</b>	<b>Same</b>
<b>Design</b>	<b>Hand-held microprocessor based</b>	<b>Hand-held manual dial</b>
<b>Force transducer</b>	<b>Solid state element</b>	<b>Spring and weight</b>
<b>Imaging transducer</b>	<b>Electronic imager</b>	<b>Operator's eye</b>
<b>IOP algorithm basis</b>	<b>Goldmann</b>	<b>Same</b>
<b>Calibration</b>	<b>None required</b>	<b>External weights</b>
<b>Contact tip</b>	<b>6mm glass coaxial</b>	<b>6 mm glass transverse</b>
<b>User interface</b>	<b>Briefly touched against eye</b>	<b>Held on patient's eye while dial is adjusted</b>
<b>Display</b>	<b>4 digit LCD</b>	<b>Scribed dial</b>
<b>Range of measurement</b>	<b>5 – 90 mmHg</b>	<b>5 –50 mmHg</b>
<b>Statistical reliability of the reading</b>	<b>The bars correspond to the goodness of fit to a linear regression of pressure against applanated diameter. The regression coefficient is 95% or better if the first bar is illuminated, 90%-95% if the second bar is illuminated, 80%-90% if the third bar is illuminated, or less than 80% if the fourth bar is illuminated.</b>	<b>No indication of reading reliability</b>
<b>Versatility</b>	<b>Patient can be in any position.</b>	<b>Same.</b>
<b>Weight (with batteries)</b>	<b>3 ounces</b>	<b>12.4 ounces</b>
<b>Dimensions</b>	<b>1" H x 1 ¼" W x 7 " L</b>	<b>1 ½" H x 1 ½" W x 11" L</b>
<b>Power source</b>	<b>2 ea 3.0 volt lithium batteries</b>	<b>4 ea. 1.5 volt AA batteries</b>

Substantial Equivalence Comparison – Tono-Pen 3 Tonometer

	EyesOPen tonometer	Tono-Pen 3
Indication	Intraocular pressure (IOP) measurement	Same
Design	Hand-held microprocessor based	Same
Force Transducer	Solid-state	Same
Imaging Transducer	Electronic CMOS imager	No imaging elements
IOP calculation basis	Goldmann	Mackay-Marg
Calibration	None required	None required
Contact tip	6mm glass coaxial	3 mm stainless steel coaxial
Measurement technique	Briefly touched against eye	Same
User control	Single push button	Same
Single use sanitary tip cover	Required	Not required
Display	4 digit LCD	Same
Range of Measurement	5 – 90 mmHg	Same
Statistical Reliability	The bars correspond to the goodness of fit to a linear regression of pressure against applanated diameter. The regression coefficient is 95% or better if the first bar is illuminated, 90%-95% if the second bar is illuminated, 80%-90% if the third bar is illuminated, or less than 80% if the fourth bar is illuminated.	The bars correspond to standard deviation expressed as a percentage of the mean pressure. The first bar is illuminated at <5%, the second bar is illuminated at 5% to 10%, the third bar is illuminated at 10% to 20%, the fourth bar is illuminated at >20%.
Versatility	Patient can be in any position.	Same.
Weight (with batteries)	3 ounces	2.3 ounces
Dimensions	1" H x 1 ¼" W x 7 " L	1" H x 1" W x 7 3/4" L
Power Source	2 ea 3 volt lithium batteries	2 ea. 1/3N 3 volt lithium batteries

**Performance Testing**

Although not required to establish substantial equivalence, the EyesOPen tonometer was tested to establish baseline performance. The EyesOPen tonometer was tested on a simulated human eye. The EyesOPen tonometer performance was consistent with the performance of the predicate tonometer devices on this simulator.



\_\_\_\_\_  
Submitter's Signature

7-28-03

\_\_\_\_\_  
Date



SEP 23 2003

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Eye-deas, LLC  
c/o Erin Sparnon  
CITECH  
5200 Butler Pike  
Plymouth Meeting, PA 19462

Re: K032769  
Trade/Device Name: EyesOPen Tonometer  
Regulation Number: 21 CFR 886.1930  
Regulation Name: Tonometer and accesories  
Regulatory Class: Class II  
Product Code: HKY  
Dated: September 5, 2003  
Received: September 8, 2003

Dear Ms. Sparnon:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4613. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

A handwritten signature in cursive script that reads "A. Ralph Rosenthal".

A. Ralph Rosenthal, M.D.  
Director  
Division of Ophthalmic and Ear,  
Nose and Throat Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

510(k) Number (if known): K032769  
Device Name: EyesOPen tonometer (manual tonometer)


Indications for Use:

The EYE-deas EyesOPen tonometer is used to measure the intraocular pressure (IOP) during routine eye examinations or when an increased IOP is suspected. This device is intended for use by Ophthalmologists, Optometrists, and other trained medical professionals.

(PLEASE DO NOT WRITE BELOW THIS LINE- CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
(Division Sign-Off)  
Division of Ophthalmic Ear,  
Nose and Throat Devices

510(k) Number K032769

Prescription Use   
(Per 21CFR801.109)

OR

Over-the-Counter Use