

K033576

# OMNI-FLOW® Medication Management System™

LifeShield® Latex-Free Primary IV Pump Set  
Distal Microbore Patient Line, Convertible Pin, 72 Inch,  
with 2 Pressure-Activated Anti-Siphon Valves, CLAVE® and  
OPTION-LOK®

DEC - 4 2003

## Special 510(K) Summary

**1. Name of Submitter:** Abbott Laboratories  
Hospital Products Division  
200 Abbott Park Road  
Abbott Park, Illinois 60064-6133  
  
Establishment Registration # 1415939

**2. Manufacturer and Establishment Registration Number:**

Manufacturer  
Abbott Hospitals, Limited  
Parque Industrial Itabo, S.A.  
Haina, San Cristobal  
Dominican Republic

Sterilization Site  
Abbott Laboratories  
Hospital Products Division  
Hwy. 301 North  
Rocky Mount, NC 27801

Establishment Registration # 9613251

Establishment Registration # 1021343

**3. Proprietary or Trade Name:** LifeShield® Latex-Free Primary IV Pump Set Distal Microbore Patient Line, Convertible Pin, 72 Inch, with 2 Pressure-Activated Anti-Siphon Valves, Inline Cassette, Collection Bag, CLAVE® and OPTION-LOK® (LN 12566)

**4. Common Name:** Set, Administration Intravascular

**5. Device Classification, Pancode and ProCode:** Class II, 80, FPA

**6. Performance Standards:** Performance standards have been established under Section 514 of the Food, Drug and Cosmetic Act for Set, Administration Intravascular. Set Administration Intravascular can be found in 21 CFR 880.5440.

**7. Intended Use:**

LifeShield® Latex-Free Primary IV Pump Set Distal Microbore Patient Line, Convertible Pin, 72 Inch, with 2 Pressure-Activated Anti-Siphon Valves, Inline Cassette, Collection Bag, CLAVE® and OPTION-LOK® is a single use device for the delivery of fluids from a container to a patient's vascular system.

The intended use and indications for use are the same as other commercially available IV administration sets designated for use with the Omni-Flow® Medication Management Systems™. This set is not intended for gravity use.

**8. Indications for Use:**

The device is used with the Omni-Flow® Medication Management System™ for intravenous infusion by or under the order of a licensed medical practitioner. The device can be used to simultaneously infuse up to four solutions/medications from both syringes and IV fluid containers.

# OMNI-FLOW<sup>®</sup> Medication Management System<sup>™</sup>

## LifeShield<sup>®</sup> Latex-Free Primary IV Pump Set Distal Microbore Patient Line, Convertible Pin, 72 Inch, with 2 Pressure-Activated Anti-Siphon Valves, CLAVE<sup>®</sup> and OPTION-LOK<sup>®</sup>

---

### 9. Device Description:

The LifeShield<sup>®</sup> Latex-Free Primary IV Pump Set Distal Microbore Patient Line, Convertible Pin, 72 Inch, with 2 Pressure-Activated Anti-Siphon Valves, Inline Cassette, Collection Bag, CLAVE<sup>®</sup> and OPTION-LOK<sup>®</sup> is designed to be used exclusively with Omni-Flow<sup>®</sup> Medication Management Systems<sup>™</sup>. The set is equipped with two pressure-activated anti-siphon valves. The two valves limit gravity flow (free flow) to 1ml per hour when the set is primed and attached to a solution bag hanging vertically at full set extension (72 inches).

### 10. Statement of Substantial Equivalence:

The LifeShield<sup>®</sup> Latex-Free Primary IV Pump Set Distal Microbore Patient Line, Convertible Pin, 72 Inch, with 2 Pressure-Activated Anti-Siphon Valves, Inline Cassette, Collection Bag, CLAVE<sup>®</sup> and OPTION-LOK<sup>®</sup> is substantially equivalent to the LifeShield<sup>®</sup> Latex-Free Primary IV Pump Set with Distal Microbore Patient Line and Gravity Flow Prevention Valve, Convertible Pin, 110 Inch, with CLAVE<sup>®</sup> and OPTION-LOK<sup>®</sup> (LN 12162).

The modified primary set (LN 12566) is substantially equivalent to currently marketed primary set (LN 12162) in that :

#### Similarities:

- 1) Both are used for intravenous infusion. with the Omni-Flow<sup>®</sup> Medication Management System<sup>™</sup>.
- 2) Both can be used with the Omni-Flow<sup>®</sup> Medication Management System<sup>™</sup> for intravenous infusion by or under the order of a licensed medical practitioner. Both devices can be used to simultaneously infuse up to four solutions/medications from both syringes and IV fluid containers.
- 3) Both are latex-free primary sets consisting of a Distal Microbore Patient Line, Convertible Pin, Inline Cassette, Collection Bag, (1) Integral Y-High Clave at the distal end of the patient line, and OPTION-LOK<sup>®</sup>.
- 4) Both sets have at least (1) Pressure Activated Anti-Siphon Valve (PAV) to provide free flow protection under specific head height conditions.
- 5) Both sets have components made from the same or similar material.
- 6) Both sets are not intended for use as a gravity set.

#### Differences:

- 1) The modified set contains two PAV's located distal to the cassette. The current set has only (1) PAV located distal to the cassette.
- 2) The modified set has a total tubing length of 72 inches (nominal). The current set has a total tubing length of 110 inches (nominal).
- 3) The PAV's in the modified set contain a comparable silicone material cleared under K790062 (IV Administration Set with Backcheck), and are currently manufactured and sold by B. Braun Medical Inc. as a bulk /non-sterile item.

# OMNI-FLOW<sup>®</sup> Medication Management System<sup>™</sup>

**LifeShield<sup>®</sup> Latex-Free Primary IV Pump Set  
Distal Microbore Patient Line, Convertible Pin, 72 Inch,  
with 2 Pressure-Activated Anti-Siphon Valves, CLAVE<sup>®</sup> and  
OPTION-LOK<sup>®</sup>**

## 11. Predicate Device Information:

Information for IV Administration sets/set components previously cleared for commercial distribution and determined to be appropriate for use as predicates is provided below.

510(k) #	Product Name	Date Submitted	Clearance Date
K790062	IV Administration Set with Backcheck Valve (B. Braun Medical)	01/10/79	04/10/79
K832755	OmniFlow Infusion Pump & IV Sets	08/09/83	11/03/83
K882594	OmniFlow Therapist Infusion System & IV Sets	06/21/88	09/20/88
K915571	Clave Connector (ICU Medical)	12/12/91	09/21/92
K944125	Abbott LTE Infusion Pump & IV Sets	08/16/94	12/06/94
K971293	Lifeshield Primary IV Set with Backcheck Valve	04/04/97	06/27/97

## 12. Comparison to Legally Marketed Device(s)

Factors	Subject Device (LN 12566)	Predicate Device(s) (LN 12162)
<b>Intended Use</b>		
	Single use device for the delivery of fluids from a container to a patient's vascular system.	Same
<b>Indications for Use</b>		
	The device is used with the Omni-Flow <sup>®</sup> Medication Management System <sup>™</sup> for intravenous infusion by or under the order of a licensed medical practitioner. The device can be used to simultaneously infuse up to four solutions/medications from both syringes and IV fluid containers.	Same
<b>Technology</b>		
➤ Basic Operating Principle	Administration/infusion of intravenous medications to a patient's vascular system with the assistance of an infusion pump.	Same
➤ Gravity Flow Protection	Provide gravity flow (free flow) protection at full vertical extension of the set (72 inches).	Provide gravity flow (free flow) protection up to a 36 inch head height.

# OMNI-FLOW<sup>®</sup> Medication Management System<sup>™</sup>

LifeShield<sup>®</sup> Latex-Free Primary IV Pump Set  
 Distal Microbore Patient Line, Convertible Pin, 72 Inch,  
 with 2 Pressure-Activated Anti-Siphon Valves, CLAVE<sup>®</sup> and  
 OPTION-LOK<sup>®</sup>

Factors	Subject Device (LN 12566)	Predicate Device(s) (LN 12162)
<b>Technology (continued)</b>		
➤ IV Line Type	Primary Patient Line	Same
➤ Infusion Pump Compatibility	Omni-Flow <sup>®</sup> Medication Management Systems <sup>™</sup>	Same
➤ Infusion Fluid Types	General Hospital	Same
<b>Design/Components</b>		
➤ Pressure Activated Anti-Siphon Valves (PAV)	(2) PAV's	(1) PAV
➤ Semi-Rigid Adapters (bonds PAV to set)	(4) Semi-Rigid Adapters	(1) Semi-Rigid Adapter
➤ Access Sites	(1) Integral Y-High Clave (distal)	Same
➤ Total Tubing Length	72 inches (Nominal)	110 inches (Nominal)
➤ Tubing Type	Distal Microbore	Same
➤ Convertible Piercing Pin	Yes	Same
➤ CAIR Clamp	Yes	Same
➤ Drip Chamber	Yes	Same
➤ In-Line Integral Cassette	(1) Patient Line Port, (1) Collection Bag Port, (4) Secondary Line Ports	Same
➤ Roller Clamp	Yes	Same
➤ OPTION-LOK <sup>®</sup> Male Adapter	Yes	Same
➤ Slide Clamp	Yes	Same
➤ Collection Bag	Yes	Same
<b>Materials</b>		
➤ Pressure Activated Anti-Siphon Valves (PAV) Housing	Polycarbonate	Same
➤ Pressure Activated Anti-Siphon Valves (PAV), Disc	Silicone	Same
➤ All Other Components	Same	Same
<b>Manufacturing Processes and Sterilization</b>		
All Components	Same	Same
<b>Packaging</b>		
➤ Assembled Set	Paperboard Carton	Same
➤ Shipping Container	Corrugated Fiberboard	Same

# OMNI-FLOW<sup>®</sup> Medication Management System<sup>™</sup>

**LifeShield<sup>®</sup> Latex-Free Primary IV Pump Set  
Distal Microbore Patient Line, Convertible Pin, 72 Inch,  
with 2 Pressure-Activated Anti-Siphon Valves, CLAVE<sup>®</sup> and  
OPTION-LOK<sup>®</sup>**

---

## 13. Summary of Safety and Effectiveness


The LifeShield<sup>®</sup> Latex-Free Primary IV Pump Set Distal Microbore Patient Line, Convertible Pin, 72 Inch, with 2 Pressure-Activated Anti-Siphon Valves, CLAVE<sup>®</sup> and OPTION-LOK<sup>®</sup> (LN 12566) as described in this submission is substantially equivalent to the predicate Lifeshield Primary Set (LN 12162), in that both sets have:

- 1) the same intended use,
- 2) the same indication for use,
- 3) the same fundamental technology and operating principle,
- 4) the same or similar materials of construction for all components,
- 5) the same manufacturing and sterilization processes, and
- 6) the same packaging.

## 14. Statement of Safety and Effectiveness

The LifeShield<sup>®</sup> Latex-Free Primary IV Pump Set Distal Microbore Patient Line, Convertible Pin, 72 Inch, with 2 Pressure-Activated Anti-Siphon Valves, CLAVE<sup>®</sup> and OPTION-LOK<sup>®</sup> (LN 12566) meets the functional claims and intended use as described in the product labeling, and is substantially equivalent to, and as safe and effective as, the LifeShield<sup>®</sup> Latex-Free Primary Pump Set with Distal Microbore Patient Line, Convertible Pin, 110 Inch, with Pressure-Activated Anti-Siphon Valve, Inline Cassette, Collection Bag, CLAVE<sup>®</sup> and OPTION-LOK<sup>®</sup> (LN 12162).

Prepared and submitted November 12, 2003 by:

  
Patricia Melerski

Manager Regulatory Affairs Device Registration  
Abbott Laboratories  
Hospital Products Division D389, J45-2N  
200 Abbott Park Road, IL 60064-6133  
Phone: 847/938-3718  
Fax: 847/938-7867



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

DEC - 4 2003

Abbott laboratories  
Ms. Patricia Melerski  
Manager, Regulatory Affairs Device Registration  
Hospital Products Division  
D-389, Bldg. J45-2N  
200 Abbott Park Road  
Abbott, Park, Illinois 60064-6133

Re: K033576

Trade/Device Name: LifeShield® Latex-Free Primary IV Pump Set Distal Microbore Patient Line, Convertible Pin, 72 Inch, with 2 Pressure-Activated Anti-Siphon Valves, Inline Cassette, Collection Bag, CLAVE® and OPTION- LOK® (LN 12566)  
Regulation Number: 880.5440  
Regulation Name: Intravascular Administration Set  
Regulatory Class: II  
Product Code: FPA  
Dated: November 12, 2003  
Received: November 13, 2003

Dear. Ms. Melerski:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4618. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

  
Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,  
Infection Control and Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

# OMNI-FLOW® Medication Management System™

LifeShield® Latex-Free Primary IV Pump Set  
Distal Microbore Patient Line, Convertible Pin, 72 Inch,  
with 2 Pressure-Activated Anti-Siphon Valves, CLAVE® and  
OPTION-LOK®

---

## Indications for Use Statement

---

510(k)  
Number  
(if known)

---

Device  
Name:

LifeShield® Latex-Free Primary IV Pump Set Distal Microbore Patient Line,  
Convertible Pin, 72 Inch, with 2 Pressure-Activated Anti-Siphon Valves, Inline  
Cassette, Collection Bag, CLAVE® and OPTION-LOK® (LN 12566)

---

Indications  
For Use:

LifeShield® Latex-Free Primary IV Pump Set (LN12566) has the following  
indications for use:

The device is used with the Omni-Flow® Medication Management System™ for  
intravenous infusion by or under the order of a licensed medical practitioner.  
The device can be used to simultaneously infuse up to four  
solutions/medications from both syringes and IV fluid containers.

PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE  
IF NEEDED

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use   
(per 21 CFR 801.109)

OR

Over-The-Counter Use

*Patricia Cucenato*

(Division Sign-Off)  
Division of Anesthesiology, General Hospital,  
Infection Control, Dental Devices

510(k) Number: 4033576