K033989 p1/2

'JAN 2 2 2004

Special 510(k) 510(k) Summary of Substantial Equivalence

SelectSite[™] Model C304 KXXXXX

510(K) SUMMARY OF SUBSTANTIAL EQUIVALENCE

Date Prepared:	22 December 2003	
Submitter:	Medtronic Inc. 7000 Central Avenue NE Minneapolis, MN 55432	
Contact:	Kristy K. Mollner, MS, RAC Sr. Regulatory Affairs Specialist (763) 514-4189	
Telephone: Fax:		
	(763) 514-6424	
Email:	kristina.mollner@medtronic.com	
Proprietary Name	Medtronic [®] SelectSite TM Model C304-S59 Deflectable Catheter System	
	Medtronic [®] SelectSite TM Model C304-L69 Deflectable Catheter System	
Common Name:	Catheter, Percutaneous	
Device Classification:	Class II, 21 CFR § 870.1250	
Product Code:	74 DQY	

Device Description

The Medtronic SelectSite Model C304 Deflectable Catheter System is designed to facilitate the introduction of transvenous devices within the chambers and coronary vasculature of the heart and the introduction of leads in the left heart via the coronary sinus.

The Model C304-S59 deflectable catheter has a single lumen for passage of transvenous devices with a maximum outer diameter of 5.6 French (1.9 mm) and a usable length of 30 cm.

The Model C304-L69 deflectable catheter has a single lumen for passage of transvenous devices with a maximum outer diameter of 5.6 French (1.9 mm) and a usable length of 40 cm.



Special 510(k) 510(k) Summary of Substantial Equivalence SelectSite[™] Model C304 KXXXXX

The deflectable catheter features a deflecting distal section controlled by the deflectable catheter handle for placement of transvenous devices within the chambers of the heart and for coronary sinus cannulation. The body of the deflectable catheter is radiopaque for visibility on fluoroscopy.

Indications for Use

The deflectable catheter system is indicated to provide a pathway through which diagnostic and therapeutic transvenous devices are introduced within the chambers and coronary vasculature of the heart, and for introducing balloon catheters into the coronary sinus or leads into vessels of the left heart via the coronary sinus.

Substantially Equivalent Devices

Predicate Device	Manufacturer	510(k) Number
Medtronic 10600 Deflectable Catheter Delivery System	Medtronic	K013517
Medtronic Attain Model 6226DEF Deflectable Catheter System	Medtronic	K032312
Introducer Valve Model 6228VAL	Cook Vascular Incorporated	K010128
Universal Slitter Model 6228SLT	Medtronic	Class 1 device exempt from pre- market notification

Summary of Studies

InVitro testing was performed to support substantial equivalence to the predicate device(s). The SelectSite Model C304 Deflectable Catheter System met all specified design and performance requirements.

Biocompatibility Information

All device materials / components were assessed for biocompatibility consistent with ISO-10993, "Biological Evaluation of Medical Devices Part-1: Evaluation and Testing". All materials were found to be biocompatible and in compliance to ISO 10993-1.

Sterilization Validation

The SelectSite Model C304 Deflectable Catheter System will be sterilized using a validated 100% Ethylene Oxide (EtO) sterilization process.

Conclusion

Through the data and information presented, as well as similarities to legally marketed devices, Medtronic, Inc. considers the SelectSite Model C304-S59 and C304-L69 Deflectable Catheter Systems to be substantially equivalent to legally marketed predicate device(s).





Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

JAN 2 2 2004

Medtronic, Inc c/o Ms. Kristy K. Mollner, M.S., RAC Senior Regulatory Affairs Specialist 7000 Central Avenue NE Mail Stop CW304 Minneapolis, MN 55432-3576

Re: K033989

Trade Name: Medtronic® SelectSite[™] Model C304-S59 Deflectable Catheter System and Medtronic® SelectSite[™] Model C304-L69 Deflectable Catheter System
Regulation Number: 21 CFR 870.1250
Regulation Name: Percutaneous Catheter
Regulatory Class: Class II (two)
Product Code: DQY
Dated: December 22, 2003
Received: December 23, 2003

Dear Ms. Mollner:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. However, you are responsible to determine that the medical devices you use as components in the kit have either been determined as substantially equivalent under the premarket notification process (Section 510(k) of the act), or were legally on the market prior to May 28, 1976, the enactment date of the Medical Device Amendments. Please note: If you purchase your device components in bulk (i.e., unfinished) and further process (e.g., sterilize) you must submit a new 510(k) before including these components in your kit. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, and labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Page 2 - Ms. Kristy K. Mollner, M.S., RAC

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on the labeling regulation, please contact the Office of Compliance at (301) 594-4646. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll free number (800) 638-2041 or (301) 443-6597, or at its Internet address <u>http://www.fda.gov/cdrh.dsma/dsmamain.html</u>.

Sincerely yours,

Bram D. Zuckerman, M.D. Director Division of Cardiovascular Devices Office of Device Evaluation Center for Devices and Radiological Health

SelectSite[™] Model C304 KXXXXX

INDICATIONS FOR USE

Page I of 1

510(k) Number (if known): <u>K03348</u>9

Device Name:

Medtronic[®] SelectSiteTM Model C304-S59 Deflectable Catheter System

Medtronic[®] SelectSiteTM Model C304-L69 Deflectable Catheter System

Indications for Use: The deflectable catheter system is indicated to provide a pathway through which diagnostic and therapeutic transvenous devices are introduced within the chambers and coronary vasculature of the heart, and for introducing balloon catheters into the coronary sinus or leads into vessels of the left heart via the coronary sinus.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use <u>×</u> (Per 21 CFR 801.109) OR

Over-the Counter Use_____

(Division Sign-Qff) Division of Cardiovascular Devices (Optional Format 1-2-96)

510(k) Number__K033991

