ARCHITECT® STAT Troponin- I Assay May 5, 2004

AUG 1 2 2004

K041192

A Fisher Scientific Company 8365 Valley Pike P.O. Box 307 Middletown, VA 22645-0307

Tel: 540-869-3200 Fax: 540-869-5249



6.3 Summary of Safety and Effectiveness



This summary of the 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990.

Applicant Name:

Josefina Infantas, MSM Sr. Regulatory Affairs Specialist Fisher Diagnostics 8365 Valley Pike P.O. Box 307 Middletown, VA 22645

Phone: 540-869-8158 Fax: 540-869-8129

Establishment Registration Number: 1181121

Identification of Device:

Device Name: ARCHITECT® STAT Troponin-I immunoassay

Proprietary/Trade Name: ARCHITECT® STAT Troponin-I immunoassay

Common Name: Troponin-I test system

Device Classification: Class II

Governing Regulation: 21 CFR 862.1215

FDA Panel: Clinical Chemistry (75)

Product Code: MMI

Identification of Predicate Device:

Access® AccuTnl Assay (K021814)

Intended Use of the Device:

ARCHITECT *STAT* Troponin-I is a Chemiluminescent Microparticle Immunoassay (CMIA) for the quantitative determination of cardiac troponin-I (cTnI) in human serum and plasma on the ARCHITECT *i* System with STAT capability. Troponin-I values are used to assist in the diagnosis of myocardial infarction (MI).

Description of the Device:

The ARCHITECT STAT Troponin-I assay is a two-step assay to determine the presence of cardiac troponin-I in human serum and plasma using CMIA technology with flexible assay protocols, referred to as Chemiflex[®].

In the first step, sample, assay diluent and anti-troponin-I antibody-coated paramagnetic microparticles are combined. After incubation and washing, anti-troponin-I acridinium labeled conjugate is added in the second step. Following another incubation and wash, pre-trigger and trigger solutions are added to the reaction mixture. The resulting chemiluminescent reaction is measured as relative light units (RLUs). A direct relationship exists between the amount of troponin-I in the sample and the RLUs detected by the ARCHITECT *i* system optics. The concentration of troponin-I is read

relative to a standard curve established with calibrators of known troponin-l concentration.

Comparison of Technological Characteristics:

The ARCHITECT® STAT Troponin-I and the ACCESS® AccuTnI assays use a chemiluminescent microparticle immunoassay (CMIA) method for the quantitative determination of cardiac troponin-I in human serum and plasma. Values obtained are used to assist in the diagnosis of myocardial infarction. Proclin® 300 is a preservative used all reagent components (microparticles, diluent and conjugate) of the ACCESS AccuTnI assay, while only the diluent and conjugate of the ARCHITECT STAT Troponin-I contain Proclin® 300 as a preservative. The ARCHITECT STAT Troponin-I microparticles contain an antimicrobial agent as a preservative. Both assays have microparticles coated with mouse monoclonal anti-troponin-I in TRIS buffer. The conjugates both contain mouse monoclonal anti-troponin-I.

Summary of Non-clinical Performance:

The ARCHITECT® STAT Troponin-I assay is substantially equivalent to the ACCESS® AccuTnI assay in terms of precison, linearity, interferences, and stability as demonstrated in non-clinical performance data in this 510(k) submission.

Summary of Clinical Performance:

The ARCHITECT *STAT* Troponin-I assay demonstrated sensitivity and specificity that is substantially equivalent to the ACCESS AccuTnI assay, using the optimal AMI "cut-off" of 0.30 ng/mL as indicated by clinical data in this 510(k) submission. The sample stability study evaluated ARCHITECT *STAT* Troponin-I assay using Lithium Heparin and Serum Separator collection tubes. There was no systematic gain or loss of the detectability of troponin-I in serum or plasma samples under any of the storage conditions evaluated in this study. A method comparison was also conducted with the ARCHITECT *STAT* Troponin-I and ACCESS AccuTnI assays and as a result, the two systems demonstrated substantial equivalence.

DEPARTMENT OF HEALTH & HUMAN SERVICES



Food and Drug Administration 2098 Gaither Road Rockville MD 20850

AUG 1 2 2004

Ms. Josefina Infantas, MSM Sr. Regulatory Affairs Specialist Fisher Diagnostics 8365 Valley Pike PO Box 307 Middletown, VA 22645

Re: k041192

Trade/Device Name: ARCHITECT® STAT Troponin-I Immunoassay

ARCHITECT STAT Troponin- I Calibrators

Regulation Number: 21 CFR 862.1215

Regulation Name: Creatine phosphokinase/creatine kinase or isoenzymes test system

Regulatory Class: Class II Product Code: MMI, JIT Dated: August 5, 2004 Received: August 5, 2004

Dear Ms. Infantas:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21

Page 2

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of *In Vitro* Diagnostic Device Evaluation and Safety at (301) 594-3084. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/dsma/dsmamain.html.

Sincerely yours,

Jean M. Corper MS, DVM. Jean M. Cooper, MS, D.V.M.

Director

Division of Chemistry and Toxicology

Office of In Vitro Diagnostic Device

Evaluation and Safety

Center for Devices and

Radiological Health

Enclosure

ARCHITECT® STAT Troponin- I Assay June 16, 2004

6.1 Indications for Use

ARCHITECT STAT Troponin-I is a Chemiluminescent Microparticle Immunoassay (CMIA) for the quantitative determination of cardiac troponin-I (cTnI) in human serum and plasma on the ARCHITECT *i* System with *STAT* protocol capability. Troponin-I values are used to assist in the diagnosis of myocardial infarction (MI).

The ARCHITECT *STAT* Troponin-I Calibrators are for calibration of the ARCHITECT *i* 2000SR System when used for the quantitative determination of cardiac Troponin-I in human serum or plasma.

510 (k) Number (if known): <u> </u>
Device Name: ARCHITECT® STAT Troponin-I Immunoassay
Over-The-Counter Use
rescription Use <u>X</u> AND/OR (21 CFR 807 Subpart C)
(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Division Sign-Off

Office of in Vitro Diagnostic
Device Evaluation and Safety

510TK) K041192