

AUG 16 2004

**Phaco Chop Needle****510(k) SUMMARY  
(per 21 CFR §807.92)**

Submitter's Name: Bausch & Lomb  
 Address: 3365 Tree Court Industrial Blvd.  
 St. Louis, MO 63122

Telephone #: (636) 226-3183  
 Fax #: (636) 226-3245

Official Correspondent: Dennis Pozzo  
 Regulatory Affairs Specialist

Date Summary Prepared: May 25, 2004

Device Name/  
 Proprietary name: Phaco Chop Needle

Classification/Common Name Phaco Needles

Class: II

Panel: Ophthalmic

Product Code: HQC

The marketed device(s) to which substantial equivalence is claimed: Bausch & Lomb  
 Phaco Needles

**PRODUCT DESCRIPTION:**

The phacoemulsification (phaco) needle is the cylindrical, metal tip which is connected to the distal end of a phaco handpiece. The needle is the component of a phaco system which, as driven by the ultrasonic handpiece, contacts and fragments the cataractous crystalline lens. Irrigation fluid flows around the external surface of the needle into the eye. The emulsified lens material and irrigant are aspirated from the eye through the phaco needle lumen. The proposed Phaco Chop Needle has the same external and internal tip design as Bausch & Lomb's currently marketed MicroFlow and MicroFlow Plus Needles. However, the internal restriction is deeper within the needle bore

**Substantial Equivalent Basis**

The proposed Phaco Chop Needle is substantially equivalent to the reusable Phaco Needles currently offered by Bausch & Lomb. They are the Standard Phaco Needle, MicroFlow Phaco Needle, MicroFlow+ Phaco Needle

**Comparison Matrix**

<b>Components</b>	<b>Currently Marketed Standard Phaco Needle</b>	<b>Currently Marketed MicroFlow Needle</b>	<b>Currently Marketed MicroFlow Plus Needle</b>	<b>Proposed Phaco Chop Needle</b>
Intended Use	Phaco needles are used for anterior segment surgery with Bausch and Lomb Microsurgical Systems	Phaco needles are used for anterior segment surgery with Bausch and Lomb Microsurgical Systems.	Phaco needles are used for anterior segment surgery with Bausch and Lomb Microsurgical Systems.	Phaco needles are used for anterior segment surgery with Bausch and Lomb Microsurgical Systems
Sterile/Non-Sterile	Sterile & Non-Sterile	Sterile & Non-Sterile	Sterile & Non-Sterile	Non-Sterile only
Single Use/Reusable	Single Use & Reusable	Single Use & Reusable	Single Use & Reusable	Reusable only
Tip Angle (degrees)	15, 30 & 45	0, 15, 30 & 45	0, 15, 30 & 45	30 & 45
Needle Tip OD (in.)	0.042	0.042	0.046	0.042
Needle Tip ID (in.)	0.036	0.036	0.036	0.036

**Statement of Indications for Use**

The Phaco Needles are intended for use with a Bausch & Lomb Millennium, Protegé or Premiere Microsurgical System for the emulsification of an opacified crystalline lens.



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

AUG 16 2004

Bauch & Lomb, Inc.  
c/o Dennis Pozzo  
3365 Tree Ct. Industrial Blvd.  
St. Louis, MO 63122-6694

Re: K041998

Trade/Device Name: Phaco Chop Needles  
Regulation Number: 21 CFR 886.4670  
Regulation Name: Phacoemulsification Needles  
Regulatory Class: Class II  
Product Code: HQC  
Dated: July 22, 2004  
Received: July 26, 2004

Dear Mr. Pozzo:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4613. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,



A. Ralph Rosenthal, M.D.  
Director  
Division of Ophthalmic and Ear,  
Nose and Throat Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

510(K) Number

Device Name: Phaco Chop Needle

**Indications for Use:**

The Phaco Needles are intended for use with a Bausch & Lomb Millennium, Protegé or Premiere Microsurgical System for the emulsification of an opacified crystalline lens.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use  OR Over-the-Counter Use

(Division Sign-Off)

*Jay L. Kaulin*

(Division Sign-Off)  
Division of Ophthalmic Ear,  
Nose and Throat Devices

510(k) Number K 041998

510(k) Number K 041998

---

Bausch & Lomb Surgical