

SEP - 7 2004

# 510(K) SUMMARY OF SUBSTANTIAL EQUIVALENCE

## Submitter

---

Medtronic, Inc.  
7000 Central Avenue N.E.  
Minneapolis, MN 55432

Contact: Paula Cordero, Regulatory Affairs Specialist  
Telephone: (763) 514-8233  
Fax: (763) 514-6424  
E-Mail: paula.cordero@medtronic.com  
Date Prepared: August 10, 2004

## Name of Device

---

Trade Name: Model 6416 Bipolar Transvenous Temporary Pacing System  
Common Name: Temporary Pacing Lead System  
Classification: Class II

## Predicate Devices

---

The predicate device for the Modified Model 6416 Bipolar Transvenous Temporary Pacing System is the current market released Model 6416 Bipolar Transvenous Temporary Pacing System.

## Device Description

---

The Model 6416 Bipolar Transvenous Temporary Pacing Lead System comprises a temporary lead and a soft-tipped Medtronic guiding catheter. The pacing lead is designed to provide pacing characteristics similar to presently marketed bipolar temporary leads. The lead is introduced percutaneously through the guiding catheter. The lead system is provided in varying lengths to allow for jugular, subclavian, or femoral approach. A removable silicone rubber torque tool is provided on the proximal portion of the lead, to facilitate rotation of the lead during implant and explant. The disposable introducer currently packaged with the system

is the Medtronic Vector X coronary guiding catheter, size 6 French. This catheter guides or steers the temporary lead during placement, and shrouds the active fixation helix from blood vessel walls and other inappropriate tissue until actual implantation in the atrium or ventricle. A hemostasis valve tightens and loosens around the lead and controls the passage of the lead through the guide catheter.

## Packaging

---

The sterile packaging for the Model 6416 Bipolar Transvenous Temporary Pacing System consists of a single pouch configuration. The pouch materials are transparent Tyvek- polyester/polyethylene laminate. The pouches are heat-sealed.

## Intended Use

---

The Model 6416 Bipolar Transvenous Temporary Pacing System features an active fixation, bipolar lead and a soft-tipped lubricated guide catheter. The system is designed for temporary intracardiac pacing and/or EGM recording. The system is disposable, for temporary single patient use with a contemplated implant duration of 7 days or less. The lead and accessories are supplied sterile. The lead is introduced transvenously using the guide catheter. Once within the appropriate chamber, the helical tip electrode of the lead is actively fixed into the endocardium. After lead placement, the guide catheter is removed by sliding it over the lead's bifurcated connector.

## Technological Characteristics

---

The technology used with the Model 6416 Bipolar Transvenous Temporary Pacing System has not changed with the modified Model 6416 Bipolar Transvenous Temporary Pacing System.

## Summary of Studies

---

Medtronic, Inc. performed system compatibility testing to support that the modified Model 6416 Bipolar Transvenous Temporary Pacing System is equivalent to the predicate device. Device testing included:

- Environmental Conditioning
- Visual Verification
- Lead/ Catheter Compatibility testing
- Preclinical confirmation of safety and performance

K042190  
p 3/3

All system tests performed have demonstrated that the modified Model 6416 Bipolar Transvenous Temporary Pacing System meets the specified requirements.

## **Sterilization**

---

The Modified Model 6416 Bipolar Transvenous Temporary Pacing System is sterilized using the same 100% Ethylene Oxide (EtO) sterilization process as the predicate device. Sterilization certification of the Modified Model 6416 Bipolar Transvenous Temporary Pacing System is based on the manufacturer's determination of substantial equivalence to the predicate device.

## **Biocompatibility**

---

Biocompatibility testing was not repeated for the modified Model 6416 Bipolar Transvenous Temporary Pacing System. No change was made to material type or material formulation from previously cleared predicate devices.

## **Conclusion**

---

Through data and information presented, numerous similarities support a determination of substantial equivalence and show the device modification does not affect the intended use of the device or alter the fundamental scientific technology of the device. Market clearance of the Modified Model 6416 Bipolar Transvenous Temporary Pacing System is supported through this Special 510(k) PreMarket Notification.



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

SEP - 7 2004

Medtronic, Inc.  
c/o Ms. Paula Cordero  
Regulatory Affairs Specialist  
7000 Central Avenue NE  
Minneapolis, MN 55432

Re: K042190

Trade Name: Medtronic® Model 6146 Bipolar Transvenous Temporary Pacing System  
Regulation Number: 21 CFR 870.3680  
Regulation Name: Cardiovascular Permanent or Temporary Pacemaker Electrode  
Regulatory Class: II (two)  
Product Code: LDF  
Dated: August 10, 2004  
Received: August 12, 2004

Dear Ms. Cordero:


We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050. This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4648. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,



Bram D. Zuckerman, M.D. *BZ*  
Director  
Division of Cardiovascular Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

# INDICATIONS FOR USE

510(k) Number (if known): N/A K042190

Device Name: Medtronic® Model 6416 Bipolar Transvenous  
Temporary Pacing System

Indications For Use: The Model 6416 Bipolar Transvenous Temporary  
Pacing System features an active fixation, bipolar  
lead and a soft-tipped lubricated guide catheter.  
The system is designed for temporary intracardiac  
pacing and/or EGM recording.

The system is disposable, for temporary single  
patient use with a contemplated implant duration of  
7 days or less. The lead and accessories are  
supplied sterile.

The lead is introduced transvenously using the  
guide catheter. Once within the appropriate  
chamber, the helical tip electrode of the lead is  
actively fixed into the endocardium. After lead  
placement, the guide catheter is removed by sliding  
it over the lead's bifurcated connector.


(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use    
(Per 21 CFR 801.109)

OR

Over-The-Counter Use

  
(Division Sign-Off)  
Division of Cardiovascular Devices

510(k) Number K042190