

K043012 1/2

NOV 17 2004

**OmniPulse™ Mini, Model 2120**

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

**I. Submitter Information:** Trimedyne, Inc.  
15091 Bake Parkway  
Irvine, CA 92618  
949-951-3800

Contact Person: Glenn Yeik  
President and COO

Summary Date: 27 September 2004

**II. Device Name**

Proprietary: OmniPulse™ Mini, Model 2120

Common: Holmium:Yttrium Aluminum Garnet (Ho:YAG) Laser System

Classification: Various Medical Lasers

**III. Predicate Device**

The predicate devices for the OmniPulse™ Mini are the Trimedyne OmniPulse™-MAX and OmniPulse Jr.™ laser systems.

**IV. Device Description**

The OmniPulse™ Mini Laser System Model 2120 is a pulsed solid-state Holmium:YAG laser system is designed to deliver infrared laser energy with a wavelength of 2.1µm and a nominal pulse width of 350 microseconds. Menu driven control options allow the operator to select the pulse energy, the pulse repetition rate or frequency, and one of three exposure modes (continuous, single burst or repetitive burst).

**V. Intended Use**

The OmniPulse™ Mini is designed for use in conjunction with Trimedyne's various arthroscopic handpieces and optical fiber delivery devices. As with the other cleared Trimedyne Ho:YAG laser systems, applications for this laser system include superficial incision, excision, resection, ablation, coagulation, hemostasis, and vaporization, with or without an endoscope, in the following indications:

- General Surgery of Soft Tissues
- Genitourinary Surgery/Urology
- Otorhinolaryngology (ENT) Surgery
- Gynecological Surgery
- Lithotripsy and Percutaneous Urinary Lithotripsy

- Orthopedic Surgery
  - Percutaneous Cervical, Lumbar, and Thoracic Disc Decompression/Discectomy
  - Therapeutic Dermatological, Therapeutic Plastic, and Aesthetic Surgical Procedures
- Gastroenterological/Gastrointestinal Surgery

When used in conjunction with the Trimedyné Side Firing Needle with Vent Sheath family of optical fiber delivery devices, this laser system may be used for interstitial incision, excision, resection, ablation, coagulation, hemostasis, and vaporization in multispecialty applications.

#### **VI. Technological Characteristics**

The laser is driven by a flash lamp and emits one pulse of laser radiation each time the flash lamp is fired. The laser pulse refers to the pulse waveform - each laser pulse consists of large numbers of overlapping spikes that have a duration of about one microsecond. These spikes compose the pulse waveform with nominal pulse duration of 350 microseconds. Laser pulses are generated at the frequency chosen by the operator. Pulse frequencies of 7, 10, 15 and 20 Hz are available.

In general, the operator may request pulse energies from 0.2 to 2.0 Joules. The range of pulse energies available depends on the pulse frequency. Not all energies are available at all pulse frequencies. The number of pulses emitted from the laser system can be controlled by the operator through the choice of one of three exposure modes.

#### **VII. Nonclinical Data**

No nonclinical data was submitted in this Premarket Notification.

#### **VIII. Clinical Data**

No clinical data was submitted in this Premarket Notification.

#### **IX. Conclusions**

All of the characteristics of the OmniPulse Mini fall with the previously cleared ranges for the predicate devices. Therefore, the OmniPulse Mini Model 2120 is substantially equivalent to the Trimedyné OmniPulse™-MAX and OmniPulse Jr.™ laser systems.



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

NOV 17 2004

Trimedyne, Inc.  
c/o Ms. Elizabeth Drew  
Medical Device Services  
Underwriters Laboratories, Inc.  
1655 Scott Boulevard  
Santa Clara, California 95050-4169

Re: K043012

Trade/Device Name: Trimedyne OnmiPulse™ Mini 20 Watt Ho:YAG Laser System,  
Model 2120

Regulation Number: 21 CFR 878.4810

Regulation Name: Laser surgical instrument for use in general and  
plastic surgery and in dermatology

Regulatory Class: II

Product Code: GEX

Dated: October 20, 2004

Received: November 2, 2004

Dear Ms. Drew:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

  
for Celia M. Witten, Ph.D., M.D.  
Director  
Division of General, Restorative  
and Neurological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number: K043012

Device Name: Trimedyne OmniPulse™ Mini 20 watt Ho:YAG Laser System, Model 2120

### Indications for Use:

Superficial incision, excision, resection, ablation, coagulation, hemostasis, and vaporization, with or without an endoscope, in the following indications:

1. Dermatology and Plastic Surgery of soft, mucosal, fatty, and cartilaginous tissues, in therapeutic plastic, therapeutic dermatological and aesthetic surgical procedures, including:
  - scars
  - tattoo removal
  - vascular lesions (including port wine stains, hemangioma, telangiectasia [facial, leg], and rosacea)
  - corns
  - papillomas
  - basal cell carcinomas
  - lesions of skin and subcutaneous tissue
  - plantar warts
  - periungual and subungual warts
  - debridement of decubitus ulcer
  - skin tag vaporization

Prescription Use   X    
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

  Miriam C. Provost    
(Division Sign-Off)  
Division of General, Restorative,  
and Neurological Devices

510(k) Number   K043012

**Indications - continued**2. Gastroenterological/Gastrointestinal Surgery, including:

- cholecystectomy
- lysis of adhesions
- appendectomy
- biopsy
- pylorostenotomy
- benign and malignant lesions
- rectal polyps of sigmoid colon
- gall bladder calculi
- biliary/bile duct calculi
- benign and malignant neoplasm
- polyps
- colitis
- ulcers
- angiodysplasia
- hemorrhoids
- varices
- esophagitis
- esophageal ulcer
- Mallory-Weiss tear
- gastric ulcer
- duodenal ulcer
- non-bleeding ulcer
- gastric erosions
- colorectal cancer
- gastritis
- bleeding tumors
- pancreatitis
- vascular malformations
- telangiectasias
- telangiectasias of the Osler-Weber-Rendu disease

3. General Surgery of soft tissue, including:

- skin incision
- tissue dissection
- excision of external tumors and lesions
- complete or partial resection of internal organs
- tumors and lesions
- tissue ablation
- mastectomy
- hepatectomy
- pancreatectomy
- splenectomy
- thyroidectomy
- parathyroidectomy
- herniorrhaphy
- tonsillectomy
- lymphadenectomy
- partial nephrectomy
- pilonidal cystectomy
- resection of lipoma
- pelvic adhesiolysis
- debridement of decubitus ulcer
- hemorrhoids
- pilonidal cyst removal and repair
- debridement of stasis ulcer
- biopsy

4. Genitourinary Surgery/Urology, including:

- superficial urinary bladder tumors
- invasive bladder carcinoma
- urethral strictures
- lesions of the external genitalia
- bladder
- urethral and ureteral tumors
- condylomas
- urethral and penile hemangioma
- bladder neck obstructions
- holmium laser incision, excision, resection, ablation, hemostasis, vaporization, and enucleation in the treatment of benign prostatic hyperplasia (BPH)

**Indications - continued**

5. Gynecological Surgery during open and endoscopic procedures, including:

- condyloma acuminata

6. Lithotripsy and Percutaneous Urinary Lithotripsy, including:

- fragmentation of urinary calculi
- fragmentation of urethral calculi
- fragmentation of kidney calculi
- treatment of distal impacted fragment of steinstrasse when guide wires cannot be passed

7. Orthopedic Surgery in pathological soft and cartilaginous tissue in small and large joints, including:

- knee meniscectomy
- knee synovectomy
- chondromalacia and tears
- loose body debridement
- lateral retinacular release
- debridement of the degenerative knee
- plica removal
- ligament and tendon release
- contouring and sculpting of articular surfaces
- debridement of inflamed synovial tissue
- capsulectomy in the knee
- chondroplasty in the knee
- chondromalacia ablation

8. Otorhinolaryngology (ENT) Surgery in soft, mucosal, cartilaginous and bony tissue, including:

- endosinus surgery
- functional endoscopic sinus surgery
- turbinate procedures (e.g. turbinectomy)
- dacryocystorhinostomy (DCR)
- ethmoidectomy
- polypectomy
- maxillary antrostomy
- frontal sinusotomy
- sphenoidotomy

**Indications - continued**

9. Percutaneous Cervical, Lumbar, and Thoracic Disc Decompression / Discectomy:  
Superficial incision, excision, resection, ablation, coagulation, hemostasis and vaporization, with or without an endoscope, in the following:
- Percutaneous Lumbar Disc Decompression/Discectomy in soft, cartilaginous, and bony tissue, including:
    - foraminoplasty
  - Percutaneous Cervical Disc Decompression/Discectomy in soft tissue, in patients with:
    - uncomplicated ruptured or herniated discs
    - neck pain with radiation down the arm
    - symptoms and signs of sensory loss, tingling, numbness, muscle weakness, and/or decreased deep tendon reflexes
    - MRI, CT, myelogram, or discogram findings of disc herniation consistent with patient signs and symptoms
    - positive electromyography and/or nerve conduction studies
    - no improvement after 12 weeks of conservative therapy (i.e., physical therapy, traction, bed rest, exercises, and medication)
  - Percutaneous Thoracic Disc Decompression/Discectomy in soft tissue, in patients with:
    - uncomplicated ruptured or herniated discs
    - thoracic and intercostal intractable pain
    - paresthesias at levels appropriate to the herniated discs visualized on MRI and CT-myelography
    - MRI, CT, myelogram, or discogram findings of disc herniation consistent with patient signs and symptoms
    - no improvement after 12 weeks of conservative therapy (i.e., physical therapy, traction, bed rest, exercises, and medication)

Interstitial incision, excision, resection, ablation, coagulation, hemostasis, and vaporization in multispecialty applications; interstitial applications should only be performed using the Trimedynce Side Firing Needle with Vent Sheath family of optical fiber delivery devices.