510(k) Summary:

JAN 1 3 2005

K643524 Prgs (42

This summary is provided as part of this Premarket Notification in compliance with 21CRF, Section 807.92.

Submitters name: B-K Medical

Address: Mileparken 34, DK2730 Herlev, Denmark

Phone: +45 44528100 Fax: +45 44528199

Contact person: Villy Braender, Regulatory Manager

Date prepared: 17 December, 2004

Trade name: Ultrasound Scanner Pro Focus 2202 Common name: Diagnostic Ultrasound System

Classification names:

Ultrasonic Pulsed Echo Imaging System (90 IYO, CFR 892.1560) Ultrasonic Pulsed Doppler Imaging System (90 IYN, CFR 892.1560) Diagnostic Ultrasonic Transducer (90 ITX, CFR 892.1570)

Identification of predicate, legally marketed device:

B-K Medical Ultrasound Scanner Type 2400, K024236 (JAN17 2003)

Device description:

Pro Focus 2202 supports the following scanning modes and combinations thereof:

B-mode, M-mode, PWD mode, CFM mode.

An optional ECG signal can be superimposed the ultrasound information in all modes and mode combinations.

The system can perform simple geometric measurements, and perform calculations in the areas of Vascular, Urology, Cardiology and OB/GYN applications.

The system can guide biopsy- and puncture needles.

An optional 3-D unit can reconstruct a series of 2-D images into a single 3-D volume and display this on the screen.

Transducers

Transducers are linear and convex arrays and mechanical sector.

The patient contact materials are biocompatible.

All transducers used together with Pro Focus 2202 are Track 3 transducers.

Acoustic output

The system controlling the Acoustic Output in Pro Focus 2202 is the same as the system in 2400. The system will assure that the acoustic output always will stay below the pre-amendments upper limits i.e. Ispta \leq 720 mW/cm² and MI \leq 1.9 (Track 3, non ophthalmic).

The Thermal Index values are maximum 6.0, i.e. $TI \le 6.0$

Clinical measurement accuracy.

Clinical measurements and calculations are described and accuracies are provided in the User Guide.

Thermal, mechanical and electrical safety.

The scanner Pro Focus 2202 has been tested by a recognized, Certified Body.

Ko 43524 Page 20f2

Acoustic Output Reporting

The Acoustic Output Reporting is made according to the standards required by "Information for Manufacturers Seeking Clearance of Diagnostic Ultrasound Systems and Transducers, FDA, CDRH, September 30, 1997"

Intended use.

2202 intended uses are contained within 2400-intended uses:

	Predicate device:	Submitted device:
	Ultrasound scanner Type 2400	Ultrasound scanner Pro
	K024236 (JAN17 2003)	Focus 2202
Modes of operation	B, M, PWD, CFM 1) and combinations.	B, M, PWD,CW, CFM 1)
	Tissue harmonic imaging	and combinations. Tissue
		and contrast harmonic
		imaging
Intended use(clinical	Abdominal	Abdominal
application)	Cardiac	Cardiac
	Fetal (incl Obstetrics)	Fetal (incl Obstetrics)
	Intraoperative	Intraoperative
	Neurosurgery	Neurosurgery
•	Pediatrics	Pediatrics
	Transrectal	Transrectal
	Small Parts (organs)	Small Parts (organs)
	Transvaginal	Transvaginal
	Peripheral vascular	Peripheral vascular
	Musculo-skeletal	Musculo-skeletal
Features	ECG (not monitoring)	ECG (not monitoring) 3D

¹⁾ CFM= Color Flow Mapping=Color Doppler and Amplitude Doppler.

Technological characteristics compared to the predicate device.

The predicate device has the same major technological characteristics as the subject device described above.

Minor differences consist: Modified processor and operating system, modified mechanical outline and 3D imaging.



Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

JAN 1 3 2005

Mr. Villy Brænder Regulatory Manager B-K Medical A/S Mileparken 34 Herlev 34, DK 2730 DENMARK

Re: K043524

Trade Name: Ultrasonic Scanner Pro Focus 2202

Regulation Number: 21 CFR 892.1550

Regulation Name: Ultrasonic pulsed doppler imaging system

Regulation Number: 21 CFR 892.1560

Regulation Name: Ultrasonic pulsed echo imaging system

Regulation Number: 21 CFR 892.1570

Regulation Name: Diagnostic ultrasonic transducer

Regulatory Class: II

Product Code: 90 IYN, IYO, and ITX

Dated: December 17, 2004 Received: December 31, 2004

Dear Mr. Brænder:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

This determination of substantial equivalence applies to the following transducers intended for use with the Ultrasonic Scanner Pro Focus 2202, as described in your premarket notification:

Transducer Model Number

1850

8661

8662

 $\frac{8803}{8811}$

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This determination of substantial equivalence is granted on the condition that prior to shipping the first device, you submit a postclearance special report. This report should contain complete information, including acoustic output measurements based on production line devices, requested in Appendix G, (enclosed) of the Center's September 30, 1997 "Information for Manufacturers Seeking Marketing Clearance of Diagnostic Ultrasound Systems and Transducers." If the special report is incomplete or contains unacceptable values (e.g., acoustic output greater than approved levels), then the 510(k) clearance may not apply to the production units which as a result may be considered adulterated or misbranded.

The special report should reference the manufacturer's 510(k) number. It should be clearly and prominently marked "ADD-TO-FILE" and should be submitted in duplicate to:

Food and Drug Administration Center for Devices and Radiological Health Document Mail Center (HFZ-401) 9200 Corporate Boulevard Rockville, Maryland 20850

This letter will allow you to begin marketing your device as described in your premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus permits your device to proceed to market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/dsma/dsmamain.html

If you have any questions regarding the content of this letter, please contact Rodrigo C. Perez at (301) 594-1212.

Sincerely yours,

Mancy C Brogdon
Nancy C. Brogdon

Director, Division of Reproductive, Abdominal and Radiological Devices

Office of Device Evaluation

Center for Devices and Radiological Health

Enclosure(s)

Diagnostic Ultrasound Indications for Use Form

System: 2202

Fill out one form for each ultrasound system and each transducer.

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

	<u> </u>	Mode of Operation										
Clinical Application	А	В	М	PWD	Tissue- and contrast narmonic imaging	Color Doppl er	Amplitude Doppler	Color Velocity Imaging	Combined (specify 1)	Continous Wave)		
Ophthalmic												
Fetal		Р	Р	Р	Р	Р	Р		Р			
Abdominal		Р	Р	Р	Р	P	Р		Р			
Intraoperative (specify)		Р	Р	Р	Р	P	P		Р			
Intraoperative Neurological		Р	Р	Р	Р	Р	Р		Р			
Pediatric		Р	Р	Р	Р	Р	Р		Р			
Small Organ (specify)		Р	Р	Р	Р	Р	P		Р	-		
Neonatal Cephalic												
Adult Cephalic												
Cardiac		Р	Р	Р	Р	Р	Р		Р			
Transesophageal												
Transrectal		Р	Р	Р	Þ	Р	Р		Р	······		
Transvaginal		Р	Р	P	Р	Р	Р		Р			
Transurethral		Р	Р	Р	Р	Р	Р		Р			
Intravascular												
Peripheral Vascular		Р	P	P	Р	Р	Р		Р			
Laparoscopic												
Musculo-skeletal Conventional		Р	Р	Р	Р	Р	Р		Р			
Musculo-skeletal Superficial		Р	Р	Р	Р	Р	Р		Р			
Other (specify)												

Other (specify)				1	1						
N= new indication; P=	previo	usly c	leare	d by F	DA; E=	added	under Appe	endix E	<u>' </u>		_
Additional Comments:	1) B+ D is F	M, B- PWD,	+D, B- C is	+C, B- Color	+D+C. B Doppler	mode in . Fetal is	cludes Tissus s often calle	e-and Contr ed Obstetri	rast Harmon CS	ic Imaging.	
								, <u>, , , , , , , , , , , , , , , , , , ,</u>		····	
	(PL	EASE D	O NOT V	VRITE BE	LOW THIS LI	NE - CONTII	NUE ON ANOTHE	R PAGE IF NEED	DED)		
		Co	ncurre	nce of	CDRH, O	ffice of D	evice Evaluat	ion (ODE)			

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)
Division of Reproductive, Abdominal, and Radiological Devices
510(k) Number + 53524

_				
С,	JST	Δŀ	n.	

2202

Transducer:

1850 (with interchangeable probes 8539,6004,6005)

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application			Mode of Operation								
General	Specific	В	М	PWD	CWD	Color	Combined	Other*			
(Track I Only)	(Tracks I & III)	-	''		CIID	Doppler	(Specify)	(Specify)			
Ophthalmic	Ophthalmic	1	!			- оррии	(0)	(Specify			
•	Fetal	<u> </u>									
	Abdominal										
	Intra-operative (Specify)	Р									
	Intra-operative (Neuro)										
	Laparoscopic										
Fetal Imaging	Pediatric					-		<u> </u>			
& Other	Small Organ (Specify)						,				
	Neonatal Cephalic										
	Adult Cephalic										
	Trans-rectal	P									
	Trans-vaginal										
	Trans-urethral										
	Trans-esoph. (non-Card.)										
	Musculo-skel. (Conventional)							-			
	Musculo-skel. (Superficial)										
	Intra-luminal										
	Other (Specify)										
	Cardiac Adult										
Cardiac	Cardiac Pediatric										
	Trans-esoph. (Cardiac)										
	Other (Specify)										
Peripheral	Peripheral vessel										
Vessel	Other (Specify)										

N = new indication; P = previously cleared by FDA(K002085); E = added under Appendix E

Additional Comments: Intraoperative: Rectum, Urethra, Urinary bladder,	
(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)	

Concurrence of Center for Devices and Radiological Health, Office of Device Evaluation

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off) Division of Reproductive, Abdominal, and Radiological Devices

the Radiosogical Devices

510(k) Number

^{*}Examples may include: A-mode, Amplitude Doppler, 3-D Imaging, Harmonic Imaging, Tissue Motion Doppler, Color Velocity Imaging

System:	2202	
Transducer:	8661	

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

CI	inical Application				Mod	de of Opera	tion	
General	Specific	В	M	PWD	CWD	Color	Combined	Amplitud
(Track I Only)	(Tracks I & III)					Doppler	(Specify)	e Doppler
Ophthalmic	Ophthalmic							
	Fetal	P	P	Р	<u> </u>	Р	P 1)	P
	Abdominal		<u> </u>					
	Intra-operative (Specify)							
	Intra-operative (Neuro)					<u> </u>	1	
	Laparoscopic	1				<u> </u>		
Fetal Imaging	Pediatric		ļ <u> </u>			<u> -</u>		
श्र Other	Small Organ (Specify)	<u>.</u>	ļ <u>.</u>					
	Neonatal Cephalic							
	Adult Cephalic		ļ					
	Trans-rectal		<u> </u>					
	Trans-vaginal	P	P	P		P	P 1)	P
	Trans-urethral		<u> </u>					
	Trans-esoph. (non-Card.)							
	Musculo-skel. (Conventional)							
	Musculo-skel. (Superficial)		<u> </u>					
	Intra-luminal							
	Other (Specify)							
	Cardiac Adult	ļ						
Cardiac	Cardiac Pediatric							
	Trans-esoph. (Cardiac)							
	Other (Specify)	ļ						
Peripheral	Peripheral vessel							
Vessel	Other (Specify)	<u> </u>						<u> </u>

N = new indication; P = previously cleared by FDA; E = added under Appendix E

Additional Com	ments:	_1)) Mode	combinati	ions:	B+M,	B+D,	B+C,	B+D+C.	(D is	PWD,	С	is
Color Flow m	apping	Doppler	including	Amplit	ude (p	oower)	Doppl	ler)				
	(PLEASE I	DO NOT WRI	TE BELOW THIS	S LINE-CO	NTINUI	E ON AI	NOTHER	R PAGE IF N	EEDED)			

Concurrence of Center for Devices and Radiological Health, Office of Device Evaluation

Prescription Use (Per 21 CFR 801.109)

Division of Reproductive, Abdominal, and Radiological Devices

510(k) Number__ X84352

^{*}Examples may include: A-mode, Amplitude Doppler, 3-D Imaging, Harmonic Imaging, Tissue Motion Doppler, Color Velocity Imaging

System:	2202
Transducer:	8662

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Cí	Clinical Application			Mode of Operation								
General	Specific	В	M	PWD	CWD	Color	Combined	Amplitud				
(Track I Only)	(Tracks I & III)					Doppler	(Specify 1)	e Doppler				
Ophthalmic	Ophthalmic											
	Fetal											
	Abdominal											
	Intra-operative (Specify 2)	P	P	P		P	P	P				
	Intra-operative (Neuro)	P	P	P		P	P	P				
	Laparoscopic											
Fetal Imaging	Pediatric	P	Р	P		-P	P	P				
& Other	Small Organ (Specify)				_							
	Neonatal Cephalic	T										
	Adult Cephalic		}									
	Trans-rectal											
	Trans-vaginal											
	Trans-urethral											
	Trans-esoph. (non-Card.)											
	Musculo-skel. (Conventional)											
	Musculo-skel. (Superficial)											
	Intra-luminal	I										
	Other (Specify)											
	Cardiac Adult											
Cardiac	Cardiac Pediatric											
	Trans-esoph. (Cardiac)											
	Other (Specify)											
Peripheral	Peripheral vessel											
Vessel	Other (Specify)											

N = new indication; P = previously cleared by FDA; E = added under Appendix E

Additional Comments:1) Mode combinations: B+M, B+D, B+C, B+D+C. (D is PWD, C is Color Flow
mapping Doppler including Amplitude(power)Doppler)
2)Intraoperative: Gall bladder
(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of Center for Devices and Radiological Health, Office of Device Evaluation

Prescription Use (Per 21 CFR 801.109)

mission of Reproductive, Abdominal, . a Radiological Devices

^{*}Examples may include: A-mode, Amplitude Doppler, 3-D Imaging, Harmonic Imaging, Tissue Motion Doppler, Color Velocity Imaging

System:	2202	
Transducer:	8803	

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

intended use: Diagnostic ditrasound imaging of									
Clinical Application		Mode of Operation							
General	Specific	В	M	PWD	Harm	Color	Combined	Amplitud	
(Track I Only)	(Tracks I & III)]	onic	Doppler	(Specify)	e Doppler	
		1			imagin				
					g				
Ophthalmic	Ophthalmic								
	Fetal	Р	Р	P	P	P	P 1)	P	
	Abdominal	P	P	P	P	P	P 1)	P	
	Intra-operative (Specify)								
	Intra-operative (Neuro)					-			
	Laparoscopic						•		
Fetal Imaging	Pediatric	P	P	P	P	P	P 1)	P	
& Other	Small Organ (Specify)								
	Neonatal Cephalic								
	Adult Cephalic								
	Trans-rectal	l					•		
	Trans-vaginal								
	Trans-urethral								
	Trans-esoph. (non-Card.)								
	Musculo-skel. (Conventional)						•	1	
	Musculo-skel. (Superficial)								
	Intra-luminal								
	Other (Specify)								
	Cardiac Adult	P	P	P	Р	P	P 1)	P	
Cardiac	Cardiac Pediatric	P	P	P	P	P	P 1)	P	
	Trans-esoph. (Cardiac)								
	Other (Specify)								
Peripheral	Peripheral vessel								
Vessel	Other (Specify)								

N = new indication; P = previously cleared by FDA; E = added under Appendix E

Additional Comments:	1)) Mode	combinations:	В+М,	B+D,	B+C,	B+D+C.	(D is	PWD,	C is
Color Flow mapping	Doppler	including Amp	litude	(power	r) Dopj	oler)			
(PLEASE	DO NOT WR	ITE BELOW THIS LINE	-CONTIN	UE ON A	ANOTH	ER PAGE IF	NEEDED)		

Concurrence of Center for Devices and Radiological Health, Office of Device Evaluation

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Offy)
Division of Reproductive, Abdominal, and Radiological Devices

510(k) Number

^{*}Examples may include: A-mode, Amplitude Doppler, 3-D Imaging, Harmonic Imaging, Tissue Motion Doppler, Color Velocity Imaging

System:	2202_	
Transducer:	8811	

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation							
General	Specific	В	M	PWD	Harm	Color	Combined	Amplitud	
(Track I Only)	(Tracks I & III)				onic	Doppler	(Specify 1)	e Doppler	
					imagin				
					g				
Ophthalmic	Ophthalmic	<u> </u>	<u> </u>				ļ		
	Fetal	1							
	Abdominal								
	Intra-operative (Specify 2)	P	P	Р	P	P	P	P	
	Intra-operative (Neuro)								
	Laparoscopic								
Fetal Imaging	Pediatric	P	Р	P	P	P	P	P	
& Other	Small Organ (Specify 3)	P	P	Р	P	P	P	P	
	Neonatal Cephalic								
	Adult Cephalic								
	Trans-rectal								
	Trans-vaginal								
ŀ	Trans-urethral								
	Trans-esoph. (non-Card.)								
	Musculo-skel. (Conventional)	Р	Р	P	P	P	P	P	
	Musculo-skel. (Superficial)	P	P	P	P	P	P	P	
	Intra-Iuminal								
	Other (Specify)		ļ						
	Cardiac Adult								
Cardiac	Cardiac Pediatric								
	Trans-esoph. (Cardiac)								
	Other (Specify)			}					
Peripheral	Peripheral vessel	P	Р	P	P	P	P	P	
Vessel	Other (Specify)]						

N = new indication; P = previously cleared by FDA; E = added under Appendix E

Additional Comments: 2)Intraoperative: Breast, liver, pancreas, biliary system

3) Small Organ: Breast, testis, penis, thyroid, parathyroid, salivary glands, lymph nodes

1) mode combinations: B, B+M, B+D, B+C, B+D+C. (D is PWD, C is Color Flow mapping Doppler including Amplitude(power)Doppler)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of Center for Devices and Radiological Health, Office of Device Evaluation

Prescription Use (Per 21 CFR 801.109)

Division of Reproductive, Abdominal, and Radiological Devices

^{*}Examples may include: A-mode, Amplitude Doppler, 3-D Imaging, Harmonic Imaging, Tissue Motion Doppler, Color Velocity Imaging