## **SECTION 6**

# 510(K) Summary CRYOcheck™ Clot S™

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

The assigned 510(k) number is:  $\frac{604357}{}$ 

**Submitters Name & Address:** 

Precision BioLogic Inc.

900 Windmill Road, Suite 100 Dartmouth, Nova Scotia B3B 1P7

Canada

**Contact Name:** 

Stephen L. Duff - Director of New Business Development

Phone: 902-468-6422 ext. 224

Fax: 902-468-6421

Email: sduff@precisionbiologic.com

**Preparation Date:** 

December 22, 2004

**Device Name & Classification:** 

CRYO*check*™ Clot S™

Common Name: Clot-based Protein S Assay

Classification Name: Test, Qualitative and Quantitative Factor Deficiency

Regulatory Class II

**Predicate Device:** 

STA® - Staclot® Protein S (K913424)

Diagnostica Stago

9, rue des Frères Chausson 92600 ASNIERES (France)

**Device Description:** 

CRYO*check™* Clot S™ consists of:

 Protein S Deficient Plasma – contains citrated pooled normal human plasma that has been depleted of protein S by

immunoadsorption, buffers and stabilizers.

 Clot S Activator – contains activated protein C, Russell's viper venom, heparin neutralizing agents, buffers and stabilizers.

Precision BioLogic Clot C & S Diluent (available separately from

Precision BioLogic).

**Device Intended Use:** 

CRYO*check*™ Clot S™ is a clot-based assay intended for the quantitative

determination of protein S activity in citrated human plasma.

## **Comparison to Predicate Device:**

Parameter	CRYO <i>check</i> ™ Clot S™	STA® - Staclot® Protein S (K913424)		
Intended Use	CRYO <i>check™</i> Clot S™ is a clot-based assay intended for use in the quantitative determination of protein S activity in citrated human plasma.	The STA® – Staclot® Protein S kit is intended for use with analyzers of the STA® brand name, for the quantitative measurement of the functional protein S level based on the principle of factor Va inhibition.		
Format	Frozen	Lyophilized		
Volume	<ul> <li>5 x 3.0 mL Protein S Deficient Plasma</li> <li>5 x 3.0 mL Clot S Activator         OR</li> <li>5 x 1.5 mL Protein S Deficient Plasma</li> <li>5 x 1.5 mL Clot S Activator</li> </ul>	<ul> <li>2 x 1 mL vials of Reagent 1 (Protein S Deficient Plasma</li> <li>2 x 1 mL vials of Reagent 2 (Human Activated Protein C)</li> <li>2 x 1 mL vials of Reagent 3 (Preparation Containing Bovine Factor Va)</li> </ul>		

#### **Correlation with Predicate Device:**

CRYO  $Check^m$   $Clot S^m$  was compared to  $STA^{\otimes}$  -  $Staclot^{\otimes}$  Protein S using 115 clinical samples from the target population for the assay. A correlation of R=0.880 was obtained.

### **Comments on Substantial Equivalence:**

It is the opinion of Precision BioLogic Inc. that CRYO*check™* Clot S™ is substantially equivalent to STA® - Staclot® Protein S, manufactured by Diagnostica Stago (France), and currently marketed in the United States by Diagnostica Stago Inc. This opinion is based on the following:

- Both products are clot-based assays.
- Both products are intended for use in the quantitative measurement of functional protein S in citrated human plasma.
- Both products provide all coagulation factors in excess by the use of protein S deficient plasma.
- Both products use exogenous activated protein C.

#### Conclusion:

CRYOcheck™ Clot S™ is substantially equivalent to STA® - Staclot® Protein S.

# DEPARTMENT OF HEALTH & HUMAN SERVICES

SERVICES . Lieu

Mr. Stephen L. Duff

MAR 1 8 2005

Food and Drug Administration 2098 Gaither Road Rockville MD 20850

Director of New Business Development Precision BioLogic Inc. 900 Windmill Road, Suite 100 Dartmouth, Nova Scotia Canada B3B 1P7

Re: k043571

Trade/Device Name: CryoCheck<sup>TM</sup> Clot S<sup>TM</sup> Regulation Number: 21 CFR § 864.7290 Regulation Name: Factor deficiency test

Regulatory Class: II Product Code: GGP Dated: February 14, 2005 Received: February 15, 2005

#### Dear Mr. Duff:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820). This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (240) 276-0484. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <a href="http://www.fda.gov/cdrh/industry/support/index.html">http://www.fda.gov/cdrh/industry/support/index.html</a>

Sincerely yours,

Robert L. Becker, Jr., MD, Pa.D

Director

Division of Immunology and Hematology

Office of In Vitro Diagnostic Device

Evaluation and Safety Center for Devices and

Radiological Health

Enclosure

# **SECTION 5**

# **Indications for Use**

	5	I0(k) Number:	K043571		
	D	evice Name:	CRYOcheck™ Ck	ot S™	
Indication	ns for Use:				
CRYO <i>ched</i> protein S	ck™ Clot S™ i activity in citr	s a clot-based as ated human plas	ssay intended for t ma.	ne quantitative determinatio	n of
indicative	ck™ Clot S™ i of an increas thrombotic e	ed risk of throm	se protein S defici poembolism. A de	ency (congenital or acquired ficiency in protein S may pro	d) which is oduce
• ty • ty a • tv	/pe I deficiend /pe II deficien ntigen levels /pe III deficie	cies correspond icies are charact of both total and ncies are defined	erized by a reduce free protein S	levels of both total and free d protein S activity but with gen level and activity of free	normal
• o • li • d • o • o • a • p	ral anticoagu ver disease lisseminated ral contracep estrogen the	lant therapy intravascular coa tives		al clinical states:	
Prescript	ion Use	✓	AND/OR	Over-The-Counter Use	
•	CFR 801 8 bpa	art D)	- ANDION	(21 CFR 807 Subpart C)	***
(PUZASE	EDO NOT W	RITE BILOW TI	HIS LINE-CONTIN	UE ON ANOTHER PAGE I	NEEDED)
n-Off	Concu	rrence of CDRH	. Office of In Vitro	Diagnostic Devices (OIVD)	

510(k) K 043571

Office of In Vitro Diagnostic Device Evaluation and Safety