

K050224

## II. 510(K) Summary of Safety and Effectiveness

(Per 21 CFR 807.92)

JUL 13 2005

### 2.1. General Information Establishment

- **Manufacturer:** NEUERO ENGINEERING INC.
- **Address:** 20F, No. 171, Cheng-Kuang Rd.,  
San-Chung City, 241, Taipei County, Taiwan
- **Registration Number:** 3004512547
- **Contact Person:** Dr. Ke-Min Jen  
Official Correspondent  
886-3-5208829 (Tel)  
886-3-5209783 (Fax)
- **Date Submitted:** January 26, 2005

### Device

- **Proprietary Name:** *3test<sup>®</sup> Glucose Monitoring System, 3TM678G*
- **Common Name:** Blood Glucose Monitoring System
- **Classification Name:** SYSTEM, TEST, BLOOD GLUCOSE, Class II,

### 2.2. Safety and Effectiveness Information

- **Predicate Device:**  
Claim of Substantial Equivalence (SE) is made to BAYER Corp. -- Glucometer Elite Blood Glucose Meter ( K964630 ).
- **Device Description:** Based on an electrochemical biosensor technology and the principle of capillary action, *3test<sup>®</sup> Glucose Monitoring System* only needs a small amount of blood. Capillary action at the end of the test strip draws the blood into

the action chamber and your blood glucose result is precisely and displayed in 5 seconds, and 100 blood glucose result memory.

- **Intended Use:**


The *3test*<sup>®</sup> glucose test strip is intended to measure the glucose in whole blood with the *3test*<sup>®</sup> Glucose Monitoring System. It is suitable for a person with diabetes to monitor their blood glucose at home by themselves. The *3test*<sup>®</sup> Glucose Monitoring System can also be used at clinical sites by nurses or professional people to test patient's glucose level in whole blood.

- **Synopsis of Test Methods and Results**

Pre-clinical and clinical data are employed upon submission of this 510(K) premarket notification according to the Guidance Document for In Vitro Diagnostic Test System; Guidance for Industry and FDA document provided by CDRH/ FDA.

- **Substantial Equivalence (SE)**

A claim of substantial equivalence is made to the predicate device BAYER Corp. -- Glucometer Elite Blood Glucose Meter (K964630). Both of them have the same working principle and technologies. The differences are control solution, dimensions of the unit and strip, the time of auto shut off, memory results, and test time. There are no safety and effectiveness aspects arising from the subject device. They are substantially equivalent.



Dr. Jen, Ke-Min  
Official Correspondent for  
NEUERO ENGINEERING INC.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

JUL 13 2005

Neuro Engineering Inc.  
c/o Dr. Ke-Min Jen  
Roc Chinese-European Industrial Research Society  
No.58, Fu-Chiun St.  
Hsin-Chu City, China  
Taiwan 300

Re: k050224  
Trade/Device Name: *3test* Glucose Monitoring System  
Regulation Number: 21 CFR 862.1345  
Regulation Name: Glucose test system  
Regulatory Class: Class II  
Product Code: NBW, CGA  
Dated: July 2, 2005  
Received: July 6, 2005

Dear Dr. Jen:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

Page 2 –

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (240) 276-0484. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Carol C. Benson, M.A.  
Acting Director  
Division of Chemistry and Toxicology  
Office of *In Vitro* Diagnostic Device  
Evaluation and Safety  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510 (K) Number ( If Known ): K 050224

Device Name: 3test<sup>®</sup> Glucose Monitoring System, 3TM678G

### Indications For Use:

The 3test<sup>®</sup> Glucose test strip is intended to measure the glucose in whole blood with the 3test<sup>®</sup> Glucose Monitoring System. It is suitable for a person with diabetes to monitor their blood glucose at home by themselves. The 3test<sup>®</sup> Glucose Monitoring System can also be used at clinical sites by nurses or professional people to test patient's glucose level in whole blood.

### NOTE:

- the 3test is to be used with capillary whole blood from the fingertip
- the 3test is not for use with neonates
- the 3test meter is to be used with the 3test Blood Glucose Test Strip, and the 3test High and Low Glucose Control Solutions

Prescription Use \_\_\_\_\_

AND/OR


Over-The-Counter Use √

(Part 21 CFR 801 Subpart D)

(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)

  
Division Sign-Off

In Vitro Diagnostic Device  
Evaluation and Safety

K050224

page 1 of 1