

K050594

Section 3
Quantia A1-AT
510(k) Summary (Summary of Safety and Effectiveness)

Submitted by:

Biokit S.A.
Can Male, Llissa d'Amunt
Barcelona
08186 Spain

MAY 13 2005

Contact Person:

Contact: Joan Guixer, Quality Assurance and Regulatory Affairs Director
Phone: 34 - 93 860 90 00

Summary Prepared:

February 24, 2005

Name of the device:

Quantia A1-AT

Classification name(s):

866.5130 Alpha-1-antitrypsin immunological test system Class II
DEM Alpha-1-Antitrypsin, Antigen, Antiserum, Control

Identification of predicate device(s):

K860894 N Antisera to Human alpha-1-Antitrypsin (Dade Behring)

Description of the device/intended use(s):

The Quantia A1-AT is intended for the *in vitro* quantitative determination of alpha-1-antitrypsin concentration in human serum or plasma (heparin with or without gel separator, EDTA) on the AEROSET[®] system as an aid in the diagnosis of juvenile and adult cirrhosis of the liver and pulmonary emphysema.

Quantia PROTEINS Control is intended for use in monitoring the quality control of results obtained with the Quantia A1-AT reagents by turbidimetry. (NOTE: This control has been also FDA 510(k) submitted for use with Quantia Beta-2 Microglobulin). For *in vitro* diagnostic use.

Quantia PROTEINS standard is intended for use in establishing the calibration curve for the Quantia A1-AT reagents by turbidimetry. For *in vitro* diagnostic use.

Statement of Technological Characteristics of the Device Compared to Predicate Device:
 Quantia A1-AT is substantially equivalent to the commercially available predicate device,
 N Antisera to Human alpha-1-Antitrypsin, in performance and intended use.

Summary of Performance Data:

In a method comparison study evaluating 111 samples with alpha-1-antitrypsin levels ranging from 42.0 to 442.5 mg/dL on the Abbott AEROSET® instrument, the slope was 1.002 and the correlation coefficient (r) was 0.9890 for Quantia A1-AT versus the predicate device.

The studies performed on the Abbott AEROSET® system with Quantia A1-AT reagents support the following precision claims:

Material	Samples /Runs	Mean (mg/dL)	CV(%) Within Run	CV(%) Between Run	CV(%) Total
Low Control (I)	2/40	77.7	1.2	0.2	1.4
Control (I + II)	2/40	156.9	0.7	0.1	0.8
High Control (II)	2/40	229.7	1.3	0.6	1.8



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
2098 Gaither Road
Rockville MD 20850

MAY 13 2005

Biokit S.A.
c/o Ms. Joan Guixer
Llissa d'Amunt
Barcelona, Spain 08186

Re: k050596

Trade/Device Name: Quantia A1-AT
Quantia Proteins Control
Quantia Proteins Standard

Regulation Number: 21 CFR 866.5130

Regulation Name: Alpha 1-antitrypsin Immunological Test System

Regulatory Class: Class II

Product Code: DEM, JJS, JJX

Dated: February 24, 2005

Received: March 8, 2005

Dear Ms. Guixer:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Ms. Joan Guixer

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0131. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,



Robert L. Becker, Jr., M.D., PhD
Director
Division of Immunology and Hematology Devices
Office of In Vitro Diagnostic Device Evaluation and Safety
Center for Devices and Radiological Health

Enclosure

Indications for Use Statement

510(k) Number (if known): k050596

Device Name: Quantia A1-AT

Indications for Use:

The Quantia A1-AT is intended for the *in vitro* quantitative determination of alpha-1-antitrypsin concentration in human serum or plasma (heparin with or without gel separator, EDTA) on the AEROSSET® system as an aid in the diagnosis of juvenile and adult cirrhosis of the liver and pulmonary emphysema.

Quantia PROTEINS Control is intended for use in monitoring the quality control of results obtained with the Quantia A1-AT reagents by turbidimetry. (NOTE: This control has been also FDA 510(k) submitted for use with Quantia Beta-2 Microglobulin). For *in vitro* diagnostic use.

Quantia PROTEINS standard is intended for use in establishing the calibration curve for the Quantia A1-AT reagents by turbidimetry. For *in vitro* diagnostic use.

Prescription Use X
(Part 21 CFR 801 Subpart D)

OR Over-The-Counter Use _____
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)

Maria M Chan
Division Sign-Off

Office of In Vitro Diagnostic
Device Evaluation and Safety

510(k) k050596