

DEC 27 2005

**510(k) Summary of Safety and Effectiveness  
Stryker® Compartmental Knee System**

**Submission Information**

Name and Address of the Sponsor  
of the 510(k) Submission: Howmedica Osteonics Corp  
325 Corporate Drive  
Mahwah, NJ 07430

For Information contact: Vivian Kelly, Regulatory Affairs Specialist  
Howmedica Osteonics Corp.  
325 Corporate Drive  
Mahwah, NJ 07430  
Phone: (201) 831-5581 Fax: (201) 831-6038

Date Summary Prepared: December 9, 2005

**Device Identification**

Proprietary Name: Stryker® Compartmental Knee System  
Common Name: Knee Prosthesis Components  
Proposed Regulatory Class: Class II

Classification Name, Reference and Product Code:

Knee Joint, Femorotibial, Polymer/Metal, Semi-constrained, Cemented Prosthesis, 21 CFR §888.3530, 87 NPJ  
Knee joint patellofemoral polymer/metal semi-constrained cemented prosthesis, 21 CFR 888.3540, 87 KRR  
Prosthesis, Knee, Femorotibial, Non-Constrained, Cemented, Metal Polymer, 21 CFR 888.3520, 87 HSX  
Prosthesis, Knee, Femorotibial, Semi-Constrained, Cemented, Metal/Polymer, 21 CFR §888.3530, 87 HRY

Description:

The Stryker® Compartmental Knee System consists of sterile, single-use components intended for replacement of the femoral side of the patellofemoral joint and/or the condyle region(s) of the femoral joint as needed. The system includes patellofemoral, femoral, and tibial components from currently marketed Howmedica Osteonics' knee systems for patellofemoral and unicompartmental arthroplasty. The system allows the physician to choose the most appropriate option to treat the patient with patellofemoral arthroplasty and/or unicompartmental arthroplasty as needed.

Indications for Use

The Stryker® Compartmental Knee System is intended to be used in cemented patellofemoral and/or unicompartmental arthroplasty in patients where conditions exist that cannot be addressed by a single device to treat the femorotibial or patellofemoral regions of the knee. The indications for the different components of the Stryker® Compartmental Knee System include conditions when the patellofemoral and/or condylar region(s) have been affected by one or more of the following conditions:

- Degenerative arthritis in the distal femur and patella,
- Patients with a history of patellar dislocation or patella fracture,
- Patients with failed previous surgery (arthroscopy, tibial tubercle elevation, lateral release) where pain, deformity or dysfunction persists,
- Moderately disabling joint disease of the knee resulting from painful osteo- or post traumatic arthritis,
- Revision of previous unsuccessful surgical procedures, either involving, or not involving, previous use of an unicompartmental knee prosthesis,
- As an alternative to tibial osteotomy in patients with unicompartmental osteoarthritis, or

- Where bone stock is of poor quality or inadequate for other reconstructive techniques as indicated by deficiencies of the femoral condyle/tibial plateau.

These components are single use only and are intended for implantation with bone cement.

**Substantial Equivalence:**

The device is substantially equivalent to its predicates for patellofemoral arthroplasty and femorotibial arthroplasty in regards to intended use, design, materials, and operational principles. The analyses demonstrate that the components from these systems are compatible when used for patellofemoral and/or femorotibial replacement.



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

DEC 27 2005

Ms. Vivian Kelly  
Regulatory Affairs Specialist  
Howmedica Osteonics Corp.  
325 Corporate Drive  
Mahwah, New Jersey 07430

Re: K052917  
Trade/Device Name: Stryker Compartmental Knee System  
Regulation Number: 21 CFR 888.3530  
Regulation Name: Knee joint femorotibial metal/polymer semi-constrained cemented prosthesis  
Regulatory Class: II  
Product Code: NPJ, KRR, HSX, HRY  
Dated: December 13, 2005  
Received: December 13, 2005

Dear Ms. Kelly:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21

Page 2 – Ms. Vivian Kelly

CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Mark N. Melkerson

Acting Director

Division of General, Restorative and Neurological Devices

Office of Device Evaluation

Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known):

Device Name: Stryker® Compartmental Knee System

Indications for Use:

The Stryker® Compartmental Knee System is intended to be used in cemented patellofemoral and/or unicompartmental arthroplasty in patients where conditions exist that cannot be addressed by a single device to treat the femorotibial or patellofemoral regions of the knee. The indications for the different components of the Stryker® Compartmental Knee System include conditions when the patellofemoral and/or condylar region(s) have been affected by one or more of the following conditions:

- Degenerative arthritis in the distal femur and patella,
- Patients with a history of patellar dislocation or patella fracture,
- Patients with failed previous surgery (arthroscopy, tibial tubercle elevation, lateral release) where pain, deformity or dysfunction persists,
- Moderately disabling joint disease of the knee resulting from painful osteo- or post traumatic arthritis,
- Revision of previous unsuccessful surgical procedures, either involving, or not involving, previous use of a unicompartamental knee prosthesis,
- As an alternative to tibial osteotomy in patients with unicompartamental osteoarthritis, or
- Where bone stock is of poor quality or inadequate for other reconstructive techniques as indicated by deficiencies of the femoral condyle/tibial plateau.

These components are single use only and are intended for implantation with bone cement.

Prescription Use   X  


(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use     

(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF  
NEEDED)

Concurrence of  Office of Device Evaluation (ODE)

**(Division Sign-Off)**

**Division of General, Restorative,  
and Neurological Devices**

**510(k) Number**   K052917