

Page 1 of 2

MAY 17 2005

**B. 510(k) SUMMARY (as required by 21 CFR 807.92)**

**Columbus Total Knee System  
(Ultra Congruent Tibial Insert)  
December 21, 2005**

**COMPANY:** Aesculap<sup>®</sup>, Inc.  
3773 Corporate Parkway  
Center Valley, PA 18034  
Establishment Registration Number: 2916714

**CONTACT:** Matthew M. Hull  
800-258-1946 (phone)  
610-791-6882 (fax)  
[matt.hull@aesculap.com](mailto:matt.hull@aesculap.com) (email)

**TRADE NAME:** Columbus Total Knee System Ultra Congruent (UC) Tibial Insert

**COMMON NAME:** Total Knee System Gliding Surface

**CLASSIFICATION NAME:** Knee joint patellofemorotibial polymer/metal/polymer semi-constrained cemented prosthesis

**REGULATION NUMBER:** 888.3560

**PRODUCT CODE:** JWH

**SUBSTANTIAL EQUIVALENCE**

Aesculap<sup>®</sup>, Inc. believes that the Ultra Congruent Tibial Insert for the Columbus Total Knee System is a line extension of the standard and deep dish tibial inserts that were cleared as components of Aesculap's Columbus Total Knee System (K022672 and K023788). It is also substantially equivalent to the ultra congruent inserts that were cleared for the VKS Knee System (K022204).

**DEVICE DESCRIPTION**

The Ultra Congruent Tibial Insert is a gliding surface for the Columbus Total Knee System made from UHMWPE. The ultra congruent insert is designed with increased contact area that provides improved stability in flexion and extension. The Columbus UC gliding surfaces range in height from 10 mm to 20 mm and in size from T0/T0+ to T5.

**INDICATIONS FOR USE**

The Columbus Total Knee System is indicated for use in reconstruction of the diseased knee joint caused by osteoarthritis, rheumatoid arthritis, post-traumatic arthritis, the need to revise failed arthroplasties or osteotomies where pain, deformity or dysfunction persist, and for patients suffering from correctable valgus or varus deformity and moderate flexion contracture.

The Columbus Knee is designed for use with bone cement.

000003

Page 2 of 2

**TECHNOLOGICAL CHARACTERISTICS(compared to Predicate(s))**

The UC tibial inserts are offered in similar in shapes and sizes as the standard and deep dish predicate devices. The same biocompatible material is used for the all of the Aesculap inserts (UC, Std, and DD) plus the UC insert for the VKS System. The Aesculap UC insert has an increased surface area over the standard and the deep dish and an area similar to the VKS UC inserts.

**PERFORMANCE DATA**

All required testing per "Draft Guidance for the Preparation of Premarket Notifications (510(k)s) Applications for Orthopedic Devices-The Basic Elements" will be done where applicable. In addition testing will be completed as applicable per the following:

- "Guidance Document for Testing Orthopedic Implants with Modified Metallic Surfaces Apposing Bone or Bone Cement",
- "Guidance for Industry on the Testing of Metallic Plasma Sprayed Coatings on Orthopedic Implants to Support Reconsideration of Postmarket Surveillance Requirements"
- "Guidance Document for Testing Non-articulating, "Mechanically Locked" Modular Implant Components"
- "Data Requirements for Ultrahigh Molecular Weight Polyethylene (UHMPE) Used in Orthopedic Devices"

000004



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

MAY 17 2006

Mr. Matthew M. Hull, RAC  
Regulatory Affairs Manager  
Aesculap, Inc.  
3773 Corporate Parkway  
Center Valley, Pennsylvania 18034

Re: K053579  
Trade/Device Name: Columbus Total Knee System Ultra Congruent (UC) Tibial Insert  
Regulation Number: 21 CFR 888.3560  
Regulation Name: Knee joint patellofemorotibial polymer/metal/polymer semi-constrained cemented prosthesis  
Regulatory Class: II  
Product Code: JWH  
Dated: May 09, 2006  
Received: May 10, 2006

Dear Mr. Hull:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21

CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Mark N. Melkerson, M.S.  
Director  
Division of General, Restorative and Neurological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

Page 1 of 1

**A. INDICATIONS FOR USE STATEMENT**

510(k) Number: K053579

Device Name: Columbus (CR) Total Knee System

**Indications for Use:**

The Columbus (CR) Total Knee System is indicated for use in reconstruction of the diseased knee joint caused by osteoarthritis, rheumatoid arthritis, post-traumatic arthritis, the need to revise failed arthroplasties or osteotomies where pain, deformity or dysfunction persist, and for patients suffering from correctable valgus or varus deformity and moderate flexion contracture.

The Columbus Knee (CR) is designed for use with bone cement.

Prescription Use   X   and/or Over-the-Counter Use                       
(per 21 CFR 801.109)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH Office of Device Evaluation (ODE)

*[Signature]*  
(Division Sign-Off)

**Division of General, Restorative,  
and Neurological Devices**

510(k) Number K053579  
000002