

FEB 21 2006

Section 3  
ACL Elite and ACL Elite Pro  
510(k) Summary

K 060162

**Submitted by:**

Instrumentation Laboratory Company  
113 Hartwell Avenue  
Lexington, MA 02421

**Contact Person:**

Carol Marble, Regulatory Affairs Director  
Phone: 781-861-4467  
Fax: 781-861-4207

**Summary Prepared:**

January 20, 2006

**Device Name:**

ACL Elite and ACL Elite Pro

**Classification Name:**

Regulation Section: Multipurpose system for *in vitro* coagulation studies (864.5425)  
Classification: Class II  
Product Code: JPA  
Panel: Hematology

**Legally Marketed Device:**

K000053 ACL 9000 (ACL 8000 and ACL 10000)

**Device Description:**

The ACL Elite and ACL Elite Pro are fully automated, high-productivity analyzers designed specifically for *in vitro* diagnostic clinical use in the hemostasis laboratory for coagulation and/or fibrinolysis testing in the assessment of thrombosis and/or hemostasis. The systems provide results for both direct hemostasis measurements and calculated parameters.

The analytical modifications to the ACL 8000/9000/10000 introduced with the new family members (ACL Elite and ACL Elite Pro) include the following:

- Modification to the test parameter (revised use of the Secondary Algorithm) for several PT and APTT assays.
- Addition of dedicated calibration for HemosIL Factor VIII and HemosIL Factor IX Deficient Plasmas, when used in conjunction with HemosIL APTT-SP.
- Introduction of automated factor parallelism, which is a manual function on the ACL 8000/9000/10000.

**Section 3 (Cont.)**  
**ACL Elite and ACL Elite Pro**  
**510(k) Summary**

**Statement of Technological Characteristics of the Device Compared to Predicate Device:**

The indications for use/intended use, operating principle and performance claims of the ACL Elite and ACL Elite Pro are substantially equivalent to the legally marketed device, the ACL 8000/9000/10000 (K000053).

**Summary of Performance Data:**

**Precision**

Within run and total precision assessed over multiple runs on the ACL Elite/ACL Elite Pro using three levels of control plasma for factors deficient plasmas, APTT and PT, and two levels of control plasmas for fibrinogen gave the following results:

| Test                                     | FVIII (% Activity)           |    |       |                 |            |
|------------------------------------------|------------------------------|----|-------|-----------------|------------|
|                                          | Material                     | n  | Mean  | Within Run % CV | Total % CV |
| HemosIL Factor VIII (K034007) w/ APTT-SP | Normal Control               | 50 | 85.3  | 4.1             | 4.5        |
|                                          | Special Test Control Level 1 | 50 | 72.4  | 2.9             | 3.6        |
|                                          | Special Test Control Level 2 | 50 | 36.3  | 4.6             | 5.4        |
|                                          |                              |    |       |                 |            |
| FIX (% Activity)                         |                              |    |       |                 |            |
| Test                                     | Material                     | n  | Mean  | Within Run % CV | Total % CV |
|                                          |                              |    |       |                 |            |
| HemosIL Factor IX (K031829) w/ APTT-SP   | Normal Control               | 50 | 115.6 | 2.7             | 5.2        |
|                                          | Special Test Control Level 1 | 50 | 76.2  | 2.3             | 4.5        |
|                                          | Special Test Control Level 2 | 50 | 41.0  | 3.3             | 5.1        |
|                                          |                              |    |       |                 |            |
| APTT (Seconds)                           |                              |    |       |                 |            |
| Test                                     | Material                     | n  | Mean  | Within Run % CV | Total % CV |
|                                          |                              |    |       |                 |            |
| HemosIL APTT-SP (K973306)                | Normal Control               | 50 | 29.3  | 1.1             | 2.0        |
|                                          | Low Abnormal Control         | 50 | 46.1  | 1.9             | 4.1        |
|                                          | High Abnormal Control        | 50 | 56.9  | 2.2             | 4.7        |
|                                          |                              |    |       |                 |            |
| APTT (Seconds)                           |                              |    |       |                 |            |
| Test                                     | Material                     | n  | Mean  | Within Run % CV | Total % CV |
|                                          |                              |    |       |                 |            |
| HemosIL SynthASil (K953981)              | Normal Control               | 50 | 29.0  | 0.8             | 1.8        |
|                                          | Low Abnormal Control         | 50 | 51.9  | 1.0             | 1.5        |
|                                          | High Abnormal Control        | 50 | 60.5  | 1.2             | 1.2        |
|                                          |                              |    |       |                 |            |

**Section 3 (Cont.)**  
**ACL Elite and ACL Elite Pro**  
**510(k) Summary**

**Summary of Performance Data (Cont.):**

**Precision (Cont.)**

| Test                                                 | Prothrombin Time (Seconds) |       |       |                 |            |
|------------------------------------------------------|----------------------------|-------|-------|-----------------|------------|
|                                                      | Material                   | n     | Mean  | Within Run % CV | Total % CV |
| HemosIL<br>PT-Fibrinogen<br>(K862301)                | Normal Control             | 50    | 12.3  | 0.6             | 1.1        |
|                                                      | Low Abnormal Control       | 50    | 20.7  | 0.9             | 1.1        |
|                                                      | High Abnormal Control      | 50    | 25.5  | 0.8             | 1.4        |
|                                                      | Fibrinogen (mg/dL)         |       |       |                 |            |
|                                                      | Material                   | n     | Mean  | Within Run % CV | Total % CV |
|                                                      | Normal Control             | 50    | 277.9 | 4.4             | 4.6        |
|                                                      | Low Fibrinogen Control     | 50    | 132.9 | 3.6             | 5.5        |
| HemosIL<br>PT-Fibrinogen<br>Recombinant<br>(K981479) | Prothrombin Time (Seconds) |       |       |                 |            |
|                                                      | Material                   | n     | Mean  | Within Run % CV | Total % CV |
|                                                      | Normal Control             | 50    | 10.7  | 1.0             | 1.1        |
|                                                      | Low Abnormal Control       | 50    | 35.0  | 1.1             | 2.0        |
|                                                      | High Abnormal Control      | 50    | 51.6  | 1.4             | 3.4        |
|                                                      | Fibrinogen (mg/dL)         |       |       |                 |            |
|                                                      | Material                   | n     | Mean  | Within Run % CV | Total % CV |
| Normal Control                                       | 50                         | 297.6 | 3.8   | 6.7             |            |
| Low Fibrinogen Control                               | 50                         | 113.1 | 5.3   | 6.4             |            |

**Method Comparison**

In comparison studies evaluating citrated plasma samples, the ACL Elite/ACL Elite Pro versus the ACL 9000 (predicate device) was shown to be statistically similar for the tests listed below:

| Test                                 | Parameter                 | n  | Slope | r      |
|--------------------------------------|---------------------------|----|-------|--------|
| HemosIL Factor VIII                  | FVIII w/ APTT-SP (% Act.) | 47 | 0.909 | 0.9938 |
| HemosIL Factor IX                    | FIX with APTT-SP (% Act.) | 49 | 0.948 | 0.9914 |
| HemosIL APTT-SP                      | APTT (Seconds)            | 54 | 1.032 | 0.9990 |
| HemosIL SynthASil                    | APTT (Seconds)            | 52 | 1.009 | 0.9981 |
| HemosIL<br>PT-Fibrinogen             | PT (sec)                  | 49 | 0.989 | 0.9985 |
|                                      | Fibrinogen (mg/dL)        | 47 | 0.996 | 0.9848 |
| HemosIL<br>PT-Fibrinogen Recombinant | PT (sec)                  | 49 | 0.972 | 0.9993 |
|                                      | Fibrinogen (mg/dL)        | 47 | 0.974 | 0.9992 |



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

Ms. Carol Marble  
Regulatory Affairs Director  
Instrumentation Laboratory Company  
113 Hartwell Avenue  
Lexington, MA 02421

FEB 21 2006

Re: k060162  
Trade/Device Name: ACL Elite and ACL Elite Pro  
Regulation Number: 21 CFR 864.5425  
Regulation Name: Multipurpose system for in vitro coagulation studies  
Regulatory Class: II  
Product Code: JPA  
Dated: January 20, 2006  
Received: January 23, 2006

Dear Ms. Marble:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820). This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

Page 2 –

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (240) 276-0484. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>

Sincerely yours,



Robert L. Becker, Jr., MD, Ph.D

Director

Division of Immunology and Hematology

Office of *In Vitro* Diagnostic Device

Evaluation and Safety

Center for Devices and

Radiological Health

Enclosure

## Indications for Use Statement

510(k) Number (if known): K060162

Device Name: ACL Elite and ACL Elite Pro

### Indications for Use:

The ACL Elite and ACL Elite Pro are fully automated, high-productivity analyzers designed specifically for *in vitro* diagnostic clinical use in the hemostasis laboratory for coagulation and/or fibrinolysis testing in the assessment of thrombosis and/or hemostasis. The systems provide results for both direct hemostasis measurements and calculated parameters.

Prescription Use    
 (Part 21 CFR 801 Subpart D)


AND/OR

Over-The-Counter Use \_\_\_\_\_   
 (21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)

  
\_\_\_\_\_  
Division Sign-Off

Office of In Vitro Diagnostic Device  
Evaluation and Safety

510(k) K060162