K060341

Section 5

Sherlock™ Tip Location System 510(k) Summary of Safety and Effectiveness 21 CFR 807.92(a)

AFT 1 COOK

5.1 General Information

Submitter Name:

Bard Access Systems, Inc. (BAS)

[Wholly owned Subsidiary of C. R. Bard, Inc.]

Address:

5425 W. Amelia Earhart Drive

Salt Lake City, UT 84116

Telephone Number:

(801) 595-0700 ext. 7136

Fax Number:

(801) 595-5425

Contact Person:

Lynn M. Kirchoff

Date of Preparation:

February 9, 2006

Registration Number:

1720496

Additional Registration Numbers:

C.R. Bard:

2212754

5.2 Subject Device Information

Device Name:

Sherlock™ Tip Location System

Trade Name:

Sherlock TM Tip Location System (TLS)

Common/Usual Name:

Vascular Access Catheter Accessories

Classification Name:

80 LJS - Accessory to Percutaneous, Implanted, Long-Term

Intravascular Catheter

21 CFR 880.5970- Class II

Classification Panel:

General Hospital

5.3 Predicate Device Information

Device Name(s):

5 Fr Dual Lumen Groshong® nXt PICC Catheter

Poly Per-Q-Cath PICC

Trade Name(s):

Groshong® nXt, Poly Per-Q-Cath®

Common/Usual Name:

Peripherally Inserted Central Catheter (PICC)

Classification Name:

80 LJS - Percutaneous, Implanted, Long-Term Intravascular Catheter

21 CFR 880.5970- Class II

Classification Panel:

General Hospital

Predicate Device Name	#510(k) **	· Clearance Date
5 Fr DL Groshong® nXt PICC Catheter	K023374	12/18/2002
Poly Per-Q-Cath® PICC	K034019	1/21/2004

Device Name:

Sensor Stylet and Sensor Guidewire

Trade Name:

Catheter Locator System

Common/Usual Name:

Vascular Access Catheter Accessories

Classification Name:

80 LJS – Accessory to Percutaneous, Implanted, Long-Term

Intravascular Catheter

21 CFR 880.5970- Class II

Classification Panel:

General Hospital

Predicate Device Name ** *** ***	510(k)	Clearance Date
Sensor Stylet and Sensor Guidewire	K935380	2/17/1995

Device Name:

Zortran Detector

Trade Name:

Zortran Detector

Common/Usual Name:

Vascular Access Catheter Accessories

Classification Name:

FOZ - Accessory to Percutaneous, Implanted, Intravascular Catheter

21 CFR 880.5970- Class II

Classification Panel:

General Hospital

Predicate Device Name	். \$510(k)∂ ்	Clearance Date
Zortran Detector	K000997	9/7/2000

5.4 Intended Use

The intended use of the subject device is the same as the intended use of the predicate devices.

Catheter stylets provide internal reinforcement to aid in catheter placement. The SherlockTM TLS Stylet contains passive magnets that generate a magnetic field. This field can be detected by the SherlockTM TLS Detector to provide the placer rapid feedback on catheter tip location.

The Sherlock TLS Detector quickly locates and confirms the position of specially designed, magnet-tipped Peripherally Inserted Central Catheters (PICCS) and Central Venous Catheters (CVCs) during and after initial placement. This device may be used by appropriate caregivers in hospitals, long-term care facilities or home-care settings. The Sherlock TLS Detector provides rapid feedback to the caregiver, but was not designed to replace conventional methods of placement verification. Users are urged to confirm correct placement according to their established institutional protocol and clinical judgment.

5.5 Indications for Use

The Indications for Use for the SherlockTM Tip Location System (TLS) is as follows:

Catheter stylets provide internal reinforcement to aid in catheter placement. The SherlockTM TLS Stylet contains passive magnets that generate a magnetic field. This field can be detected by the SherlockTM TLS Detector to provide the placer rapid feedback on catheter tip location.

5.6 Summary of Changes

The following are modifications made to the currently marketed predicate PICC stylet and Sensor Stylet:

• The addition of encapsulated passive magnetic material in the distal portion of the stylet to allow for detection when used with a magnetic sensor.

5.7 Device Description

The Sherlock Tip Location System (TLS) consists of the Sherlock TLS Detector and Sherlock TLS Stylet.

The Sherlock TLS Stylet has been developed to aid in the placement of Bard Access Systems catheters using current placement techniques. The stylets are designed to give the catheters added support and stiffness while traversing the patient's venotomy. Also, should the clinician choose to do so, the stylets have been designed to be used in conjunction with Sherlock TLS Detector to allow for rapid feedback of catheter tip placement. The information provided by the Sherlock TLS is not meant to replace conventional methods of catheter placement verification. Clinicians are urged to confirm correct catheter placement according to their established institutional protocol and clinical judgment.

5.8 Technological Comparison to Predicate Devices

The technological characteristics of the Sherlock TLS are substantially equivalent to those of the predicate devices in terms of intended use, application, user population, basic design, performance, labeling, packaging and sterilization method.

5.9 510(k) Substantial Equivalence Decision Tree

New device is compared to Marketed Device? Yes.

Does the new device have the same indication statement and intended use as the predicate? Yes. However, there are minor modifications to the indication verbiage.

Does the new device have the same technological characteristics, e.g. design, materials, etc.? Not in all regards. The principles of operation and basic design are the same as the predicate devices. The main changes in design are the addition of passive magnetic material encapsulated in the distal portion of the stylet.

Could the new characteristics affect safety or effectiveness?

Yes. The design changes may affect safety or effectiveness of the device.

Do the new characteristics raise new types of safety and effectiveness questions? No. Safety and effectiveness questions are the same as for the predicate devices.

Do accepted scientific methods exist for assessing effects of the new characteristics? Yes. Testing was based on FDA guidance document and standards to evaluate the devices' performance:

- Coronary and Cerebrovascular Guidewire FDA Guidance, dated 1/95
- ISO 11070:1998, Sterile, single-use intravascular catheter introducers
- AAMI/ANSI/ISO 11135:1994, Medical Devices Validation and Routine Control of Ethylene Oxide Sterilization
- AAMI/ANSI/ISO 10993-1:1997, Biological Evaluation of Medical Devices Part 1: Evaluation and Testing, and the FDA Modified ISO 10993 Test Profile

Are performance data available to assess effects of new characteristics?

Yes. Bench testing was based on the above referenced guidance document and standards. All test results confirm that the modified device is substantially equivalent to the predicate devices.

5.10 Conclusion

The Sherlock™ Tip Location System (TLS) met all the performance criteria of the tests performed and, based on FDA's decision tree, is substantially equivalent to its predicate devices covered by: K023374, K034019, K935380, and K000997.





Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

APR 1 4 2006

C.R. Bard, Incorporated C/O Ms. Lynn M. Kirchoff Regulatory Affairs Specialist Bard Access Systems, Incorporated 5425 West Amelia Earhart Drive Salt Lake City, Utah 84116

Re: K060341

Trade/Device Name: Sherlock™ Tip Location System (TLS)

Regulation Number: 880.5200

Regulation Name: Intravascular Catheter

Regulatory Class: II Product Code: FOZ Dated: February 9, 2006 Received: February 10, 2006

Dear Ms. Kirchoff:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device. Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food. Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations. Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act): 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/industry/support/index.html.

Sincerely yours,

Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital, Infection Control and Dental Devices Office of Device Evaluation Center for Devices and Radiological Health

Enclosure

Section 4 Indications for Use

10(k) Number (if known): 19060341
Device Name: Sherlock™ Tip Location System (TLS)
ndications for Use:
Catheter stylets provide internal reinforcement to aid in catheter placement. The Sherlock TM TLS Stylet ontains passive magnets that generate a magnetic field. This field can be detected by the Sherlock TM TLS Detector to provide the placer rapid feedback on catheter tip location.
Prescription Use X Over-The-Counter Use (Part 21 CFR 801 Subpart D) AND/OR (21 CFR 801 Subpart C)
(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF NEEDED)
Concurrence of CDRH, Office of Device Evaluation (ODE)