

MAY 18 2006

510(k) Summary
Dimension® TACR Flex® reagent cartridge

This summary of 510(k) safety and effectiveness information is submitted in accordance with the requirements of SMDA 1990 and 21 CFR §807.92.

The assigned 510(k) number is: K060502

1. Manufacturer's Name, Address, Telephone, and Contact Person, Date of Preparation

Manufacturer: Dade Behring Inc.
P.O. Box 6101
Newark, DE 19714

Contact Information: Dade Behring Inc.
P.O. Box 6101
Newark, DE 19714
Attn: Yuk-Ting Lewis
Tel: 302-631-7626

Date of Preparation: May 2, 2006

2. Device Name / Classification

Dimension® TACR Flex® reagent cartridge / Class II

3. Identification of the Predicate Device

Abbott IMx® Tacrolimus II Assay, P970007
(Note: Tacrolimus test systems have been reclassified into Class II since the predicate was approved.)

4. Device Description

The Dimension® TACR Flex® reagent cartridge is an in vitro diagnostic device that consists of prepackaged reagents in a plastic eight-well cartridge for use on the Dade Behring Dimension® clinical chemistry system.

5. Device Intended Use

The Dimension® TACR Flex® reagent cartridge is an in vitro diagnostic test intended to quantitatively measure tacrolimus in human whole blood on the Dimension® clinical chemistry system. Measurements of tacrolimus are used as an aid in the management of tacrolimus therapy in kidney and liver transplant patients.

6. Medical device to which equivalence is claimed and comparison information

The Dimension® TACR Flex® reagent cartridge is substantially equivalent in intended use and technological characteristics to the Abbott IMx® Tacrolimus II Assay. Both devices are immunoassays intended for use in the quantitative measurement of tacrolimus in human whole blood. The Dimension® TACR Flex® reagent cartridge has an assay range of 1.2-30 ng/mL. The Abbott IMx® Tacrolimus II Assay has an assay range of 1.5-30 ng/mL

Comparison Information

Method comparison studies were conducted at two external sites comparing the Dimension® TACR Flex® reagent cartridge against two predicates:

- liquid chromatography / tandem mass spectrometry (LC/MS/MS) and
- the Abbott IMx® Tacrolimus II Assay.

Samples from 2 transplant patient groups (liver and kidney) were included in the studies. The data from both sites were pooled and analyzed by Passing Bablok linear regression.

Comparative Method

LC/MS/MS	Slope	Intercept	r	n
All samples	1.13	-0.27	0.88	184
Kidney	1.16	-0.60	0.88	103
Liver	1.06	0.49	0.88	81

Abbott IMx® Tacrolimus II Assay

All samples	0.92	0.10	0.85	175
Kidney	0.90	-0.06	0.87	97
Liver	0.93	0.45	0.82	78



Food and Drug Administration
2098 Gaither Road
Rockville MD 20850

MAY 18 2006

Ms. Yuk-Ting Lewis
Regulatory Affairs & Compliance Manager
Dade Behring, Inc.
PO Box 6101, M/S 514
Newark, DE 19714

Re: k060502
Trade/Device Name: Dimension® TACR Flex® reagent cartridge
Regulation Number: 21 CFR§862.1678
Regulation Name: Tacrolimus test system
Regulatory Class: Class II
Product Code: MLM
Dated: April 20, 2006
Received: April 21, 2006

Dear Ms. Lewis:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

Page 2 –

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (240) 276-0484. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Alberto Gutierrez, Ph.D.

Director

Division of Chemistry and Toxicology

Office of In Vitro Diagnostic Device

Evaluation and Safety

Center for Devices and

Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K060502

Device Name: Dimension® TACR Flex® reagent cartridge

Indications For Use:

The Dimension® TACR Flex® reagent cartridge is an in vitro diagnostic test intended to quantitatively measure tacrolimus in human whole blood on the Dimension® clinical chemistry system. Measurements of tacrolimus are used as an aid in the management of tacrolimus therapy in kidney and liver transplant patients.


Prescription Use x
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 801)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)



Division Sign-Off

Office of In Vitro Diagnostic Device
Evaluation and Safety

510(k) K060502