K061048

510(k) Summary – COBAS INTEGRA Glucose HK Gen. 3

1 1**7** - 2 15 2

Introduction	According to the requirements of 21 CFR 807.92, the following information provides sufficient detail to understand the basis for a determination of substantial equivalence
Submitter name, address, contact	Roche Diagnostics 9115 Hague Rd Indianapolis IN 46250 (317) 521-3723
	Contact person: Corina Harper
	Date prepared: April 11, 2006
Device Name	Proprietary name: COBAS INTEGRA Glucose HK Gen. 3 test
	Common name: Glucose HK Gen. 3
	Classification name: Glucose Test System
Device Description	The cassette COBAS INTEGRA Glucose HK Gen. 3 contains an in vitro diagnostic reagent system intended for use on COBAS INTEGRA SYSTEMS for the quantitative determination of glucose in serum, plasma, urine, and cerebrospinal fluid (CSF).
	The test principle is an enzymatic reference method with hexokinase.
Intended use	In vitro test for the quantitative determination of glucose in serum, plasma, urine and cerebrospinal fluid (CSF) on COBAS INTEGRA systems.
Predicate Device	We claim substantial equivalence to the COBAS INTEGRA Glucose HK Liquid cleared as K972250.
Substantial equivalency – Similarities	The table below indicates the similarities between the modified COBAS INTEGRA Glucose HK Gen. 3 test and its predicate device (COBAS INTEGRA Glucose HK Liquid, K972250).

Feature	Predicate device: Glucose HK Liquid (K972250)	Modified device: Glucose HK Gen. 3	
General			
Intended Use/	The cassette COBAS INTEGRA	In vitro test for the quantitative	
Indications for	Glucose HK Liquid contains an in	determination of glucose in serum,	
Use	vitro diagnostic reagent system	plasma, urine and cerebrospinal	
	intended for use on COBAS	fluid (CSF) on COBAS INTEGRA	
	INTEGRA systems for the	systems	
	quantitative determination of		
	glucose concentration in serum,		
	fluid (CSF)		
Specimen type	Serum, plasma, urine, CSF	Same	
Test principle			
Reference	Enzymatic reference method with	Same	
method	hexokinase.		
Reagent informati	on		
Stability - shelf	2-8 °C until expiration date	Same	
life and on-board			
	COBAS INTEGRA 400		
	o weeks at 10 to 15 C		
	COBAS INTEGRA 700/800		
	8 weeks at 8°C		
Calibrator	Calibrator f.a.s.	Same	
Or 11 to a sector 1	Interval: each lot	Sama .	
Quality control	Serum and plasma:	Same	
	Precipath and Precipath U Plus		
	Urine:		
	Quantitative urine controls		
	<u>CSF:</u>		
	Quantitative CSF controls		
Traceability	Standardized against Isotone	Same	
Traccability	Dilution Mass Spectrometry	Same	
Performance char	acteristics		
Measuring range	0-40 mmol/L (0-720 mg/dL)	0.12-40 mmol/L (0.12-720 mg/dL)	
	Extended measuring range with	Extended measuring range with	
	recommended post dilution factor of $10: 0.400 \text{ mmol/L} (0.7200 \text{ mmol/L})$	recommended post dilution factor $a = 12, 0, 12, 400 \text{ mmol/L}$ (2.16.7200	
	10. 0.400 mmore (0-7200 mg/dE)	mg/dL)	

Expected values (literature reference) Additional values are referenced in the method sheet	Plasma (fasting): 3.88-6.38 mmol/L Urine: 1 st morning urine 0.3-1.1 mmol/L 24 h urine 0.11-0.50 mmol/24h Serum/plasma: Adults 4.11-5.89 mmol/L	Plasma (fasting): 3.88-6.38 mmol/L Urine: 1 st morning urine 0.3-1.1 mmol/L 24 h urine 0.3-0.96 mmol/L Same
		· · · · · · · · · · · · · · · · · · ·

Substantial	The table below indicates the similarities between the modified COBAS
equivalency –	INTEGRA Glucose HK Gen. 3 test and its predicate device (COBAS
Differences	INTEGRA Glucose HK Liquid, K972250).

.

Feature	Predicate device: Glucose HK Liquid (K972250)	Modified device: Glucose HK Gen.3
Reagent informat	ion	
R1 R2	Mono reagent in vial A and B (liquid) 100 mmol/L MOPS, 12 mmol/L ATP, 6 mmol/L NAD ⁺ , 10 mmol/L	R1: 5.0 mmol/L MES, ≥4.5 mmol/L ATP, 24 mmol/L Mg ⁺⁺ , ≥7.0 mmol/L NADP, pH 6.0
	Mg ⁺⁺ , ≥50 µcat/L HK(yeast), ≥50 µcat/L G6PDH (microbial), 0.09% Sodium azide, pH 7.1	R2: 4.0 mmol/L Mg ⁺⁺ , 200 mmol/L HEPES, \geq 300 µcat/L HK (yeast), \geq 300 µcat/L G6PDH (microbial), pH 8.0
Performance char	acteristics	

Precision	<u>Serum and plasma:</u> Within run: 1.7% @ 5.3 mmol/L 0.72% @ 33.2 mmol/L Between run: 2.6% @ 5.3 mmol/L 1.5% @ 33.2 mmol/L	Serum and plasma: Within run CV%: 0.41% @ 4.48 mmol/L 0.47% @ 12.48 mmol/L Between day: 1.09% @ 4.44 mmol/L 0.90% @ 12.46 mmol/L
	Urine application Within run: 1.7% @ 1.7 mmol/L 1.8% @ 37.1 mmol/L Between run: 4.3% @ 1.7 mmol/L 2.9% @ 37.1 mmol/L	Urine application Within run: 1.35% @ 0.83 mmol/L 0.64% @ 2.42 mmol/L Between day: 0.75% @ 0.84 mmol/L 0.83% @ 2.43 mmol/L
	<u>CSF application</u> Within run: 1.6% @ 1.7 mmol/L 1.8% @ 3.3 mmol/L Between run: 2.3% @ 1.7 mmol/L 1.9% @ 3.3 mmol/L	CSF application Within run: 1.13% @ 3.20 mmol/L 1.49% @ 9.31 mmol/L
Linearity	0-40 mmol/L (before dilution)	0.12-40 mmol/L (before dilution)
Lower detection limit	Serum and plasma: 0.033 mmol/L <u>Urine application:</u> 0.22 mmol/L <u>CSF application</u> : 0.023 mmol/L	Serum, plasma, urine and CSF: 0.12 mmol/L
Endogenous interferences	Hemolysis no significant interferences Icterus no significant interferences	Hemolysis: up to 1200 H Index Icterus: up to 60 I Index
	Lipemia no significant interferences	Lipemia: up to 1900 L Index

Ex Int	ogenous terferences	Falsely low results may be caused by elevated pyruvates levels	Tetracyclin at therapeutic concentration gives falsely low results in urine samples	
		Gammopathy, in particular IgM, may cause unreliable results in rare cases	Same	
Proposed Propo Labeling direct of the require		oposed labeling sufficient to describe the rections for use can be found in Section V the device labeling presented contains al quired per 21 CFR 809.10.	device, its intended use, and the /. We believe the proposed version l of the technical information	
Validatio Design C	on and D Control pr of pr	Development activities were conducted under appropriate design control procedures and the overall product specifications were met. The Declaration of Conformity with Design Controls and Results of Risk Analysis are provided in Section 5.1. Analytical Performance.		
Confide	entiality R or be	Roche Diagnostics Corporation requests that the FDA not disclose the nature or existence of this submission until the substantial equivalence decision has been reached.		
Closing	M th la Tl su G	Modification of the COBAS INTEGRA Glucose HK Gen.3 does not affect the intended use or indications for use of the device as described in the labeling, nor does it alter the fundamental scientific technology of the device. Therefore, we trust the information provided in this Special 510(k) will support a decision of substantial equivalence of the COBAS INTEGRA Glucose HK Gen.3 to the predicate.		
	If he	If you have any questions or require further information, please do not hesitate to contact this office.		
	•	Phone: (317) 521-3831 FAX: (317) 521-2324 email: corina.harper@roche.com		

• • •

DEPARTMENT OF HEALTH & HUMAN SERVICES



Public Health Service

Food and Drug Administration 2098 Gaither Road Rockville MD 20850

MAY 1 7 2006

Ms. Corina Harper Regulatory Affairs Consultant Roche Diagnostics 9115 Hague Road PO Box 50416 Indianapolis, IN 46250-0416

Re: k061048 Trade/Device Name: COBAS INTEGRA Glucose HK Gen 3 Regulation Number: 21 CFR§ 862.1345 Regulation Name: Glucose test system Regulatory Class: Class II Product Code: CFR Dated: April 14, 2006 Received: April 17, 2006

Dear Ms. Harper:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

Page 2 –

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (240) 276-0484. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/industry/support/index.html.

Sincerely yours,

Alberto Gutierrez, Ph.D. Director Division of Chemistry and Toxicology Office of In Vitro Diagnostic Device Evaluation and Safety Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K06/048

Device Name: COBAS INTEGRA Glucose HK Gen 3

Indications For Use:

In vitro test for the quantitative determination of glucose in serum, plasma, urine and cerebrospinal fluid (CSF) on COBAS INTEGRA systems.

Glucose measurements are used in diagnosis and treatment of carbohydrate metabolism disorders including diabetes mellitus, neonatal hypoglycemia, and idiopathic hypoglycemia, and pancreatic islet cell tumors.

Prescription Use XXXX (Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)

Division \$ign-Off

Office of In Vitro Diagnostic Device **Evaluation and Safety**

61048