

K061496

**510(k) SUMMARY OF SAFETY AND EFFECTIVENESS  
SAS™ Legionella Test**

JUN 28 2007

The 510(k) summary of safety and effectiveness submission is in accordance with the requirements of SMDA 1990 and 21 CFR 807.92

Submitted by: SA Scientific, Ltd.  
4919 Golden Quail  
San Antonio, TX 78240

Establishment Reg. No. 1645225

Contact Person: Veronica DeLeon

Date Prepared: May 30, 2006

Proprietary Name: SAS™ Legionella Test

Common Name: SAS™ Legionella Test, SAS™ *Legionella pneumophila* Test

Classification Name: Legionella, SPP., Elisa

Device Classification: 21 CFR 866.3300

Regulatory Class: Class II

Product Code: MJH

Substantial Equivalence: Binax™ Now® Legionella Urinary Antigen Test manufactured by Binax™ Inc., Portland, Maine

Device Description: The SAS™ Legionella test utilizes a combination of polyclonal antibodies against the antigens of *Legionella pneumophila*. The SAS™ Legionella test begins with the addition of urine to the test device. The specimen is absorbed by the sample pad and then moves through the conjugate pad which contains dried gold conjugated antibodies which are specific for *Legionella pneumophila* antigens; if the Legionella antigens are present in the urine sample, a "half-sandwich" immunocomplex is formed. This immuno-complex then migrates via capillary action along a nitrocellulose membrane containing immobilized antibodies to *Legionella pneumophila* antigens. The immobilized antibodies bind the "half-sandwich" immuno-complex to

form a “whole sandwich” immuno-complex. Thus, when the “whole sandwich” is formed, a visible, pink colored line develops in the specimen zone on the test device. In the absence of a Legionella antigen, a “sandwich” immuno-complex is not formed and a negative result is indicated. To serve as a procedural control, a pink colored control line will always appear in the control zone regardless of the presence or absence of Legionella antigen. The test is available in cassette format.

**Intended Use:** The SAS™ Legionella Test is a visually read, in vitro immunochromatographic rapid assay for the presumptive qualitative detection of *Legionella pneumophila* serogroup 1 antigens in human urine. This test is intended to aid in the presumptive diagnosis of Legionnaires' disease in conjunction with culture and other methods for patients with signs and symptoms of pneumonia. This test is for prescription use only.

**Quality Controls:** The SAS™ Legionella Test provides two (2) internal procedural controls. It is recommended that external quality controls should be performed on each new test kit box opened. Positive and negative external controls are supplied separately.

**Device Comparison:** The SAS™ Legionella Test and Binax™ Now® Legionella Urinary Antigen Test are rapid immunoassays tests utilizing immunochromatographic technology for the visualization of *Legionella pneumophila* antigen. Each utilizes an antibody conjugated to colored particles and an antibody printed onto a membrane.

**Performance Summary:** The SAS™ Legionella Test performed substantially equivalent to the predicate device, Binax™ Now® Legionella Urinary Antigen Test and to culture. This was verified by comparison of frozen and fresh urine specimens.

Cross reactivity and interference studies were performed on viral and bacterial strains commonly found in human urine. None of the organisms interfered or cross-reacted with the performance of the SAS™ Legionella Test.

Prepared by: Veronica DeLeon  
Veronica DeLeon  
Regulatory Affairs

Date: 6/21/07



Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

JUN 28 2007

Veronica DeLeon  
Regulatory Affairs  
SA Scientific, Ltd.  
4919 Golden Quail  
San Antonio, TX 78240

Re: k061496

Trade/Device Name: SAS™ Legionella Test  
Regulation Number: 21CFR 866. 3300  
Regulation Name: Haemophilus spp. Serological Reagents  
Regulatory Class: Class II  
Product Code: MJH  
Dated: May 30, 2006  
Received: May 31, 2006

Dear Ms. DeLeon:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice

Page 2 –

requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820). This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (301) 594-3084. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150, or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,



Sally A. Hojvat, M.Sc., Ph.D.  
Director  
Division of Microbiology Devices  
Office of *In Vitro* Diagnostic Device  
Evaluation and Safety  
Center for Devices and Radiological Health

# INDICATIONS FOR USE

Applicant: SA Scientific, Inc.  
4919 Golden Quail  
San Antonio, TX. 78240  
Ph: (210) 699-8800 Fax: (210) 699-6545

Establishment Reg. No.: 1645225

501(k) Number: 061496

Device Names: SAS™ Legionella Test

## Indications for Use:

The SAS™ Legionella Test is a visually read, in vitro immunochromatographic rapid assay for the presumptive qualitative detection of *Legionella pneumophila* serogroup 1 antigens in human urine. This test is intended to aid in the presumptive diagnosis of Legionnaires' disease in conjunction with culture and other methods for patients with signs and symptoms of pneumonia. This test is for prescription use only.


Prescription Use X  
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)

  
Division Sign-Off

Page 1 of 1

Office of In Vitro Diagnostic  
Device Evaluation and Safety

vi

510(k) K061496