

**510(k) Summary for the
Dimension Vista™ System Total Triiodothyronine Calibrator
(T3 CAL – KC250)**

A. 510(k) Number: K061885

AUG 14 2006

B. Analyte: Total Triiodothyronine (T3)

C. Type of Test: Calibrator Material

D. Applicant: Dade Behring Inc., P.O. Box 6101, Newark, DE 19714-6101
Victor M. Carrio, Regulatory Affairs and Compliance Manager
Office: (302) 631-0376 Fax: (302) 631-6299

E. Proprietary and Established Names:

Dimension Vista™ System Total Triiodothyronine Calibrator (T3
CAL – KC250)

F. Regulatory Information:

1. Regulation section: 21 CFR § 862-1150 – Calibrator
2. Classification: Class II
3. Product Code: JIT – Calibrator, Secondary
4. Panel: Clinical Chemistry

G. Intended Use: The T3 CAL is an *in vitro* diagnostic product for the calibration of Total Triiodothyronine (T3) method on the Dimension Vista™ System.

H. Device Description:

The T3 CAL is a liquid, human serum based product containing L-triiodothyronine. The kit consists of six vials, three vials of Calibrator A (2.0 mL per vial) and three vials of Calibrator B (1.5 mL per vial). T3 CAL is ready for use, no preparation is required.

I. Substantial Equivalence Information:

1. Predicate Device: K032697 – Dimension® Total Triiodothyronine Calibrator.
2. Comparison with Predicate:

	Device	Predicate
Item	Dimension Vista™ System Total Triiodothyronine Calibrator	Dimension® Total Triiodothyronine Calibrator
Intended Use	The T3 CAL is an <i>in vitro</i> diagnostic product for the calibration of Total Triiodothyronine (T3) method on the Dimension Vista™ System.	The T3 Calibrator is intended for use in the calibration of Total Triiodothyronine (T3) method on the Dimension® clinical chemistry system with the Heterogeneous Immunoassay Module.
Analytes	Total Triiodothyronine (T3).	Total Triiodothyronine (T3).
Form	Liquid	Liquid
Traceability	USP ¹ L-triiodothyronine (USP Catalog # 36800)	USP ¹ L-triiodothyronine (USP Catalog # 36800)
Matrix	Stripped human serum.base.	Stripped human serum.base.
Levels	Two levels.	Five levels.

¹ United States Pharmacopeia

J. Standard/Guidance Document Referenced:

1. Guidance: Guidance for Industry - Abbreviated 510(k) Submissions for In Vitro Diagnostic Calibrators; Final, 02/22/1999
Guidance for Industry and FDA Staff - Use of Symbols on Labels and in Labeling of In Vitro Diagnostic Devices Intended for Professional Use, 11/30/2004
2. Standards: CEN 13640 Stability testing of In-Vitro Diagnostic Devices
ISO 14971:2000 Medical devices -Application of risk management to medical devices

K. Performance Characteristics:

1. Stability: Target shelf life for the Dimension Vista™ System Total Triiodothyronine (T3) Calibrator is 24 months. Calibrator shelf life is determined by comparing results of the product stored at 4°C with control stored at -20°C. The method is calibrated from this stored material. The 4°C material values are recovered versus the calibration. Recovery versus time is monitored and percent change over time is determined. Percent change should be less than or equal to 8 %. Shelf-life stability (expiration) dating assignment at commercialization reflects the real-time data on file at Dade Behring, Inc.
A vial punctured by the instrument and stored on board is stable for one day.
An open vial not stored on board of the instrument, but recapped and stored in a refrigerator is stable for 31 days.

For testing, vials are opened /punctured on day zero. For opened vials, a quantity sufficient for multiple calibrations is removed and the vials are recapped and stored at 2 – 8 °C. Opened vials are tested on days 0, 2, 8, 15, 22, and 32 versus freshly opened vials. Punctured vials are tested once after at least 25 hours versus freshly opened vials.

2. Traceability: The assigned values of the Total Triiodothyronine (T3) Calibrator are traceable to the United States Pharmacopeia (USP) L-triiodothyronine (USP Catalog # 36800).

3. Value Assignment:

Master Pool is manufactured by weighing USP L-triiodothyronine (USP Catalog # 36800) into stripped human serum at five levels. The Master Pool is stored frozen, -10 to -20 °C. The bottle value for the Master Pool is assigned for each level by testing N = 45 replicates per level on multiple instruments, using a previous Master Pool as a control.

A Commercial Lot is manufactured by weighing Purified T3 into stripped human serum at two levels. The bottle value assignment of each level is tested using multiple instruments calibrated with Master Pool for N = 45 total replicates per level. A previous lot of commercial calibrator is used as control.



Food and Drug Administration
2098 Gaither Road
Rockville MD 20850

Mr. Victor M. Carrio
RA/QS Compliance Manager
Dade Behring, Inc.
PO Box 6101
Newark, DE 19714-6101

AUG 14 2006

Re: k061885
Trade/Device Name: Dimension Vista™ Total Triiodothyronine Calibrator
(T3 CAL, KC250)
Regulation Number: 21 CFR§862.1150
Regulation Name: Calibrator
Regulatory Class: Class II
Product Code: JIT
Dated: June 30, 2006
Received: July 3, 2006

Dear Mr. Carrio:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

Page 2 –

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (240) 276-0484. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Alberto Gutierrez, Ph.D.
Director
Division of Chemistry and Toxicology
Office of In Vitro Diagnostic Device
Evaluation and Safety
Center for Devices and
Radiological Health

Enclosure

Indications For Use Statement

510(k) Number (if known):

Device Name: **K061885**

Dimension Vista™ Total Triiodothyronine Calibrator (T3 CAL, KC250)

Indications for Use:

The T3 CAL is an *in vitro* diagnostic product for the calibration of Total Triiodothyronine (T3) method on the Dimension Vista™ System.

Prescription Use X
(Per 21 CFR 801 Subpart D)

AND/OR

Over-the-counter Use _____
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of -In Vitro Diagnostic Devices (OIVD)


Division Signatory
Office of In Vitro Diagnostic Device
Evaluation and Safety
510(k) K061885