

~~DEC 29 2006~~

## 9. 510(K) SUMMARY OF SAFETY AND EFFECTIVENESS

This summary of safety and effectiveness information is being submitted in accordance with the requirements of The Safety Medical Devices Act of 1990 (SMDA 1990) and 21 CFR Part 807.92.

Assigned 510(k) Number:

K062787

Date of Summary Preparation: September 15, 2006

Manufacturer:

Phadia AB  
Rapsgatan 7  
SE-751 37 Uppsala, Sweden

510 (k) Contact Person:

Martin Mann  
Regulatory Affairs Manager  
Phadia US Inc.  
4169 Commercial Avenue  
Portage, Mi 49002 , USA  
+1 (-269-492) -1957 (Phone)  
+1 (-269-492) -7541 (Fax)  
[martin.mann@phadia.com](mailto:martin.mann@phadia.com)

Device Name:

EliA™ Celikey IgA Well  
EliA™ Celiac Control

Common Name:

Tissue transglutaminase autoantibodies  
immunological test system, Antigen and Control

### Classification

<u>Product Name</u>	<u>Product Code</u>	<u>Class</u>	<u>CFR</u>
EliA™ Celikey IgA	MVM	II	866.5660
EliA™ Celiac Control	MVM	II	866.5660

### Substantial Equivalence to

Celikey Tissue Transglutaminase (human, recombinant) IgA Antibody Assay  
510(k) number: K041174

### **Intended Use Statement of the New Device**

EliA Celikey IgA is intended for the in vitro semi-quantitative measurement of IgA antibodies directed to tissue transglutaminase (tTG) in human serum and plasma. EliA Celikey IgA is based on recombinant human tissue transglutaminase as antigen and is useful as an aid in the clinical diagnosis of patients with celiac disease. EliA Celikey IgA uses the EliA IgA method on the instrument ImmunoCAP 100 and ImmunoCAP 250.

EliA Celiac Control is intended for laboratory use in monitoring the performance of in vitro measurement of antibodies to tissue transglutaminase (tTG) and gliadin with ImmunoCAP 100 or 250 using the EliA IgG or IgA method.

#### Special condition for use statement

The device is for prescription use only.

#### Special instrument requirements

ImmunoCAP100/ImmunoCAP250 are fully automated immunoassay analyzers, which include software for evaluation of test results.

### **General Description of the New Device**

The new device belongs to a fully integrated and automated system for immunodiagnostic testing. It comprises a Fluorescence-Immunoassay test system using EliA single wells as the solid phase and is intended to be performed on the instruments ImmunoCAP 100 and ImmunoCAP 250. The conjugate for the EliA IgA method is mouse anti-human IgA beta-galactosidase, which uses 4-Methylumbelliferyl- $\beta$ D-Galactoside as substrate. The total IgA calibration is based on a set of six WHO-standardized IgA Calibrators derived from human serum. They are used to establish an initial calibration curve, which may be used for up to 28 days on additional assays and can be stored by the instrument. Each additional assay includes calibrator (curve) controls that have to recover in defined ranges to ensure that the stored calibration curve is still valid. The Fluorescence-Immunoassay test system includes test-, method specific and general reagents that are packaged as separate units.

### **Test Principle of the New Device**

The EliA Celikey IgA Wells are coated with human recombinant tTG. If present in the patient's specimen, antibodies to tTG bind to their specific antigen. After washing away non-bound antibodies, enzyme-labeled antibodies against human IgA antibodies (EliA IgA Conjugate) are added to form an antibody-conjugate complex. After incubation, non-bound conjugate is washed away and the bound complex is incubated with a Development Solution. After stopping the reaction, the

**EliA™ Celikey IgA – New Device**  
**510(k) Submission**  
**Section 9. Summary of Safety and Effectiveness**

---

fluorescence in the reaction mixture is measured. The higher the response value, the more specific IgA is present in the specimen. To evaluate test results, the response for patient samples is compared directly to the response for calibrators.

**Device Comparison**

The new and the predicate device both represent non-competitive solid phase EIAs. Both IVDs are used as an aid in the diagnosis of Coeliac Disease.

**Laboratory equivalence**

The comparability of predicate device and new device is supported by a data set including

- results obtained within a comparison study between new and predicate device
- results obtained for clinically defined sera
- results obtained for samples from apparently healthy subjects (normal population).

In summary, all available data support that the new device is substantially equivalent to the predicate device.



Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

Mr. Martin R. Mann  
Phadia US Inc.  
4169 Commercial Ave.  
Portage, MI 49002

DEC 29 2006

Re: k062787  
Trade/Device Name: EliA™ Celikey IgA Immunoassay, EliA™ Celiac Control  
Regulation Number: 21 CFR 866.5660  
Regulation Name: Tissue transglutaminase autoantibodies immunological test system,  
Antigen and Control  
Regulatory Class: Class II  
Product Code: MVM  
Dated: September 15, 2006  
Received: September 18, 2006

Dear Mr. Mann:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

Page 2 –

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (240) 276-0490. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,



Robert L. Becker, Jr., MD, Ph.D  
Director

Division of Immunology and Hematology  
Office of In Vitro Diagnostic Device Evaluation & Safety  
Center for Devices and Radiological Health

Enclosure

EliA™ Celikey IgA – New Device  
510(k) Submission  
Section 1. Indications for Use Statement

---

**Indications for Use**

510(k) Number:           K062 787          

Device Name: **EliA™ Celiac Control**

**Indications For Use:**

EliA Celiac Control is intended for laboratory use in monitoring the performance of in vitro measurement of antibodies to tissue transglutaminase (tTG) with ImmunoCAP 100 or 250 using the EliA IgG or IgA method.

  
Division Sign-Off

Office of In Vitro Diagnostic  
Device Evaluation and Safety

510(k)           K062 787          

Prescription Use              
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use              
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

**EliA™ Celikey IgA – New Device**  
**510(k) Submission**  
**Section 1. Indications for Use Statement**

---

**Indications for Use**

510(k) Number:           K062787          

Device Name: **EliA™ Celikey IgA Immunoassay**

**Indications For Use:**

EliA Celikey IgA Immunoassay is intended for the in vitro semi-quantitative measurement of IgA antibodies directed to tissue transglutaminase (tTG) in human serum and plasma (EDTA, citrate). EliA Celikey IgA Immunoassay is based on recombinant human tissue transglutaminase as antigen and is useful as an aid in the clinical diagnosis of patients with celiac disease, in conjunction with other laboratory and clinical findings. EliA Celikey IgA Immunoassay uses the EliA IgA method on the instrument ImmunoCAP 100 and ImmunoCAP 250.

  
Division Sign-Off

Office of In Vitro Diagnostic  
Device Evaluation and Safety

510(k)           K062787          

Prescription Use              
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use              
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)