



K070223

## 510(k) Summary

### A. Submitter

MAR 13 2007

Aalto Scientific, Ltd.  
1959 Kellogg Ave.  
Carlsbad, CA 92008  
Telephone: (760) 431-7922  
Fax: (760) 431-6824

### B. Contact Person

Alan Vekich  
Regulatory Affairs Specialist  
Telephone: (760) 431-7922 ext. ~~21~~ 120  
E-mail: [avekich@aaltoscientific.com](mailto:avekich@aaltoscientific.com)

### C. Date of Summary Preparation

December 1, 2006

### D. Device Identification

Product Trade Name:	Audit™ MicroCV™ Lipids Linearity Set (Low)
Common Name:	Lipids Linearity
Classification Name:	Assay QC Material
Device Classification:	Class I
Regulation Number:	21 CFR 862.1660
Panel:	75
Product Code:	JJY

### E. Device to Which Substantial Equivalence is Claimed

Audit™ MicroCV™ General Chemistry Linearity Set  
Aalto Scientific, Ltd., Carlsbad, CA  
K042318



## F. Description of the Device

The Audit™ MicroCV™ Lipids Linearity Set (Low) is a human based, lyophilized, five level set of QC material, with each level containing 6 analytes. It is used to confirm the proper calibration, linear operating range, and reportable range of Lipids methods for the analytes listed. Level A is near the lower limit level and Level E has concentrations near the upper limit of instruments. Levels B – D are related by linear dilution of Level A and Level E.

## G. Statement of Intended Use

The Lipids Linearity Set is intended to simulate human patient serum samples for the purpose of verifying and validating the Analytical Measurement Range for non-waived Lipids testing methods as identified in the package insert.

## I. Summary of Performance Data

Stability studies have been performed to determine the reconstituted stability and shelf life for the Audit™ MicroCV™ Lipids Linearity Set (Low). All supporting data is retained on file at Aalto Scientific, Ltd. Product claims are as follows:

*Reconstituted Stability:* Once a vial has been reconstituted, all analytes will be stable for 24 hours when stored tightly capped at 2 - 8° C.

*Shelf Life:* Two years, when stored unopened at 2 - 8° C.

*Note:* Real time studies are ongoing to support the shelf life of this product.

## H. Technical Characteristics Compared to Predicate Device

Characteristics	Audit™ MicroCV™ Lipids Linearity Set (Low) (New Device)	Audit™ MicroCV™ General Chemistry Linearity Set (K042318)
Intended Use	Audit™ MicroCV™ Lipids Linearity Set (Low) is assayed quality control material consisting of Delipidized Human Serum. It is intended to simulate human patient serum samples for the purpose of monitoring the precision and to detect systematic analytical deviations of laboratory testing procedures. This product may also be used as unassayed quality control material for these same analytes and may be used for proficiency testing in inter-laboratory surveys. In addition, this product may also be used to perform CLIA directed calibration verification for these same analytes with similar reagents on similar instrumentation in accordance with current CLIA-88 guidelines and regulations.	Audit™ MicroCV™ General Chemistry Linearity Set is assayed quality control material consisting of human based serum. It is intended to simulate human patient serum samples for the purpose of monitoring the precision and to detect systematic analytical deviations of laboratory testing procedures. This product may also be used as unassayed quality control material for these same analytes and may be used for proficiency testing in inter-laboratory surveys. In addition, this product may also be used to perform CLIA directed calibration verification for these same analytes with similar reagents on similar instrumentation in accordance with current CLIA-88 guidelines and regulations.
Number of Analytes per vial	6	30
Number of levels per set	5	5
Contents	5 x 2 mls	5 x 5 mls
Matrix	Delipidized Human Serum	Human Based Serum
Type of Analytes	Lipids	General Chemistry
Form	Lyophilized	Lyophilized
Stabilizers	Sucrose	None
Preservatives	Sorbitol	Sorbitol
Storage	2 to 8° C Until expiration date	2 to 8° C Until expiration date
Reconstituted Stability	24 hours at 2 to 8° C	24 hours at 2 to 8° C



## J. Conclusions

Based upon the purpose of the device, the descriptions and labeling of the predicate device, the safety and efficacy, and the stability data generated, the product is substantially equivalent to the predicate device.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

Aalto Scientific, Ltd.  
c/o Mr. Alan Vekich  
Regulatory Affairs Specialist  
1959 Kellogg Avenue  
Carlsbad, CA 92008

**MAR 13 2007**

Re: k070223  
Trade/Device Name: Audit™ Microcv™ Lipids Linearity Set (Low)  
Regulation Number: 21 CFR§862.1660  
Regulation Name: Quality control material (assayed and unassayed).  
Regulatory Class: Class I  
Product Code: JJY  
Dated: January 17, 2007  
Received: January 24, 2007

Dear Mr. Vekich:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

Page 2 –

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (240) 276-0490. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address at <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

*Jean M. Cooper, M.S., D.V.M.*

Jean M. Cooper, M.S., D.V.M.

Director

Division of Chemistry and Toxicology

Office of *In Vitro* Diagnostic Device

Evaluation and Safety

Center for Devices and

Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): unknown

Device Name: Audit™ MicroCV™ Lipids Linearity Set (Low)

Indications For Use:

The Audit™ MicroCV™ Lipids Linearity Set (Low) consists of five levels. Each level contains the following analytes: Apolipoprotein A1, Apolipoprotein B, Cholesterol, HDL Cholesterol, LDL Cholesterol and Triglyceride and may be used for proficiency testing in interlaboratory surveys and to perform CLIA directed calibration verification for these same analytes with similar reagents on similar instrumentation in accordance with current CLIA-88 guidelines and regulations.

In addition, Level A – E of this product may be used as unassayed quality control material for these analytes or as an assayed quality control material for the analyzer systems specified in the package insert. It is not intended to be used as an assayed quality control material for any other analyzer systems.

Prescription Use   X    
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)

Carol Benson  
Division Sign-Off

Office of In Vitro Diagnostic Device  
Evaluation and Safety

K070223