

## Section 5: 510(k) Summary

K071047

JUN - 8 2007

**Submitted by:** Masimo Corporation  
40 Parker  
Irvine, CA 92618  
PHONE (949) 297-7000  
FAX (949) 297-7001

**Official Correspondent:** James Cronin, Vice President, Regulatory Affairs

**Company Contact for this Submission:** Marguerite Thomlinson, Manager, Regulatory Affairs

**Date Summary Prepared:** April 12, 2007

**Trade Name:** Masimo Patient Safety Net

**Common Name:** System, Network and Communication, Physiological Monitors

**Device Class:** Class II

**Product Code:** MSX

**Classification Name:** System, Network and Communication, Physiological Monitors

**Substantially Equivalent Devices:** Bernoulli Management System  
510(k) Number - K061932

### Description of Masimo Patient Safety Net

The Masimo Patient Safety Net (PSN) is a supplemental alarm system for physiological monitoring devices. It is not intended to replace any part of the patient monitoring procedures already existing for the patient monitoring devices.

The PSN communicates with multiple patient monitoring devices and distributes physiological monitoring information remotely. Skilled clinicians remotely receive wireless transmission of patient physiological monitoring information, from patient monitoring device(s) to their pagers. The transmitted information includes alarm information and physiological parameters.

### Intended use/ Indications for Use

The Masimo Patient Safety Net (PSN) is intended to be used as a supplemental alarm system, communicating with multiple patient monitoring devices. The PSN provides secondary display of physiological monitoring parameters. It enables the viewing and monitoring of patient physiological conditions. The PSN is used in hospitals or hospital-type environments.

## Section 5: 510(k) Summary

K071047

### Technology

The Masimo Patient Safety Net (PSN) includes standard telecommunication and IT hardware. Masimo develops the application software which is used on the system server and PC computer (central monitoring station).

The PSN functions and transmits data similarly to the predicate device. The PSN and the predicate device use similar data transmission and communication technologies.

### Test Summary

The Masimo Patient Safety Net (PSN) complies with the voluntary standards as detailed in Section 9 of this submission. The following quality assurance measures were applied to the development of the PSN:

- Risk Analysis
- Design Reviews
- Component Level Testing
- System Level Testing
- Performance Testing
- Safety Testing
- Environmental Testing

### Conclusions

The information in this 510(k) submission demonstrates that the Masimo Patient Safety Net is substantially equivalent to the predicate device as a supplementary alarm system, with respect to safety, effectiveness, and performance.



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

JUN - 8 2007

Masimo Corporation  
c/o Marguerite Thomlinson  
Manager, Regulatory Affairs  
40 Parker  
Irvine, California 92618

Re: K071047

Masimo Patient Safety Net  
Regulation Number: 21 CFR 870.2300  
Regulation Name: Cardiac monitor (including cardiometer and rate alarm)  
Regulatory Class: Class II  
Product Code: MSX  
Dated: April 12, 2007  
Received: April 13, 2007

Dear Ms. Thomlinson:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Page 2 – Ms. Marguerite Thomlinson

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Bram D. Zuckerman, M.D.  
Director  
Division of Cardiovascular Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Section 4 - Indications for Use

510(k) Number (if known): K071047

Device Name: Masimo Patient Safety Net

### Indications For Use:

The Masimo Patient Safety Net (PSN) is intended to be used as a supplemental alarm system, communicating with multiple patient monitoring devices. The PSN provides secondary display of physiological monitoring parameters. It enables the viewing and monitoring of patient physiological conditions. The PSN is used in hospitals or hospital-type environments.

Prescription Use   X  

AND/OR

Over-The-Counter Use

(Per 21 CFR 801 Subpart D)

(Per 21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF  
NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

*B. Munn*  
(Division Sign-Off)  
Division of Cardiovascular Devices  
510(k) Number K071047

0011